



Remodeling or new setup?

Whether you're completely remodeling your office or just want to freshen up an exam room, Midmark has what you need to make your office run more efficiently... such as our trusted Ritter products.

Midmark carries a full line of quality healthcare equipment including power and manual examination tables, sterilizers, casework, lighting, seating, digital diagnostic products and more.

Simplicity, clean lines and ease of use, which are hallmarks of all our products, have earned us both national and international recognition, as well as our position as the industry leader.

Contact your sales Representative for details on how we can help you provide efficient patient care in your practice, or visit midmark.com.



We Bring More to the Table.TM

Ritter[®]
by MIDMARK

Alko Enterprises

Intelligent weighing and measuring for your practice.



seca 703

Functional, user-friendly, precise: seca scales and measuring devices are reliable aids to help you survive the strain of everyday medical life. seca products are not only extremely robust and durable but are tailor-made to meet the requirements of the user. Whether baby, column or floor scales or measuring devices for children or adults – seca products all have useful details such as transport castors, cable-linked remote displays, TARE, HOLD and BMI functions. So you can calmly look forward to your day's work at the practice.



seca 374



seca 416



BREWER ACCESS HIGH-LOW EXAM TABLE

Chair-to-table design effortlessly takes patient from consultation to examination. 18" low height for easy and safe patient transfers. Return to chair feature easily takes table back to chair position in one, simple step (on power back model only). Standard grab bars on all models provide added safety and comfort. Provides patient a stable and secure grab point to use when getting on and off table. Conveniently swings out for lateral wheelchair transfers. Bars are easily removable. Abundant storage capacity. Patented pass-through side drawers allow left or right access. Additional spacious front drawer. Convenient covered paper roll storage holds up to six rolls and keeps paper clean. Pass-through work surface pulls out on right or left side. Saves space by eliminating the need for additional stands or trays. Perfect for charts, instruments, procedure preparation, or portable computer. Intuitive foot control, on/off switch (located on back of table) allows table's foot control to be disabled and prevents patients and/or children from moving table. Comfortable, ergonomic upholstery with easy to clean surface; seamless design. Extendable, adjustable ergonomic stirrups. Available in Clamshell, Gunmetal, Adriatic, Burgundy, Emerald, Mauve, Azure Blue, Royal Plum, Plum, Deep Tea, Petal, Blue Fog, Forest Green, Spanish Tile, Chocolate, Black, Putty, Graphite, Arctic Plush, Charcoal Plush, Shetland Plush, Mocha Plush, Berry Plush, Admiral Plush, CAL 133, and Special Color.

6000	6000 Model w/pneumatic back
6001	6001 Model w/pneumatic back, pelvic tilt & drawer warmer
6115-5X	Upholstery Set (Top & Toe Pad), Specify Color
6115-SP	Upholstery Set (Top & Toe Pad), Special Colors
6115-XX	Upholstery Set (Top & Toe Pad), Standard Colors, Specify Color
6500	6500 Model w/power back & return to chair
6501	6501 Model w/power back, Return to Chair, Pelvic Tilt & Drawer Warmer



7000-00



7000-04


BREWER ASSIST POWER PROCEDURE TABLE

Exclusive, easy-to-adjust three-point pivot headrest system. Convenient table paper storage accepts 18" and 21" rolls. Comfortable, ergonomic upholstery that is easy to clean. Easy wheelchair patient transfer at 19". Supports 450 lbs. Exclusive footrest extends an additional 5". Intuitive foot and handheld controls. Height: 19"-40". Length: 69" with footrest and headrest extends to 90". Width: 28". Leg rest: 16.5"W x 11.25"H. Headrest: 13"W x 11"H. 115 VAC nominal, 60Hz. Safety Extra-Low Voltage (SELV). UL listed and CUL listed. Three year warranty. Available in Clamshell, Gunmetal, Adriatic, Burgundy, Emerald, Mauve, Azure Blue, Royal Plum, Plum, Deep Teal, Custom Color, Petal, Blue Fog, Forest Green, Spanish Tile, Chocolate, Black, Putty, Graphite, and CAL 133.

7000-00	Power Table w/o upholstery
7000	Power Table w/upholstery

Accessories


7115	Upholstery Only
99500	Hand Pendant
99501	Side Rails Including Hardware (Set of 2)
99502	Universal Socket Clamp (needed for all accessories)
99503	Arm Board
99504	Deluxe Articulating Knee Crutches (right & left)
99505	Stainless Steel Drain Pan




6100




6115




6160



6170



6173



6190



6192



6196

CLINTON BARIATRICS SERIES TABLES

2" high-density urethane foam padding. Easy to clean laminate surfaces. All tops are plywood. All tops are attached with bolt through the top tee nut fastening system. All stress points reinforced with steel angles. All cabinets are assembled using steel fasteners and metal supports. 600 lbs capacity. Heavy duty construction throughout.

6100	Bariatrics Power Table w/power backrest, Hands-Free Foot Control, 72"L x 20"-30"H x 33"W
6115	Bariatrics Echo-Image Cabinet Style Table w/pillow, hinged drop leaf & storage compartments, 78"L x 28"H x 33"W
6160	Bariatrics Full Cabinet Treatment Table w/3 storage areas w/adjustable shelf, 78"L x 31"H x 33"W
6170	Bariatrics Echo-Image Power Table w/pillow & hands-free foot pedal control, 76"L x 20"-30"H x 33"W
6173	Bariatrics Echo-Image Power Table w/pillow, built-in slip-down casters & hands-free foot pedal control, 76"L x 20"-30"H x 33"W
6190	Bariatrics Treatment Table w/laminate full length shelf w/H-braces & pillow, 6-Leg Design, Built-In Middle Leg Levelers, 78"L x 31"H x 30"W
6192	Bariatrics Flat Top Power Table w/hands-free foot control, 76"L x 20"-30"H x 33"W
6196	Low Height Bariatrics Treatment Table w/laminate full length shelf w/H-braces & pillow, 6-Leg Design, Built-In Middle Leg Levelers, 78"L x 26"H x 30"W



5190



5192

CLINTON POWER TABLES

2" thick, firm, foam padding. All welded steel scissors lifting mechanism. Easy-clean, attractive plastic shroud hides and protects lifting mechanism. Sturdy, welded steel base with easily adjusted foot levelers. Overload protection. Hospital-grade plug. 400 lbs lifting capacity.

5190	Flat Top Power Table w/foot operated controls, 72"L x 22 1/2"-35"H x 27"W
5192	Power Table w/adjustable backrest & hands-free foot control, 72"L x 22 1/2"-35"H x 27"W

Tables

HAUSMANN ECONO-LINE™ ELECTRIC HI-LO 5-SECTION TREATMENT TABLE

Motorized height adjustment from 37" high to 18" low for wheelchair transfer. Five sections include foot section, head section, seat section and two armrests. Gas spring rear section raises to 60° and headrest raises up to 30° and down to 22° below horizontal. Armrests adjust up and down for localized treatment. 2" thick urethane foam top. Breathing slot and plug in head section. Fully retractable casters. Heavy duty steel frame. Powder coated epoxy in neutral cream finish. Includes pneumatic foot control. Limit switches. Hospital grade plug.

6057 74"L x 27"W x 18"-37"H Table

**HAUSMANN ECONO-LINE™ POWER CLINIC EXAMINATION TABLE**

Dual direct-drive electrical lift actuators adjust height from 24" to 35". Convenient hand control operation. Two door and two drawer storage cabinet. Includes 36" long (10) position adjustable backrest. Four position adjustable 18" long legrest. Vinyl-clad stirrups adjust in length. Includes paper dispenser and cutter.

4430 Examination Table, 56"-74"L x 27"W x 24"-35"H
 4435 Examination Table, 56"-74"L x 27"W x 24"-35"H Left Hand Table
 4731 Examination Table w/o stirrups, 56"-74"L x 27"W x 24"-35"H

**HAUSMANN ELECTRIC HI-LO TREATMENT TABLE**

Adjustable height hydraulic cylinders are built into each four corner steel legs. Hands free foot control activation. Apron and shelf are almond laminate. Black powder coated legs with levelers. Upholstery comes in 14 standard vinyl colors.

4705 Crank Hi-Lo Treatment Table, Same as #4705 Except Hand Crank Operation w/fold-away handle
 4707 Electric Hi-Lo Treatment Table, 72"L x 30"W x 27"-37"H



6040-709



6056-709

HAUSMANN ELECTRIC TILT TABLES

Gray vinyl upholstery. 2" thick urethane foam top. Tilt angle indicator from 0-90°. Three extra-wide safety straps. Includes foot board. Four 3" locking casters. Heavy duty steel frame is powder coated in neutral cream finish. 300 lb capacity.

6040-709 78"L x 28"W x 33"H Table, Heavy Duty, Pneumatic Hand Control, Removable Footboard
 6056-709 78"L x 28"W x 25-39"H Table, Two Heavy Duty Motors, Hi-Lo, Pneumatic Hand Control, Two Support Handles, Stainless Steel Footboard (removable w/tools)

HAUSMANN ELECTRIC WORK TABLE

Direct-drive electric lift mechanism elevates top from 27" to 43"H. Ideal height range for sitting, standing or wheelchair use. Sleek, compact pedestal design offers more leg room for patients. 36" x 60" Natural Oak high pressure laminate top with black protective bullnose molding on perimeter edge. Black laminate base plate with leveling legs. Easily accommodates four ambulatory patients or two wheelchairs and therapists. Hand control panel is conveniently mounted under work surface. 100V, 60Hz, 1.8 amps. Hospital grade plug.

4375 Electric Work Table

**HAUSMANN MULTI-PURPOSE ECONO-LINE ELECTRIC TABLE**

Electric motorized height adjustment from 20" to 28". Air spring pneumatic headrest raises to 80° and to 20° below horizontal. Hands-free foot control operation. One adjustable foot for floor leveling. Includes paper dispenser and cutter. Heavy duty steel frame with cream powder coated epoxy finish.

4710 74"L x 25"W x 20"-38"H Table



4795



595

HAUSMANN POWERMATIC® BARIATRIC HI-LO TREATMENT TABLE

Designed specifically for the treatment of bariatric patients. Three direct drive motors - 2 synchronized for height adjustment and 1 for the backrest to facilitate patient positioning. 600 lb lifting capacity. Spacious size and 2" thick urethane foam padding for patient comfort. Elevates from 20" to 29" for wheelchair accessibility. Remote (wired) control for easy operation. Ships fully assembled. Black laminate base with satin aluminum trim and levelers.

4795 Treatment Table, 78"L x 32"W, Height Adjusts From 27"-29"

Accessories

595 Digital Scale



HAUSMANN POWERMATIC® EXAMINATION TABLES

Hi-lo, smooth, powerful motorized lift system. 400 lb capacity. Hands-free pneumatic foot control. Vinyl upholstered pad with built-in pillow. Pneumatic backrest with dual adjustment from 0 to 75°. Sliding 21"W x 16"L leg rest extension with removable pad. Adjustable concealed stirrups made of heavy duty steel with vinyl coated heel rests. Laminate drawers on nylon rollers. Black laminate base plate with leveling legs. Paper dispenser and paper cutter.

- 4440 55" to 71" L x 26" D Table, Motorized Lift From 26"-36", Two Front Drawers, Two Side Drawers
 4445 55" to 71" L x 26" D, Left Hand Table, Motorized Lift From 26"-36", Two Front Drawers, Two Side Drawers

Accessories

- 089 Motorized Power Backrest
 091 Pelvic Tilt
 094 Drawer Warmer
 095 Armboard
 097 Irrigation Pan
 098 "Pop-Up" Elevating Leg Rest
 17 Side Mounting Rails For Surgical Boot Stirrups For Models 4450 & 4455



4750



4766

HAUSMANN POWERMATIC® PROCEDURE TABLES

2" thick urethane foam top. 400 lb capacity. Includes paper dispenser/cutter. Foot control for hands-free operation. Thermal overload protection. Limit switches. Hospital grade plug. Formica base plate with leveling legs. 76"L x 27"W.

- 4750 Procedure Table Height Adjusts From 19"-31", Wheelchair Accessible
 4756 Procedure Table Height Adjusts From 19"-31", Wheelchair Accessible, Pneumatic 30" Backrest Elevates to 75° w/side controls
 4757 Procedure Table Height Adjusts From 19"-31", Wheelchair Accessible, Power 30" Backrest Elevates to 75° w/foot control
 4766 Procedure Table Height Adjusts From 26"-38", Pneumatic 30" Backrest Elevates to 75° w/side lever controls
 4767 Procedure Table Height Adjusts From 26"-38", Power 30" Backrest Elevates to 75° w/foot control



4746



4747

HAUSMANN POWERMATIC® TREATMENT TABLE

All Formica power table with ample storage consisting of (2) drawers and (1) door cabinet with adjustable shelf. 76" x 27" with power adjustable height range of 26"-38". Available in your choice of 114 Formica colors and 12 vinyl upholstery colors at no additional charge. Paper dispenser and paper cutter included. Available with 99L Left Hand Table option.

- 4746 Table w/airlift backrest
 4747 Table w/electric power backrest

HAUSMANN POWERMATIC® WHEELCHAIR ACCESSIBLE ADA EXAM TABLE

Custom designed for patients with special needs, including bariatric, geriatric and ambulatory patients. Hi-Lo motorized lift from 20" to 32". Hands-free foot control. Vinyl upholstered removable pad, with contoured padded head section is attached with Velcro® brand fasteners. Air-Spring pneumatic backrest with dual controls adjusts from 0 to 75°. Concealed stirrups adjust in length and lateral positions. Vinyl coated heelrests for patient comfort. Includes slide-out legrest with pad. Includes paper dispenser and cutter. Hospital grade duplex outlet. 100V, 60Hz, 3.5 amps. UL Listed. 927 Folkstone Gray or your choice of any high pressure laminate solid color or woodgrain. Upholstery comes in 14 standard vinyl colors.

- 4460 ADA Exam Table, 55"-71" L x 26" D x 20"-32" H
 4465 ADA Exam Table, Left Hand Table, 55"-71" L x 26" D x 20"-32" H



MIDMARK 622 PREMIUM UPHOLSTERY

The Midmark 622 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 622 powers-up and down for unlimited patient positioning. Designed especially for the exam room, the 622 features new technology and a new look. It complies with ADA guidelines.

- 622-001-235 Premium Upholstery, Blueberry
 622-001-236 Premium Upholstery, Silver Sage
 622-001-237 Premium Upholstery, Tea Green
 622-001-238 Premium Upholstery, Fossil Gray
 622-001-239 Premium Upholstery, Cashmere Blue
 622-001-240 Premium Upholstery, Terra Cotta
 622-001-247 Premium Upholstery, Belagio
 622-001-248 Premium Upholstery, Siena
 622-001-249 Premium Upholstery, Tuscany
 622-001-250 Premium Upholstery, Firenze
 622-001-251 Premium Upholstery, Capri
 622-001-252 Premium Upholstery, Treviso



MIDMARK 622 ULTRA PREMIUM UPHOLSTERY

The Midmark 622 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 622 powers-up and down for unlimited patient positioning. Designed especially for the exam room, the 622 features new technology and a new look. It complies with ADA guidelines.

- 622-001-241 Ultra Premium Upholstery, Pacific Blue
 622-001-242 Ultra Premium Upholstery, Spruce
 622-001-243 Ultra Premium Upholstery, Olivine
 622-001-244 Ultra Premium Upholstery, Artic
 622-001-245 Ultra Premium Upholstery, Sky
 622-001-246 Ultra Premium Upholstery, Canyon

MIDMARK 622 WITH HEATER AND PELVIC TILT - PREMIUM UPHOLSTERY

The Midmark 622 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 622 powers-up and down for unlimited patient positioning. Designed especially for the exam room, the 622 features new technology and a new look. It complies with ADA guidelines.

622-002-235	Premium Upholstery, Blueberry
622-002-236	Premium Upholstery, Silver Sage
622-002-237	Premium Upholstery, Tea Green
622-002-238	Premium Upholstery, Fossil Gray
622-002-239	Premium Upholstery, Cashmere Blue
622-002-240	Premium Upholstery, Terra Cotta
622-002-247	Premium Upholstery, Belagio
622-002-248	Premium Upholstery, Siena
622-002-249	Premium Upholstery, Tuscany
622-002-250	Premium Upholstery, Firenze
622-002-251	Premium Upholstery, Capri
622-002-252	Premium Upholstery, Treviso

MIDMARK 622 WITH HEATER AND PELVIC TILT - ULTRA PREMIUM UPHOLSTERY

The Midmark 622 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 622 powers-up and down for unlimited patient positioning. Designed especially for the exam room, the 622 features new technology and a new look. It complies with ADA guidelines.

622-002-241	Ultra Premium Upholstery, Pacific Blue
622-002-242	Ultra Premium Upholstery, Spruce
622-002-243	Ultra Premium Upholstery, Olivine
622-002-244	Ultra Premium Upholstery, Artic
622-002-245	Ultra Premium Upholstery, Sky
622-002-246	Ultra Premium Upholstery, Canyon
623-002-241	Ultra Premium Upholstery, Pacific Blue
623-002-242	Ultra Premium Upholstery, Spruce
623-002-243	Ultra Premium Upholstery, Olivine
623-002-244	Ultra Premium Upholstery, Artic
623-002-245	Ultra Premium Upholstery, Sky
623-002-246	Ultra Premium Upholstery, Canyon



MIDMARK 623 PREMIUM UPHOLSTERY

The Midmark 623 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 623 powers-up and down for unlimited patient positioning. It features a power back section and an auto return control, which brings the table to the home position with a touch of a button. Designed especially for the exam room, the 623 features new technology and a new look. It complies with ADA guidelines.

623-001-235	Premium Upholstery, Blueberry
623-001-236	Premium Upholstery, Silver Sage
623-001-237	Premium Upholstery, Tea Green
623-001-238	Premium Upholstery, Fossil Gray
623-001-239	Premium Upholstery, Cashmere Blue
623-001-240	Premium Upholstery, Terra Cotta
623-001-247	Premium Upholstery, Belagio
623-001-248	Premium Upholstery, Siena
623-001-249	Premium Upholstery, Tuscany
623-001-250	Premium Upholstery, Firenze
623-001-251	Premium Upholstery, Capri
623-001-252	Premium Upholstery, Treviso



MIDMARK 623 ULTRA PREMIUM UPHOLSTERY

The Midmark 623 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 623 powers-up and down for unlimited patient positioning. It features a power back section and an auto return control, which brings the table to the home position with a touch of a button. Designed especially for the exam room, the 623 features new technology and a new look. It complies with ADA guidelines.

623-001-241	Ultra Premium Upholstery, Pacific Blue
623-001-242	Ultra Premium Upholstery, Spruce
623-001-243	Ultra Premium Upholstery, Olivine
623-001-244	Ultra Premium Upholstery, Artic
623-001-245	Ultra Premium Upholstery, Sky
623-001-246	Ultra Premium Upholstery, Canyon

MIDMARK 623 WITH HEATER AND PELVIC TILT - PREMIUM UPHOLSTERY

The Midmark 623 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 623 powers-up and down for unlimited patient positioning. It features a power back section and an auto return control, which brings the table to the home position with a touch of a button. Designed especially for the exam room, the 623 features new technology and a new look. It complies with ADA guidelines.

623-002-235	Premium Upholstery, Blueberry
623-002-236	Premium Upholstery, Silver Sage
623-002-237	Premium Upholstery, Tea Green
623-002-238	Premium Upholstery, Fossil Gray
623-002-239	Premium Upholstery, Cashmere Blue
623-002-240	Premium Upholstery, Terra Cotta
623-002-247	Premium Upholstery, Belagio
623-002-248	Premium Upholstery, Siena
623-002-249	Premium Upholstery, Tuscany
623-002-250	Premium Upholstery, Firenze
623-002-251	Premium Upholstery, Capri
623-002-252	Premium Upholstery, Treviso

MIDMARK 623 WITH HEATER AND PELVIC TILT - ULTRA PREMIUM UPHOLSTERY

The Midmark 623 High-Low Examination Table with a low height of 18 inches and an extraordinary comfort seating system is comfortable, easy to access, and reduces patient anxiety. The Midmark 623 powers-up and down for unlimited patient positioning. It features a power back section and an auto return control, which brings the table to the home position with a touch of a button. Designed especially for the exam room, the 623 features new technology and a new look. It complies with ADA guidelines.

623-002-241	Ultra Premium Upholstery, Pacific Blue
623-002-242	Ultra Premium Upholstery, Spruce
623-002-243	Ultra Premium Upholstery, Olivine
623-002-244	Ultra Premium Upholstery, Artic
623-002-245	Ultra Premium Upholstery, Sky
623-002-246	Ultra Premium Upholstery, Canyon

RITTER 222 AND 223 BARRIER-FREE™ EXAM TABLE

Getting on an exam table can be an anxious, insecure moment, especially for elderly and special needs patients. The Ritter 222 and 223 Barrier-Free tables lower to 18 inches from the floor, making it easier for patients to get on the table and eliminating excess strain on their backs, as well as the backs of physicians and staff. Patients can now get on the table with no or minimal staff assistance, helping them maintain their dignity. Measuring 37 inches in highest position, these exam tables have the absolute lowest and highest height combination on the market, making exams more comfortable. This maximum low and high height combination also eliminates bending and stooping to help your patients onto the table surface. Once your patient is on the table, with a touch of the foot control, the Ritter Barrier-Free tables allow precise and smooth positioning with minimal effort on everyone's part, helping your patients and you feel relaxed. With unprecedented low and high positioning heights, the Ritter 222 has a gas spring assisted back section with an easy access backrest mechanism, while the Ritter 223 has a power back system. Both tables include: Built-in pillow; Hideaway stirrups that lock in one of four lateral positions; Foot extension shelf to expand the table length; Paper roll holder hidden under headrest; Storage for four 21" x 3.5" diameter paper rolls; Polystyrene treatment pan; Foot control; Exam Assistant Drawer System™ and Leveling screws. Available in Sand Gray, Slate Gray, Iris Blue, Nile Green, Terra Rosa, Royal Lavender, Pebble Gray, Perfect Plum, Moss, Navy, Shadow, Dusty Blue, and Clay.

222-015	Receptacle, Soft Touch Top, Base Only
222-015-00	Receptacle, Soft Touch Top
222-015-01	Receptacle, Seamless Top
222-016	Receptacle, Drawer, Heater, Pelvic Tilt, Soft Touch Top, Base Only
222-016-00	Receptacle, Drawer, Heater, Pelvic Tilt, Soft Touch Top
222-016-01	Receptacle, Drawer, Heater, Pelvic Tilt, Seamless Top
223-015	Receptacle, Soft Touch Top, Base Only
223-015-00	Receptacle, Soft Touch Top
223-015-01	Receptacle, Seamless Top
223-016	Receptacle, Drawer, Heater, Pelvic Tilt, Soft Touch Top, Base Only
223-016-00	Receptacle, Drawer, Heater, Pelvic Tilt, Soft Touch Top
223-016-01	Receptacle, Drawer, Heater, Pelvic Tilt, Seamless Top



MIDMARK 630 UNIVERSAL POWER PROCEDURE TABLES

With a low height of 19" and an exclusive Premium Comfort System™, it provides easy access for patients who are pregnant, elderly, or have disabilities. Full eight-way power, an incredible height range of 19" to 40", 450 lb weight capacity, a hand and foot control, unprecedented patient comfort, and multitude of available accessories make this procedure table the deal choice for improving the look, feel and efficiency of the practice. #630-002 offers the added feature of programmable positions over 630-001. #630-003 offers the added feature of auxiliary receptacles over 630-001. #630-004 offers the added feature of four programmable positions and auxiliary receptacles over 630-001.

630-001	Midmark 630 Table, 4-Function, Non-Programmable w/out receptacle & w/out upholstered top
630-002	Midmark 630 Table, 4-Function, Programmable w/out receptacle & w/out upholstered top
630-003	Midmark 630 Table, 4-Function, Non-Programmable w/receptacle & w/out upholstered top
630-004	Midmark 630 Table, 4-Function, Programmable w/receptacle & w/out upholstered top



MIDMARK 75L POWER PROCEDURE TABLES

With a low height of 19" and an exclusive Premium Comfort System™, it provides easy access for patients who are pregnant, elderly, or have disabilities. Full eight-way power, an incredible height range of 19" to 40", 450 lb weight capacity, a hand control, auxiliary receptacles, unprecedented patient comfort, and multitude of available accessories make this procedure table the ideal choice for improving the look, feel and efficiency of the practice.

75-020	Midmark 75L Table, 4-Function, Non-Programmable, Hand Control Only w/receptacle & w/out upholstered top
--------	---



MIDMARK 75LP POWER PROCEDURE TABLE

With a low height of 19" and an exclusive Premium Comfort System™, it provides easy access for patients who are pregnant, elderly, or have disabilities. Full eight-way power, an incredible height range of 19" to 40", 450 lb weight capacity, a hand control, 4 programmable positions, auxiliary receptacles, unprecedented patient comfort, and multitude of available accessories make this procedure table the ideal choice for improving the look, feel and efficiency of the practice. The Midmark 75LP (75-021) offers the added feature of four programmable positions over the Midmark 75L (75-020).

75-021	Midmark 75LP Table, 4 Function, Programmable, Hand Control Only w/receptacles & w/out upholstered top
--------	---



RITTER 244 BARIATRIC TREATMENT TABLE WITHOUT TOP

244-001	Ritter 244 Bariatric Treatment Table w/out top
244-001-201	Ritter 244 Bariatric Treatment Table, Slate Grey
244-001-202	Ritter 244 Bariatric Treatment Table, Iris Blue
244-001-203	Ritter 244 Bariatric Treatment Table, Nile Green
244-001-204	Ritter 244 Bariatric Treatment Table, Terra Rosa
244-001-205	Ritter 244 Bariatric Treatment Table, Royal Lavender
244-001-216	Ritter 244 Bariatric Treatment Table, Pebble Grey
244-001-229	Ritter 244 Bariatric Treatment Table, Perfect Plum
244-001-230	Ritter 244 Bariatric Treatment Table, Moss
244-001-231	Ritter 244 Bariatric Treatment Table, Navy
244-001-232	Ritter 244 Bariatric Treatment Table, Shadow
244-001-233	Ritter 244 Bariatric Treatment Table, Dusty Blue
244-001-234	Ritter 244 Bariatric Treatment Table, Clay

Tables

RITTER 319 POWER PROCEDURE TABLE

A truly universal power table with eight-way power, allows you to precisely and effortlessly position patients for a multitude of universal procedures. The power height, back, foot and tilt functions are all operated by an easy to use foot control. The low height of 26" eliminates you and your staff bending and stooping. Features and accessories like a removable foot section, articulating headrest, knee break and adjustable stirrups are durable and functional. A paper roll holder is hidden under the headrest and comes with paper holding tear straps to keep paper in place. A sculptured, form-fitting cushion made of seamless upholstery, makes clean-up faster and more effective. Replace XXX in item number to specify upholstery colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender.

319-003-XXX Power Examination Table w/top

**RITTER 391 ENT CHAIR**

Designed specifically with the otolaryngologist in mind. That's because it offers a unique, thin back section allowing extra close patient access and a magnetic headrest, which adjusts into position precisely and effortlessly. For added convenience, the membrane touch pad switches, which move the base and back sections of the chair, are located on both sides of the chair in a convenient location. Other features of this product include: one programmable memory position, a home function which with the single touch of a switch returns the chair to a seated position for patient exit, full function foot control, flat table positioning, swing out armrests, 360° base rotation with foot-lock on both sides, and seamless vacuum-formed upholstery. Replace XXX in item number to specify upholstery colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender.

391-002 ENT Chair w/out upholstered top
391-002-XXX ENT Chair w/upholstered top

RITTER 391 ENT CHAIR ACCESSORIES

154083100 Nile Green, Upholstery Set Only
154083101 Royal Lavender, Upholstery Set Only
154083102 Terra Rosa, Upholstery Set Only
154083103 Slate Grey, Upholstery Set Only
15408395 Sand Grey, Upholstery Set Only
15408396 Iris Blue, Upholstery Set Only
154087100 Articulating Headrest For 391, Nile Green
154087101 Articulating Headrest For 391, Royal Lavender
154087102 Articulating Headrest For 391, Terra Rosa
154087103 Articulating Headrest For 391, Slate Grey
15408795 Articulating Headrest For 391, Sand Grey
15408796 Articulating Headrest For 391, Iris Blue
154088100 IV Armboard, Nile Green
154088101 IV Armboard, Royal Lavender
154088102 IV Armboard, Terra Rosa
154088103 IV Armboard, Slate Grey
15408895 IV Armboard, Sand Grey
15408896 IV Armboard, Iris Blue
9A240001 Sunnex Light For 391 (Field Installed Only)

**RITTER 75 EVOLUTION POWER EXAMINATION TABLE**

A truly universal power table with eight-way power, allows you to precisely and effortlessly position patients for a multitude of universal procedures. The power height, back, foot and tilt functions are all operated by an easy to use foot control. The low height of 26" eliminates you and your staff from bending and stooping. Features and accessories like a removable foot section, articulating headrest, knee break and adjustable stirrups are durable and functional. A paper roll holder is hidden under the headrest and comes with paper holding tear straps to keep paper in place. A sculptured, form-fitting cushion made of seamless upholstery, makes clean-up faster and more effective. Replace XXX in item number to specify upholstery colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender.

75-015 Ritter 75 Evolution Table, 4-Function, Foot Control, w/out upholstered top
75-015-XXX Ritter 75 Evolution Table, 4-Function, Foot Control, w/upholstered top

**RITTER 75E UNIVERSAL PROCEDURE TABLE**

Non-programmable w/foot control only

75-019 Ritter 75e Universal Procedures Table w/out top
75-019-201 Ritter 75e Universal Procedures Table, Slate Grey
75-019-202 Ritter 75e Universal Procedures Table, Iris Blue
75-019-203 Ritter 75e Universal Procedures Table, Nile Green
75-019-204 Ritter 75e Universal Procedures Table, Terra Rosa
75-019-205 Ritter 75e Universal Procedures Table, Royal Lavender
75-019-216 Ritter 75e Universal Procedures Table, Pebble Grey



BREWER ACCESS™ EXAM TABLE

Spacious, secure and adjustable patient step. Large abundant storage for supplies. Versatile pass-thru back drawers that open from both sides of the table. Reversible outlets that can be easily positioned on the right or left side. Comfortable, easy to clean upholstery. Smooth, quiet pneumatic cylinder for effortless positioning. Specify upholstery color code when ordering: (1) Clamshell, (2) Gunmetal, (3) Adriatic, (4) Burgundy, (5) Emerald, (6) Mauve, (7) Azure Blue, (8) Royal Plum, (9) Plum, (10) Deep Teal, (11) Petal, (12) Blue Fog, (13) Forest Green, (14) Spanish Tile, (15) Chocolate, (16) Black, (17) Putty, (18) Graphite, (SP) Special Custom Color (additional charges may apply).

5000-XX	Access Exam Table, Specify Color
5001-XX	Access Exam Table w/pelvic tilt & drawer heater, Specify Color
5000-CAL133	Access Exam Table, 5000 Model, Cal 133
5001-CAL133	Access Exam Table w/pelvic tilt and drawer heater, Cal 133

Accessories

5115-CAL133	Upholstered Replacement Top, Cal 133
5115-XX	Upholstered Replacement Top (Standard Colors Only) Specify Color



BREWER BASIC EXAM TABLE

Designed with an emphasis on comfort and practicality. Comes standard with a pneumatic cylinder which makes patient positioning a smooth, easy and relaxing experience for you and your patient. Has two convenient front drawers (19.2" x 19.7" x 4") and two large side drawers (22.5" x 20" x 7") that are available on the right or left side of the table. Features a sturdy, retractable patient step with a durable mat. The upholstered top is made of multiple layers of resilient foam, contoured and curved to the body's natural shape. Seamless upholstery, void of stitching and embossing provides a smooth surface for easy cleaning. Equipped with an easy to change paper roll which accepts 21" x 3.5" and 18" x 3.5" paper. Comes with a 3 year warranty. Weight capacity: 325lbs; Height: 32"; Table Length: 57.5" extends to 71"; Table Width: 26.5"; Leg Extension: 15.2". Available in Clamshell-01, Gunmetal-02, Adriatic-03, Burgundy-04, Emerald-05, Mauve-06, Azure Blue-07, Royal Plum-08, Plum-09, Deep Teal-10, Petal-11, Blue Fog-12, Forest Green-13, Spanish Tile-14, Chocolate-15, Black-16, Putty-17, and Graphite-18.

400-XX-L-E	4001 Model w/electrical outlets, Pelvic Tilt & Drawer Warmer, Left Handed Configuration
4000-XX-L	4000 Model, Left Handed Configuration
4000-XX-L-E	4000 Model, Left Handed Configuration, Elec. Outlet
4000-XX-R	4000 Model, Right Handed Configuration
4000-XX-R-E	4000 Model, Right Handed Configuration, Elec. Outlet
4001-XX-L	4001 Model w/pelvic tilt and drawer warmer, Left Handed Configuration
4001-XX-L-E	4001 Model w/electrical outlets, Pelvic Tilt & Drawer Warmer, Left Handed Configuration
4001-XX-R	4001 Model w/pelvic tilt and drawer warmer, Right Handed Configuration
4001-XX-R-E	4001 Model w/electrical outlets, Pelvic Tilt & Drawer Warmer, Right Handed Configuration



HAUSMANN MANUAL EXAMINATION TABLE

Two drawers on front and three drawers on side. Recessed, full width black plastic pulls are standard. Satin aluminum pulls are optional. Paper dispenser and cutter. Pneumatic backrest. Removable Velcro® vinyl top. Four-position laterally adjustable stirrups. Slide-out leg rest with Velcro® pad. Extra wide pull-out footstep with rounded corners. Specify color when ordering. 55"-71"L x 26"W x 33"H.

4400	Right Hand Model Table
4405	Left Hand Model Table



4420



4423

HAUSMANN MAXIMUM VALUE™ EXAMINATION TABLES

Vinyl upholstered removable Velcro top pad. Pneumatic backrest with dual adjustment controls from 0-75°. Sliding 21"W x 16"L legrest extension with removable pad. Heavy duty concealed stirrups adjust in length and are laterally adjustable in four positions. Built-in rubber tread foot step 21" wide. Paper dispenser and cutter. Hospital grade duplex outlet. Black formica base with leveling legs. 55"-71"L x 26"W. Laminate folkstone gray with selection of trim colors.

4420	Table, w/2 front drawers & 1 side door w/shelf
4421	Table, Left Hand, w/2 front drawers & 1 side door w/shelf
4423	Table, w/2 front drawers & 3 side doors
4424	Table, Left Hand, w/2 front drawers & 3 side doors



HAUSMANN X-L POWER-BACK EXAM TABLE

Highest patient weight capacity at 600 lbs. 30" wide one-piece upholstered top. Strongest slide-out footstep. Backrest is motorized to facilitate patient positioning. Heavy duty concealed stirrups. Equipped with five drawers with laminate drawer fronts and 3mm protective edging. Welded heavy-duty steel chassis with durable powder coated finish.

4416	X-L Power-Back Exam Table
------	---------------------------

Tables



1020 Series



1010 Series



1030 Series



1013 Series



533-38



1520-38



1570-38

CLINTON CLASSIC LINE TREATMENT TABLES

2" high-density urethane foam padding. Dark Cherry or Natural finish. Solid hardwood legs. Secure H-brace construction. Scratch and stain resistant surfaces. 350 lbs capacity. Matching laminate shelving. Backrest standard on most models. Tops secured to frames with steel angles. Premium, stain-resistant, woven, knit-backed vinyl upholstery. Welt strips conceal upholstery strips. Three-year limited warranty. 16 top colors available.

- 1010-24 Straight Line Treatment Table w/hardwood legs & adjustable backrest, 72"L x 31"H x 24"W
- 1010-27 Straight Line Treatment Table w/hardwood legs & adjustable backrest, 72"L x 31"H x 27"W
- 1010-30 Straight Line Treatment Table w/hardwood legs & adjustable backrest, 72"L x 31"H x 30"W
- 1013-27 Straight Line Treatment Table w/hardwood legs, laminate shelf, two drawers & adjustable backrest, 72"L x 31"H x 27"W
- 1013-30 Straight Line Treatment Table w/hardwood legs, laminate shelf, two drawers & adjustable backrest, 72"L x 31"H x 30"W
- 1020-24 Straight Line Treatment Table w/hardwood legs, laminate shelf & adjustable backrest, 72"L x 31"H x 24"W
- 1020-27 Straight Line Treatment Table w/hardwood legs, laminate shelf & adjustable backrest, 72"L x 31"H x 27"W
- 1020-30 Straight Line Treatment Table w/hardwood legs, laminate shelf & adjustable backrest, 72"L x 31"H x 30"W
- 1030-24 Straight Line Treatment Table w/hardwood legs, laminate shelf, storage compartments & adjustable backrest, 72"L x 31"H x 24"W
- 1030-27 Straight Line Treatment Table w/hardwood legs, laminate shelf, storage compartments & adjustable backrest, 72"L x 31"H x 27"W
- 1030-30 Straight Line Treatment Table w/hardwood legs, laminate shelf, storage compartments & adjustable backrest, 72"L x 31"H x 30"W
- 1520-38 Mat Top Treatment Table w/drawer & full length laminate shelf, 78"L x 31"H x 30"W
- 1533-38 Mat Top Treatment Table w/laminate storage compartment, drawer & full laminate shelf, 78"L x 31"H x 30"W
- 1570-38 Cabinet-Style Mat Top Treatment Table w/2 laminate sliding doors & adjustable laminate shelf, 78"L x 31"H x 30"W



7114



7122



7133



7144




7177

CLINTON SELECT SERIES TREATMENT TABLES

Easy to clean laminate surfaces. Vertical spanner adds support to panel legs. 2" thick foam. Welt strip conceals upholstery staples. Scratch and mar resistant surfaces. 350 lbs capacity. Chip resistant Slate Gray with PVC edge trim. All stress points reinforced by steel. Three-year limited warranty. Fifteen laminate colors. Sixteen vinyl colors.

- 7114 Infant Treatment/Changing Table w/full shelf, rimmed top on 3 sides, removable pad & 2 round stainless steel top rails, 46"L x 36"H x 24"W
- 7120 Straight Line Treatment Table w/full length storage shelf, 72"L x 31"H x 27"W
- 7122 Straight Line Treatment Table w/full length storage shelf & adjustable backrest, 72"L x 31"H x 27"W
- 7130 Treatment Table w/storage shelf & 2 Euro-style drawer, 72"L x 31"H x 27"W
- 7133 Treatment Table w/adjustable backrest, storage shelf & 2 Euro-style drawer, 72"L x 31"H x 27"W
- 7140 Cabinet Style Laminate Treatment Table w/2 storage compartments w/sliding doors, 72"L x 31"H x 27"W
- 7144 Cabinet Style Laminate Treatment Table w/adjustable backrest, 2 storage compartments w/sliding doors, 72"L x 31"H x 27"W
- 7170 Cabinet Style Treatment Table w/2 sliding doors front & back, 72"L x 31"H x 27"W
- 7174 Cabinet Style Laminate Treatment Table w/4 doors, adjustable backrest & 2 storage compartments, 72"L x 31"H x 27"W
- 7177 Cabinet Style Treatment Table w/2 sliding doors front & back & adjustable backrest, 72"L x 31"H x 27"W
- 7178 Cabinet Style Laminate Treatment Table w/4 doors & 2 storage compartments, 72"L x 31"H x 27"W

Call for current pricing!



2020 Series




2021 Series



2010 Series



2030 Series



2040



2070

CLINTON SIGNATURE SERIES OAK TABLES

2" high-density urethane foam padding. Matching laminate shelves and cabinets. Solid oak legs with leg levelers for stability. Scratch and stain resistant finish. Tops are reinforced and securely attached with steel angles. Backrest standard on all models. 350 lbs capacity. Premium, stain-resistant, woven, knit-backed vinyl upholstery. Welt strips conceal upholstery strips. Natural and Dark Cherry finish with choice of 16 top colors. Three-year limited warranty.

- 2010-24 Oak Table w/solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 24"W
- 2010-27 Oak Table w/solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 27"W
- 2010-30 Oak Table w/solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 30"W
- 2020-24 Oak Table w/laminated full shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 24"W
- 2020-27 Oak Table w/laminated full shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 27"W
- 2020-30 Oak Table w/laminated full shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 30"W
- 2021-24U Oak Table w/vinyl covered full shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 24"W
- 2021-27U Oak Table w/vinyl covered full shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 27"W
- 2021-30U Oak Table w/vinyl covered full shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 30"W
- 2030-27 Oak Table w/solid hardwood legs, H-braces, adjustable backrest, laminated full shelf & storage unit, 72"L x 31"H x 27"W
- 2030-30 Oak Table w/solid hardwood legs, H-braces, adjustable backrest, laminated full shelf & storage unit, 72"L x 31"H x 30"W
- 2040-27 Oak Table w/laminated partitioned full length shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 27"W
- 2040-30 Oak Table w/laminated partitioned full length shelf, solid hardwood legs, H-braces & adjustable backrest, 72"L x 31"H x 30"W
- 2070-27 All-Laminate Table w/solid oak legs, 2 sliding doors, adjustable shelf & adjustable backrest, 72"L x 31"H x 27"W
- 2070-30 All-Laminate Table w/solid oak legs, 2 sliding doors, adjustable shelf & adjustable backrest, 72"L x 31"H x 30"W

Accessories

- 079 Optional Triple Bolted Legs For Extra Strength & Stability



9013 Series



9044 Series



9140 Series




9133 Series

CLINTON STYLE LINE TREATMENT TABLES

2" high-density urethane foam padding. Easy to clean laminate surfaces. Vertical spanner adds support to panel legs. Euro-style drawers with all steel and laminate construction. Scratch and mar resistant surfaces. Chip resistant PVC edge trim. All stress points reinforced with steel angles. 350 lbs capacity. Backrest standard on most models. Premium, stain-resistant, woven, knit-backed vinyl upholstery. Welt strips conceal upholstery strips. Three-year limited warranty. 16 top colors.

- 9013-27 Straight Line Treatment Table w/2 drawers, shelf & adjustable backrest, 72"L x 31"H x 27"W
- 9013-30 Straight Line Treatment Table w/2 drawers, shelf & adjustable backrest, 72"L x 31"H x 30"W
- 9020-27 Straight Line Treatment Table w/full length storage shelf & adjustable backrest, 72"L x 31"H x 27"W
- 9020-30 Straight Line Treatment Table w/full length storage shelf & adjustable backrest, 72"L x 31"H x 30"W
- 9044-27 Cabinet Style Laminate Treatment Table w/adjustable backrest, 2 storage compartments w/sliding doors & adjustable self in each compartment, 72"L x 31"H x 27"W
- 9044-30 Cabinet Style Laminate Treatment Table w/adjustable backrest, 2 storage compartments w/sliding doors & adjustable self in each compartment, 72"L x 31"H x 30"W
- 9070-27 Cabinet Style Treatment Table w/2 sliding doors front & back plus adjustable backrest, 72"L x 31"H x 27"W
- 9070-30 Cabinet Style Treatment Table w/2 sliding doors front & back plus adjustable backrest, 72"L x 31"H x 30"W
- 9074-24 Cabinet Style Laminate Treatment Table w/4 doors, adjustable backrest & 2 storage compartments each w/adjustable shelf, 72"L x 31"H x 24"W
- 9074-27 Cabinet Style Laminate Treatment Table w/4 doors, adjustable backrest & 2 storage compartments each w/adjustable shelf, 72"L x 31"H x 27"W
- 9074-30 Cabinet Style Laminate Treatment Table w/4 doors, adjustable backrest & 2 storage compartments each w/adjustable shelf, 72"L x 31"H x 30"W
- 9102-30 Straight Line Treatment Table w/14" wide double pull-out shelves & full length shelf, 72"L x 31"H x 30"W
- 9104-38 Cabinet Style Laminate Treatment Table w/2 enclosed storage compartments w/adjustable shelves, 22"W pull-out shelf w/locking length positions & 2 adjustable backrests, 78"L x 31"H x 30"W
- 9105-30 Cabinet Style Laminate Treatment Table w/3 drawers, adjustable backrest & storage compartment w/adjustable shelf, 72"L x 31"H x 30"W
- 9106-38 Laminate Cast Table w/adjustable backrest & 3 large storage sections w/adjustable shelves, 78"L x 31"H x 30"W
- 9107-38 Cabinet Style Laminate Treatment Table w/3 drawers, 2 storage compartments w/2 adjustable shelves & adjustable backrest, 78"L x 31"H x 30"W
- 9110-24 Space Saver Table w/adjustable shelf, 2 doors, adjustable footrest & adjustable backrest, 72" extended/56" folded L x 31"H x 24"W
- 9110-27 Space Saver Table w/adjustable shelf, 2 doors, adjustable footrest & adjustable backrest, 72" extended/56" folded L x 31"H x 27"W
- 9133-27 Laminate Treatment Table w/3 drawers, storage shelf, adjustable backrest & matching stool, 72"L x 31"H x 27"W
- 9133-30 Laminate Treatment Table w/3 drawers, storage shelf, adjustable backrest & matching stool, 72"L x 31"H x 30"W

Tables



CLINTON STYLELINE LAMINATED TREATMENT TABLE ACCESSORIES

003	Slide-Out Shelf, 9"W
004	Drawer, 14 1/2"D x 16 1/4"W x 3 1/4"H
005	Table Extra Height, 36" Maximum
006	Table Extra Length, 78" Maximum
008	Armboard, Adjustable
009	Nosehole, 9"L x 2"W
010	Stirrups, Adjustable, Chrome-Plated, Heels Covered w/plastic coating
020	Small Pillow, 14" x 12" x 3"
030	Paper Dispenser, Accommodates Paper 18" to 21"
034	Paper Dispenser & Cutter Combo
040	Paper Cutter, Clear Plastic
060	Full-Size Head Pillow, 14" x 21" x 3"
077	Safety Strap, 2" x 92"



HAUSMANN ALL-PURPOSE TREATMENT TABLE

Three all laminate drawers on nylon rollers. Two fully-enclosed door cabinets. Open center section is 24" wide, 2" thick urethane foam top. Vinyl upholstery finish. Black formica covered plywood base. 400 lb capacity. Available with 99L Left Hand Table option.

4834	78"L x 30"W x 31"H Table
------	--------------------------



HAUSMANN DELUXE TREATMENT TABLE

Two drawers on nylon roller bearings. Hinged door storage section with one shelf. Includes backrest, paper dispenser and paper cutter. Laminate finish. Urethane foam top. Vinyl upholstery finish. 400 lb capacity. Available with 99L Left Hand Table option.

4796	76"L x 27"W x 31"H Table
------	--------------------------



HAUSMANN ENCLOSED CABINET TREATMENT TABLE

Table has one adjustable shelf and two sliding doors.

4043	24"W x 72"L x 31"H Table
4043-027	27"W x 72"L x 31"H Table
4043-030	30"W x 72"L x 31"H Table
4043-030-078	30"W x 78"L x 31"H Table



4602



4624



4643

HAUSMANN EXAM/TREATMENT TABLES

Adjustable head and foot rests. Measures 24" x 72" x 31" (41" with rests in down position). Stirrups are standard. Available with 99L Left Hand Table option.

4602	Table w/H-brace
4624	Table w/H-brace & Full Shelf
4643	Table w/fully enclosed cabinet, 2 Hinged Doors & Adj. Shelf



HAUSMANN MD SERIES™ TREATMENT TABLES

Available with 99L Left Hand Table option.

4259	72" x 27" x 31" Table, Two Base Cabinets w/6 laminate drawers on nylon rollers, 25" Open Center Section
------	---



HAUSMANN MODEL 4298 TREATMENT TABLE

Unique 8" high pull-out, non-slip foot stool. Foot stool slides into base for storage. Three large all laminate drawers on nylon rollers. Full width shelf. 2" thick foam top. 350 lb weight capacity. Available with 99L Left Hand Table option.

4298 72"L x 27"W x 32"H Treatment Table

HAUSMANN MULTI-SECTIONAL TREATMENT TABLE

Unique multi-functional design for therapy, sports or medical use. Includes manually adjustable elevating 0° to 45° headrest with nose cutout. Includes pneumatic gas-spring 0° to 75° backrest with one lever control. Manually adjustable 18" x 18" legrest raises to horizontal position. Fully enclosed Oak laminate cabinet with sliding doors and shelf. 2" high-density urethane foam vinyl top. Upholstery comes in 14 standard vinyl colors.

4044 Multi-Sectional Treatment Table, 72"-90"L x 27"W x 31"H



HAUSMANN OVERSIZE TREATMENT TABLE

Six laminate legs and a sturdy plywood sub-top. 600 lb capacity. H-Brace with full laminate shelf. 2" thick urethane foam top. Vinyl upholstery finish.

4086 78"L x 30"W x 26"H Table
4087 78"L x 30"W x 31"H Table



HAUSMANN PORTABLE MASSAGE TABLE

73" L x 30"W. Height adjustable from 23" to 33". All-hardwood construction. Reinforced bracing. 400 lb capacity. Foam top. Carrying handles. Folds easily for compact storage.

7600-706 Slate Blue Massage Table

Accessories

081-706 Cushioned Face Cradle, Slate Blue
082 Carrying Case w/shoulder strap



HAUSMANN PORTABLE TREATMENT TABLES

72"L x 24"W x 28"H. Sturdy hardwood frame with eight tubular satin aluminum legs. 1" thick urethane foam top. Carrying grip handles, corner guards and bottom bumpers. 300 lb capacity. Folds to a compact 24" x 36" x 6". Not for outdoor use.

7500-079-701 Oak Brown Table w/nose cutout
7500-079-706 Slate Blue Table w/nose cutout
7500-079-711 Grotto Green Table w/nose cutout
7500-701 Oak Brown Table
7500-706 Slate Blue Table
7500-711 Grotto Green Table



4524



4541



4543

HAUSMANN PROFESSIONAL TREATMENT TABLES

20" wide laminate drawer on nylon rollers. 2" thick urethane foam velcro-lock removable top pad. Rimless, black laminate sub-top design eliminates uncomfortable rim edges. Vinyl upholstery finish in 12 colors. 350 lb capacity. 78"L x 30"W x 31"H. Available with 99L Left Hand Table option.

4524 Table, H-Brace w/full width velcro-lock laminate shelf
4541 Table, H-Brace w/laminate shelf & open shelf storage cabinet
4543 Table, Fully-Enclosed Storage Cabinet w/sliding doors & one shelf



4102



4143

HAUSMANN SPACE-SAVER TREATMENT TABLES

Full width 20" long foot end folds down to reduce closed table size to only 52". Positive action locking devices allows the foot end to be easily elevated into three positions. 2" thick urethane foam top. Specify color when ordering. Available with 99L Left Hand Table option.

4102 Table w/recessed H-Brace in Table Leg
4124 Table, H-Brace w/full width laminate shelf
4143 Table, Fully Enclosed Cabinet w/doors & adjustable shelf



4002



4024



4041

HAUSMANN TREATMENT TABLES WITH H-BRACE

Table tops generously padded with 2" thick urethane foam and covered with durable vinyl. Choice of 10 upholstery colors. High-pressure formica finish resists stains and scratches. H-brace assembly with lock-tite increases strength and durability. Available in gray, cherry and natural oak. 31"H. Available with 99L Left Hand Table option.


4002	24"W x 72"L Table w/H-brace
4002-027	27"W x 72"L Table w/H-brace
4002-030	30"W x 72"L Table w/H-brace
4002-078	24"W x 78"L Table w/H-brace
4024	24"W x 72"L Table w/H-brace & Full Shelf
4024-027	27"W x 72"L Table w/H-brace & Full Shelf
4024-030	30"W x 72"L Table w/H-brace & Full Shelf
4024-078	24"W x 78"L Table w/H-brace & Full Shelf
4041	24"W x 72"L Table w/H-brace, Shelf & Open-Shelf Storage Cabinet
4041-027	27"W x 72"L Table w/H-brace, Shelf & Open-Shelf Storage Cabinet
4041-030	30"W x 72"L Table w/H-brace, Shelf & Open-Shelf Storage Cabinet
4002-030-078	30"W x 78"L Table w/H-brace
4024-030-078	30"W x 78"L Table w/H-brace & Full Shelf
4041-030-078	30"W x 78"L Table w/H-brace, Shelf & Open-Shelf Storage Cabinet



ITTER 95 CLASSIC SERIES EXAMINATION TABLE

Ideal for new medical practices, industrial first aid areas, school nurses' stations or as an additional table for any practice. Its hardwood construction and h-brace frame give it strength and stability, while its simple design allows easy access to patients without obstructions and is your most economical option available in many colors.

95-100-200	27"W x 72"L x 30"H Exam Table, Sand Grey
95-100-201	27"W x 72"L x 30"H Exam Table, Slate Grey
95-100-202	27"W x 72"L x 30"H Exam Table, Iris Blue
95-100-203	27"W x 72"L x 30"H Exam Table, Nile Green
95-100-204	27"W x 72"L x 30"H Exam Table, Terra Rosa
95-100-205	27"W x 72"L x 30"H Exam Table, Royal Lavender
95-100-206	27"W x 72"L x 30"H Exam Table, Wedgewood
95-100-207	27"W x 72"L x 30"H Exam Table, Monaco
95-100-208	27"W x 72"L x 30"H Exam Table, Royal Carriage
95-100-209	27"W x 72"L x 30"H Exam Table, Monet
95-100-210	27"W x 72"L x 30"H Exam Table, Soft Pearl
95-100-211	27"W x 72"L x 30"H Exam Table, Waterford
95-100-212	27"W x 72"L x 30"H Exam Table, English Gardens
95-100-213	27"W x 72"L x 30"H Exam Table, Cabernet
95-100-214	27"W x 72"L x 30"H Exam Table, Tapestry
95-100-215	27"W x 72"L x 30"H Exam Table, Cathedral



7370-71

CLINTON FIND-4 TREATMENT TABLES

A few minutes can seem like hours to a child waiting for the doctor or nurse to come into the exam room. Entertain the young patient with a "Find-4" table while he or she waits. It's designed to be easy for staff, yet engaging for children. Ask the child to find one of the fun objects or animals hidden in the picture. There are: 4 brown eggs, 4 cats, 4 birds and 4 mice in each scene. Keep the number of the objects a secret. Just ask the child how many he or she found. Whether it is all or a few, every child's count is correct to you. Next visit ask them to find one of the other 3 remaining objects. Keep young patients amused for many visits to come. Each table is 72"L x 31"H x 27"W and features 4 doors, adjustable backrest and 2 storage compartments each with an adjustable shelf.

7370-71	Sweet Dreams Candy Factory Scene Treatment Table
7370-72	Wings 'N Things Scene Treatment Table
7370-73	Secret Garden Scene Treatment Table
7370-74	Fairy Tale Dale Scene Treatment Table



7020



7030



7040



7822



7833



7844

CLINTON FUN SERIES TABLES

Easy to clean laminate surfaces. Graphics are embedded in the laminate. Graphics on three sides. Durable Perma Block top coated vinyl. Fire-rated 2" thick foam padding. Doors mounted with Euro-style hinges. Large storage capacity.

7020	Zoo Bus Table
7030	Engine K-9 Table
7040	Outback Buggy Table
7050	Cool Camper
7060	Medicine Show Wagon
7822	Zoo Bus Table w/scale
7833	Engine K-9 Table w/scale
7844	Outback Buggy Table w/scale

CLINTON SELECT SERIES GRAPHICS

Fun accents that brighten up any Select Series table, pediatric cabinets or anywhere a touch of fun is desired. Printed on tough, adhesive-backed, durable vinyl. They are washable and resist tearing and peeling.

9610	Baseball Table Graphics
9640	Play Blocks Table Graphics
9650	Under The Sea Table Graphics
9660	Sports Table Graphics
9680	Alphabet Table Graphics
9721	Zebra Graphics, Approx. 14"H
9722	Fish Graphics, Approx. 14"H
9723	Airplane Graphics, Approx. 10"H
9724	Tennis Ball Graphics, Approx. 11"H
9725	Firefighter Graphics, Approx. 14"H
9726	Koala Graphics, Approx. 14"H
9727	Basketball Graphics, Approx. 12"H
9729	Friendly Bear Graphics, Approx. 14"H
9731	T-Depressor Eater Graphics, Approx. 6"H
9732	Cottontail Graphics, Approx. 11"H
9733	Tiger Buddies Graphics, Approx. 14"H
9735	Bandage Tiger Graphics, Approx. 7"H
9736	Gauze Raccoon Graphics, Approx. 4"H
9737	Smiling Seals Graphics, Approx. 14"H
9741	Swab Beaver Graphics, Approx. 6"H
9742	Baseball Graphics, Approx. 12"H
9743	Froggie Graphics, Approx. 10"H

CLINTON SELECT SERIES SCALE TABLES

Comfortable work height of 36". Top is rimmed on three sides for safety. The front and back rim edges are beveled and the corners rounded for safety. Easy-clean, replaceable center pad. Built-in infantometer with folding ends for measuring height. Easy access to weight, measurement and storage. Built-in concealed paper dispenser and cutter. Recessed front toe kick. Adjustable built-in levelers provide a stable surface. Easy-clean all laminate base. Ample storage with various door and drawer configurations. Scale features include: dependable accurate Tanita scale, easy-read digital display, autocalibration on start-up, auto shut-off, electric or battery operation (6 AA batteries required), weights in increments of lbs., ozs. or kilos, maximum weight 40 lbs.

7810	Pediatric Scale/Treatment Table w/1 drawer & 2 open shelf storage areas, 43"L x 36"H x 24"W
7820	Pediatric Scale/Treatment Table w/2 doors & 1 adjustable shelf, 43"L x 36"H x 24"W
7830	Pediatric Scale/Treatment Table w/2 doors & 1 drawer storage, 43"L x 36"H x 24"W
7840	Pediatric Scale/Treatment Table w/2 doors & 2 drawer storage, 43"L x 36"H x 24"W



7836



7837

CLINTON THEME SERIES SCALE TABLES

A picture is worth a thousand words and these tables are a picture of fun that needs no words. Brightly adorned with color, these whimsical scenes will enchant and delight young patients and parents alike. All tables are 24"D x 43"L x 36"H and come with one drawer and two doors.

7831	Cool Pals Scale Table
7832	Rainforest Follies Scale Table
7834	Dino Day's Scale Table
7835	Space Place Scale Table
7836	Ocean Commotion Scale Table
7837	Alley Cats & Dogs Scale Table



7931



7932

CLINTON THEME SERIES TREATMENT TABLES

Easy to clean laminate surfaces. All graphics are permanently embedded in the high pressure laminate and will not peel off. 2" of firm foam padding. Four door storage with two adjustable shelves. Slat Gray ends and backs. Scratch and mar resistant surfaces. Chip resistant PVC edges. Premium, stain-resistant, woven, knit-backed vinyl upholstery.

7931	Cool Pals Treatment Table w/large storage capacity & adjustable backrest, 68"L x 33"H x 27"W
7931-1	Cool Pals Treatment Table w/large storage capacity, 68"L x 33"H x 27"W
7932	Rainforest Follies Treatment Table w/large storage capacity & adjustable backrest, 68"L x 33"H x 27"W
7932-1	Rainforest Follies Treatment Table w/large storage capacity, 68"L x 33"H x 27"W
7934	Dino Days Treatment Table w/large storage capacity & adjustable backrest, 68"L x 33"H x 27"W
7934-1	Dino Days Treatment Table w/large storage capacity, 68"L x 33"H x 27"W
7935	Space Place Treatment Table w/large storage capacity & adjustable backrest, 68"L x 33"H x 27"W
7935-1	Space Place Treatment Table w/large storage capacity, 68"L x 33"H x 27"W
7936	Ocean Commotion Treatment Table w/large storage capacity & adjustable backrest, 68"L x 33"H x 27"W
7936-1	Ocean Commotion Treatment Table w/large storage capacity, 68"L x 33"H x 27"W
7937	Alley Cats & Dogs Treatment Table w/large storage capacity & adjustable backrest, 68"L x 33"H x 27"W
7937-1	Alley Cats & Dogs Treatment Table w/large storage capacity, 68"L x 33"H x 27"W



HAUSMANN DALMATION TREATMENT TABLE

Sturdy H-brace design. All-laminate construction is easy to clean and maintain. Choice of 10 Resist-All vinyl upholstery colors.

4904	58"L x 24"W x 34"H Treatment Table
------	------------------------------------

HAUSMANN PEDIATRIC DIGITAL SCALE TABLE

Digital scale with lb-oz readout. Reads in 1 oz increments. Lb-oz/kg switch. Auto power shut-off. Electric and battery operation. 39" long infantometer to measure infant. 12"H slide-out drawer for desk charger. Two 6"H drawers and one 18"W door. Paper dispenser and cutter. Padded Navy Blue vinyl top with rim on three sides. If ordered without Wall Panel 4992, scale head attaches to table.

4947	Digital Scale Table, 43"L x 24"D x 36"H
------	---

Accessories

4992	Triangular Wall Panel For Use w/#4947
------	---------------------------------------

Tables

HAUSMANN PEDIATRIC TREATMENT TABLE

Finished in child friendly decorative laminate colors. Three 18"W yellow laminate drawers on nylon rollers. Three 18"W blue laminate doors with shelf. Paper dispenser and paper cutter. Removable 70" x 24" navy blue top pad on Velcro® brand fasteners.

4990 Treatment Table, 72"L x 24"D x 34"H

Accessories

4991 Rectangular Wall Panel For Use w/#4990 Table



4960



4961



4966



4967



4968, 4969

HAUSMANN PEDIATRIC TREATMENT TABLES

The colorful characters are permanently embedded into the high-pressure laminate surfaces. They won't peel off or fade away.

- 4960 Tonka™ Treatment Table, Yellow w/black trim, Black base w/four 10" wheels, Front End Loader Serves as Paper Roll Dispenser, Plexi-Mirror Side Window & Windshield
- 4961 Candy Land™ Treatment Table w/o Plate, w/39" ruler, 3 Self-Closing Doors & 3 Drawers on Nylon Roller Self-Closing Slides
- 4962 Candy Land™ Treatment Table with Plate, w/39" ruler, 3 Self-Closing Doors & 3 Drawers on Nylon Roller Self-Closing Slides
- 4966 Mr. Potato Head™ Treatment Table w/39" ruler, 3 Self-Closing Doors & 3 Drawers on Nylon Roller Self-Closing Slides
- 4967 Mr. Potato Head™ Cabinet
- 4968 Mr. Potato Head™ Hamper
- 4969 Mr. Potato Head™ Smart Cart



HAUSMANN PETE'S PICKUP™ PEDIATRIC TABLE

Child friendly molded blue plastic car chassis. Built-in toddler car seat with dual steering wheels. Two fixed wheels and black cabinet base trim. One blue laminate drawer on nylon rollers. Two blue laminate doors on concealed, self-closing euro-hinges. Yellow vinyl 52" x 24" padded upholstered top with 18" x 18" legrest extension. Legrest elevates up to 90° and extends top length to 70". Includes paper dispenser and cutter. 2" high-density urethane foam tops.

4975 Pediatric Table, 60"-78"L x 24"W x 34"H



HAUSMANN TRAIN CAR DESIGN PEDIATRIC TABLE

Colorful exam/treatment table minimizes young patient's anxiety while enhancing the room. Padded top for comfort. Storage drawer plus three roomy cabinets provide generous space for storage. 24"W x 72"L x 34"H. Upholstery: #704 Nordic Blue or #707 Black. Specify upholstery color.

4970 Tables for Tois™

Accessories

2216A Matching Blue Footstep, 20" x 16" x 8"H



HAUSMANN WALL FOLDING TREATMENT/CHANGING TABLE

Fold-away design saves floor space. 2" thick foam upholstered top. Sturdy, floor-supported wall mounting plate for added strength. Two folding Formica® natural oak laminate legs. Upholstery available in 10 vinyl colors.

4075 72"L x 24"W x 31"H Table



309-005



309-006

RITTER 309 PEDIATRIC EXAMINATION TABLE

Everything needed to treat infants in one convenient unit. Includes a standard, built-in measurement arm. Arm pivots out to measure and folds back in, out of physician's way. Head end has an upholstered panel to make measuring easier and more accurate. Hidden paper roll holder and padded one-piece top with a protective finish make cleaning easy. Top features beveled sides and upholstered end panel which helps keep baby from rolling off. Flat sides and no sharp edges increase safety. Table's 36" height minimizes bending and stooping. Extra large, easy to clean storage drawer and bulk storage underneath provide ample room for diapers and equipment. Built-in digital scale weighs babies right on counter. 309-005 Pediatric Exam Table with digital scale (.2 lbs to 40 lbs, 100 grams to 90 kilograms range in .2 increments; measures given in both imperial or metric). 309-006 Pediatric Exam Table without digital scale. Replace XXX in item number with specify upholstery colors: 200-Sand Grey/Sand Grey, 201 Slate Grey/Slate Grey, 201-1 Slate Grey/Sand Grey, 202 Iris Blue/Iris Blue, 202-1 Iris Blue/Sand Grey, 203 Nile Green/Nile Green, 203-1 Nile Green/Sand Grey, 204 Terra Rosa/Terra Rosa, 204-1 Terra Rosa/Sand Grey, 205 R. Lavender/R. Lavender, 205-1 R. Lavender/Sand Grey, 206-1 Wedgewood/Sand Grey, 207-1 Monaco/Sand Grey, 208-1 Royal Carriage/Sand Grey, 209-1 Monet/Sand Grey, 210-1 Soft Pearl/Sand Grey, 211-1 Waterford/Sand Grey, 212-1 English Gardens/Sand Grey, 213-1 Cabernet/Sand Grey, 214-1 Tapestry/Sand Grey, 215-1 Cathedral/Sand Grey, 999 Special Colors.

309-005 Pediatric Examination Table w/scale

309-006 Pediatric Examination Table

CLINTON ECHO IMAGE SERIES TABLES

2" high-density urethane foam padding. Stationary models feature easy to clean laminate surfaces available in several colors. The top and frame are upholstered together to prevent dirt from collecting in seams. All metal slide bar supports window drop. 11 position adjustable headrest. 400 lbs capacity. Power units feature: all welded steel frame, single drive motor system, hospital-grade 100V plug, fully enclosed lift system, and hands-free foot control.

- 4615 Echo Image Table w/full triple leg support, 76"L x 31"H x 27"W, 10"W x 12"D Drop Leaf Window
- 4620 Power Echo Image Table, 76"L x 20"-30"H x 33"W, 8"W x 11"D Drop Leaf Window
- 4630 Power Echo Image Table w/built-in flip-down casters, 76"L x 20"-30"H x 33"W, 8"W x 11"D Drop Leaf Window
- 4690V Vascular Table w/tilting table top, 71"L x 24"-34"H x 32"W, 8"W x 11"D Fold-Down Section



1702-30



1703-30

CLINTON SPORTS TRAINING TABLES

- 1702-30 Heavy Duty Taping Table w/solid wood legs, H-braces & laminated full shelf, 350 lbs Capacity, 48"L x 36"H x 30"W
- 1703-30 Heavy Duty Taping Table w/solid wood legs, H-braces & laminated 2 level storage compartments & shelf, 350 lbs Capacity, 48"L x 36"H x 30"W
- 1701-27 Heavy Duty Taping Table w/solid wood legs & H-braces, 350 lbs Capacity, 48"L x 36"H x 27"W
- 1701-30 Heavy Duty Taping Table w/solid wood legs & H-braces, 350 lbs Capacity, 48"L x 36"H x 30"W
- 1702-27 Heavy Duty Taping Table w/solid wood legs, H-braces & laminated full shelf, 350 lbs Capacity, 48"L x 36"H x 27"W
- 1703-27 Heavy Duty Taping Table w/solid wood legs, H-braces & laminated 2 level storage compartments & shelf, 350 lbs Capacity, 48"L x 36"H x 27"W
- 1763 Knee Gatch Table w/panel legs, 2 level laminated storage compartments & open shelf, Adjustable Backrest & Adjustable Foot Section, 350 lbs Capacity, 78"L x 31"H x 30"W
- 1764 Double Leg Lift Table w/panel legs, 2 level laminated storage compartments & shelf, Each Leg Rest Can Be Raised & Lowered Independently, 78"L x 31"H x 30"W
- 1772 2 Station Taping Table w/3 storage sections w/adjustable shelves, Self-Closing Doors, Center Drawer, 76"L x 36"H x 30"D
- 1780 3 Station Taping Table/Treatment Table Combo w/3 storage sections w/adjustable shelves, Enclosed Cabinet w/2 doors & adjustable shelf, 78"L x 36"H x 30"W



HAUSMANN BARIATRIC ELECTRIC TILT TABLE

Electric tilt mechanism with heavy-duty motor for up to 500 lb patients. Tilt indicator displays angle of tilt 0 to 90°. Extra-wide 34" x 77" long top unholstered in 709 Gray vinyl with 2" thick high-density urethane foam. 23" extra-wide x 18" deep fixed foot board. Three extra-long patient safety straps. Fixed table height of 33". Pneumatic hand control operable from either side of the table. Set of four 3" locking casters. Heavy-duty gauge steel frame is powder coated in neutral cream finish. 3.5" clearance under frame. 110V, 60Hz, 3.5 amps. Hospital grade plug, UL listed.

- 6045-709 Bariatric Electric Tilt Table



4325



4326



4327

HAUSMANN CRANK HEIGHT HI-LO WORK TABLES

Adjustable from 26" to 38"H. Removable handle. Accommodates wheelchair or ambulatory patients. Heavy gauge steel frame, powder coated in neutral cream finish.

- 4324 60"L x 28"W Table w/natural oak laminate top
- 4325 60"L x 28"W Table w/butcher block top
- 4326 66"L x 48"W Table w/butcher block top
- 4327 66"L x 48"W Table w/natural oak laminate top

HAUSMANN ECHO-SCAN™ TABLE

Two laminated solid panel legs. Convenient laminate storage shelf. Folkstone Gray laminate finish. 3" high-density urethane foam top. 10"W x 12"D drop-away section on durable piano hinge. Sliding latch activates drop-away section.

- 4892 Echo-Scan™ Table, 76"L x 27"W x 31"H



HAUSMANN ECONOMY STAND-IN TABLE

Designed to promote vertical weight bearing activities. Sturdy laminated wood construction. Fixed work surface height of 45". Height of standing platform can be adjusted in 10" vertical range to accommodate different size patients. Standing platform has anti-slip treads. Laminate top is 30" x 30" with 15"W x 12"D cutout. Adjustable height padded knee rest. Heavy duty hinges and positive latch. Folkstone Gray laminate finish.

- 6125 Economy Stand-In Table

Tables



6175



6180

HAUSMANN ELECTRIC STAND-IN TABLES

Motorized height adjustment of top. Top is 33"W x 24"D. Cutout size is 19"W x 10"D. Hand control operates electric tabletop height and motorized patient lift. Includes lift support harness. Padded knee, chest and back supports. Heavy gauge steel frame.

- 6175 Econo-Line Stand-In Table w/2 locking & 2 non-locking heavy-duty casters, Can Lift a 300 lb Patient
6180 Stand-In Table w/retractable locking swivel casters, Can Lift a 500 lb Patient



HAUSMANN EXTRA-WIDE BARIATRIC TREATMENT TABLE

3" extra-thick high-density urethane foam top for bariatric use. Extra-wide 34" x 78" long upholstered top. Reinforced oak laminated wood construction can accommodate patients up to 800 lbs. Three storage cabinets, each with laminate interiors and one adjustable self. Doors have concealed, self-closing, Euro-hinges. Natural Oak laminate. Choice upholstery in 14 standard vinyl colors.

- 4074 Bariatric Treatment Table, 78"L x 34"W x 30"H

HAUSMANN HI-LINE TREATMENT TABLES

All purpose medical table for treatment, examination or therapy. Sturdy panel legs with durable laminate finish. Folkstone Gray edges on legs and shelf. 2" high-density urethane foam top. Resist-All™ vinyl upholstery finish.

- 4206 Hi-Line Treatment Table w/storage shelf, 72"L x 27"W x 31"H
4246 Hi-Line Treatment Table w/4 sliding laminate doors (2 on each side), 72"L x 27"W x 31"H



HAUSMANN HI-LO WHEELCHAIR ACCESSIBLE TILT TABLE

Medium gray vinyl upholstery. Comfortable 74"L x 24"W wide cushioned top. Equipped with three extra-wide safety straps with hook and loop locks. Electric hi-lo lift mechanism adjusts from 19" low wheelchair height to 34". Tilts from horizontal to 90° to fixed rate of speed of 4 degrees per second. Removable rubber tread metal foot board. Convenient hand control on coiled cord operates height and tilt adjustments. Includes four 3" diameter locking casters. Heavy duty steel frame is powder coated in neutral cream finish. 350 lb weight capacity.

- 6058 Hi-Lo Wheelchair Accessible Tilt Table

Accessories

- 6027-58 Patient Activity Tray For 6058



HAUSMANN HORSE-SHOE THERAPY TABLE

Accommodates up to four wheelchair patients. Durable, cloud wash formica top. Seven adjustable height metal legs. Hand control knob on each leg.

- 6674 72"L x 48"W x 26"-34"H Therapy Table



4300



4320

HAUSMANN HYDRAULIC WORK TABLES

Convenient foot pedal-activated hydraulic lift adjust table height from 29" to 44". Natural oak formica top. Oak laminate lift enclosure and base plate with four leveling legs.

- 4300 60"L x 36"W Table
4320 66"L x 48"W Table



HAUSMANN MANUAL TILT TABLE

Operates with a threaded worm gear mechanism that ensures automatic stop and locks in any position from 0 to 90°. Heavy duty steel frame, 4" locking casters. Removable footboard; restrainer straps. Hand crank folds flat when not in use. 78"L x 28"W x 33"H. Gray upholstery.

- 6020-709 Tilt Table



HAUSMANN ORTHOPEDIC/CAST TABLE

Choice of standard vinyl or nylon-reinforced vinyl top. Urethane foam top. Three fully encased door storage cabinet sections with adjustable shelf. 400 lb capacity. Available with 99L Left Hand Table option.

- 4073 78"L x 30"W x 31"H Table



HAUSMANN ORTHOPEDIC/HAND SURGERY TABLE

All-laminate drawer on nylon rollers. Two fully enclosed door cabinets. Slide out 24" wide shelf in center section. 2" thick urethane foam top. Vinyl upholstery. Black formica covered plywood base. 400 lb capacity. Available with 99L Left Hand Table option.

4839 78"L x 30"W x 31"H Table



HAUSMANN POWERMATIC® ECHO-SCAN™ TABLE

3" thick urethane foam top. Top has 10"W x 12"D drop-away section on durable piano hinge. Sliding latch activates drop-away action. Two direct-drive synchronized stain aluminum life mechanisms. Power height adjustment from 20" wheelchair height to 29"H. Hands-free foot control. Hospital grade plus. 100V, 60Hz, 2.8 amps.

4790 Echo-Scan™ Table, 76"L x 27"W x 20"-29"H

Accessories

074 Four 4" Locking Casters For 4790 Only



HAUSMANN POWERMATIC® WORK TABLE

Electrical power height adjustment from 27" to 39". Hand control. Formica top 66" x 48". Black bumper Tee-molding around top. Rounded corners. Black laminate lift enclosure and sturdy chrome-trimmed base plate with four leveling legs. Limit switches. Thermal overload. Hospital grade plug.

4380 Work Table



HAUSMANN REVERSABLE™ ACTIVITY/HAND THERAPY TABLE

360° rotating insert with pegboard. Can be locked in 2 positions. Pegboard has 100 peg holes with plastic pegs in assorted colors. 1 1/2" Cloud Wash Formica top with black edges. 26" to 34" adjustable height metal legs.

6272 Quarter Round Design Therapy Table, 94" Front Radius x 28"D,
(2) 16" x 16" Pegboards



HAUSMANN ROUND HYDRAULIC WORK TABLE

Hydraulic foot activated mechanism raises height from 27 1/2" to 34"H. 48" diameter 1 1/8" oak formica top swivels 360° and locks in place.

4335 Round Work Table



HAUSMANN SECTIONAL ELECTRIC TREATMENT TABLES

Gray vinyl upholstery. 2" thick urethane foam. Breathing slot in head section. Four 3" locking casters (except 6050). Heavy duty steel frame is powder coated in neutral crème finish. 350 lb capacity. High/Low adjustment.

6050-709 76"L x 26"W x 22"-40"H, 3-Section Table, Manual Headrest, Pneumatic Foot Control

Accessories

072 Four 3" Locking Casters, For 6050-709, 4/set

HAUSMANN SPLIT LEG LIFT TABLE

Designed for leg lift exercises for training and rehabilitation. Each padded leg rest can be adjusted and locked into desired angle up to 45°. Open shelf storage, accessible from both sides of the table. Natural Oak laminate. Upholstery comes in 14 standard vinyl colors.

4065 Split Leg Lift Table, 78"L x 30"W x 31"H

Tables

**HAUSMANN STAND-IN HYDRAULIC TABLE**

Features include parallel bars, non-slip ramp for easy access, and an oversized backrest that slides and locks into place. Once standing, adjustable padded back and knee supports hold patient comfortably and securely. Front mounted hydraulic control raises and lowers the 31" x 41" work surface from 41" to 58" high.

6100 Stand-In Hydraulic Table

**HAUSMANN THERAPY TABLE**

Height of spring-powered table can be easily adjusted with dual hand control levers. 14" height range allows easy use sitting, standing or in wheelchair. Folkstone Gray Formica® 3/4" tops. Black epoxy heavy-duty dual steel uprights and U base.

- 6282 Hand Therapy Table w/4 leveling legs, 32"L x 18"D x 26 1/2"-40 1/2"H
 6283 Hand Therapy Table w/4 leveling legs & drawer, 32"L x 18"D x 26 1/2"-40 1/2"H
 6284 Mobile Therapy Table w/four 2" twin-wheel nylon casters, 32"L x 24"D x 28 1/2" - 42 1/2" w/19"W x 6"D cutout

Accessories

- 075 Set of (4) 2" Nylon Locking Casters

**HAUSMANN WORK TABLES**

- 4331-048 48"L x 30"W x 26"-34"H Table, No Apron, Cloud Wash Formica Top, 4 Adjustable Metal Legs, Hand Control Knob
 4337-048 48"L x 36"W x 26"-34"H Table, No Apron, Cloud Wash Formica Top, 4 Adjustable Metal Legs, Hand Control Knob
 4337-060 60"L x 36"W x 26"-34"H Table, No Apron, Cloud Wash Formica Top, 4 Adjustable Metal Legs, Hand Control Knob

**LUMEX THREE POSITION RECLINERS**

Blow-molded table either slide on and off or can be locked in one of five positions. Three easy-to-attain positions: sitting, elevated leg or full recline. Contoured seat back and articulating headrest. Zinc plated caster horns for improved rust resistance.

- | | |
|-----------|---|
| 574G406 | Pecan Recliner |
| 574G407 | Terra Cotta Recliner |
| 574G413 | Dove Recliner |
| 574G427 | Blue Ridge Recliner |
| 574G454 | Royal Blue Recliner |
| 574G851 | Doe Skin Recliner |
| 574G857 | Jade Recliner |
| 574G858 | Grey Taupe Recliner |
| 574G863 | Rosewood Recliner |
| FR574G427 | Blue Ridge Recliner, Fire Rated Version (CA133) |
| FR574G454 | Royal Blue Recliner, Fire Rated Version (CA133) |
| FR574G851 | Doe Skin Recliner, Fire Rated Version (CA133) |
| FR574G857 | Jade Recliner, Fire Rated Version (CA133) |
| FR574G858 | Grey Taupe Recliner, Fire Rated Version (CA133) |
| FR574G863 | Rosewood Recliner, Fire Rated Version (CA133) |

MJM SHOWER GURNEYS AND STRETCHERS "900" SERIES

- 910 Shower Gurney w/3-position elevating headrest, Close Cell Waterproof Foam Pad, Drain Pan, Drain Hose, 6-5" x 1 1/4" Heavy Duty Threaded Stem Casters
 910-DB Bariatric Shower Gurney w/three position elevating headrest, Close cell waterproof foam pad, Drain pan and Drain hose / 900lb capacity


**A complete line of
stretchers is available.
Contact us for
information.**




3713-16



3710-16



3600 Series



3620 Series

CLINTON RECOVERY COUCHS

2" ultra-firm foam padding. 1" thick, solid plywood, frame construction. All cabinets are assembled using steel fasteners for strength and stability. Easy clean laminate models or classic chrome or wood leg models. Color coordinated laminate and vinyl. 300 lbs capacity. Premium, stain-resistant, woven, knit-backed vinyl upholstery. Three-year limited warranty.

- 3600-24 Chrome Leg Couch w/non-adjustable wedge headrest, 72"L x 18"H x 24"W
- 3600-27 Chrome Leg Couch w/non-adjustable wedge headrest, 72"L x 18"H x 27"W
- 3610-24 Panel Leg Couch w/non-adjustable wedge headrest, 72"L x 18"H x 24"W
- 3610-27 Panel Leg Couch w/non-adjustable wedge headrest, 72"L x 18"H x 27"W
- 3620-24 Hardwood Leg Couch w/non-adjustable wedge headrest, 72"L x 18"H x 24"W
- 3620-27 Hardwood Leg Couch w/non-adjustable wedge headrest, 72"L x 18"H x 27"W
- 3710-10 Panel Leg Couch w/full shelf & non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3710-15 Panel Leg Couch w/full shelf & adjustable flat foam headrest, 72"L x 20"H x 27"W
- 3710-16 Panel Leg Couch w/full shelf & adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3713-10 Upholstered Apron Couch w/double drawer storage & non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3713-15 Upholstered Apron Couch w/double drawer storage & adjustable flat foam headrest, 72"L x 20"H x 27"W
- 3713-16 Upholstered Apron Couch w/double drawer storage & adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3770-10 Upholstered Apron Couch w/sliding doors & non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3770-15 Upholstered Apron Couch w/sliding doors & adjustable flat foam headrest, 72"L x 20"H x 27"W
- 3770-16 Upholstered Apron Couch w/sliding doors & adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3790-10 Solid Base Couch w/non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3790-15 Solid Base Couch w/adjustable flat foam headrest, 72"L x 20"H x 27"W
- 3790-16 Solid Base Couch w/adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3810-10 Panel Leg Couch w/non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3810-15 Panel Leg Couch w/adjustable headrest, 72"L x 20"H x 27"W
- 3810-16 Panel Leg Couch w/adjustable pillow wedge, 72"L x 20"H x 27"W
- 3813-10 Panel Leg Couch w/double drawer storage & non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3813-15 Panel Leg Couch w/double drawer storage & adjustable headrest, 72"L x 20"H x 27"W
- 3813-16 Panel Leg Couch w/double drawer storage & adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3870-10 Panel Leg Couch w/sliding doors & non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3870-15 Panel Leg Couch w/sliding doors & adjustable headrest, 72"L x 20"H x 27"W
- 3870-16 Panel Leg Couch w/sliding doors & adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3890-10 Solid Base Couch w/non-adjustable pillow wedge headrest, 72"L x 20"H x 27"W
- 3890-15 Solid Base Couch w/adjustable headrest, 72"L x 20"H x 27"W
- 3890-16 Solid Base Couch w/adjustable pillow wedge, 72"L x 20"H x 27"W

Accessories

- 016 Adjustable Wedge Headrest



7010



7030

HAUSMANN ECONO-LINE™ RECOVERY COUCH

Sturdy end panel legs and apron. 2" thick urethane foam vinyl top. Folkstone Gray laminate finish.

- 7010 72" x 27" x 22" Couch
- 7030 72" x 27" x 22" Couch w/cabinet & 2 sliding doors

HAUSMANN FLEETWOOD COUCH

2" thick urethane foam vinyl top. Ten position, 0-45° adjustable headrest. Natural Oak laminate legs. Upholstery comes in 14 standard vinyl colors.

- 7155 Fleetwood Couch, 72"L x 27"W x 18"H

HAUSMANN INWOOD COUCH

3" thick urethane foam vinyl top. Ten position, 0-45° adjustable headrest. Natural Oak laminate legs. Upholstery comes in 14 standard vinyl colors.

- 7905 Inwood Couch, 72"L x 27"W x 18"H

CLINTON MANAGED CARE QUICK CABINETS

Pre-assembled units save time. Easy-clean all laminate construction. Tough, chip resistant PVC edges. Assembled with steel fasteners. Four adjustable leveling feet on the base units. Three-way adjustable 110° hinges. Euro-style drawer fronts with metal sides that hold up to 75 lbs. Recessed toe kick. Easy hook and rail wall installation. Edgebanding on all edges. Door/drawer pulls available in three colors. Large color selection/match or complement. Easily removes for relocation. Drawers can be added to any base unit. Doors and drawer can be moved to accommodate sinks.

- 8024 Single Base Cabinet w/2 doors & 1 drawer, 24"L x 35"H x 18"D
- 8036 Base Cabinet w/2 doors & 2 drawers, 36"L x 35"H x 18"D
- 8042 Base Cabinet w/3 doors & 2 drawers, 42"L x 35"H x 18"D
- 8048 Base Cabinet w/4 doors & 2 drawers, 48"L x 35"H x 18"D
- 8054 Base Cabinet w/3 doors & 3 drawers, 54"L x 35"H x 18"D
- 8066 Base Cabinet w/5 doors & 3 drawers, 66"L x 35"H x 18"D
- 8072 Base Cabinet w/6 doors & 3 drawers, 72"L x 35"H x 18"D
- 8212 Short Wall Cabinet, 24"L x 12"H x 12"D
- 8224 Single Wall Cabinet w/2 doors, 24"L x 24"H x 12"D
- 8236 Wall Cabinet w/2 doors, 36"L x 24"H x 12"D
- 8242 Wall Cabinet w/2 doors, 42"L x 24"H x 12"D
- 8248 Wall Cabinet w/2 doors, 48"L x 24"H x 12"D
- 8254 Wall Cabinet w/3 doors, 54"L x 24"H x 12"D
- 8266 Wall Cabinet w/3 doors, 66"L x 24"H x 12"D
- 8272 Wall Cabinet w/3 doors, 72"L x 24"H x 12"D

Accessories

- 011 Stainless Steel Sink & Faucet
- 033 Waste Drop Door
- 044 Add a Drawer
- 055 Door & Drawer Locks
- 066 Door Lock & Inside Latch Combo
- 076 Wall Mount Desk w/1 leg & 1 drawer, 24"L x 30"H x 18"D

Cabinets

CLINTON THEME SERIES QUICK CABINETS

Pre-assembled units save time. Easy-clean all laminate construction. All graphics are permanently embedded in the high pressure laminate and will not peel off. Tough, clip-resistant PVC edges. Assembled with steel fasteners. Four adjustable leveling feet on each unit. Three-way adjustable 110° hinges. Euro-style drawer fronts with metal sides that hold up to 75 lbs. Recessed toe kick. Easy hook and rail wall installation. Slate Gray edgebanding on all edges. Black chrome door/drawer pulls. Easily removes for relocation. Drawers can be added to any base unit. Doors and drawer can be moved to accommodate sinks.

- 6120-B School House Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6120-W School House Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6130-B Firehouse Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6130-W Firehouse Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6131-B Cool Pals Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6131-W Cool Pals Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6132-B Rainforest Follies Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6132-W Rainforest Follies Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6134-B Dino Days Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6134-W Dino Days Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6135-B Space Place Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6135-W Space Place Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6136-B Ocean Commotion Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6136-W Ocean Commotion Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D
- 6137-B Alley Cats & Dogs Base Cabinet, 4 Doors w/adjustable shelves & 3 drawers, Shipped in 1 Unit, 61"L x 35"H x 18"D
- 6137-W Alley Cats & Dogs Wall Cabinet, 4 Doors w/adjustable shelves, Shipped in 1 Unit, 61"L x 24"H x 12"D

CLINTON TREATMENT CABINETS

All laminate construction. Euro-style drawers with metal sides.

- 8805 Treatment Cabinet w/5 drawers, 24"L x 35"H x 18"D
- 8820 Treatment Cabinet w/2 doors
- 8821 Treatment Cabinet w/2 doors & 1 drawer, 24"L x 35"H x 18"D
- 8822 Treatment Cabinet w/2 doors & 2 drawers, 24"L x 35"H x 18"D
- 8900 Mobile Equipment Cabinet w/1 flip-down door, 18 1/2"L x 30"H x 19 1/2"D, Mounted on Dual Wheel Swivel Casters
- 8910 Mobile Cabinet w/1 drawer & open storage compartment, 18 1/2"L x 30"H x 19 1/2"D, Mounted on Dual Wheel Swivel Casters
- 8921 Mobile Treatment Cabinet w/2 doors & 1 drawer, 24"L x 34"H x 18"D, Mounted on Dual Wheel Swivel Casters
- 8922 Mobile Treatment Cabinet w/2 doors & 2 drawers, 24"L x 34"H x 18"D, Mounted on Dual Wheel Swivel Casters
- 8930 Mobile Cabinet w/3 drawers, 18 1/2"L x 30"H x 19 1/2"D, Mounted on Dual Wheel Swivel Casters
- 8950 Mobile Treatment Cabinet w/5 drawers, 24"L x 34"H x 18"D, Mounted on Dual Wheel Swivel Casters

Accessories

- 089 Optional Wire Basket, 17 1/2"L x 5"H x 6 3/4"D

**HAUSMANN ACCENT™ CABINETS**

Attractive gray cabinets are available with matching handles or with handles and trim in a choice of five accent colors/pull colors (966 Burgundy/Burgundy pulls, 969 Navy Blue/Navy Blue pulls, 318 Shadow Blue/Gray pulls, 956 Victorian Teal/Gray pulls or 927 Folkstone Gray/Gray pulls). Four lower cabinet options and two types of upper cabinets provide high efficiency and versatility. A coordinated desk (see GLR-DK) module may be attached to either end of a base cabinet or fastened to the wall (optional). Cabinets can be free-standing or wall-attached. End panels and counter tops have square front edges in gray laminate. Available in 16 different combinations. Sinks are standard on some units. Locks are available. Specify left or right hand - optional sink on order.

- GLR-A1 All laminate free-standing cabinet system 50"W x 77"H. 18"D base section, 2 Drawer/Door Cabinet w/sink cabinet, stainless steel sink, faucet & towel drop, Upper section is 10"D w/4 framed plexiglas doors
- GLR-B4 All laminate free standing cabinet system 50"W x 77"H. 18"D base section w/2 doors & 5 drawers w/o sink cabinet, on nylon rollers, Upper Section is 10"D w/4 solid doors

HAUSMANN MOBILE MEDICAL CABINET

High pressure laminate top slides to left and to right to reveal 6" deep storage well. Three drawers with full width black pulls. Easy rolling twin-wheel casters. Available in choice of any high pressure laminate color or woodgrain.

- GM-18-38 18"W x 32"H x 18"D Medical Cabinet

**HAUSMANN PRO-LINE™ PROFESSIONAL CABINETS**

Modular laminate, dowel and glue cabinetry. Neutral color interiors. Square-edged counter tops. Satin chrome wire pulls or nylon wire pulls in selected colors. Self-closing laminate drawer slides with nylon rollers and positive stops. Doors have concealed, self-closing 3-way adjustable hinges. Gray laminate. Choice of counter top colors. Other models available.

- B-18-5D 18"W x 36"H x 18 1/2"D Base Cabinet w/5 drawers
- B-18-L 18"W x 36"H x 18 1/2"D Base Cabinet, 1 Door Hinged Left
- B-18-R 18"W x 36"H x 18 1/2"D Base Cabinet, 1 Door Hinged Right
- B-24 24"W x 18"D x 36"H Cabinet, 2 Doors
- B-24-2D 24"W x 18"D x 36"H Cabinet, 2 Doors & 2 Drawers
- B-24-5D 24"W x 36"H x 18 1/2"D Base Cabinet w/5 drawers
- B-401 Proline Cabinet Combo Consists of: B-24-2D & B-24
- B-402 Proline Cabinet Combo Consists of: B-18-L, B-18-2D, W-18-24-L & W-18-24-R
- B-403 Proline Cabinet Combo Consists of: B-24, B-24-2D & W-24-24
- B-404 Proline Cabinet Combo Consists of: B-24, B-24-1D, D-24, D-EP & W-24-24
- W-18-24-L 18"W x 24"H x 12"D Wall Cabinet, 1 Door, Hinged Left
- W-18-24-R 18"W x 24"H x 12"D Wall Cabinet, 1 Door, Hinged Right
- W-24-24 24"W x 24"H x 12"D Wall Cabinet, 2 Doors


Accessories

- LOCK Lock
- LOCK-BOLT Lock & Bolt
- S Stainless Steel Sink w/faucet

HAUSMANN TREATMENT CABINETS

One laminate drawer on nylon rollers. All laminate interiors with one shelf. Cloud wash formica top.

- 8728-346 30"L x 16"W x 30"H Cabinet, Natural Oak
8728-927 30"L x 16"W x 30"H Cabinet, Folkstone Gray



344-005

344-006

RITTER DRESSING NOOK

When opened, it's a mirrored dressing area with storage for personal items and clothing; closed, it appears as a color-coordinated wall panel. 24"W x 7"D x 63"H. Available panel colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

- 344-005 RH Door Dressing Nook
344-006 LH Door Dressing Nook



340-003

341-005

347-005

RITTER TREATMENT CABINETS

Keeps equipment and disposables at your fingertips. Specify door and drawer front color. Available colors in: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

- 340-003 24"W x 18"D x 36"H Cabinet, One 4" Drawer, 2 Doors w/shelf, Sand Grey Kydex Top
340-004 24"W x 18"D x 36"H Cabinet, One 4" Drawer, 2 Doors w/shelf, Laminate Top (specify color)
341-005 24"W x 18"D x 36"H Cabinet, Two 4" Drawers, 2 Doors, Sand Grey Kydex Top
341-006 24"W x 18"D x 36"H Cabinet, Two 4" Drawers, 2 Doors, Laminate Top (specify color)
341-007 24"W x 24"D x 36"H Cabinet, Two 4" Drawers, 2 Doors, Sand Grey Kydex Top
341-008 24"W x 24"D x 36"H Cabinet, Two 4" Drawers, 2 Doors, Laminate Top (specify color)
347-005 24"W x 18"D x 36"H Cabinet, Five 4" Drawers, Sand Grey Kydex Top
347-006 24"W x 18"D x 36"H Cabinet, Five 4" Drawers, Laminate Top (specify color)
347-007 24"W x 24"D x 36"H Cabinet, Five 4" Drawers, Sand Grey Kydex Top
347-008 24"W x 24"D x 36"H Cabinet, Five 4" Drawers, Laminate Top (specify color)

**HAUSMANN MOBILE CABINET FOR SPLINTING & SUPPLIES**

Gray laminate finish on all surfaces. Heavy-duty twin wheel swiveling locking casters. Grab-bar handle on left side. Locking door cabinet on right side. Drawer with rounded inside corners and glides.

- 6695 42"L x 24"D x 31"H Cabinet



85161H

8550LG

KENDALL IN-ROOM SYSTEM WALL ENCLOSURES AND GLOVE BOXES

Sturdy, locked sharps container mounting allows for secure point of use disposal of sharps. For convenience, all cabinets are universally keyed. Sturdy, durable, long lasting ABS plastic provides a secondary barrier of protection. Easily monitored fill level through viewing windows on front and sides. Keep gloves readily available with optional glove dispenser. Beige, unobtrusive style fits any room décor.

- 85161H Wall Enclosure For 5qt Container, 1/cs
85165H 5qt Enclosure w/window, 5/cs
8516HDL 85161H Replacement Door w/window & lock, 1/cs
85301H 3 gal, Enclosure w/window, 1/cs
85401H 3 & 4 gal, Low Profile Wall Enclosure, 1/cs
8546 Glove Box Prevention Adapter, 5/cs
8550B Glove Boxes w/brackets, 5/cs
8550LG Glove Dispenser, Top Loading, Attaches to 5 & 12 qt Wall Enclosures, 6.75"H x 4.5"D x 12"W, 10/cs
8552B Accessory Tray w/bracket, 5/cs
8553B Accessory Tray, 10/cs
8554B Double-Up Wall Bracket, 5/cs
8555SA Multi-Glove Dispenser, Opaque Beige, 10/cs

**OMNIMED BEAM® PASS-THRU CABINET**

Designed for all in wall pass-thru convenience. Each door is hinged with spring closure. Adjustable trim flange allows mounting to wall thickness of 3 1/2" - 5 1/2". Specify key lock or thumb lock. Constructed of 22 gauge stainless steel.

- 181785 Pass-Thru Cabinet, 12"H x 11 1/2"W x 6"D



8154-346

8248-84-346

HAUSMANN TALL STORAGE CABINETS

All laminate dowel and glue wood construction. Easy to clean laminate on all exterior and interior surfaces. Laminate shelves. Durable, 5 knuckle institutional hinges. Satin chrome wire pulls.

- 8154-346 Storage Cabinet, 1 Door, 5 Adjustable Shelves, 24"W x 23"D x 76"H, Natural Oak
8154-927 Storage Cabinet, 1 Door, 5 Adjustable Shelves, 24"W x 23"D x 76"H, Folkstone Gray
8248-84-346 Storage Cabinet, 2 Doors, 8 Adjustable Shelves, 48"W x 23"D x 76"H, Natural Oak, Includes 4 Drawers
8248-84-927 Storage Cabinet, 2 Doors, 8 Adjustable Shelves, 48"W x 23"D x 76"H, Folkstone Gray, Includes 4 Drawers

Cabinets



3826MAW



3902MAW

LAKESIDE HOMED MEDICAL STORAGE CABINETS

- 3824MAW Cabinet, 2 Shelves - 1 Fixed & 1 Adjustable, 12"D x 30"W x 30"H, 2 Clear Acrylic Windows, Cut Key Lock, Autumn White
- 3826MAW Cabinet, 2 Shelves - 1 Fixed & 1 Adjustable, 12"D x 24"W x 26"H, 2 Clear Acrylic Windows, Cut Key Lock, Autumn White
- 3902MAW Cabinet, 3 Shelves - 1 Fixed & 2 Adjustable, 30"H x 30"W x 10"D, 22 Gauge, Spring Clip Door Latch, Wall-Mount Option, High Security Cam Lock, Autumn White
- 6030MAW Cabinet, 3 Shelves - 1 Fixed & 2 Adjustable, 30"W x 18"D x 60"H, Peg Holes on Doors, Inside Hangers, Securely Lockable, Baked Enamel Finish, Autumn White
- 6035MAW Cabinet, 4 Adjustable Shelves, 30"W x 18"D x 60"H, Peg Holes on Doors, Inside Hangers, Securely Lockable, Baked Enamel Finish, Autumn White

OMNIMED BEAM® WALL STORAGE CABINETS

Store all regularly used supplied or equipment at your fingertips. Designed for easy mounting on most walls, doors or vertical surfaces. Smaller sizes can be mounted horizontally in drawers, under shelves or inside other cabinetry. The ambi design allows you the convenience of mounting cabinets for left or right opening. Security can be provided with a variety of key lock styles or allow unlimited access with latch cabinet closure and handle. Constructed of heavy gauge coated steel, and are extremely strong and durable, yet lightweight. Easy to clean and requires no maintenance for years of trouble free use. Specify color: beige or gray.

- 291619 Wall Storage Cabinet, 26 3/4"H x 16"W" x 4"D w/2 shelves

Accessories

- 291618-1 Shelf, 16"W x 4"D w/4 brackets
- 291618-2 Shelf, 16"W x 8"D w/4 brackets

HARLOFF NARCOTICS CABINET

Two high security pick-resistant tubular locks for regulatory compliance. Narcotics box is delivered ready to use (door handle must be customer installed). Powder coat finish (13 colors available 2710/2720/2730). Also available in stainless steel (2712/2722/2732). 20 gauge unibody steel construction. Stainless steel pinned door hinge.

- 2710 Narcotics Box, Compact, Single Door, Double Lock, ea
- 2712 Stainless Steel Compact Narcotics Box, Single Door w/dual tubular locks separately keyed, 10"H x 7"W x 7"D
- 2720 Medium Narcotics Box, Double Door w/dual tubular locks separately keyed, 16"H x 12"W x 9"D
- 2722 Stainless Steel Medium Narcotics Box, Double Door w/dual tubular locks separately keyed, 16"H x 12"W x 9"D
- 2730 Large Narcotics Box, Single Door w/dual tubular locks separately keyed, 29.5"H x 23.5"W x 10.5"D
- 2732 Stainless Steel Large Narcotics Box, Single Door w/dual tubular locks separately keyed, 29.5"H x 23.5"W x 10.5"D
- 2750 Dual Lock Narcotic Cabinet, Large, for Gastroenterology
- 40401 Additional/Replacement Key - Key number must be specified, ea
- AL75PC Desk/Wall Mount Cabinet for Punchcards including 3 Exchange baskets, capacity of 75 punchcards, ea

Accessories

- 40400RPL Replacement Lock Kit, Specify Old Lock & Key#, Includes 1 Key For The Existing Lock, 1 New Lock & 2 New Keys



3801MAW



3802MAW



3812MAW



3822MAW



3901MAW



3901MAW-2

LAKESIDE HOMED NARCOTIC CABINETS

- 3801MAW Single Door/Single Lock Wall Cabinet, 1 Adjustable Shelf, 3 1/2" D x 10"H x 14"W, 16 Gauge, High Security Cam Lock, Autumn White
- 3802MAW Single Door/Double Lock Box, 7 1/4"H x 8 3/4"W x 6"D, 20 Gauge, Two High Security Cam Lock, Autumn White
- 3812MAW Single Door/Double Lock Cabinet, 1 Fixed Shelf, 14"W x 8"D x 12"H, 20 Gauge, One Position Key Release, Autumn White
- 3821MAW Double Door/Double Lock Cabinet, 2 Adjustable Shelves, 18"W x 8"D x 30"H, Full Length Piano Hinges, Two High Security Cam Lock, Autumn White
- 3822MAW Double Door/Double Lock Cabinet, 1 Adjustable Shelf, 12"W x 8"D x 15"H, 20 Gauge, Two High Security Cam Lock, Autumn White
- 3901MAW Single Door/Single Lock Cabinet, 3 Shelves - 1 Fixed & 2 Adjustable, 24"W x 30"H x 10"D, 18 Gauge, Mounting Hardware Included, 1 High Security Cam Lock, Autumn White
- 3901MAW-2 Single Door/Double Lock Cabinet, 2 Shelves - 1 Fixed & 1 Adjustable, 24"W x 10"D x 30"H, 18 Gauge, Mounting Hardware Included, 2 High Security Cam Lock, Autumn White

**OMNIMED BEAM® CORNER NARCOTIC CABINETS**

Creates additional storage space in an otherwise "lost" area. Four aluminum, adjustable and removable shelves and can be right or left hinged. Cabinet constructed of 22 gauge stainless steel. Piano hinge suspension.

- 181765 Corner Narcotic Cabinet, Single Door w/double lock, 24"H x 24"W x 8"D



OMNIMED BEAM® DOUBLE DOOR NARCOTIC CABINET ACCESSORIES

181601-1	Replacement Shelf w/4 brackets For #181601, 4"D
181602	Replacement Lock Kit For #181601, 181611, 181651
181603	Replacement Key/Wafer Key, Min 2 Keys
181604	Replacement Wafer Lock & 3 Keys
181605	Replacement Push Button Combination Lock
181681-1	Replacement Shelf For #181681, Large w/4 brackets
181682	Replacement Lock Kit For #181681
181683-1	Maxi Shelf For #181683 w/4 brackets
181691	Replacement Lock w/3 Flat Keys
181691-1	Replacement Flat Key Blanks



181501



181601



181680

OMNIMED BEAM® DOUBLE DOOR NARCOTIC CABINETS

Double door concept makes each cabinet a controlled system for storing and dispensing small but adequate quantities of narcotics from the limited space at nurse's station. Double door security. Two key/two person safety system. Triple bolt wafer lock on outer door. Inner door lock has flat key for easy identification. Wafer lock is virtually pick proof for extra security and protection. Easily mounts on wall or suitable for ambulance use. Meets NY State Schedule 2 Drug Storage requirements. Complies with Dept. of Justice Drug Enforcement and Joint Commission requirements. Check your state codes. Welded corners, easy to clean. Shipped complete with locks and 3 keys for each lock. Cabinet constructed of 22 gauge stainless steel. Shelves are aluminum, adjustable and removable. Piano hinge suspension. Can be right or left hinged.

181501	Mini Narcotic Cabinet, 9"H x 8"W x 5 5/8"D
181502	Narcotic Cabinet, 8"H x 16"W x 8"D
181601	Narcotic Cabinet w/2 shelves, 15"H x 11"W x 4"D
181650	Narcotic Cabinet, Push Button Keyless Combination Lock, 15"H x 11"W x 8"D
181651	Narcotic Cabinet w/2 shelves, 15"H x 11"W x 8"D
181680	Narcotic Cabinet, Push Button Keyless Combination Lock, 24"H x 16"W x 8"D
181681	Narcotic Cabinet w/4 shelves, 24"H x 16"W x 8"D
181683	Maxi Narcotic Cabinet, 30"H x 22"W x 12"D



182100



182150



182175

OMNIMED BEAM® ECONOMY NARCOTIC CABINETS

Made of heavy gauge beige steel and meets all regulatory compliance requirements. Features an ambi top for left or right hinge, double door - double lock security and adjustable shelves. Ideal for wall mount point of care convenience. Call for custom size and color considerations.

182100	Economy Mini Narcotic Cabinet, 9"H x 8"W x 5 5/8"D
182125	Economy Narcotic Cabinet w/2 shelves, 15"H x 11"W x 4"D
182150	Economy Narcotic Cabinet w/2 shelves, 15"H x 11"W x 8"D
182175	Economy Narcotic Cabinet w/2 shelves, 24"H x 16"W x 8"D



181401



181451

OMNIMED BEAM® SINGLE DOOR NARCOTIC CABINETS

Features two wafer locks spaced 3" apart to meet schedule two drug storage requirements. Each cabinet is supplied with 3 keys per lock and adjustable shelves. Cabinet is constructed of 22 gauge stainless steel. Shelves are aluminum, adjustable and removable. Piano hinge suspension. Can be left or right hinged.

181401	Narcotic Cabinet, 2 Shelves, 15"H x 11"W x 4"D
181451	Narcotic Cabinet, 2 Shelves, 15"H x 11"W x 8"D
181481	Narcotic Cabinet, 4 Shelves, 24"H x 16"W x 8"D



181801




181851

OMNIMED BEAM® TWIN NARCOTIC CABINETS

Assembly in 5 easy steps. Same specifications as our standard Narcotic Cabinets, but mounted side-by-side. Double door-double lock. Two lock/two key system. Includes a fascia and cap to cover top and center connectors. Single door models available.

181801	Twin Narcotic Cabinet, 15"H x 22"W x 4"D
181851	Twin Narcotic Cabinet, 15"H x 22"W x 8"D
181881	Twin Narcotic Cabinet, 24"H x 32"W x 8"D



AKRO ATTACHED LID CONTAINERS

Injection molded from high density polyethylene. Containers are completely modular and stack one on top of the other. Lids are attached so they can never be out of reach. Available in three different sizes. Nest for compact storage when they are not in use. Ergonomically designed molded-in handle grips. Anti-jamming ribs. Padlock eyes for security. Full length stainless steel hinge pins for durability. Textured bottoms increase container conveyability. Grey and Blue.

39085	8.5 Gallon Container, 21 1/2"L x 15"W x 9"H, 6/ctn
39120	12 Gallon Container, 21 1/2"L x 15"W x 12 1/2"H, 6/ctn
39170	17 Gallon Container, 21 1/2"L x 15"W x 17"H, 3/ctn
39160	16 1/2 Gallon Container, Grey Only
39175	17 1/2 Gallon Container, Grey Only
39280	28 1/2 Gallon Container, Grey Only

Accessories

RU843HR1420	Dolly For 39085, 39120 & 39170 on 3" Poly Casters
-------------	---

Cabinets

**AKRO BINS**

Polypropylene storage bins that hang or stack. Extra-wide stacking ledge assures secure stacking. Anti-slide stop prevents loaded bins from shifting forward when stacked. Large label slot for clear identification. Reinforcing ribs prevent bins from spreading when loaded. Finger grips allow easy handling. Available in red, yellow, blue or stone (specify when ordering).

30006	1 Sided, 12 Rail Assembly, Accommodated 96 each #30-210 or 30-220; 24 each 30-240; 12 each 30-250; 24 each 30-239
30008	1 Sided, 16 Rail Assembly, Accommodated 48 each #30-230; 24 each 30-235
30012	2 Sided, 12 Rail Assembly, Accommodated 192 each #30-210 or 30-220; 48 each 30-240; 24 each 30-250; 48 each 30-239
30016	2 Sided, 16 Rail Assembly, Accommodated 96 each #30-230; 48 each 30-235
30148	Steel Rail, 48"L x 3"H, Wall-Mounted Unit
30210	Bins, 5 3/8"L x 4 1/8"W x 3"H, 24/cs
30220	Bins, 7 3/8"L x 4 1/8"W x 3"H, 24/cs
30224	Bins, 10 7/8"L x 4 1/8"W x 4"H, 12/cs
30230	Bins, 10 7/8"L x 5 1/2"W x 5"H, 12/cs
30234	Bins, 14 3/4"L x 5 1/2"W x 5"H, 12/cs
30235	Bins, 10 7/8"L x 11"W x 5"H, 6/cs
30239	Bins, 10 3/4"L x 8 1/4"W x 7"H, 6/cs
30240	Bins, 14 3/4"L x 8 1/4"W x 7"H, 12/cs
30250	Bins, 14 3/4"L x 16 1/2"W x 7"H, 6/cs
30255	Bins, 10 7/8"L x 16 1/2"W x 5"H, 6/cs
30265	Bins, 18"L x 8 1/4"W x 9"H, 6/cs
30270	Bins, 18"L x 16 1/2"W x 11"H, 3/cs
30120ESD	ESD Bins, 11 5/8"L x 4 1/8"W x 4"H, 24/ctn
30124ESD	ESD Bins, 23 5/8"L x 4 1/8"W x 4"H, 12/ctn
30128ESD	ESD Bins, 17 7/8"L x 4 1/8"W x 4"H, 12/ctn
30210ESD	ESD Bins, 5 3/8"L x 4 1/8"W x 3"H, Black Only, 24/ctn
30220ESD	ESD Bins, 7 3/8"L x 4 1/8"W x 3"H, Black Only, 24/ctn
30230ESD	ESD Bins, 10 7/8"L x 5 1/2"W x 5"H, Black Only, 12/ctn
30235ESD	ESD Bins, 10 7/8"L x 11"W x 5"H, Black Only, 6/ctn
30240ESD	ESD Bins, 14 3/4"L x 8 1/2"W x 7"H, Black Only, 12/ctn
30280	Bins, 20"L x 12 3/8"W x 61"H, 4/cs
30281	Bins, 20"L x 12 3/8"W x 8"H, 3/cs
30282	Bins, 20"L x 12 3/8"W x 12"H, 2/cs
30283	Bins, 20"L x 18 3/8"W x 12"H
30290	Bins, 29 1/4"L x 18 3/8"W x 12"H

Accessories

RU842HR2018	Dolly For 30283 w/2" Polyo Casters, 20"L x 18 3/8"W x 3 1/2"H
RU842HR2918	Dolly For 30290 w/2" Polyo Casters, 29 1/4"L x 18 3/8"W x 3 1/2"H
RU843HR1818	Dolly For 30283 w/3" Polyo Casters
RU843HR1827	Dolly For 30290 w/3" Polyo Casters

AKRO BINS ACCESSORIES

30136	Louvered Panel, 35 3/4"L x 19"H
30211	Lids For #30-210, 24/ctn
30221	Lids For #30-220, 24/ctn
30231	Lids For #30-230, 12/ctn
30236	Lids For #30-235, 6/ctn
30299	Security Clips
30424	Mobile Kit - Casters & Handle
30618	Louvered Panel, 18"L x 19"H, 4/ctn
30626	Snap Fastener Louvered Rack
30636	Louvered Panel, 35 3/4"L x 19"H, 4/ctn
30653	Two Sided Louvered Rack
30812	Bin Cart, 38 1/2"L x 24"W x 36 1/2"H
40210	Dividers For #30-210, 6/pkg
40220	Dividers For #30-220, 6/pkg
40230	Dividers For #30-230 & 30-235, 6/pkg
40239	Dividers For #30-239, 6/pkg
40245	Dividers For #30-240 & 30-250, 6/pkg
40265	Dividers For #30-265, 6/pkg
40270	Dividers For #30-270, 6/pkg
98600	Bench Rack, 27 3/4"L x 8"W x 21"H
98636	Bench Rack, 36"L x 8"W x 21"H

AKRO LID STORAGE ORGANIZERS

Separate lid compartment holds small items. When lid is shut a unique lid design traps contents in place to keep items from gravitating to other compartments. The main compartment holds larger items. Carrying handle. Adjustable and removable dividers. Blue with clear lids and black latches.

06115	Lid Storage Organizer, 15"L x 11 5/16" x 3 1/4"H
06118	Lid Storage Organizer, 18 1/4"L x 13 3/8"W x 3 5/8"H
06217	Dual Lid Storage Organizer, Flat Bottom, 6/ctn

AKRO MODULAR CABINET ACCESSORIES

20228	Replacement Drawers For #19-228, 56/ctn
20320	Replacement Drawers For #19-320, 40/ctn
20416	Replacement Drawers For #19-416, 32/ctn
20715	Replacement Drawers For #19-715, 30/ctn
20909	Replacement Drawers For #19-909 & 19-109, 15/ctn
40228	Dividers For #19-228, 25/pkg
40320	Dividers For #19-320, 25/pkg
40416	Dividers For #19-416, 25/pkg
40715	Dividers For #19-715, 25/pkg
40909	Dividers For #19-909 & 19-109, 25/pkg

**AKRO MODULAR CABINETS**


Perfect for small part storage. Cabinets feature steel frames and drawer supports, a charcoal gray baked enamel finish, and five drawer sizes. Drawers are molded of gray styrene with transparent windows and a retaining tab to prevent spilling. Finger-grip drawer fronts for easy opening. Windowed drawers with molded-in label holders. Two dividers and one label furnished per drawer. Key-hole slots for easy mounting. Built-in nesting for safe stacking. 17"W x 11"H x 11"D.

19109	Cabinet w/9 drawers & locking door
19228	Cabinet w/28 drawers
19320	Cabinet w/20 drawers
19416	Cabinet w/16 drawers
19715	Cabinet w/15 drawers
19909	Cabinet w/9 drawers

AKRO NSBEES™ NEST & STACK BINS


Injection molded co-polymer polypropylene for economy, strength and durability. They are waterproof, resistant to grease, oil, solvents, chemicals and autoclavable up to 250°F. Reinforced side ribs. Built-in stacking posts and postholes. Generous hopper front gives easy access to stored items and the large label area allows quick content identification. Specify color when ordering: blue or yellow.

30716	Bins, 11 7/8" L x 6 5/8" W x 5" H, 10/ctn
30718	Bins, 11 7/8" L x 8 3/8" W x 5" H, 8/ctn
30776	Bins, 17 7/8" L x 6 5/8" W x 7" H, 10/ctn
30778	Bins, 17 7/8" L x 8 3/8" W x 7" H, 8/ctn
30796	Bins, 8 7/8" L x 6 5/8" W x 5" H, 10/ctn

**AKRO SHELF BINS**

Durable polypropylene shelf bins for storage of parts and components. Designed for use on 12", 18", and 24" deep shelving or vertical storage and retrieval units. Built-in hanglock catches underside of shelf above, allowing each bin to be tilted out for complete access. Reinforced edges for additional strength. Available in red, yellow or blue (specify when ordering).

30110	Shelf Bin, 11 5/8" L x 2 3/4" W x 4" H, 24/ctn
30120	Shelf Bin, 11 5/8" L x 4 1/8" W x 4" H, 24/ctn
30124	Shelf Bin, 23 5/8" L x 4 1/8" W x 4" H, 12/ctn
30128	Shelf Bin, 17 7/8" L x 4 1/8" W x 4" H, 12/ctn
30130	Shelf Bin, 11 5/8" L x 6 5/8" W x 4" H, 12/ctn
30138	Shelf Bin, 17 7/8" L x 6 5/8" W x 4" H, 12/ctn
30150	Shelf Bin, 11 5/8" L x 8 3/8" W x 4" H, 12/ctn
30158	Shelf Bin, 17 7/8" L x 8 3/8" W x 4" H, 12/ctn
30164	Shelf Bin, 23 5/8" L x 6 5/8" W x 4" H, 6/ctn
30170	Shelf Bin, 11 5/8" L x 11 1/8" W x 4" H, 12/ctn
30174	Shelf Bin, 23 5/8" L x 11 1/8" W x 4" H, 6/ctn
30178	Shelf Bin, 17 7/8" L x 11 1/8" W x 4" H, 12/ctn
30184	Shelf Bin, 23 5/8" L x 8 3/8" W x 4" H, 6/ctn

**AKRO SHELF BINS ACCESSORIES**

30101	Bin Cup, Small, 3 1/4" x 2" x 3", 48/ctn
30102	Bin Cup, Large, 5" x 2 3/4" x 3", 24/ctn
40120	Bin Divider For #30-120, 30-124 & 30-128, 24/pkg
40130	Bin Divider For #30-130, 30-138 & 30-164, 24/pkg
40150	Bin Divider For #30-150, 30-158 & 30-184, 24/pkg
40170	Bin Divider For #30-170, 30-174 & 30-178, 24/pkg
40410	Extended Label Holder, 10° Angle, 24/pkg
40445	Extended Label Holder, 45° Angle, 24/pkg

AKRO-NEST & STACK TOTE CONTAINERS

Injection molded from high density polyethylene. Resistant to solvents and chemicals. Won't rust, corrode or bend out of shape. Complies with FDA regulations for use in meat and poultry plants. Three colors: red, blue, gray. Optional lids and label holders available.

35180	Container, 18" L x 11" W x 6" H, 6/ctn
35190	Container, 19 1/2" L x 15 1/2" W x 10" H, 6/ctn
35200	Container, 19 1/2" L x 13 1/2" W x 8" H, 6/ctn
35230	Container, 23 1/2" L x 19 1/2" W x 13" H, 3/ctn
35240	Container, 23 1/2" L x 15 1/2" W x 12" H, 3/ctn
35300	Container, 29 1/2" L x 19 1/2" W x 15" H, 3/ctn
35185	Container, 18" L x 11" W x 9" H, 6/ctn
35195	Container, 19 1/2" L x 15 1/2" W x 13" H, 6/ctn
35225	Container, 23 1/2" L x 19 1/2" W x 10" H, 3/ctn

AKRO-STRAIGHT WALL CONTAINER SYSTEM

Available in mesh or solid walls and bases. Open mesh sides allow high visibility and good air flow. Mesh sided lighter than solid sided. Solid sided containers have smooth inner walls. Stackable. High density polyethylene. Complies with all FDA regulations.

37208	Mesh Side/Base Container, 15 3/4" L x 11 3/4" W x 8 1/4" H, 12/ctn
37672	Mesh Side/Solid Base Container, 23 3/4" L x 15 3/4" W x 12 1/4" H, 3/ctn
37278	Mesh Side/Solid Base Container, 15 3/4" L x 11 3/4" W x 8 1/4" H, 12/ctn
37288	Solid Side/Base Container, 15 3/4" L x 11 3/4" W x 8 1/4" H, 12/ctn
37608	Mesh Side/Base Container, 23 3/4" L x 15 3/4" W x 8 1/4" H, 4/ctn
37612	Mesh Side/Base Container, 23 3/4" L x 15 3/4" W x 12 1/4" H, 3/ctn
37616	Mesh Side/Base Container, 23 3/4" L x 15 3/4" W x 16 1/8" H, 2/ctn
37676	Mesh Side/Solid Base Container, 23 3/4" L x 15 3/4" W x 16 1/8" H, 2/ctn
37678	Mesh Side/Solid Base Container, 23 3/4" L x 15 3/4" W x 8 1/4" H, 4/ctn
37682	Solid Side/Base Container, 23 3/4" L x 15 3/4" W x 12 1/4" H, 3/ctn
37686	Solid Side/Base Container, 23 3/4" L x 15 3/4" W x 16 1/8" H, 2/ctn
37688	Solid Side/Base Container, 23 3/4" L x 15 3/4" W x 8 1/4" H, 4/ctn

Accessories

RU843HR1624	Dolly For 37608, 37678, 37688, 37612, 37682, 37616, 37676, 37686 on 3" Poly Casters
-------------	---

AKRO-TUB

Infection molded from polyethylene for durability and uniformity. Autoclavable up to 250°F. Complies with FDA regulations for use in meat and poultry plants. Rounded corners. Rolled top rim.


34240	Storage Bin, 24 1/2" L x 19" W x 9 1/2" H, 6/ctn
-------	--



HAUSMANN COMBINATION TREATMENT/WORK TABLE/DESK

Dual-purpose that is ideal for schools or clinics. 2" thick removable reversible urethane foam pad with nylon reinforced vinyl cover. Top pad has straps and can hang vertically when not in use. Folkstone Gray high pressure laminate sub-top with black tee-mold edging. Upholstery pad is blue on one side and gray on reverse side with 2" gray border.

4082 72"L x 24"W x 32"H Work Table/Desk



3810MAW-A



3810MAW-M

LAKESIDE HOMED WALL DESK

Heavy duty 16 gauge construction. Top mounted key lock. Single pocket across desk. Hardware is not included (due to variations in mounting surfaces).

3810MAW-A Automatic Close Wall Desk, 27"W x 18"H x 4"D, Writing Surface:

27"W x 18"D, Top Mounted Key Lock, Autumn White

3810MAW-M Manual Close Wall Desk, 27"W x 18"H x 4"D, Writing Surface:

27"W x 18"D, Top Mounted Key Lock, Autumn White



OMNIMED BEAM® DELUXE WALL DESK ACCESSORY

291550 Replacement Lock Assembly w/3 keys

OMNIMED BEAM® DRY ERASE WALL DESK KIT

The Dry Erase Wall Desk Kit is designed as an add-on accessory to our Elite, Deluxe and Slim Line Wall Desks. A magnetic dry erase pad adheres to the writing surface and a clear acrylic panel with a "living hinge" attaches on to the surface as well, to protect any messages or information from being smeared while charting etc. Each kit is supplied with an eraser and a black dry erase marker.

741304 Portable IV Pole

OMNIMED BEAM® ELITE WALL DESKS

Features solid wood frames with smooth rounded corners in natural oak or rich cherry finish. Increase patient/resident chart confidentiality and security with your choice of lock systems. Push button keyless entry lock provides endless combination options and the standard flat key lock offers durability and economy. Desk conveniently holds two charts and the large writing surface provides sample space for point-of-care record keeping or medication administration. Available with self-close and manual close door styles.

291533CHR	Elite Wall Desk Single Size Self Close Cherry
291533OAK	Elite Wall Desk Single Size Self Close Oak
291534HCR	Elite Wall Desk Single Size Manual Close Cherry
291534OAK	Elite Wall Desk Single Size Manual Close Oak
291535CHR	Elite Wall Desk Full Size Self Close Cherry
291535OAK	Elite Wall Desk Full Size Self Close
291536CHR	Elite Wall Desk Full Size Cherry
291536OAK	Elite Wall Desk Full Size
291537CHR	Elite Wall Desk Full Size Self Close
291537OAK	Elite Wall Desk Full Size Manual Close Oak
291538CHR	Elite Wall Desk Full Size Manual Close Cherry
291538OAK	Elite Wall Desk Full Size Manual Close Cherry
291539CHR	Elite Wall Desk Single Size Self Close
291539OAK	Elite Wall Desk Single Size Self Close
291540CHR	Elite Wall Desk Single Size Manual Close Cherry
291540OAK	Elite Wall Desk Single Size Manual Close Oak

OMNIMED BEAM® LAPTOP WALL DESK

291559-BG Laptop Wall Desk



OMNIMED BEAM® STANDARD WALL DESK

Writing surface easily pulls down for use. Features two large pocket compartments. Two clips conveniently mounted to organize memos and notes. 20 3/4"W x 15"H x 3"D.

291501 Wall Desk, Aluminum

291503 Wall Desk, Stainless Steel



OMNIMED BEAM® WALL POCKET

Can be mounted on wall, door, bed end or file cabinet. Holds charts, notes, records, computer printouts, magazines, etc. Easily attaches with the supplied pressure sensitive tape or screws (not included). Easy access for records to be filed or used. Tapered to 1" at bottom.

255713 Aluminum Wall Pocket, 14"W x 7"H

255715 Aluminum Wall Pocket, 17"W x 9 3/4"H



RITTER MOBILE DESK

24"x 18"D x 30"H. Two locking casters in front with two carpet casters in back. Available in standard colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

3740-002 Mobile Desk, 1 - 1 1/2" Drawer, Laminate (Specify Top Color & Drawer Panel Color)



WELCH ALLYN DIRECT FOCUSING HEADLIGHT

Lightweight and comfortable to wear. Focusing sleeve allows you to adjust the spot size according to the procedure. Light spot can adjust from 1/2" spot at 8" or 6 1/2" spot at 13".

46003 Headlight Complete w/headband & 6V transformer, 110-130V, 60 Hz

WELCH ALLYN DIRECT FOCUSING HEADLIGHT ACCESSORIES

460152 Lens Assembly Only (lamp not included)
 460167-501 Polyethylene Headband Only
 46080 Polyethylene Headband Only, w/lamp housing & quick-detach connector
 791003-502 Transformer Only



WELCH ALLYN HALOGEN HEADLIGHT

Allows practitioner to use both hands, while providing a focal source of illumination. Even light intensity across all spot sizes. Open ribbed construction prevents excess heat build-up.

49003 Complete headlight w/headband, wall transformer & mounting plate, 110-130V, 50/60 Hz
 490074 Complete Headlight w/o wall transformer

Accessories

460167-503 Headband Only
 790030-501 Transformer Only



BREWER HALOGEN 35 EXAM LIGHT

Dichroic filters out infrared (hot light). Color-corrected halogen light source. Cool housing. Reliable cool light output. 35 watt halogen bulb. Long bulb life (4000 hours). Crisp white light, reducing eye strain and fatigue. Low voltage. Made in the U.S.A.

16100 Halogen 35 Exam Light



BREWER INCANDESCENT BASIC EXAMINATION LIGHT

6 1/2" reflector system allows for reduced shadows during examinations. Our non-metallic head eliminates the dangers of a hot metal lampshade. The power cord is hidden to avoid entanglement. Chrome dipped bulb will protect your eyes from unwanted filament images causing pupil dilation and irritation. Manufactured in USA.

15100 Basic Examination Lamp



BREWER REPLACEMENT BULBS & CASTER BASES FOR LIGHTS

001525 Replacement Bulb For #15100
 001635 Replacement Bulb For #16100
 016150 Replacement Bulb For 163001, 163002, 163003, 163004
 100165 Caster Base For #16100 and #15100 Lamp

Lighting

BURTON SUPER BRIGHT SPOT™ EXAMINATION LIGHT

3,000 footcandles. Color temperature of 3200°K. Dichroic reflector filters out infrared to keep the beam cool. 20-watt halogen bulb with 3,000 hour life provides 5"-9" (13-23cm) spot size. Floorstand includes 24" rugged gooseneck and telescoping upright (with a range of 35"-62"). Unique space-saving base with 4 built-in casters.

0152042 Wall Mount Light w/extended spring arm
0152113 115V Light, Floorstand Model (casters included)

Accessories

0007000PK Set of 4 ESX Replacement Bulbs

**BURTON SUPER NOVA® EXAM LIGHT**

2,600 footcandles. 3100°K color temperature. 20-watt halogen bulb with 4,000 hour average bulb life; clear safety cover over bulb. 24" sturdy gooseneck for easy positioning. Cool housing temperature. Smooth surface for easy cleaning. Easy-grip handles for positioning. Easy bulb changes (from top) without tools. Unique space-saving base with 4 built-in casters. Superior arm system - will not sag.

0990020 115V Light, Floorstand Model (casters included)

Accessories

0008004PK Set of 4 Nova Replacement Bulbs

**CLINTON GOOSENECK LAMPS**

T-10 Lamp w/chrome finish

**RITTER 152 HALOGEN EXAM LIGHT**

A halogen light that uses a 35 watt quartz halogen bulb. It provides 1,500 footcandles at 24" for an optimum amount of light for general use without high cost. The 3,000K color temperature provides a cool white light for optimum color recognition. Uses dichroic coating to eliminate infrared light, which permits only cool light to the doctor's work area. Uses a unique heat sink, which effectively removes heat from the bulb and its socket. This feature increases the life of the bulb for minimal downtime. A small head section features high quality lighting without being burdensome. Listed by ETL to UL2601, the medical and dental equipment standard, and assures the customer of their safety. Designed in accordance to the strict aesthetic standards of the Ritter line of medical products. Cord is located at the lowest point of the light base, for optimum cord management. Heavy gauge, vinyl covered spring steel in gooseneck to prevent sagging.

152-001 Halogen Exam Light w/35 watt quartz halogen bulb



48210



52641

WELCH ALLYN HALOGEN EXAM LIGHT™ III ACCESSORIES

48200 Fiber Optic Light Pipe, 49 1/2"
48210 Fiber Optic Neonate Light Pipe For 48606
48605 Focusing Sleeve
48740 Exam Light III Light Box
48850 Mobile Stand w/locking casters for Exam Light III
48859 Universal Table Bracket For Exam Light III
52640 Disposable Sheaths, 52", 25/bx, 5 bx/cs
52641 Disposable Sheath Dispenser, 5/cs

SummitDoppler

Listening to Life™

....in Every Practice

LifeDop ABI for PAD (Peripheral Arterial Disease)

- Meets ABI reimbursement criteria
- Portable bi-directional doppler with compact printer
- Bi-directional probe
- Label paper
- Forms
- ABI video
- Examination guide
- Quick reference guide
- ABI glossary of terms



Obstetrical



Vascular



L150 Series: Non-Display Dopplers

- Superior sound quality
- Battery indicator light
- Interchangeable probes for both obstetrical and vascular applications
- Includes your choice of one probe: waterproof, 2 or 3 MHz for OB, and 4, 5, or 8 MHz for vascular



L250 Series: Digital Display Dopplers

- Display, which shows the heart rate, battery life and signal quality
- Superior sound quality
- Interchangeable probes for both obstetrical and vascular applications
- Includes your choice of one probe: waterproof, 2 or 3 MHz for OB, and 4, 5, or 8 MHz for vascular





Break Down The Barriers™



It's the Right Thing to Do



The Ritter 222 is one of a full line of Barrier-Free™ exam tables, offering the highest level of efficiency, safety, and comfort for both physicians and their patients.

"I took my mother to her doctor appointment today. It's hard watching her age and become more and more dependent. She had so much trouble just getting on the exam table... as she held to the doctor's arm with all her might, I could sense her insecurity and discomfort."

Prepare for an increase in patients, especially the elderly, who have difficulty getting onto an exam table. Equip your exam rooms with Barrier-Free™ tables that will allow all types of patients to get on the exam table effortlessly. And, Barrier-Free™ tables help you and your staff by decreasing the chance of back injury, extra labor and unproductive time.


Contact your sales Representative for details on how we can help you provide efficient patient care in your practice, or visit midmark.com.

Alko Enterprises

Ritter®
by MIDMARK



0960112



0960120

BURTON GLEAMER® TASK LIGHT

3,100 footcandles (Standard model). Color temperature of 3300°K. Cool housing temperature and low infrared; light can be in continuous use for hours without causing discomfort to physician or patient. Standard Model produces a spot of 6"; Wide Beam model produces a spot of 16". Easy grip ridges. Optimal form and function; 17" (43cm) rugged gooseneck and 16" (41cm) extension arm. Trouble-free design; 4,000 hour 35-watt halogen bulb.

0960112	Gleamer Light Floorstand
0960120	Gleamer Light w/extension arm & universal mount
0980112	Gleamer Wide Beam Spot Floorstand
0980120	Gleamer Wide Beam Spot w/extension arm & universal mount

Accessories


0009600PK	Set of 4 Gleamer 35 Watt Replacement Bulbs
0009608PK	Wide Beam Spot Replacement Bulbs, 4/pk

**WELCH ALLYN HALOGEN EXAM LIGHT™ III**

Intense halogen lamp provides brilliant 3200°K light, 750 footcandles of a 4" spot at 16". Yields smooth, even, cool illumination at distal end (48"). Extremely small distal output with easy-to-use focusing system for adjusting and concentrating the light spot. Includes fiber optic light pipe, focusing sleeve, and light box. Mounts on wall, table, desk, or mobile stand. 110-130V, 50/60 Hz.

48700	Exam Light III (stand not included)
48760	Exam Light III w/neonate light pipe

CIDA®



44300



44310

WELCH ALLYN LS135™ HALOGEN EXAMINATION LIGHT

Computer-designed reflector and front cover bezel provide homogeneous 6" spot of light with greater depth of field than conventional gooseneck lights. Halogen illumination's whiter light allows for truer flesh tones. Flexible 24" smooth encased arm. Available as floorstand, wall mount, table mount, or mobile examination light. UL listed, CSA certified, and IEC 601-1 approved.

44300	Halogen Examination Light w/integral floor stand
44310	Halogen Examination Light w/wall or table mount

**WELCH ALLYN LS150™ HALOGEN EXAM LIGHT**

50 watt lamp delivers superior light quality across the variable spectrum for more accurate color rendition. Specially redesigned light head keeps luminaire cool to the touch-even during multiple exams. Lamp has been offset at a 12" angle relative to the luminaire housing, allowing an unobstructed view of the field of interest. Reduced base profile frees up valuable exam room floor space. Encased, flexible 24" arm enables accurate one-hand positioning. Multiple installation options, including table-mounted, mobile stand and wall-mounted, enable customization to be suit needs.

44350	Caster Base For LS135 & LS150
44500	Halogen Exam Light w/floor stand
44510	Halogen Exam Light w/wall mount

Lighting



0008100PK



0112012



0112022



1007101

BURTON FOCUSABLE REFLECTOR LIGHTS

More than 3,000 footcandles. Includes a 100-watt halogen bulb with 2,000 hour life. Bulb housed in a heat absorbing, color-correcting glass cylinder. Design reduces shadows in field of vision. Choice of 12" (.30m) or 14" (.36m) reflectors. Battery option. Two SteriHandles and Adapters (available free with new units; can be purchased for older units).


- 0112012 12" Diameter Floorstand Model
0112022 12" Diameter Single Head Ceiling Mount Model

Accessories

- 0006126PK Halogen Reflector 100W Replacement Bulbs, 115V, 4/pk
0008100PK Sterile Disposable Light Covers, 25/pk
1003016 29" Drop Tube Rod
1006760 Heat Glass Cylinder Replacement
1007101 Reflector SteriHandle Assembly Replacement Model, 2/pk



44200



44210



44250

WELCH ALLYN LS200™ HALOGEN PROCEDURE LIGHT

Three lamp design provides shadow-free illumination and assures that if one lamp fails, two continue to operate. Allows working distance of 24" with 7" spot and 12" depth of field. 2500 footcandles, 3500-K color temperature. Drift-free maneuverability with a wide range of positioning options. Meets standards for UL 544, CSA 22.2-125, and IEC 601-1.

- 44200 Procedure Light w/arm & luminaire section

Accessories

- 44210 Wall Mount Accessory
44250 Mobile Stand Accessory



163001



163002



163003

BREWER FEATHERLIGHT™ 3000 SERIES

Featherlight™ movement is ease of positioning with counterbalance precision. No gas cylinders or tension springs to wear out. Multi light source provides shadow free characteristics with cool illumination. Remarkable 540° rotation at the arm and yoke provides flexibility of positioning. 5000 foot candle. Made in the U.S.A.

- 163001 Wall Model Exam Light
163002 Floor Model Exam Light
163003 Ceiling Model Exam Light

BURTON COOLSPOT™ II MINOR SURGICAL LIGHT

Up to 10,000 footcandles at 24" (61cm) and 4,000 footcandles at 1 meter. Up to 4700°K color temperature. Variable intensity wheel and high intensity switch. New, improved housing design that's easy to clean and includes a convenient removable, autoclavable SteriHandle™. Friction knob allows user to position self-balancing arm without a tool. Floorstand includes footswitch, locking casters and cord wrap.


- 0204370 Single Light w/Fastrac & Single Trolley
0204372 Double Light w/Fastrac & Double Trolley
0244312 Light w/Fleximount Floor Stand
0244322 Light w/Fleximount Single Ceiling Mount Kit
0244332 Light w/Fleximount Double Ceiling Mount Kit
0244342 Light w/Fleximount Wall Mount

Accessories

- 0007006PK Set of 4 EKE Replacement Bulbs

BURTON FASTRAC™ & FLEXMOUNT™ ASSEMBLIES

- 1008060 Outpatient II/Coolspot II 58" Ceiling Mount Drop Tube
0100770 Outpatient II/Coolspot II Fastrac Mount for 1 Trolley, 6' Length
0102110 Outpatient II Fleximount Floor Stand Base & 5' Upright (casters included)
0102160 Outpatient II Fleximount Wall Mount Kit
0102180 Outpatient II Single Ceiling Trolley (for one lighthouse)
0102190 Outpatient II Double Ceiling Trolley (for two lighthouses)
0102540 Outpatient II Fleximount Single Ceiling Mount Kit
0102550 Outpatient II Fleximount Double Ceiling Mount Kit



0204470



0224112



0224122



0224132



0224142

BURTON OUTPATIENT® II MINOR SURGICAL LIGHT

Enhanced illumination, plus focusability. 7,000 footcandles (40% brighter). Focuses with central removable, autoclavable SteriHandle™ II. Friction knob allows user to position self-balancing arm without a tool. Uses three “optically improved” 50-watt halogen bulbs - each with 2,000 hour life. Floorstand includes footswitch, locking casters and cord wrap.

0204470	Single Light w/Fastrac & Single Trolley
0204472	Double Light w/Fastrac & Double Trolley
0224112	Light w/Fleximount Floor Stand
0224122	Light w/Fleximount Single Ceiling Mount Kit
0224132	Light w/Fleximount Double Ceiling Mount Kit
0224142	Light w/Fleximount Wall Mount

Accessories

0006130PK	Set of 3 Hi-Intensity 50 Watt FTY Bulbs
1008428	Outpatient II SteriHandle Replacement Model

BURTON OUTPATIENT® II/COOLSPOT™ II SERIES COMBOS

A double mount with one Outpatient II and one CoolSpot II light to provide greater flexibility in high intensity lighting in a single fixed or overhead track system.

0204572	Outpatient II/Coolspot II Light w/Fastrac & Double Trolley
0234552	Outpatient II/Coolspot II Light w/Fleximount Double Ceiling Mount Kit

BURTON OUTPATIENT® PLUS MINOR SURGICAL LIGHT ACCESSORIES

0006123PK	Set of 3 EXT Watt Glass Covered Replacement Bulbs
1007696	Universal SteriHandle® Replacement Model



355-022



355-025



355-031


RITTER MINOR SURGERY LIGHTS

It's counterbalanced to give you precise positioning without drifting. Less than 2 1/2 lbs of force allow you to easily maneuver the light. Features include a large pattern size with an even distribution of light throughout, unsurpassed shadow control and superior foot-candle illumination. Single ceiling lights include one 100W quartz halogen bulb, dual ceiling lights include two 100W quartz halogen bulbs.

355-022	Mobile Exam Light
355-025	8 ft Single Ceiling Light
355-028	9 ft Single Ceiling Light
355-031	8 ft Dual Ceiling Light
355-034	9 ft Dual Ceiling Light

Accessories

9A231001	Conversion Kit w/cover, Convert Existing Single Light System To Dual System
----------	---



354-022



354-025



354-031

RITTER MINOR SURGERY SPOTLIGHT

The Ritter 354 has a color temperature of 4,200K to provide a whiter light and a truer representation of tissue color. Easy to use levers allow the amount of light and spot size to be adjusted. Thanks to the counterbalanced system, it takes just one hand to maneuver the spotlight allowing precise positioning every time with no need to worry about drifting. In addition to being durable and dependable, it has smooth, rounded, soft lines to make the exam room even more appealing. Single spotlight models include one 150W halogen bulb while double spotlight models include two 150W halogen bulbs.

354-022	Mobile Minor Surgery Light
354-025	Single Spotlight, 8 ft Ceiling
354-028	Single Spotlight, 9 ft Ceiling
354-031	Double Spotlight, 8 ft Ceiling
354-034	Double Spotlight, 9 ft Ceiling



RITTER MINOR SURGERY/SPOTLIGHT COMBO

The incredible features of both the Ritter 354 Spotlight and the Ritter 355 Minor Surgery Light can now be realized in a single product. It offers a unique solution for the physician who performs multiple procedures. Includes one 150W halogen bulb for the spotlight and one 100W halogen bulb for the surgery light.

354-037	354/355 Combo Light, 8 ft. Ceiling
354-040	354/355 Combo Light, 9 ft. Ceiling

BURTON GENESIS® PLUS & GENIE II ACCESSORIES

0002004PK	Genesis® Plus 50 Watt Halogen Bulbs, 6/set
0002008PK	Genie™ II 50W Halogen Replacement Bulbs, 4/pk
2100288	Genesis/Genie™/Genie™ II Replacement Spare Focus Handle
2100374	Steel Plate Assembly For Special Height Ceilings



BURTON GENESIS® PLUS MAJOR O.R. LIGHT

30" (76cm) diameter lighthouse. 14,800 footcandles (160,000 lux) at 1 meter. 4000°K color temperature. 6"-10" (15-25cm) spot size. Six 50-watt halogen bulbs wired in parallel to assure continuous illumination; each bulb rated for 2,000 hour life. Focus adjustments on center and side positioning handles. Continuous dimming. Completely smooth front surface for easy cleaning. Easy bulb change. Mounting system for 10' (3.03m) standard; special ceiling plates available for lower and higher ceilings.

0100322	Genesis® Plus Light, Single Ceiling
0100332	Genesis® Plus Light, Double Ceiling

BURTON GENIE™ II/GENIE™ PLUS MAJOR O.R. LIGHT

22" (56cm) diameter lighthouse. Incredible 12,000 footcandles (130,000 lux) at 1 meter for Genie II; 9,300 footcandles (100,000 lux) at 1 meter for Genie Plus. 4000°K color temperature. 6"-10" (15-25cm) spot size. Four 50-watt halogen bulbs wired in parallel to assure continuous illumination; each bulb rated for 2,000 hour life. Focus adjustment on center positioning handle. Continuous dimming. Completely smooth front surface for easy cleaning. Easy bulb change. Mounting system for 10' (3.03m) standard; special ceiling plates available for lower and higher ceilings.

0300522	Genie™ II Light Single Ceiling
0300532	Genie™ II Light Double Ceiling
0600392	Genesis/Genie Combo Unit



BURTON VISIONARY® O.R. LIGHT

20" (51cm) diameter lighthouse. 7,500 footcandles (81,000 lux) at 1 meter. 4000°K color temperature. 4"-8" (94-206mm) spot size. 360° limitless rotation around vertical axes. Focuses with autoclavable handle. Cool operation assured with special patented heat filter and reflector system. Four 35-watt quartz halogen bulbs - each with 2,000 hour bulb life. 115V standard; 110V, 230V/240V and battery (APS) models also available.

VS10DC	Double Ceiling Light, 115V
VS10SC	Single Ceiling Light, 115V
VS12DC	Double Ceiling Light, 230V
VS12SC	Single Ceiling Light, 230V

Accessories

0002000PK	Set of 4 Halogen Bulb
1019014	Visionary® Replacement/Spare Focus Handle



BURTON ULTRAVIOLET EXAMINATION LIGHTS

Portable, hand-held lights for maximum flexibility. Optically ground and polished glass magnifier lens with 3-power magnification and focal length of 8". Available with magnifier and white bulbs for application need. Unlike acrylic, Burton's polished glass lens will not scratch. 3-year warranty. Made in U.S.A.

0315010	Light w/two UV bulbs (no magnifier) - 115V
0316020	Light w/3X magnifier & Two UV Bulbs - 115V
0316030	Light w/3X magnifier, Two White & Two UV Bulbs - 115V

Accessories

0001127PK	Set of 4 Replacement White Bulbs
1007090	Wall Bracket For 31602, 31622, 31603 & 31623


BURTON EPIC MAGNIFIER

Will command attention based on aesthetics - but its real value is tied to its world class optics, performance and precision movement. The patented arm system stops exactly where positioned and will not require constant tightening or lock mechanism adjustments. Features a 3.5 diopter polished glass lens, 6.5" (165mm) diameter circular viewing area, powerful 22-watt bulb, quick start electronic ballast and integral handles and cover for lens protection.

- 0194040 115V Magnifier w/floorstand (casters included)
 0194042 115V Magnifier w/wall or table mount



0194040



0194042

BURTON WAVE™ PLUS FLUORESCENT MAGNIFIER

(2) 13-watt bulbs for shadow reduction. 3.5 diopter fluorescent magnifier. Arm allows complete freedom of movement. "Quick-start" switch and ballast to activate fluorescent bulbs. 10,000 hour bulb life. 3-year warranty. UL listed.

- 0184040 Magnifier, Floor Stand Model (casters included)
 0184042 Magnifier, Wall Mount Model

DURO-MED PAGE MAGNIFIER

Magnify entire page at once. Framed for ease of use.

- 599-8193-0255 Page Magnifier, 7 1/2" x 11 3/4"



6000X-F



6010-F

CLINTON BARIATRIC SERIES BLOOD DRAWING CHAIRS

Seat has additional depth and an extra center support for added strength. Extra dimensions make it easier for patient and technician. Adjustable height and depth flip arms. Armrests lock in place when in use, then flips out of the way for easy exit. Comfortable 2" padding throughout. Padded and upholstered seat, back and armrests. Easy-clean, antimicrobial, heavy weight, vinyl upholstery. Neutral Grey powder-coated frame. 2 knobs lock chrome armrest into position. Adjustable foot levelers. Chairs may be special ordered with straight armrests instead of flip arms (677 Option). Accommodates larger patient or child with parent up to 700 lbs. combined weight.

- 6000B-F Bariatric Chair w/2 flip arms, 35" x 21" Seat, 20"H x 48"W x 32"D
 6000B-FD Bariatric Chair w/flip arm & drawer, 35" x 21" Seat, 20"H x 58"W x 32"D

**CLINTON EXTRA WIDE BLOOD DRAWING CHAIRS**

- 6099X-F Extra Wide & Tall Chair w/2 flip arms, 35" x 18" Seat, 26"H x 48"W x 28"D
 6099X-FD Extra Wide & Tall Chair w/flip arm & drawer, 35" x 18" Seat, 26"H x 58"W x 28"D



6000X-F



6010-F



6020-F



6022-F

CLINTON LAB-X SERIES BLOOD DRAWING CHAIRS

1 1/4" square, heavy duty, all-welded, tubular steel frame. Adjustable height and depth flip arm is padded and upholstered. Armrest locks in place when in use, then flips out of the way for easy exit. Comfortable 2" padding throughout. Neutral Grey powder-coated frame. 2 knobs lock chrome armrest into position. Easy-clean durable antimicrobial vinyl. Adjustable foot levelers. Complete series of chairs to choose from. 400 lb Capacity.

- 6000X-F Extra Wide Chair w/2 flip arms, 35" x 18" Seat, 20"H x 48"W x 28"D
 6000X-FD Extra Wide Chair w/flip arm & drawer, 35" x 18" Seat, 20"H x 58"W x 28"D
 6010-F Chair w/flip arm, 17" x 16" Seat, 20"H x 34"W x 26"D
 6011-F Extra Tall Chair w/flip arm, 17" x 16" Seat, 26"H x 34"W x 26"D
 6020-F Chair w/flip arm & drawer, 17" x 16" Seat, 20"H x 43"W x 26"D
 6022-F Extra Tall Chair w/flip arm & drawer, 17" x 16" Seat, 26"H x 43"W x 26"D

HAUSMANN BLOOD DRAWING CHAIRS

Gray powder-coated steel frame. One-piece contoured black polypropylene seat. Armrests are interchangeable and can be set up for left or right hand use. One armrest pivots 360°. The other armrest is non-pivoting. Both armrests have 10" height range. Easy-to-clean gray laminate finish. Shipped ready-to-assemble.

- 2191 Blood Chair w/armrests that adjust from 26" to 36"H w/control knobs, 29"L x 18"D x 19"H
 2192 Tall Blood Chair w/footrest, Armrests Adjust From 34" to 44"H, 29"L x 26"D x 26"H
 2196 Blood Chair w/supply drawer & armrests that adjust from 26" to 36"H w/control knobs, 37 1/2"L x 18"D x 19"H
 2197 Tall Blood Chair w/footrest & supply drawer, Armrests Adjust From 34" to 44"H, 37 1/2"L x 26"D x 26"H

Chairs



165-005-200



165-006-200

RITTER BLOOD DRAWING CHAIR

Wide and comfortable with padded arm support that rotates 360°. Arm rest adjusts to accommodate any size/height patient. Storage drawer optional. Stain-resistant upholstery. Upholstery colors: Wedgewood, Monaco, Royal Carriage, Monet, Soft Pearl or coordinating designer patterns.

- 165-005 Blood Drawing Chair
165-006 Blood Drawing Chair w/drawer



1200 Series, 1250 Series

BREWER SIDE CHAIRS

- | | |
|---------|--|
| 1200-01 | Side Chair with Arms, Clamshell |
| 1200-02 | Side Chair with Arms, Gunmetal |
| 1200-03 | Side Chair with Arms, Adriatic |
| 1200-04 | Side Chair with Arms, Burgundy |
| 1200-05 | Side Chair with Arms, Emerald |
| 1200-06 | Side Chair with Arms, Mauve |
| 1200-07 | Side Chair with Arms, Azure Blue |
| 1200-08 | Side Chair with Arms, Royal Plum |
| 1200-09 | Side Chair with Arms, Plum |
| 1200-10 | Side Chair with Arms, Deep Teal |
| 1200-SP | Side Chair with Arms, Special Color |
| 1250-01 | Side Chair without Arms, Clamshell |
| 1250-02 | Side Chair without Arms, Gunmetal |
| 1250-03 | Side Chair without Arms, Adriatic |
| 1250-04 | Side Chair without Arms, Burgundy |
| 1250-05 | Side Chair without Arms, Emerald |
| 1250-06 | Side Chair without Arms, Mauve |
| 1250-07 | Side Chair without Arms, Azure Blue |
| 1250-08 | Side Chair without Arms, Royal Plum |
| 1250-09 | Side Chair without Arms, Plum |
| 1250-10 | Side Chair without Arms, Deep Teal |
| 1250-SP | Side Chair without Arms, Special Color |



C-10



C-20



C-40



C-50

CLINTON SIDE CHAIRS

- C-10 Plastic Seat Chair, Black Powder Coat Frame Finish, Nylon Glides, Gray
C-20 2" Thick Padded Seat, White Powder Coat Frame, Nylon Glides, 16 Colors of Upholstery Available
C-40 2" Thick Firm Foam Padded Seat w/out arms, Black Powder Coated Frame Finish, Nylon Glides, Contoured Seat & Back, 16 Colors of Upholstery Available
C-50 2" Thick Firm Foam Padded Seat w/arms, Black Powder Coated Frame Finish, Nylon Glides, Contoured Seat & Back, 16 Colors of Upholstery Available

HAUSMANN SIDE & ARM CHAIRS

Black epoxy steel frame. Upholstered seat and backrest. Specify color: Gray, Slate Blue, Rose, Grotto Green, Black.

- 2158 Arm Chair
2159 Side Chair



273008



273009

CIDA PNEUMATIC STOOLS

Features a 5-leg, die cast aluminum 22" diameter base. The base has a satin finish top with complementary gray sides and matching gray casters. Equipped with a smooth operating air cylinder lift mechanism, operated by a convenient semi-ring under seat. Seat cushion is 4" thick, 16" diameter. Available in upholstery colors to match any exam room.

- 273008 Stool 17.5"-24.5"
273009 Stool w/upholstered backrest



273015



273017



273019

CIDA SCREW-ADJUSTABLE REVOLVING STOOLS

Value priced adjustable examination stools feature all-welded construction and chrome-plated, tubular steel legs. Height adjusts from 20" to 28". Stools have a 14" seat diameter and 4" of padding. All models feature rubber, ball bearing casters.

273015	Revolving Stool
273017	Revolving Stool w/round foot ring
273019	Revolving Stool w/square foot ring



CIDA VALUE PLUS PNEUMATIC STOOL

Pneumatic height adjustment of 17 1/2" - 22.5". Features 5-leg, 21" diameter composite base. Poly foam seat cushion is 4" thick, 15" diameter. Available in upholstery colors to match any exam room

273006	Stool
--------	-------



21340



21340B

BREWER FOOT OPERATED PNEUMATIC ADJUSTABLE STOOL

All new ergonomic design. Hands free foot control, pneumatic height adjustment 18.5"-24". Powder coat black 23" aluminum base. 18" wide contoured polyfoam seat cushion. High impact black plastic backrest shroud. Seamless upholstery standard on this model. Available without backrest (features 16" wide, 4" thick round seat cushion).

21340	Foot Operated Stool
21340B	Foot Operated Stool w/backrest



BREWER LABORATORY STOOL

Designed for hospital, clinical and reference laboratories. Easy-to-clean, seamless upholstery and exceptional comfort with superb durability. 18" diameter seat cushion with 3" foam padding. Pneumatic height. Adjustable heel ring. Nylon glides allow for easy movement. Stable 23" five point base.

21520B	Laboratory Stool w/backrest, 18" Contour Seat, Adjusts 27" to 33"
21521B	Laboratory Stool w/backrest, Adjusts 24" to 29.5"



BREWER MEDICAL TASK CHAIR


21435B	Chair w/two-way seat/back adjustable
21435BA	Chair w/two-way seat/back adjustable & adjustable armrests

Stools

**BREWER PNEUMATIC STOOL**

Features a 5-leg, die cast aluminum 23" diameter base. The base has a satin finish top with complementary gray sides and matching gray casters. Equipped with a smooth operating air cylinder lift mechanism, operated by a convenient semi-ring under seat. Seat cushion is 4" thick, 16" diameter. 50 standard colors available, plus a variety of special order colors.

11001	Stool 17.5"-24.5"
11001B	Stool w/upholstered backrest
11001BFR	Stool Same as 11001B w/adjustable round footring
11001BVFR	Stool Same as 1101BFR w/seamless upholstery
11001FR	Stool w/adjustable round footring 20" - 27.5"
11001V	Stool w/seamless upholstery
11001VFR	11001FR Stool w/seamless upholstery
11001BV	Stool w/backrest & seamless upholstery

**BREWER SCREW-ADJUSTABLE REVOLVING STOOLS**

Value priced adjustable examination stools feature all-welded construction and chrome-plated, tubular steel legs. Stools have a 13 1/2" seat diameter and 4" of padding. All models feature rubber, ball bearing casters.

1010	Operator's Stool, 15" Round Seat, Range 17.75"-26" w/black accent base
1010B	1010 Stool w/backrest
23050	Revolving Stool 20"-28"
23051	Revolving Stool w/round foot ring 20"-28"
23052	Revolving Stool w/square foot ring 20"-28"

**BREWER THREADLESS SCREW ADJUSTABLE STOOL**

Threadless precision height adjustment 17.5-22.5". Comfortable 15" poly-foam seat cushion. 5 leg, 21" diameter black ABS plastic composite base. Available with 9.5" x 15.75" contoured backrest with durable black plastic shroud. Optional seamless/vacuum formed upholstery available.

22400	Stool
22400B	Stool w/backrest
22400BV	Stool w/backrest & seamless upholstery
22400V	Stool w/seamless upholstery

**BREWER VALUE PLUS PNEUMATIC STOOL**

Pneumatic height adjustment. 15" diameter. 4" thick seat cushion. 21" diameter black plastic composite base.

22500	Stool 17.5"-22.5"
22500B	Stool w/backrest
22500BV	Stool w/backrest & seamless upholstery
22500V	Stool w/seamless upholstery

**CLINTON CLASSIC SERIES STOOLS**

Full chrome frame. 1" diameter, S-style all-welded frame. Exclusive "Soft Feel" poly foam padding provides two layers of comfort. 4" thick seat, 2" rubber or nylon dual wheel casters. Solid steel machine screw spin height adjustment. All models have chrome foot ring. Cover cap conceals top welds. Seat is mounted to a 1/4" thick x 6" square plate. Premium, stain-resistant, woven, knit-backed, vinyl upholstery. 16 standard color choices. Seat diameter: 14"; Height range: 18 1/2" x 25 1/2".

2100	Stool w/2" rubber casters & square foot ring
2100-CC	Stool w/2" dual wheel hooded casters & square foot ring
2102	Stool w/2" rubber wheel, ball bearing casters & round foot ring
2102-CC	Stool w/2" dual wheel hooded casters & round foot ring
2111	Non-Adjustable Stool, 15" x 15"

CLINTON EPIC SERIES SPECIALTY STOOLS

Contemporary styling. Five-leg cast aluminum base. Soft Feel poly foam padding. Soft roll, dual wheel, hooded casters. Pneumatic cylinder height adjustment. Premium, stain-resistant, woven, knit-backed vinyl upholstery. Sixteen standard color choices.

- 2145-W Hands-Free Stool w/adjustable height back, 23" Black Powder Coated Base, 18" Contour Seat, Foot Activated Height Control
- 2145 Hands-Free Stool, 23" Black Powder Coated Aluminum Base, 4" Thick Round Seat, Foot Activated Height Control
- 2145-SP Hands-Free Stool w/adjustable height armrest w/extra firm padding, 4" Thick Round Seat, Foot Activated Height Control, 23" Black Powder Coated Aluminum Base
- 2155-SP Stool w/adjustable height armrest w/extra firm padding, Satin Cast Aluminum Base, 4" Thick Round Seat

CLINTON EPIC SERIES STOOLS

Contemporary styling. Cast aluminum base, 23". "Soft Feel" laminated seat cushion. Contour seat for extra comfort. Pneumatic cylinder height adjustment. Fully adjustable tilting seat. Premium, stain-resistant, woven, knit-backed, vinyl upholstery. Soft roll, dual wheel, hooded casters. Adjustable contour backrest with high impact shroud. 16 standard color choices.

- 2165-W Ergonomic Task Chair w/arms, Powder Coated Black Base
- 2166-W Lab Stool w/out arms, Brushed Base
- 2167-W Ergonomic Task Chair w/out arms, Powder Coated Black Base
- 2175-W Office Chair w/out arms, Brushed Base
- 2188-W Lab Stool w/arms, Brushed Base



CLINTON SELECT SERIES STOOLS

Premium, die cast, satin finish, 23" aluminum base. 5-leg design for greater stability. Exclusive "Soft Feel" poly foam padding is engineered for comfort. 4" thick seat. Hooded, dual wheel, nylon casters. Single lever, pneumatic, or easy spin lift height adjustment. 3/4" thick plywood top. Premium, stain-resistant, woven, knit-backed, vinyl upholstery. 16 standard color choices. Customizable with options.

- 2150 Spin Lift Stool w/17 1/4" - 21 1/2" Height Range
- 2155 Pneumatic Lift Stool w/18" - 23 1/2" Height Range
- 2165 Pneumatic Laboratory Stool w/22 1/2" - 30" Height Range, Adjustable Foot Ring & Black Nylon Glides
- 2166 Pneumatic Laboratory Stool w/22 1/2" - 30" Height Range, Adjustable Foot Ring & Dual Wheel Casters

CLINTON STAINLESS STEEL STOOL

Welded stainless steel frame. Smooth machine screw height adjustment. Height ranges from 22" up to 30". 2" ball-bearing casters. Round foot ring. 14" diameter seat. 4" thick seat for comfort.

- SS-2142 Stool
- SS-2142-21 Stool w/backrest



CLINTON STYLELINE SERIES STOOLS

Contemporary style. One-piece BIFMA tested nylon base with 3 color options. 5-leg design for greater stability. Exclusive "Soft Feel" poly foam padding provides super comfort. 4" thick seat, 2" nylon dual wheel casters. Single lever, pneumatic, or easy spin lift height adjustment. 3/4" thick plywood top. Premium, stain-resistant, woven, knit-backed, vinyl upholstery. 16 standard color choices.

- 2130 Spin Lift Stool w/18"-22 1/4" Height Range, Black Base
- 2135 Pneumatic Stool w/18"-23" Height Range, Black Base
- 2135-P Pneumatic Stool w/18"-23" Height Range, Putty Base
- 2135-S Pneumatic Stool w/18"-23" Height Range, Slate Base
- 2135-M Pneumatic Stool w/18"-23" Height Range



RITTER 125 AIR-LIFT STOOL

Equipped with a built-in shock absorber, an extra wide seat made of dense foam, and a beadless upholstery seam that eliminates the seat edge from digging into your legs. An air cylinder quickly and effortlessly raises and lowers the stool to the height that is ideal for you. Five-caster base. Height adjusts pneumatically from 18 1/2" to 24 1/2" with the handi-reach ring or lever. Many colors available.

- 125-001 Air-Lift Stool

RITTER 126 AIR-LIFT STOOL

Features a chrome five-caster, star-shaped base allowing physicians to get in close curing exams. It comes with a backrest to promote better posture and leave physicians feeling less fatigued at the day's end. The thick, but firm, foam padding enhances comfort, while the dual hooded silent casters allow the physician to sit on the edge of the seat, roll around the room and never tip. Height adjusts from 16 1/2" to 24 1/2". Many colors available.

- 126-001 Air-Lift Stool w/back



RITTER 127 SURGEON STOOL

Similar to the Air-Lift Stool, but with adjustable, padded body support. A foot-operated height adjustment provides precise positioning. From 18 1/2" to 24 1/2". Many color options

- 127-001 Air-Lift Surgeon Stool w/foot release

Stools

**RITTER 128 SURGEON STOOL**

All of the features of #127 stool, but with a hand operated height control. Many colors available.

128-001 Air-Lift Surgeon Stool w/hand release



195-001-200



196-001-200

RITTER 195/196 AIR-LIFT STOOLS

Features an amply padded, wide seat cushion and five-caster, black composite base with pneumatic height adjustment from 18 1/2" to 24 3/4". The thick, but firm, foam padding enhances comfort, while the dual hooded silent casters allow the physician to sit on the edge of the seat, roll around the room and never tip.

195-001 Air-Lift Five Leg Stool w/black base

196-001 Air-Lift Five Leg Stool w/black base & back



123-001



124-001



193-001



194-001

RITTER PHYSICIAN STOOLS

Height adjusts with a fixed height rotation and easy-to-use collar. Set at any height between 16 1/2" and 23 1/2". Base made of polished cast aluminum. Available in colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 R. Lavender, 206 Wedgewood, 207 Monaco, 208 Royal Carriage, 209 Monet, 210 Soft Pearl, 211 Waterford, 212 English Gardens, 213 Cabernet, 214 Tapestry, 215 Cathedral, 216 Pebble Grey, 235 Blueberry, 236 Silver Sage, 237 Tea Green, 238 Fossil, 239 Cashmere Blue, 240 Terra Cotta, 247 Belagio, 248 Siena, 249 Tuscany, 250 Firenze, 251 Capri, 252 Treviso, 241 Special Color Pacific Blue, 242 Special Color Spruce, 243 Special Color Olivine, 244 Special Color Arctic, 245 Special Color Sky, 246 Special Color Canyon, 999 Special Color.

123-001 Five Caster Stool

124-001 Five Caster Stool w/back

193-001 Five Caster Stool w/black base

194-001 Five Caster Stool w/back & w/black base

**HAUSMANN AIR-LIFT STOOLS**

16" diameter seat with 3 1/2" thick foam. Durable 5-leg, black composite base. (5) 2" twin wheel black nylon casters with urethane tread.

2113 Stool, Hand Control Pneumatic Height Adjustment

2123 Stool, Hand Control Pneumatic Height Adjustment, w/padded, spring backrest 10"H x 15"W

2133 Stool, Foot Control Pneumatic Height Adjustment

2143 Stool, Foot Control Pneumatic Height Adjustment, w/padded spring backrest 10"H x 15"W

**HAUSMANN CHOICE LINE™ AIR-LIFT STOOLS**

Hand operated air-spring pneumatic height adjustment (except 2150). 15" diameter seat, 3" thick foam. Durable 5-leg, black composite base. Five 2" twin wheel nylon casters.

2150 Revolving Stool w/adjustable worm gear height range of 7", Fixed Height Rotation

2153 Air-Lift Stool w/single lever height control handle, 6" Height Range

2157 Air-Lift Stool w/360° circular ring control handle, 6" Height Range

HAUSMANN P.T. "SCOOTER" STOOL

14" x 18" padded vinyl seat. Fixed height of 18". Four chrome legs with swivel casters.

2300 Stool

HAUSMANN PNEUMATIC LAB STOOL WITH BACKREST

16" diameter seat with 3 1/2" thick foam. Hand control pneumatic height adjustable 28 1/4" to 35" with padded backrest. 5 leg stationary base. Equipped with foot ring.

2160 Stool w/backrest

**BREWER STAINLESS STEEL STOOLS**

11200SS Step Stool

Accessories

11220SS Handrail For Step Stool 11200SS



11200-1




11220

BREWER STEP STOOL

Ideally suited for exam rooms, X-ray, patient rooms or any area where sturdy support is required. This stool has a chrome-plated finish with steel legs that are welded to the steel top. The top has a skid-resistant mat. The rubber leg-tips also provide skid-resistance. The chrome-steel handrail model reaches 35" in height from the floor to the top of the handrail.

- 11200-1 Step Stool
- 11200-2 Step Stool, 2/ctn
- 11220 Step Stool w/handrail
- 21200 Bariatric Step Stool, 600 lb Capacity
- 21220 Bariatric Step Stool w/handrail, 600 lb Capacity
- 31200 Heavy Duty Two-Step, Step Stool
- 31220 Heavy Duty Two-Step, Step Stool w/handrail



6110



6142

CLINTON BARIATRICS FOOT STOOLS

Heavy duty construction. Rubber tread top for safety. 600 lbs capacity.

- 6110 X-Large Step Stool, 23"L x 8"H x 23"W, Laminated Plywood Base, Sand
- 6142 Bariatrics Foot Stool, 14 1/4" x 11 1/4" x 9", Chrome-Plated Base



10-H



10-K

CLINTON PEDIATRIC STEP STOOL

- 10-H Step Stool w/hyenas (laminated)
- 10-K Step Stool w/koalas (laminated)



CIDA STEP STOOL

Ideally suited for exam rooms, X-ray, patient rooms or any area where sturdy support is required. This stool has a chrome-plated finish with steel legs that are welded to the steel top. The top has a skid-resistant mat. The rubber leg-tips also provide skid-resistance. The chrome-steel handrail model reaches 35" in height from the floor to the top of the handrail.

- 272022 Step Stool w/handrail

TECH-MED MAYO STAND

Secure height adjustments from 28" to 48". Mobile base includes locking casters. Base width 20". Removable stainless steel tray 12 1/2" x 18". 11 lb weight capacity.

- 4352 Mayo Stand w/U-Shaped Base, Height Adjusts From 31" - 50", Stainless Steel
- 4365 Mayo Stand w/mobile base

ALLIED "E" TRANSPORT SYSTEM

- 65361-E A-Straight Post Transport System, Includes: Cart, E Cylinder (Empty), 0 - 8 Regulator, Cannula w/7" tubing
- 65363-E B-Post Transport System w/on-off lever, Includes: Cart, E Cylinder (Empty), 0 - 8 Regulator, Cannula w/7" tubing
- 65364-E C-Post Transport System w/on-off lever & gauge, Includes: Cart, E Cylinder (Empty), 0 - 8 Regulator, Cannula w/7" tubing



BREWER I.V. STANDS

Heavy bar steel base to achieve the lowest possible center of gravity. Base is epoxy coated to resist corrosion. Ram's horn type hooks prevent accidental spilling. Height adjusts 51 3/4" - 94". No loss knob.

- 43403 Two-Hook Model I.V. Stand
- 43406 Four-Hook Model I.V. Stand
- 11290 I.V. Pole w/2-Ram's Horn Welded Hooks, Chrome Plated Steel Tubing, Four-Leg w/23" steel base, Height 47" - 81"
- 11320 I.V. Pole w/2-Ram's Horn Hooks, Chrome Plated Steel Tubing, Five-Leg w/21" plastic base, Height 51.75" - 94"
- 11321 I.V. Pole w/4-Ram's Horn Hooks, Chrome Plated Steel Tubing, Five-Leg w/21" plastic base, Height 51.75" - 94"



43428



43429



43430



43431



43433


BREWER IV STAND ACCESSORIES

Easy-to-change accessories make equipping an IV stand or infusion pump stand quick and easy. For Operating Rooms, Intensive Care, Emergency Room and Patient Care Units.

- 43428 Steering Handle
- 43429 Oxygen Cylinder Holder (holds D or E size)
- 43430 Infusion Pump Frame (holds multiple pumps)
- 43431 Personal Belongings Tray
- 43433 Universal Clamp (Required for all accessories)



11300



11350

BREWER PROFESSIONAL QUALITY I.V. STANDS

Professional quality IV stands for hospital or home use. Conventional use or as an infusion pump stand. Height adjusts from 47" - 84 1/2". All units have twin wheel casters.

- 11300 2-Rams Horn Hooks, 4-Leg Base I.V. Stand 47"-84 1/2"
- 11350 2-Rams Horn Hooks, Aluminum 5-Leg Base I.V. Stand 52.5"-93.75"
- 11360 4-Rams Horn Hooks, Aluminum 5-Leg Base I.V. Stand 52.5"-93.75"



BREWER STAINLESS STEEL I.V. STAND

Stainless steel. Height adjusts from 60" to 104". Comes with 2 ram's horn hooks. 22" diameter base has four legs and rolls on four 2" ball bearing style, conductive rubber casters.

- 48403 I.V. Stand w/4 casters

CLINTON IV STANDS

- IV-30 Value IV Pole, Chrome Plated Steel w/2 hooks & 4 legs
- IV-31 Value IV Pole, Heavy Base w/2 hooks & 4 leg base, Chrome Plated Steel w/knob height adjustment
- IV-314 Value IV Pole, Heavy Base w/4 hooks & 4 leg base, Chrome Plated Steel w/knob height
- IV-33 Select IV Pole, 6 Leg Heavy Base w/2 hooks
- IV-35 Select IV Pole, 5 Leg Base w/2 hooks
- IV-354 Select IV Pole, 5 Leg Base w/4 hooks
- IV-37 Select IV Pole, 5 Leg Nylon Base w/2 hooks
- IV-374 Select IV Pole, 5 Leg Nylon Base w/4 hooks
- IV334 Select IV Pole, 6 Leg Heavy Base w/4 hooks
- SS-31 IV Pole w/2 hooks, Stainless Steel

HAUSMANN 4-HOOK IV POLE

4-hook detachable top cap with push button release. Four leg anti-tip heavy steel base equipped with four 2" diameter twin-wheel nylon casters. Pole height adjusts from 49" to 85"H.

- 2188 IV Pole

LAKESIDE I.V. STAND ACCESSORIES

- 4870 Steering Handle, 10"W x 9"L x 1"H, Black Vinyl Non-Slip
- 4871 Universal Clamp, Attaches To Poles Up To 1" Diameter

Everything you need
for your
state-of-the-art
facility.



272011



272013

CIDA I.V. STANDS

Professional quality IV stands for hospital or home use. Conventional use or as an infusion pump stand. Height adjusts from 47" - 84 1/2". Epoxy coated steel base and no loss knob.

- 272011 2-Rams Horn Hooks I.V. Stand, Epoxy Coated Steel Base
- 272013 2-Rams Horn Hooks I.V. Stand, Aluminum Base
- 272015 4-Rams Horn Hooks I.V. Stand, Aluminum Base



272007



272017

CIDA HEAVY DUTY I.V. STANDS

Conventional use or as an infusion pump stand. The 1 1/4" lower column will accommodate all clamp-on infusion pumps. Extra heavy, wide stance base has six 3" rubber wheel casters for superior stability and easy rolling. Bases are epoxy coated to resist corrosion. Height adjusts from 58"-102".

- 272007 2-Hook Model I.V. Stand, 6-Leg Heavy Base
- 272009 4-Hook Model I.V. Stand, 6-Leg Heavy Base
- 272017 2-Hook Model I.V. Stand, 4-Leg Heavy Base
- 272019 4-Hook Model I.V. Stand, 4-Leg Heavy Base



272000



272003

CIDA MAYO INSTRUMENT STANDS

Available in two designs, both of durable heavy gauge chrome-plated steel tubing. Screw lock adjustment provides secure lock at desired height. Large stainless steel instrument tray is removable for easy cleaning. Two rubber wheel, ball bearing casters.


- 272000 Mayo Instrument Stand
- 272003 California Style Base Instrument Stand, 2 Pole Support



BREWER ELECTRIC NEEDLE STAND

Mounting panel measures 6 3/4" x 10" x 9/16". Height 38 1/2". Support post is 1" steel tubing. Stable, 5 point cast aluminum base with gray nylon casters. Strong round edge wire hook.

- 43577 Electric Needle Stand



43408



43416



43417

BREWER HEAVY DUTY INFUSION PUMP/I.V. STAND

Conventional use or as an infusion pump stand. The 1 1/4" lower column will accommodate all clamp-on infusion pumps. Extra heavy, wide stance base has six 3" rubber wheel casters for superior stability and easy rolling. Bases are epoxy coated to resist corrosion. No loss knob.

- 43408 I.V. Stand w/2 hooks
- 43409 I.V. Stand w/4 hooks
- 43416 I.V. Small Wheelbase w/twin ram's horn hooks
- 43417 I.V. Small Wheelbase w/4 ram's horn hooks

Stands

BREWER MRI ANTI-MAGNETIC PRODUCTS

Stainless steel.

43704 Stainless Steel Bowl Only



48450



48451

BREWER SOLUTION BOWL STANDS

48450 Single Solution Bowl Stand w/bowl

48451 Double Solution Bowl Stand w/bowl

**OMNIMED BEAM® ARTHROSCOPY STAND**

Featuring a large 16" x 16" stainless steel tray and security bracket to hold pumps up to 8 1/2" high. Hanger assembly holds up to four 5000 ml containers. Weighted base and rubber casters provide stability and smooth transport.

741302 Arthroscopy Stand

**OMNIMED BEAM® DOPPLER STAND**

Offers protection as well as convenience. Expandable coil cord secures doppler to stand, reducing disappearance or misplacement. One size pocket fits all. Optional basket to store monitor, earpiece or other items. Cast base for stability. Smooth rolling casters or mobility. 41" high.

350010 Doppler Instrument Stand

**OMNIMED BEAM® HEAVY WEIGHT™ INFUSION STAND**

The beefy column and weighted five leg base provide stability and holds up to four 5000 ml containers. Manually raises to a maximum height of nine feet. Cast aluminum lock securely holds at the desired height. Air-glide system controls descent of pole when unlocked.

741300 Heavy Weight™ Infusion Stand



350002



350005

OMNIMED BEAM® INSTRUMENT STAND ACCESSORIES

350002 Transport Basket, 10" x 13" x 4"

350005 Wire Basket

**OMNIMED BEAM® LAPTOP TRANSPORT STAND**

The large tray features dual handles for transport and height adjustment for standing or seated use. Simply depress the foot ring to adjust the surface height from 31" to 43". Star leg base and 3" rubber casters provide a smooth, quiet ride.

350305 Laptop Transport Stand



OMNIMED BEAM® POWER LIFTER® AUTOMATED IRRIGATION STAND

Heavy duty pneumatic irrigation stand featuring a rechargeable battery and weighted base. The touch of a foot pedal automatically raises irrigation solution containers to a maximum height of 9 feet. Depress the foot pedal and the pole lowers quickly and smoothly. Holds up to four 5000 ml irrigation containers. UL approved external battery charger. Can be operated by one person.

741314 Power Lifter Irrigation Stand

Accessories

741316 Low Voltage Light, Factory Installed

741317 Low Voltage Light Kit



OMNIMED BEAM® STAND

Features a fixed, stainless steel shelf (16 1/4" x 12 3/4") with an adjustable padded security clamp to hold the instrument in place. A sturdy 5 leg cast aluminum base provides excellent stability and the 3" silent rolling rubber casters assure smooth mobility. 42" floor to shelf height.

350000 Beam® Stand



RITTER INSTRUMENT PANEL

Houses physicians' most commonly used diagnostic equipment in one convenient location. Available in standard colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

9067-002 30"W x 19"H Instrument Panel

9067-003 34"W x 19"H Instrument Panel



BREWER KICK BUCKETS

Stainless steel bucket and frame.

48711 Kick Bucket w/frame

48712 Kick Bucket Pail Only

48714 Kick Bucket Frame Only

POLAR WARE KICK BUCKETS

Pails are tapered and ears for handles are an integral part of pail, resisting breakoff and allowing for easy cleaning. Stainless steel.

1014 Kick Bucket Only

PROFEX ANTI-STATIC KICK BUCKET LINERS

60 Dispenser Box Liners, 50/bx, 10 bx/cs

61 Dispenser Box Liners, 100/bx, 5 bx/cs

65 Liners, Bulk Pack, 500/cs

65B Economy Liners, Blue, 500/cs

65R3 Red Liners, Standard, 3 mil, 250/cs

75 Liners For Oval Kick Buckets, Bulk Packed, 500/cs

75B Liners For Oval Kick Buckets, Economy, Blue, 500/cs

75R3 Liners For Oval Kick Buckets, Red, 3 mil, 250/cs



6077



6088

CLINTON SPECIALTY BLOOD DRAWING CHAIRS

6077 Infant Station, 20"W x 32"L Top w/2 rims, 2 Safety Straps, Plastic Work Tray & Removable Middle Padded Section

6088 Wheelchair Station, Wall Mounted, Adjustable Height & Depth, 6"W x 24"L Armrest

Carts



6220



6221

HARLOFF LONG TERM CARE TREATMENT CART

Comes with 4 shelves with 4 standard bins, 1 shelf with 4 deep bins, and 1 storage tub on roll-out shelf. Dual flush-folding utility doors with Velcro® strips to hold doors open. Key lock with slam-lock feature, supplied with two keys. Ball-bearing full extension drawer slides on bottom tub shelf. Unibody steel construction. 5" sealed ball bearing premium casters, all full swivel, 2 with brake. Removable plastic top. Corner bumpers. Powder coat finish available in 12 colors.

6220 41.5"H x 32"W x 22"D Treatment Cart

6221 44.5"H x 38.5"W x 22"D Treatment Cart w/stainless steel top rail/push handle, 2 gallon waste container & 3 gallon sharps container & bracket



9029-8-346



9029-8-927



9038-20-346

HAUSMANN SMART CARTS

All laminate medical therapy carts 18 1/2"W x 18 1/2"D x 30"H; top surface 17"W x 17 1/2"D. Choice of natural oak or folkstone gray laminate finish. Back of cabinet is finished in matching laminate. Ideal for equipment, supplies or modalities. All cabinet case edges are high impact oak edge on oak case and black edge on gray case. High impact polystyrene drawer interiors with rounded inside corners for ease of cleaning. Equipped with easy-sliding, self-lubricating polystyrene track. Easy-to-use drawer pulls in satin aluminum with oak cabinet, or black nylon with gray cabinet. Casters are durable, 2 1/4" hooded ball casters for ease of mobility. Durable and maintenance-free laminate surfaces on all cabinet exteriors and inside cupboards. Shipped fully assembled.

9018-20-346 Drawer/Door Cart w/one adjustable shelf, Natural Oak (UPS-able)

9018-20-927 Drawer/Door Cart w/one adjustable shelf, Folkstone Gray (UPS-able)

9029-24-927 Drop Front Door Cart For Surface Mounted or Recessed EKG, Folkstone Gray (specify EKG) (UPS-able)

9029-8-346 One Drawer Cart, Natural Oak (UPS-able)

9029-8-927 One Drawer Cart, Folkstone Gray (UPS-able)

9038-20-346 Three Drawers Cart, Natural Oak (UPS-able)

9038-20-927 Three Drawers Cart, Folkstone Gray (UPS-able)

9018-10-346 Drawer/Door Cart w/cabinet, Natural Oak (UPS-able)

9018-10-927 Drawer/Door Cart w/cabinet, Folkstone Gray (UPS-able)

**ITTER 345 SPECIAL PROCEDURES CARTS**

Primarily used in examination rooms for general examinations and minor procedures. There are two storage drawers which accommodate supplies used during examinations. Replace XXX in item number to specify upholstery colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

345-004 Cart

345-005 Cart w/heater & Electrical Outlet

345-006 Cart w/gas lift

**ITTER 3702 MOBILE TREATMENT CARTS**

Offers the durability physicians need at a price they can afford. It depreciates much faster than built-ins, and you have a spectacular treatment cart. Comes with two 4" drawers and 2 doors and is 24"W x 18"D x 34"H or 24"W x 24"D x 34"H. Ritter door/drawer front and laminate countertop colors: replace XXX in item number to specify upholstery colors: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

3702-009 24"W x 18"D x 34"H Cart w/Sand Grey Kydex Top (Specify Door/Drawer Front Color)

3702-010 24"W x 18"D x 34"H Cart w/Sand Grey Kydex Top & Handles (Specify Door/Drawer Front Color)

3702-011 24"W x 18"D x 34"H Cart w/laminate top (Specify Color For Laminate Top & Door/Drawer Front)

3702-012 24"W x 18"D x 34"H Cart w/laminate top & handles (Specify Color For Laminate Top & Door/Drawer Front)

3702-013 24"W x 24"D x 34"H Cart w/Sand Grey Kydex Top (Specify Color For Door/Drawer Front)

3702-014 24"W x 24"D x 34"H Cart w/Sand Grey Kydex Top & Handles (Specify Color For Door/Drawer Front)

3702-016 24"W x 24"D x 34"H Cart w/laminate top & handles (Specify Color For Laminate Top & Door/Drawer Front)



BITTER 3705 MOBILE TREATMENT CART

Offers the durability physicians need at a price they can afford. Comes with five 4" drawers and is 24"W x 18"D x 34"H or 24"W x 24"D x 34"H. Door/drawer front and laminate countertop colors, replace XXX in item number to specify upholstery: 200 Sand Grey, 201 Slate Grey, 202 Iris Blue, 203 Nile Green, 204 Terra Rosa, 205 Royal Lavender, 999 Special Colors.

3705-009	Cart w/Sand Grey Kydex Top, 24"W x 18"D x 34"H w/five 4" drawers (Specify Door/Drawer Front Color)
3705-010	Cart w/Sand Grey Kydex Top & Handles, 24"W x 18"D x 34"H w/five 4" drawers (Specify Door/Drawer Front Color)
3705-011	Cart w/laminate top, 24"W x 18"D x 34"H w/five 4" drawers (Specify Laminate Top & Door/Drawer Front Color)
3705-012	Cart w/laminate top & handles, 24"W x 18"D x 34"H w/five 4" drawers (Specify Laminate Top & Door/Drawer Front Color)
3705-013	Cart w/Sand Grey Kydex Top, 24"W x 24"D x 34"H w/five 4" drawers (Specify Door/Drawer Front Color)
3705-014	Cart w/Sand Grey Kydex Top w/handles, 24"W x 24"D x 34"H w/five 4" drawers (Specify Door/Drawer Front Color)
3705-015	Cart w/laminate top, 24"W x 24"D x 34"H w/five 4" drawers (Specify Laminate Top & Door/Drawer Front Color)
3705-016	Cart w/laminate top & handles, 24"W x 24"D x 34"H w/five 4" drawers (Specify Laminate Top & Door/Drawer Front Color)

WATERLOO DUAL-SIDED MEDICATION CARTS

Best lock key system. Dual access drawers. Different drawer sizes. Lockable portable asset storage with bins. Bin sizes available are 3", 4.5" or 6". Four 5" full swivel locking casters with locking tow brake. Additional storage with the accessory module with built in waste receptacle. Economical cup dispenser that fits in on top of cart.. Wrap-around rubber bumper base provides perimeter protection. Pull-out shelf.

DMC-3	Dual-Sided Medication Cart Shell Only (order for custom carts)
DMC-3W	Dual-Sided Medication Cart w/o side module, Shell Only (order for custom carts)
DMC3-4S10	Dual-Sided Med Cart Combo Includes 4-Tier Cassettes & Bins, Must Identify 3", 4.5" or 6" Bins Desired
DMC3-55T10	Dual-Sided Med Cart Combo Bins Sold Separately, 5 Bin Drawers & Supply Drawers
DMC3-LTC1	Dual-Sided Med Cart Combo Includes EB-2B, Long Term Care Cart, No Patient Bins
DMC3-LTC3	Dual-Sided Med Cart Combo, 5 Bin Drawers, No Patient Bins Included Must Order Separately
DMC3-LTC6	Dual-Sided Med Cart Combo Includes EB-2B, Long Term Care Cart, No Patient Bins
DMC3-LTC7	Dual-Sided Med Cart Combo Includes EB-2B, Long Term Care Cart, No Patient Bins

WATERLOO UNI-CART® MED JUNIOR SERIES MEDICAL CARTS

Casters are 4" quality with 2 brake and 2 swivel. Options: two shell sizes, two shell colors eight drawer front colors, configurable drawer sizes, open lower storage space available. No standard lock. Optional gate lock available.

JSG-4369-LTG	4-Drawer Cart, All Light Gray, No Lock
JSG-436S-LTG	4-Drawer Cart, All Light Gray, No Lock
JSW-433336-PEDS	6-Drawer Cart, White/Multi-Colored Drawers, No Lock
JSW-4369-WHT	4-Drawer Cart, All White, No Lock
JSW-436S-WHT	3-Drawer Cart w/storage, All White, No Lock
JTG-433369-LTG	6-Drawer Cart, All White, No Lock, Lt. Grey
JTG-43669-LTG	5-Drawer Cart, All White, No Lock, Lt. Grey
JTW-433369-WHT	6-Drawer Cart, All White, No Lock
JTW-43669-WHT	5-Drawer Cart, All White, No Lock

WATERLOO UNI-CART® TREATMENT & PROCEDURE CARTS

Wrap around bumper provides perimeter protection. Ships assembled. Pre-plated zinc coated steel painted with a superior process of powder coat. Full extension drawers. White ABS plastic top and push handle and key lock standard. Optional push button lock with key override available. Contact Waterloo for more information and part numbers. Twenty brilliant color choices for drawer front colors. Eight color choices for the 9 drawer split drawer cart.

USGKU-4369-GRA	4-Drawer Cart, Grey Shell, Charcoal Drawers
UTGKU-433369-GRA	6-Drawer Cart, Grey Shell, Charcoal Drawers
UTGKU-9SR-GRA	9-Drawer Cart, Top Drawers Are Split (Drawer 1 is 4"H, Drawer 2 is 3"H), Grey Shell, Charcoal Drawers
UTGKU-BT3R	2-Drawer Cart w/doors & pull-out shelf, Grey Shell, Charcoal Drawers
UTGSU-9SR-GRA	9-Drawer Cart, Top Drawers Are Split (Drawer 1 is 4"H, Drawer 2 is 3"H), Grey Shell, Charcoal Drawers, Push Button Lock

AKRO UTILITY/PRO CART

Quick change from flat top to box top. I-Beam construction yielding the most accessibility for the lower shelf contents. Capability of adding rails to hand Akro-Bins or modular cabinets. Heavy duty structural foam for commercial usage. Heavy duty 5" non-marking casters. Two locking. 400 lb capacity limit (200 lb load limit per shelf). Unassembled.

30930	40"L x 19.5"W x 32"H Cart
30936	45"L x 24"W x 32"H Cart
Accessories	
30900	Rail Kit For 30930, Small Cart
30906	Rail Kit For 30936, Large Cart
30936SHELF	Middle Shelf For Large ProCart



BREWER MEDIUM/HEAVY DUTY STAINLESS STEEL UTILITY CART

Three shelves. 3" rubber wheel casters. Sturdy structural stainless steel handle. 300 lb capacity.

63500	Utility Cart
-------	--------------



Cart shown with basket

CIDA MULTIPURPOSE CART

Suited for almost any task. Carries two 6" wire baskets or one 12" basket (not included). Folds flat for storage. Chrome plated steel frame. No assembly required.

275015 Cart

Accessories

275017 6" Basket, 2/ctn

**CIDA UTILITY CART**

Plastic, all purpose and sturdy to accommodate almost all sensitive equipment and disposable supplies. 100 lb capacity. Three ABS plastic shelves and 3" swivel casters.

275010 Cart, 27 1/4" x 17" Shelf Size, 35 1/2"H, 11" Between Shelves

CLINTON MOBILE CARTS

TC-233 Chrome Cart Frame w/o baskets, 19 1/2"H x 30"D x 39"H

TC-50 Plastic Utility Cart, 22 1/2"W x 16"D x 32"H

8920 Mobile Equipment Cart w/2 drawers, Mounted On 5-Leg Black Nylon Base, Laminated

8940 Mobile Equipment Cart w/1 drawer & 1 shelf, Mounted On 5-Leg Black Nylon Base, Laminated

Accessories

TB-206 Wire Basket, 6" Deep For TC-233 Cart

TB-212 Wire Basket, 12" Deep For TC-233 Cart

ENRICHMENT LARGE SUPPLY CART

Specially designed shelves for safe distribution of equipment and supplies. 3 removable shelves with a load capacity of 150 pounds per shelf. Usable shelf space of 42" x 20". Includes 5" casters (2 locking, 2 non-locking).

40400 Large Supply Cart, 45" x 55" x 24"

ENRICHMENT UTILITY SUPPLY CART

Specially designed shelves for safe distribution of equipment and supplies. 3 removable shelves with a load capacity of 150 pounds per shelf. Usable shelf space of 24" x 16". Includes 3" casters (2 locking, 2 non-locking).

40450 Small Supply Cart, 33"H x 32"W x 20"D

GRAHAM FIELD UTILITY CART

Chrome-plated tubular steel frame with 4" casters. Three stainless steel shelves. Pushes from either end. Shipped knocked down; assembles in minutes. 400 lb load capacity.

8146 35 3/4"H x 28"L x 16"W Cart

**LAKESIDE ERGO-ONE™ UTILITY CARTS**

Ergonomic benefits are designed into every cart. All stainless steel construction and electronically welded. 20 gauge stainless steel shelves with reinforced edges. Leg bumpers protect walls and furnishings. Comfortable vinyl handle grips are available. Easy-to-clean and sanitize, simply wipe down or steam clean. Lake-Glide Casters: 5" No-Mark® cushion tread. Model 6600 caster assembly required.

- 6600 Standard Duty Cart w/2 shelves, 300lb Capacity, 16 1/2"W x 28"L x 44 3/8"H, 4-4" Swivel Casters
- 6610 Standard Duty Cart w/2 shelves, 300lb Capacity, 19"W x 31 1/8"L x 44 3/8"H, 4-4" Swivel Casters
- 6800 Medium Duty Cart w/2 shelves, 500lb Capacity, 16 1/2"W x 28"L x 44 3/8"H, 4-5" Swivel Casters
- 6810 Medium Duty Cart w/2 shelves, 500lb Capacity, 19"W x 31 1/8"L x 44 3/8"H, 4-5" Swivel Casters
- 6820 Medium Duty Cart w/2 shelves, 500lb Capacity, 22"W x 39 1/8"L x 44 3/8"H, 4-5" Swivel Casters
- 6830 Medium Duty Cart w/2 shelves, 500lb Capacity, 22"W x 53 1/8"L x 44 3/8"H, 2-5" Swivel & 2-5" Fixed Casters
- 7010 Heavy Duty Cart w/2 shelves, 700lb Capacity, 18 5/8"W x 35 3/8"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7015 Heavy Duty Cart w/3 shelves, 700lb Capacity, 18 5/8"W x 35 3/8"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7020 Heavy Duty Cart w/2 shelves, 700lb Capacity, 21 5/8"W x 41 3/8"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7025 Heavy Duty Cart w/3 shelves, 700lb Capacity, 21 5/8"W x 41 3/8"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7030 Heavy Duty Cart w/2 shelves, 700lb Capacity, 21 5/8"W x 57 3/16"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7120 Tough Transport Cart w/2 shelves, 1000lb Capacity, 21 5/8"W x 41 3/8"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7130 Tough Transport Cart w/2 shelves, 1000lb Capacity, 21 5/8"W x 57 7/16"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters
- 7140 Tough Transport Cart w/2 shelves, 1000lb Capacity, 24 5/8"W x 50 7/16"L x 46 3/4"H, 2-5" Swivel & 2-8" Fixed Casters



LAKESIDE GUARD RAIL CARTS

Ideal for any situation where you are transporting items which could fall off the cart easily. Protective guard rail helps eliminate spills and breakage of materials. Integrated push handle allows for easy maneuverability of your cart. All welded stainless steel construction allows for durability and easy sanitation. Standard handle and leg bumpers protect cart, walls and doorways from damage.

- | | |
|-----|---|
| 316 | Cart, 16 1/4"W x 27 1/2"L x 33 3/8"H, 300 lb Capacity, 3 Shelves, 3 1/2" Swivel Casters |
| 526 | Cart, 19"W x 31"L x 33 3/4"H, 500 lb Capacity, 3 Shelves, 4" Swivel Casters |
| 725 | Cart, 19 3/8"W x 32 5/8"L x 34 1/2"H, 700 lb Capacity, 2 Shelves, 5" Swivel Casters |
| 726 | Cart, 19 3/8"W x 32 5/8"L x 34 1/2"H, 700 lb Capacity, 3 Shelves, 5" Swivel Casters |
| 729 | Cart, 22 3/8"W x 38 5/8"L x 37 1/4"H, 700 lb Capacity, 2 Shelves, 5" Swivel Casters |
| 730 | Cart, 22 3/8"W x 38 5/8"L x 37 1/4"H, 700 lb Capacity, 3 Shelves, 5" Swivel Casters |



LAKESIDE UTILITY CART - PLASTIC

Specifically engineered for healthcare's quick and easy sanitation. Patented leg design for permanent snap together assembly without tools. Shelves precision molded from high-density polyethylene to resist stains and odors. Durable, versatile and resists dents.

- | | |
|------|---|
| 2000 | Cart w/3 Shelves, 17 1/8" x 32 1/2" x 34 7/8", 4" Swivel Casters, Beige |
| 2001 | Cart w/4 Shelves, 16" x 24", 4" Swivel Casters, Beige |
| 2300 | Cart w/6 Shelves, 17 1/4" x 26 5/8" x 35 1/2", 3" Swivel Casters, Beige |
| 2500 | Cart w/3 Shelves, 17 1/4" x 31 3/4" x 38", 4" Swivel Casters, Black |
| 2505 | Cart w/3 Shelves, 17 1/4" x 43" x 33", 4" Swivel Casters, Black |



LAKESIDE UTILITY CART - STAINLESS STEEL "U" FRAME

Heavy Duty 700 lb Capacity. Angled "U" shaped steel frame and electronic welding provide the strength needed for heavy duty jobs. Rugged 18 gauge reinforced stainless steel shelves are stain and rust resistant. Extra load casters with cushion tread non-marking wheels. Easy to clean and sanitize, simply wipe down or steam clean. Leg and handle bumpers protect walls and furniture. Caster Type and Size: 2 Swivel/2 Fixed, 5" dia./8" dia.

- | | |
|-----|---|
| 710 | Heavy Duty Cart w/2 shelves, 16 1/4"W x 30"L x 34 1/4"H |
| 711 | Heavy Duty Cart w/3 shelves, 16 1/4"W x 30"L x 34 1/4"H |
| 758 | Heavy Duty Cart w/2 shelves, 22 3/8"W x 54 5/8"L x 37"H |
| 759 | Heavy Duty Cart w/3 shelves, 22 3/8"W x 54 5/8"L x 37"H |
| 721 | Heavy Duty Cart w/2 Shelves, 19 3/8" x 32 5/8" x 34 1/2"H |
| 722 | Heavy Duty Cart w/3 Shelves, 19 3/8" x 32 5/8" x 34 1/2"H |
| 743 | Heavy Duty Cart w/2 Shelves, 22 3/8" x 38 5/8" x 37 1/8"H |
| 744 | Heavy Duty Cart w/3 Shelves, 22 3/8" x 38 5/8" x 37 1/8"H |



LAKESIDE UTILITY CARTS STAINLESS STEEL ANGLE LEG

All welded stainless steel construction allows for increased durability and easy sanitation. Handle and leg bumpers help protect cart, walls and doorways. Integrated push handle allows for easy maneuverability. Utility Carts (except NSF models) feature new and improved sound deadening material, which allows for quieter operation.

- | | |
|-----|--|
| 311 | Standard Duty Cart w/3 Shelves, 16 1/4"W x 27 1/2"L x 32 1/8"H, 300 lb Capacity, 3 1/2" Swivel Casters |
| 322 | Standard Duty Cart w/3 Shelves, 18 3/8"W x 30 3/4"L x 33"H, 300 lb Capacity, 3 1/2" Swivel Casters |
| 411 | Medium Duty Cart w/3 Shelves, 16 3/4"W x 27 5/8"L x 33 3/8"H, 500 lb Capacity, 3 1/2" Swivel Casters |
| 422 | Medium Duty Cart w/3 Shelves, 19"W x 31"L x 33 3/8"H, 500 lb Capacity, 3 1/2" Swivel Casters |
| 444 | Medium Duty Cart w/3 shelves, 22 3/8"W x 39 1/4"L x 37 1/4"H, 500 lb Capacity, 5" Swivel Casters |
| 459 | Medium Duty Cart w/2 shelves, 22 3/8"W x 54 1/8"L x 37 1/4"H, 500 lb Capacity, 5" Casters-2 Swivel/2 Fixed |

Carts



LAKE SIDE UTILITY CARTS - STAINLESS STEEL TUBULAR MODEL

All welded stainless steel construction allows for increased durability and easy sanitation. Handle and leg bumpers help protect cart, walls and doorways. Integrated push handle allows for easy maneuverability. Utility Carts (except NSF models) feature new and improved sound deadening material, which allows for quieter operation.

- 210 Medium Duty Cart w/2 shelves, 17 1/2"W x 27"L x 35 3/4"H, 500 lb Capacity, 4" Swivel Casters
- 211 Medium Duty Cart w/3 shelves, 17 1/2"W x 27"L x 35 3/4"H, 500 lb Capacity, 4" Swivel Casters
- 221 Medium Duty Cart w/2 Shelves, 20"W x 30"L x 35 3/4"H, 500 lb Capacity, 4" Swivel Casters
- 222 Medium Duty Cart w/3 Shelves, 20"W x 30"L x 35 3/4"H, 500 lb Capacity, 4" Swivel Casters
- 243 Medium Duty Cart w/2 shelves, 22"W x 36"L x 40 5/8"H, 500 lb Capacity, 5" Swivel Casters
- 244 Medium Duty Cart w/3 shelves, 22"W x 36"L x 40 5/8"H, 500 lb Capacity, 5" Swivel Casters
- 472 Medium Duty Cart w/2 shelves, 17 1/2"W x 27"L x 35 3/4"H, 500 lb Capacity, 4" Swivel Casters, Chrome Plated Legs/Frame
- 473 Medium Duty Cart w/3 shelves, 17 1/2"W x 27"L x 35 3/4"H, 500 lb Capacity, 4" Swivel Casters, Chrome Plated Legs/Frame
- 479 Medium Duty Cart w/2 shelves, 17 1/2"W x 27"L x 35"H, 500 lb Capacity, 4" Swivel Casters, Chrome Plated Legs/Frame
- 489 Medium Duty Cart w/3 Shelves, 20"W x 27"L x 35"H, 500 lb Capacity, 3 1/2" Swivel Casters, Chrome Plated Legs/Frame
- 492 Medium Duty Cart w/2 shelves, 23"W x 36"L x 40 1/8"H, 500 lb Capacity, 4" Swivel Casters, Chrome Plated Legs/Frame
- 493 Medium Duty Cart w/3 shelves, 23"W x 36"L x 40 1/8"H, 500 lb Capacity, 4" Swivel Casters, Chrome Plated Legs/Frame
- 499 Medium Duty Cart w/2 shelves, 23"W x 36"L x 37"H, 500 lb Capacity, 4" Swivel Casters, Chrome Plated Legs/Frame


LAKE SIDE UTILITY CARTS - TOUGH TRANSPORT

1000 lb Capacity - "U" Frame Leg Models. Angled "U" shaped steel frame and electronic welding provide the strength needed for heavy duty jobs. All joints are continuous heliarc welded to withstand heavy loads and continuous use. Rugged 14 gauge reinforced stainless steel shelves. Extra load casters with cushion tread non-marking wheels. Stainless steel design resists staining and rusting and is easy to maintain, clean and sanitize. Steam cleans for your convenience. Caster Type and Size: 2 Swivel/2 Fixed, 5" dia./8" dia.

- 938 U Frame Cart w/2 shelves, 19 3/4"W x 33"L x 34 1/4"H, 1000 lb Capacity
- 939 U Frame Cart w/3 shelves, 19 3/4"W x 33"L x 34 1/4"H, 1000 lb Capacity
- 943 U Frame Cart w/2 shelves, 22 3/4"W x 39"L x 37 3/8"H, 1000 lb Capacity
- 944 U Frame Cart w/3 shelves, 22 3/4"W x 39"L x 37 3/8"H, 1000 lb Capacity
- 947 U Frame Cart w/2 shelves, 25 7/8"W x 42"L x 37 3/8"H, 1000 lb Capacity
- 949 U Frame Cart w/3 shelves, 25 7/8"W x 42"L x 37 3/8"H, 1000 lb Capacity
- 953 U Frame Cart w/2 shelves, 25 3/4"W x 48"L x 37 3/8"H, 1000 lb Capacity
- 954 U Frame Cart w/3 shelves, 25 3/4"W x 48"L x 37 3/8"H, 1000 lb Capacity
- 958 U Frame Cart w/2 shelves, 22 3/4"W x 55"L x 37"H, 1000 lb Capacity
- 959 U Frame Cart w/3 shelves, 22 3/4"W x 55"L x 37"H, 1000 lb Capacity



264620




264626



264627



264629



264630


OMNIMED BEAM® UTILITY CARTS

Serves as an all-purpose cart. Use with 6" or 12" gray powder coated baskets. Features 1" metal tubing on 3" silent rubber casters. Easily folds flat for storage and shipping. Table-top storage racks mount securely to transport patient charts. Add a basket to accommodate completed charts or other necessary items.

- 264620 Utility Cart, 38"H x 19"W x 30"L

Accessories

- 264626 6" Deep Basket, 6"H x 14 1/4"W x 24"L
- 264627 12" Deep Basket, 12"H x 14 1/4"W x 24"L
- 264630 Table Top Storage Rack for 6 Capacity Cart
- 264629 Table Top Storage Rack for 8 Capacity Cart



OMNIMED BEAM® WIRE SHELF UTILITY CART

Perfect for transporting records, supplies and equipment throughout your facility. Consider for use in Medical Records, Central Supply, Housekeeping, Pharmacy Warehouse, etc. Features include three heavy-duty wire shelves and 3" rubber casters for smooth, quiet transport over any surface. 200 lb recommended maximum capacity.

- 264650 Wire Shelf Utility Cart, 38"H x 23"W x 17"D

WATERLOO UNI-CART® GENERAL & UTILITY CART ACCESSORIES

Color options for drawers and panel (replace XXX with color code): LTG-Light Gray, GRA-Charcoal, ELB-Electric Blue, DKB-Dark Blue, LTB-Light Blue, TLG-Teal Green, FWG-Fairway Green, SPG-Spring Green, FRG-Forest Green, YEL-Yellow, ORG-Orange, PNK-Pink, BRY-Berry, VIL-Violet, WNE-Wine, RED-Red, ELR-Electric Red, WHT-White, NTL-Natural, TPE-Taupe.

12D-XXX	12" Drawer (Select Color)
12DF-XXX	12" Drawer Front (Select Color)
3D-XXX	3" Drawer (Select Color)
3DF-XXX	3" Drawer Front (Select Color)
4D-XXX	4" Drawer (Select Color)
4DF-XXX	4" Drawer Front (Select Color)
6D-XXX	6" Drawer (Select Color)
6DF-XXX	6" Drawer Front (Select Color)
9D-XXX	9" Drawer (Select Color)
9DF-XXX	9" Drawer Front (Select Color)
ABH-2	Ambu Bag Hooks For ER2000 (2)
ANES-PKG	Anesthesia Cart Accessory Pkg Includes: SR-1 Side Rail, BSK-1 Waste Receptacle, BRS-2 Base Rail, TB-6 Tilt Bin, TD-1 Tape Dispenser, ORG-1S Small Clear Organizer
ANES-PKGDIX	Deluxe Pkg Includes: SR-1 Side Rail, BSK-1 Waste Receptacle, TB-6 Tilt Bin, TD-1 Tape Dispenser, ORG-1S Small Clear Organizer, SHLF-210 Base Rail System, Stabilizer Frame Bumper & Shelf, GMH-1 Glove/Mask Holder
BAG-3	Sealable Bags (20) For ER Tray 3
BAG-6	Sealable Bags (20) For ER Tray 6
BRK-ER	Bracket Attaches Shelf to UC
BRS-2	Base Rail System
BSK-1	Waste Receptacle w/clip
CB-1	Cross Bar
CBB-6	Cardiac Board & Basket
CBB-ER	Cardiac Board For ER2000
CBB-PC345	Cardiac Board & Brackets For PC345 & PC343HC Carts
CH-3	Catheter Holder
CHT-1	Chart Holder
CMB-3	Cassette Mount Bracket
CP-DIV-TRAY	Clear Plastic Dividers Trays
CR-1	Coat Rack
DIV-10L	Large Dividers, 10/ctn
DIV-10S	Small Dividers, 10/ctn
DIV-1S	Simple Divider Set, 10/ctn
DIV-5C	Divider 5" w/clip
DIV-CC2	2" Tall, Cut & Clip Dividers
DIV-CC4	4" Tall, Cut & Clip Dividers
DIV-DRW12	Divider Set For Crash Cart Drawers 1 & 2
DIV-DRW34	Divider Set For Crash Cart Drawers 3 & 4
DIV-DRW5	Divider Set For Crash Cart Drawer 5
DIV-DRW6	Divider Set For Crash Cart Drawer 6
EMB-1	Portable Emergency Med Box
EMB-1C	Portable Emergency Med Box w/caddy
EMER-BOX	3" Removable Med Box w/clear lid
EMER-BOX4	4" Removable Med Box w/clear lid
EMER-BOX4WDIV	4" Removable Med Box w/clear lid & dividers
EMER-BOXWDIV	3" Removable Med Box w/clear lid & dividers
EMER-DV	Dividers For EMER-BOX
EMER-PKGM	Emergency Cart Accessory Pkg Includes: CBB-6 Cardiac Board, DIV-CC2 Drawer Dividers, CP-DIV-TRAY Drawer Trays, IV-2 IV Pole, OTB-4 Oxygen Tank Bracket, PS-100 Seals, R-3 Railing, SSG-2C Side Shelf, UH-2 Hooks
ER-BIN-L	Large Door Bin ER Series
ER-BIN-S	Small Door Bin ER Series
ER-LBL	Labels For ER2000 Drawers
ER-SHLF	Swivel Monitor Shelf ER Series
ER-SHLFM	Swivel Monitor Shelf Uni-Cart
ER-TRAY3	3" Transparent Tray
ER-TRAY6	6" Transparent Tray
GLSGK-13B	Gate Lock For Short Cart, All Light Gray, Key Lock
GLSGK-22	Gate Lock For Short Cart, All Light Gray, Key Lock
GLSWB-433336	Gate Lock For Short Cart, All White, Key Lock

WATERLOO UNI-CART® GENERAL & UTILITY CART ACCESSORIES (Continued)

GLSWK-13B	Gate Lock For Short Cart, All White, Key Lock
GLSWK-22	Gate Lock For Short Cart, All White, Key Lock
GLTGG-28	Gate Lock For Tall Cart, All Light Gray, Key Lock
GLTWK-28	Gate Lock For Tall Cart, All White, Key Lock
GMH-1	Glove & Mask Holder
INF-PKG	Infection Control Cart Accessory Pkg Includes: CR-1 Coat Rack, R-3 Railing, UH-2 Hooks
IV-2	IV Pole For Uni-Cart
IV-3	IV Pole For ER2000
LT-1	Magnetic Flex Light
MINI-BX	Portable Mini Box
O-2	Oxygen Tank Bracket For ER2000
ORG-IL	Large Clear Organizer Box
ORG-IS	Small Clear Organizer Box
OS-1	Outlet Strip Hospital Grade
OTB-4	Oxygen Tank Bracket For Uni-Cart
OTB-PC345	Oxygen Tank Brackets
PC1000	Silent Knight Replacement Bags, 1000/pkg
PNL-XXX	9" Factory Installed Panel (Select Color)
PRO-PKG	Procedure Cart Accessory Pkg Includes: DIV-1 Dividers, DT-6 & DT-7 Drug Trays, R-3 Railing, SSG-2 Shelf, UH-2 Hooks
PS-100	Plastic Breakaway Seals, 100/ctn
R-3	Railing
SB-1	Sharps Bracket
SB-2L	Sharps Bracket, Lockable
SD-3	Security Drawer, Lockable
SF-4	Stabilizer Frame Bumper
SH-LTG	Side Handle, Gray
SHLF-2	Shelf, Uni-Cart Single Tier
SR-1	Side Rail
SS-M	Suction Shelf
SSG-2C	Side Shelf, Gray
T-2	Replacement Top
TB-3	Tilt Bin w/3 large compartments
TB-4	Tilt Bin w/4 large compartments
TB-5	Tilt Bin w/5 large compartments
TB-6	Tilt Bin w/6 large compartments
TD-1	Tape Dispenser
TT-1	Transfer Tray w/o dividers
TT-1A	Transfer Tray w/bins
TT-1B	Transfer Tray w/bins
UH-2	Utility Hooks
WALLY	Waste Receptacle
WALLY-ATT	Waste Receptacle w/attached bracket
WALLY-REM	Waste Receptacle w/removable bracket w/side rail
WB-1L	Large Wire Basket
WB-1S	Small Wire Basket
WDC-2025	Plastic Tray w/dividers Fits in 3" Drawers or Larger
WDC-2050	Plastic Tray w/dividers Fits in 6" Drawers or larger
WNH-1	Wire Notebook Holder
WSB-1	Wire Storage Basket

AKRO CARTS

Polyethylene container. 16-gauge zinc-plated steel tubing frame. Latch bars that extend around the frame for added strength. Joints are rived and heavy duty 90 gal. joints are welded. All are available with or without a piano hinged lid which props open during use. Specify color when ordering: green, gray or blue.

76460	Medium Duty 60 Gallon Capacity Cart w/o lid
76461	Medium Duty 60 Gallon Capacity Cart w/lid
76490	Medium Duty 90 Gallon Capacity Cart w/o lid
76491	Medium Duty 90 Gallon Capacity Cart w/lid
76590	Heavy Duty 90 Gallon Capacity Cart w/o lid
76591	Heavy Duty 90 Gallon Capacity Cart w/lid


**ALLIED SMALL CYLINDER CART**

Chrome-plated steel with 6" heavy duty wheels and rubber handgrip.

- 65061 Step Saver 4 Small Cylinder Carrier
- 65070 Small Cylinder Cart, Holds Single "E" Sized Cylinder, Fixed Handle, 1/cs
- 65076 Small Cylinder Cart, Holds Single "E" Sized Cylinder, Adjustable Handle
- 65090 Twin Cylinder Cart, Small Cylinders
- 65096 Multi Cylinder Cart, Small Cylinders, Holds 12

**BURDICK CART WITH NYLON CASTERS**

- 008046 Universal Cart

**BURDICK CLINIC ECG CARTS**

Designed to accommodate all Burdick ECG products. Electrocardiographs mount securely to the carts. Their sturdy design make these carts a functional and attractive addition to the primary care office. Shipped fully assembled. Mar-resistant finish. Large bottom tray for supply storage. Specify the Universal Clinic Cart for the following electrocardiographs: EK 10, Elite II, Eclipse LE, Eclipse LEi, E350, E350i, E560. Specify the Eclipse Clinic Cart for the following electrocardiographs: Eclipse 4, 400, 8, 800 and 850.

- 008022 Universal Clinic ECG Cart, 34"H x 16"W x 26"D

BURDICK STRESS ECG CART

Accommodates the Burdick electrocardiographs and the associated stress systems. Shipped fully assembled. Features: swing out/lock away handle, 5" tracking locked swivel casters, accessory drawer and defibrillator tray.

- 008018 Stress ECG Cart, 32 3/4"H x 19"W x 25"D

ENRICHMENT EMERGENCY CART

Versatile cart for secure transfer of selected emergency equipment. Support racks positioned for optimum protection of 2 E-cylinder tanks. Three 2" deep shelves for additional supplies. Electrical cord wrap brackets. Front hook for CPR board. 5" casters (2 locking, 2 non-locking). Sure-Chek dust cover standard; colors available are: 33-Tan; 37-White; 39-Light Blue.

- 40650 Emergency Cart w/dust cover, 35"H x 37"W x 20"D

ENRICHMENT ICE CART

Ideal for beverage distribution. 36 qt. Cooler, ice scoop with pouch, drainage pail and Sure-Chek skirt included. 3 shelves for easy accessible serving and storage. Skirt colors available are: 33-Tan; 37-White; 39-Light Blue. 3" casters (2 locking, 2 non-locking).

- 40200 Ice Cart w/skirt, 41"H x 28"W x 20"D

ENRICHMENT LASER/MONITOR CART

Specifically designed for use in laser surgery. Standard with 5 shelves, monitor tie-down strap, E-cylinder holder and electrical cord wrap brackets. 4 locking 3" casters.

- 40100 Laser/Monitor Cart, 57"H x 36"W x 20"D

ENRICHMENT UTILITY/EKG CART

Features a full length drawer for supplies and equipment. Two full shelves. 3" casters (2 locking, 2 non-locking).

- 40700 Utility/EKG Cart w/drawer, 33"H x 32"W x 20"D

HARLOFF ANESTHESIA CART

Mechanical combination pushbutton lock (7451, 7351, 7350, 7456). Key lock with two keys (6451, 6450, 6351, 6456, 6350). Hospital grade 6 socket electrical outlet with 15 ft cord and cord wrap (6456, 7456). Two 3" exchange trays and two 6" exchange trays (6456, 7456). Full wrap-around vinyl bumper on separate steel frame. 5" sealed ball-bearing premium casters, all full swivel, two with brake. Removable plastic top. Raised back rail system with two horizontal rails (7451, 6451, 6351, 7351, 6456, 7456). Four aluminum push handle/side rails (6350 & 7350 have two). Ball-bearing full extension drawer slides. Cart is delivered with wheels installed; some accessories require user installation. Unibody steel construction. Powder coat finish (12 colors available).

- 6350 4 Drawer Keylock Workstation Cart w/premium rubber wheels, 35.5"H x 32"W x 22"D
- 6351 4 Drawer Specialty Pkg Keylocked Workstation w/premium rubber wheels, 51.5"H x 32"W x 22"D
- 6450 6 Drawer Workstation Cart w/key lock, Standard Pkg, 41.5"H x 32"W x 22"D
- 6451 6 Drawer Workstation Cart w/key lock, Specialty Pkg, 57.75"H x 32"W x 22"D
- 6456 Deluxe Anesthesia Cart, Keylock w/premium rubber wheels, 57.5"H x 32"W x 22"D
- 7350 4 Drawer Combination Lock Workstation Cart w/premium rubber wheels, 36.5"H x 32"W x 22"D
- 7351 4 Drawer Specialty Pkg Combination Locked Workstation w/premium rubber wheels, 52"H x 32"W x 22"D
- 7450 6 Drawer Combination Lock Workstation w/premium rubber wheels, 42.5"H x 32"W x 22"D
- 7451 6 Drawer Workstation Cart w/mechanical combination pushbutton lock, Specialty Pkg, 58.25"H x 32"W x 22"D
- 7455 5 Drawer Anesthesia Cart, Special Accessory Pkg, Combination Lock
- 7456 Deluxe Anesthesia Cart, Combination Lock w/premium rubber wheels, 58"H x 32"W x 22"D
- 7650 Standard Monitor Cart Combination Lock w/premium rubber wheels
- 7651 Specialty Pkg Combination Locked Monitor Cart w/premium rubber wheels

HARLOFF CAST CART ACCESSORIES

602001	Stainless Steel Folding Shelf
602002	Stainless Steel Side Storage Cabinet
602003	Tourniquet Pump Pole
602003R	Retrofit Bucket Bracket w/swivel tourniquet pump pole, For Field Installation
602004	Scissors Security Cable (scissors not included)
602005	Clear Label Holders, 10/pk
602007	Tourniquet Pump Pole Only
602012	Black Vinyl Drawer Mat
602013	3" Drawer, Stainless Steel, Cast Cart
602014	3" Drawer Divider Kit, Cast Cart
602016	6" Drawer, Stainless Steel, Cast Cart
602017	6" Drawer Divider Kit, Cast Cart
602019	9" Drawer, Stainless Steel, Cast Cart
602106	Top Compartment/Organizer, Cast Cart - NEW - Holds BSN and 3M Supplies
603013	3" Drawer, Painted, Cast Cart
603016	6" Drawer, Painted, Cast Cart

HARLOFF E-SERIES CART

E24-3B	E-Series Cart 24" Cabinet, 3 Drawers (1-6", 2-9"), Breakaway Lock
E24-3K	E-Series Cart 24" Cabinet, 3 Drawers (1-6", 2-9"), Key Lock w/2 keys
E24-4B	E-Series Cart 24" Cabinet, 4 Drawers (2-3", 1-6", 1-12"), Breakaway Lock
E24-4EMS	E-Series Emergency Cart, 24" Cabinet, 4 Drawers, Breakaway Lock & Accessory Pkg (accessories must be customer installed)
E24-4K	E-Series Cart 24" Cabinet, 4 Drawers (2-3", 1-6", 1-12"), Key Lock w/2 keys
E24-5B	E-Series Cart 24" Cabinet, 5 Drawers (3-3", 1-6", 1-9"), Breakaway Lock
E24-5K	E-Series Cart 24" Cabinet, 5 Drawers (3-3", 1-6", 1-9"), Key Lock w/2 keys
E24-6B	E-Series Cart 24" Cabinet, 6 Drawers (4-3", 2-6"), Breakaway Lock
E24-6K	E-Series Cart 24" Cabinet, 6 Drawers (4-3", 2-6"), Key Lock w/2 keys
E24-7B	E-Series Cart 24" Cabinet, 7 Drawers (6-3", 1-6"), Breakaway Lock
E24-7K	E-Series Cart 24" Cabinet, 7 Drawers (6-3", 1-6"), Key Lock w/2 keys
E24-8B	E-Series Cart 24" Cabinet, 8 Drawers (8-3"), Breakaway Lock
E24-8K	E-Series Cart 24" Cabinet, 8 Drawers (8-3"), Key Lock w/2 keys
E30-10B	E-Series Cart 30" Cabinet, 10 Drawers (10-3"), Breakaway Lock
E30-10K	E-Series Cart 30" Cabinet, 10 Drawers (10-3"), Key Lock w/2 keys
E30-2K	E-Series 30" Cabinet, 2 Drawers (1-3", 1-6"), Locking Storage Compartment w/2 doors (21"H w/shelf), Key Lock w/4 keys
E30-3B	E-Series Cart 30" Cabinet, 3 Drawers (1-6", 2-12"), Breakaway Lock
E30-3K	E-Series Cart 30" Cabinet, 3 Drawers (1-6", 2-12"), Key Lock w/2 keys
E30-4B	E-Series Cart 30" Cabinet, 4 Drawers (1-3", 1-6", 1-9", 1-12"), Breakaway Lock
E30-4K	E-Series Cart 30" Cabinet, 4 Drawers (1-3", 1-6", 1-9", 1-12"), Key Lock w/2 keys
E30-5B	E-Series Cart 30" Cabinet, 5 Drawers (2-3", 1-6", 2-9"), Breakaway Lock
E30-5K	E-Series Cart 30" Cabinet, 5 Drawers (2-3", 1-6", 2-9"), Key Lock w/2 keys
E30-6B	E-Series Cart 30" Cabinet, 6 Drawers (4-3", 1-6", 1-12"), Breakaway Lock
E30-6EMS	E-Series Emergency Cart, 30" Cabinet, 6 Drawers, Breakaway Lock & Accessory Pkg (accessories must be customer installed)
E30-6K	E-Series Cart 30" Cabinet, 6 Drawers (4-3", 1-6", 1-12"), Key Lock w/2 keys
E30-7B	E-Series Cart 30" Cabinet, 7 Drawers (5-3", 1-6", 1-9"), Breakaway Lock
E30-7K	E-Series Cart 30" Cabinet, 7 Drawers (5-3", 1-6", 1-9"), Key Lock w/2 keys
E30-8B	E-Series Cart 30" Cabinet, 8 Drawers (6-3", 2-6"), Breakaway Lock
E30-8K	E-Series Cart 30" Cabinet, 8 Drawers (6-3", 2-6"), Key Lock w/2 keys
E30-9B	E-Series Cart 30" Cabinet, 9 Drawers (8-3", 1-6"), Breakaway Lock
E30-9K	E-Series Cart 30" Cabinet, 9 Drawers (8-3", 1-6"), Key Lock w/2 keys
EMS1	Accessory Pkg: Stainless Steel Top Rail/Push Handle, Side-Mounted Drop Shelf, IV Pole, Cardiac Board, Oxygen Bottle Holder, (2) 65830-P1 Drawer Trays & (4) Utility Hooks

HARLOFF ECONOMY TREATMENT CART

Key lock with two keys. Removable plastic top. Drawer fronts have integrated pull handles. Roller drawer slides. Unibody steel construction. 3" steam casters, two with brake. Taupe powder coat finish.

ETC-3	27.25"H x 23.38"W x 15.13"D Three Drawer Treatment Cart w/bottom compartment
ETC-5	27.25"H x 23.38"W x 15.13"D Five Drawer Treatment Cart

HARLOFF INFECTION CONTROL/ISOLATION CART

Key lock with two keys (6510, 6511). Mechanical combination pushbutton lock (7511). Stainless steel breakaway locking system (6500). Full wrap-around vinyl bumper on separate steel frame. 5" sealed ball-bearing premium casters, all full swivel, two with brake. Removable plastic top. Aluminum push handle (6500, 6510, 7510). Stainless steel top rail/push handle (7511, 6511, 6501). Ball-bearing full extension drawer slides. Cart is delivered with wheels installed; some accessories require user installation. Unibody steel construction. Powder coat finish (12 colors available).

6500	Infection Control Cart, Breakaway Lock, Standard Pkg, 41.5"H x 32"W x 22"D
6501	Infection Control Cart, Breakaway Lock, Specialty Pkg, 62.5"H x 32"W x 22"D
6510	Infection Control Cart, Key Lock, Standard Pkg, 41.5"H x 32"W x 22"D
6511	Infection Control Cart, Key Lock, Specialty Pkg, 62.5"H x 32"W x 22"D
6520	4 Drawer Standard Package Infection Control Cart, Key Lock
6521	4 Drawer Specialty Package Infection Control Cart, Key Lock
7510	Infection Control Cart, Mechanical Pushbutton Combination Lock, Standard Pkg, 42.5"H x 32"W x 22"D
7511	Infection Control Cart, Mechanical Pushbutton Combination Lock, Specialty Pkg, 62.5"H x 32"W x 22"D

HARLOFF INSTRUMENT CART

Key lock with two keys. Stainless steel top rail/push handle (6726/6736/6746). Removable plastic top (6726/6736/6746). Aluminum push handle (6725/6735/6745). Ball-bearing full extension drawer slides. 3" stem caster all full swivel, two with brake. Cart is delivered with wheels and push handle installed. Powder coat finish (12 colors available).

6725	6" Instrument Cart w/key lock, Standard Pkg, 35.75"H x 26"W x 20.5"D
6726	6" Instrument Cart w/key lock, Specialty Pkg, 38.5"H x 26"W x 21.5"D
6735	12" Instrument Cart w/key lock, Standard Pkg, 35.75"H x 23.5"W x 21.5"D
6736	12" Instrument Cart w/key lock, Specialty Pkg, 38.5"H x 23.5"W x 21.5"D
6745	18" Instrument Cart w/key lock, Standard Pkg, 35.75"H x 26"W x 20.5"D
6746	18" Instrument Cart w/key lock, Specialty Pkg, 38.5"H x 26"W x 21.5"D

HARLOFF MAGNETIC RESONANCE SAFE ANESTHESIA CART

Aluminum construction for cabinet and all drawers. Replaceable plastic top with integrated push handles for easy to maneuverability. Stainless steel pull-out shelf. Key and breakaway lock models available. Standard drawer configurations. Ball-bearing full extension drawer slides. White cabinet color with the choice of 12 colors for the drawer

MR6K	6 Drawer MR Safe Cart w/key lock, Standard Pkg, 45.75"H x 30.5"W x 23.75"D
MR6K-MAN	6 Drawer MR Safe Cart w/key lock, Anesthesia Pkg, 59.75"H x 32"W x 25"D
MR7K	7 Drawer MR Safe Cart w/key lock, Standard Pkg
MR7K-MAN	7 Drawer MR Safe Cart w/key lock, Anesthesia Pkg

HARLOFF MAGNETIC RESONANCE SAFE CART ACCESSORIES

MR-4TILT	4-Compartment Tilt Bin
MR-5TILT	5-Compartment Tilt Bin
MR-6TILT	6-Compartment Tilt Bin
MR-CARDBRD	Cardiac Board & Mounting Brackets
MR-IV2	IV Pole
MR-O2	O2 Bracket
MR-POSELF	Replacement Pullout Shelf
MR-RBR2	2-Tier Raised Back Rail
MR-RBR3	3-Tier Raised Back Rail
MR-RLHK	Rail Hook
MR-SR	Side Rail
MR-UHOOKS4	Utility Hooks, 4/pk
MR-3GWASTE	Waste Container, 3 gallon, with mounting bracket, MR-Safe
MR-8QWASTE	Waste Container, 8 quart, with mounting bracket, MR-Safe

Carts

HARLOFF MAGNETIC RESONANCE SAFE EMERGENCY CART

Aluminum construction for cabinet and all drawers. Replaceable plastic top with integrated push handles for easy to maneuverability. Stainless steel pull-out shelf. Key and breakaway lock models available. Standard drawer configurations. Ball-bearing full extension drawer slides. White cabinet color with the choice of 12 colors for the drawers.

MR6B	6 Drawer MR Safe Cart w/breakaway lock, Standard Pkg, 45.75"H x 30.5"W x 23.75"D
MR6B-EMG	6 Drawer MR Safe Cart w/breakaway lock, Emergency Pkg, 45.75"H x 34.5"W x 25.5"D
MR7B	7 Drawer MR Safe Cart w/breakaway lock, Standard Pkg
MR7B-EMG	7 Drawer MR Safe Cart w/breakaway lock, Emergency Pkg

HARLOFF MALIGNANT HYPERTHERMIA CART

Demonstrates JCAHO compliance and/or best practice for all or suites, Level 1 Trauma Units, Surgery Centers, and Ambulatory Facilities. Cart includes: integrated 1.6 cu ft refrigerator/freezer with hospital grade plug, 10 glycol packs, LED-illuminated digital thermometer with cart power safety feature, 3" and 6" storage drawers, stainless steel breakaway locking system on both the refrigerator and the two storage drawers, full wrap-around vinyl bumper on separate steel frame, ball-bearing full extension drawer slides, removable plastic top, 5-inch premium wheels with one brake and one tracking caster, powder coat finish in 12 colors. The cart increases patient safety and departmental efficiency by providing a convenient, portable and economic solution for storage and transport of cooled glycol packs, cold saline IV solutions, and dantrolene (Dantrium) used to stabilize MH. Cart can hold critical MH emergency supplies at the proper temperature in close proximity to areas where they are most likely needed.

6640	MH Cart, Standard Package, 1-Compartment w/pre-installed refrigerator, 41.5" H x 32" W x 22" D
6640K	Malignant Hyperthermia Cart Standard Package, Keylocking Drawers, Breakaway lock on refrigerator, ea
6641	MH Cart, Specialty Package, 1-Compartment w/pre-installed refrigerator, 44.5" H x 32" W x 22" D, Accessories Included: IV Pole (extends to 85"/216cm), Full Drawer Tray Organizer, Side-Mounted Drop Shelf & Stainless Steel Top Rail/Handle

HARLOFF MEDICATION CART LOCKS AND KEYS

40406-XX	Lock, Cabinet, Economy Medication Cart
40680	BEST® Lock, Replacement, Narcotics Drawer
40681	BEST® Lock, Replacement, Cabinet
40682-KEY	BEST® Lock Replacement Key, Cabinet or Narc. Drawer
40693-XX	Key for 40406-XX Lock
40695	Master Key for 40406 Locks
40696-XX	Lock, Narcotics Drawer, 5 Wide/Econ Med Carts
40723-XX	Key for 40696-XX Lock



3145B



3145B-EMG



3156B

HARLOFF MINI-LINE EMERGENCY CARTS

Narrow design with reduced footprint. Stainless steel breakaway locking system. Ball-bearing full extension drawer slides. All steel construction. 3" stem casters, all full swivel, 2 with brake. Removable plastic top. Powder coat finish available in 12 colors.

3145B	32"H x 18"W x 18"D Short Five Drawer Emergency Cart
3145B-EMG	58.5"H x 25.5"W x 18"D Short Five Drawer Emergency Cart w/IV pole, oxygen tank holder, cardiac board & aluminum push handle/side rail
3156B	38.25"H x 18"W x 18"D Tall Six Drawer Emergency Cart



3145K



3145K-ANS



3156B-EMG



3156K

HARLOFF MINI-LINE WORKSTATION CARTS

Narrow design with reduced footprint. Key lock with two keys. Ball-bearing full extension drawer slides. All steel construction. 3" steam casters, all full swivel, 2 with brake. Removable plastic top. Powder coat finish in 12 colors.

3145K	32"H x 18"W x 18"D Five Drawer Short Workstation Cart
3145K-ANS	45.63"H x 18"W x 18"D Five Drawer Short Workstation Cart w/2-tier raised back rail system, 2 side rails & a 6-tilt bin organizer
3156B-EMG	64.75"H x 25.5"W x 18"D Tall Six Drawer Emergency Cart w/IV pole, oxygen tank holder, cardiac board & aluminum push handle/side rail
3156K	38.25"H x 18"W x 18"D Six Drawer Tall Workstation Cart
3156K-ANS	54.75"H x 18"W x 18"D Six Drawer Tall Workstation Cart w/2-tier raised back rail system, 2 side rails & a 6-tilt bin organizer

HARLOFF MONITOR ANESTHESIA CART

Short overhead shelf system with two shelves (6650). Tall overhead shelf system with two shelves offering multiple mounting locations (6651). Key lock with two keys. Full wrap-around vinyl bumper on extended width steel frame. 5" sealed ball-bearing premium casters, all full swivel, two with brake. Removable plastic top. Aluminum push handle/side rail. Hospital grade 6 socket electrical outlet with 15 ft cord and cord wrap. Ball-bearing full extension drawer slides. Cart is delivered with wheels installed; accessories require user installation. Unibody steel construction. Powder coat finish (12 colors available).

6650	Standard Monitor Cart Keylocked w/premium rubber wheels, 66.5"H x 40"W x 25"D
6651	Specialty Pkg Keylocked Monitor Cart w/premium rubber wheels, 75"H x 40"W x 25"D

Ask about our
new products.

HARLOFF PHYSICIAN CART

Key lock with four keys (6150/6151). Key lock with two keys (6132/6142/6146). Stainless steel breakaway locking system (6160/6161/6130/6140/6145). 18" storage compartment, locking door, and adjustable shelf (6150/6160/6151/6161). Full wrap-around vinyl bumper on separate steel frame. 5" casters, all full swivel, two with brake. Removable plastic top (6151/6161). Stainless Steel top rail (6151/6161). Aluminum push handle (6150/6160/6130/6132/6142/6140/6145/6146). Ball-bearing full extension drawer slides. Cart is delivered with wheels installed. Unibody steel construction. Powder coat finish (12 colors available).

6130	GP Line 5 Drawer Cart w/breakaway lock, Standard Pkg, 40.25"H x 28.5"W x 22.25"D
6132	GP Line 5 Drawer Cart w/key lock, Standard Pkg, 40.25"H x 28.5"W x 22.25"D
6133	GP Line 5 Drawer Cart w/o lock, Standard Pkg
6140	GP Line 6 Drawer Cart w/breakaway lock, Standard Pkg, 40.25"H x 28.5"W x 22.25"D
6142	GP Line 6 Drawer Cart w/key lock, Standard Pkg, 40.25"H x 28.5"W x 22.25"D
6145	GP Line 4 Drawer Cart w/breakaway lock, Standard Pkg, 40.25"H x 28.5"W x 22.5"D
6146	GP Line 4 Drawer Cart w/key lock, Standard Pkg, 40.25"H x 28.5"W x 22.5"D
6150	GP Line 3 Drawer & Storage Compartment Cart w/key lock, Standard Pkg, 40.25"H x 28.5"W x 22.25"D
6151	GP Line 3 Drawer & Storage Compartment Cart w/key lock, Specialty Pkg, 43"H x 28.5"W x 22.25"D
6160	GP Line 3 Drawer & Storage Compartment Cart w/breakaway lock, Standard Pkg, 40.25"H x 28.5"W x 22.25"D
6161	GP Line 3 Drawer & Storage Compartment Cart w/breakaway lock, Specialty Pkg, 43"H x 28.5"W x 22.25"D

HARLOFF TREATMENT/STORAGE CART

21" storage compartment, locking doors and adjustable shelf. Key lock with four keys. Full wrap-around vinyl bumper on separate steel frame. 5" sealed ball-bearing premium casters, all full swivel, two with brake. Removable plastic top. Stainless steel top rail/push handle (6201). Aluminum push handle (6200). Ball-bearing full extension drawer slides. Cart is delivered with wheels installed; some accessories require user installation. Unibody steel construction. Powder coat finish (12 colors available).

6200	Multi-Treatment Cart, Standard Pkg, 41.5"H x 32"W x 22"D
6201	Multi-Treatment Cart, Specialty Pkg, 44.5"H x 32"W x 22"D

HARLOFF WOOD MEDICATION FURNITURE

WD270BOX	Wood Med Cart, 8 Drawer, Up to 270 Unit Dose Bed Box Capacity, No Narcotics Box, Choice of 6 finishes and 3 top colors
WD270BOXN	Wood Med Cart, 8 Drawer, Up to 270 Unit Dose Bed Box Capacity, with internal Narcotics Box included, Choice of 6 finishes and 3 top colors
WD480PC	Wood Med Cart, 4 Drawer, Up to 480 Card Capacity, No Narcotics Box, Choice of 6 finishes and 3 top colors
WD480PCN	Wood Med Cart, 4 Drawer, Up to 480 Card Capacity with internal Narcotics Box included, Choice of 6 finishes and 3 top colors
WD7952	Waste Basket, aluminum, laminated to match wood cabinet, Choice of 6 finishes
WD7953	Waste Basket with top cover, laminated to match wood cabinet, Choice of 6 finishes



LAKESIDE HOMED CLASSIC EMERGENCY CARTS

18 and 20 gauge steel with electronically-welded seams. Epoxy-polyester powder coating - resists marring, chipping and peeling. Rubber corner bumpers extend out 1 1/4" for wall and cart protection. Removable ABS plastic top allows for thorough cleaning. Premium casters - donut shaped, all swivel, two lock. Roller bearing slides or compound glides on drawers. Pre-drilled holes allow for easy accessory installation. Push handle standard.

3903C-XX	3-Drawer Cart, Short Shell, 26 5/8"W x 18"D x 37"H, Drawers: (2) 3", (1) 6", 13" Bottom Compartment, Locking Bar on Drawers, Butterfly Lock on Door (Specify color: RD-Red, BL-Blue)
3914C-XX	4-Drawer Cart, Short Shell, 37 1/2"H x 18"D x 30 1/2"W, Drawers: (2) 3", (1) 6", (1) 12", Key Lock (Specify color: RD-Red, BL-Blue)
3926C-XX	6-Drawer Cart, Tall Shell, 44"H x 18"D x 30 1/2"W, Drawers: (4) 3", (1) 6", (1) 12", Security Coded Dimple Lock (Specify color: RD-Red, BL-Blue)



23954M-XX



23965M-XX



23966M-XX

LAKESIDE HOMED PREFERRED EMERGENCY CART

Heavy duty steel construction with electronically welded seams. Premium 5" (12.7cm) donut-shaped casters make this cart easy to maneuver (2 locking, 2 swivel). Redesigned caster system makes casters the easiest to install. Double-wall construction for added stability and strength. ABS plastic top removes for thorough cleaning. Two-thirds bumper system for extra protection. Built-in push handle maximizes cart mobility. Built-in side shelf slides in and out for quick access. Built-in stabilizing base prevents accidents from happening. Full length drawer pulls simplify opening the drawers. Pre-drilled holes speed up accessorizing. Hybrid powder coat finish endures and looks good years down the hall. Electronic lock available. Security-coded dimpled key and lock system is master key-able. "Easy Turn" gang locking mechanism for quick locking. Roller bearing glides make drawers easy to open.

23954M-XX	Cart, Short Shell, 35"H x 22"D x 26.5"W, Drawers: (2) 3", (1) 6", (1) 12", Key Lock (Specify color: RD-Red, BL-Blue, AW-Autumn White)
23954MEL-XX	Cart, Short Shell, 35"H x 22"D x 26.5"W, Drawers: (2) 3", (1) 6", (1) 12", Electronic Lock (Specify color: RD-Red, BL-Blue, AW-Autumn White)
23965M-XX	Cart, Tall Shell, 43"H x 22"D x 26.5"W, Drawers: (2) 3", (2) 6", (1) 12", Key Lock (Specify Color: RD-Red, BL-Blue, AW-Autumn White)
23965MEL-XX	Cart, Tall Shell, 43"H x 22"D x 26.5"W, Drawers: (2) 3", (2) 6", (1) 12", Electronic Lock (Specify color: RD-Red, BL-Blue, AW-Autumn White)
23966M-XX	Cart, Tall Shell, 43"H x 22"D x 26.5"W, Drawers: (4) 3", (1) 6", (1) 12", Key Lock (Specify color: RD-Red, BL-Blue, AW-Autumn White)
23966MEL-XX	Cart, Tall Shell, 43"H x 22"D x 26.5"W, Drawers: (4) 3", (1) 6", (1) 12", Electronic Lock (Specify Color: RD-Red, BL-Blue, AW-Autumn White)

**LAKESIDE HOMED PREFERRED MEDICATION CARTS**

- 23944MAW 4-Drawer, Punch Card Medication Cart, 43"H x 22"D x 32"W w/side mount storage cabinet, Drawers: (1) 3", (3) 9", 9 Punch Card Dividers, 29000NC Locking Narcotic Box, Includes Waste Receptacle & Storage Cabinet, Key Lock, Autumn White
- 23970MAW Unit Dose Medication Cart, 43"H x 22"D x 34"W, Drawers: (1) 3", (1) 6", (1) 12", 29000NC Locking Narcotic Box, Includes Waste Receptacle & Storage Cabinet, Security Coded Dimple Key & Lock, Autumn White

**LAKESIDE HOMED PREFERRED PROCEDURE/ TREATMENT CARTS**

Heavy duty steel construction with electronically welded seams. Premium 5" (12.7cm) donut-shaped casters make this cart easy to maneuver (2 locking, 2 swivel). Redesigned caster system makes casters the easiest to install. Double-wall construction for added stability and strength. ABS plastic top removes for thorough cleaning. Two-thirds bumper system for extra protection. Built-in push handle maximizes cart mobility. Built-in side shelf slides in and out for quick access. Built-in stabilizing base prevents accidents from happening. Full length drawer pulls simplify opening the drawers. Pre-drilled holes speed up accessorizing. Hybrid powder coat finish endures and looks good years down the hall. Electronic lock available. Security-coded dimpled key and lock system is master key-able. "Easy Turn" gang locking mechanism for quick locking. Roller bearing glides make drawers easy to open.

- 23922M-XX Procedure Cart, 2-Drawer/2-Door, 43"H x 22"D x 26.5"W, Drawers: (1) 3", (1) 6", 20" Doors w/1 adjustable pull-out shelf, Key Lock (Specify color: RD-Red, YE-Yellow, AW-Autumn White)
- 23968M Treatment Cart, Tall Shell, Latex Free, 43"H x 22"D x 26.5"W, Drawers: (7) 3", (1) 9", Key Lock, Autumn White w/orange drawers
- 23969M-XX Treatment Cart, Tall Shell, 43"H x 22"D x 26.5"W, Drawers: (8) 3", (1) 6", Key Lock (Specify Color: RD-Red, BL-Blue, YE-Yellow, AW-Autumn White)

**MIDMARK DIAGNOSTIC CART**

Includes: ECG Acquisition Module Gooseneck Accessory, printer shelf and storage basket.

- 4746-001 Diagnostic Cart

**OMNIMED BEAM® RUBBER CASTER BUMPERS**

An excellent way to protect your furniture, walls and chart racks against accidents. 5" diameter rubber bumpers fit over casters on all Omnicart® Vertical and Flat Storage Racks.

- 264004 Rubber Bumpers, 4/pkg

WATERLOO MEDTRX & HALFTRX LONG TERM CARE CARTS

Removable, replaceable ABS top with large unobstructed work area. Styled push handles. Full size and half size drawers. Optional double lock narcotics storage available. Roller bearing slides with quick release levers. Wrap around bumper. Four 5" swivel casters, one brake, one tracking. Best® key lock system is standard. Options are push-button lock system or electronic lock with key override and narcotic locks for drawers for additional charge.

- | | |
|--------------------|------------------------------------|
| HALFTRX-1.0 | Half Trax Cart Shell Only |
| MEDTRX-1.5 | Full MedTrx Cart Shell Only |
| HTK1-33311 | Half Trax Combo w/BLS-1GC |
| HTK1-33TTT6 | Combo Includes Cassettes/Bins, BLS |
| HTK1-6TTTT | Combo Includes Cassettes/Bins, BLS |
| HTK1-999 | Half Trax Combo w/BLS-1GC |
| MTK1-4411/33311 | Med Trax Combo w/BLS-1GC |
| MTK1-4411/3611 | Med Trax Combo w/BLS-1GC |
| MTK1-4411/6TTTT | Combo Includes Cassettes/Bins, BLS |
| MTK1-4411/999 | Med Trax Combo w/BLS-1GC |
| MTK1-4411/TFTF6 | Combo Includes Cassettes/Bins, BLS |
| MTK1-4771/33FTF6 | Combo Includes Cassettes/Bins, BLS |
| MTK1-4771/33TTTT | Combo Includes Cassettes/Bins, BLS |
| MTK1-4771/999 | Med Trax Combo w/BLS-1GC |
| MTK1-4771/TFTFT | Combo Includes Cassettes/Bins, BLS |
| MTK1-8A8A8A/8A8A8A | Med Trax Combo w/BLS-1GC |
| MTK1-999/33311 | Med Trax Combo w/BLS-1GC |

WATERLOO UNI-CART® ANESTHESIA CARTS

Wrap around bumper provides perimeter protection. Ships assembled. Carts are pre-plated zinc coated steel painted with a superior process of powder coat. Full extension drawers. Roller bearing slides with quick release levers. Key lock standard. Optional push button lock with key override available, as well as other drawer front colors.

- | | |
|------------------|---|
| USGKU-4369-DKB | 4-Drawer Anesthesia Cart, Grey Shell, Dark Blue Drawers, Key Lock |
| UTGKU-433369-DKB | 6-Drawer Anesthesia Cart, Grey Shell, Dark Blue Drawers |
| UTGSU-46612-FWG | 4-Drawer Anesthesia Cart, Grey Shell, Fairway Green Drawers, Push Button Lock |

WATERLOO UNI-CART® EMERGENCY DEPARTMENT CARTS

Wrap around bumper provides perimeter protection. Ships assembled. Carts are pre-plated zinc coated steel painted with a superior process of powder coat. Full extension drawers. Lever lock standard. Standard white ABS plastic top and push handle. 5" casters (2 swivel, 1 brake, 1 tracking). Optional push button lock with key override and gate locks available.

Special Order Cart	Special Order Cart
USBLU-436P-ELB	3-Drawer Cart w/panel, All Blue
USGKU-4369-ELB	4-Drawer Cart, Grey Shell, Blue Drawers
USGKU-4369-RED	4-Drawer Cart, Grey Shell, Red Drawers, Key Lock
USRLU-4369-RED	4-Drawer Cart, All Red
USRLU-436P-RED	3-Drawer Cart w/panel, All Red
UTBKU-433369-ELB	6-Drawer Cart, All Blue
UTBLU-433369-ELB	6-Drawer Cart, All Blue
UTBLU-43336P-ELB	5-Drawer Cart w/panel, All Blue
UTGLU-8PEDS	8-Drawer Peds Cart
UTGLU-9PEDS	9-Drawer Peds Cart
UTRKU-433369-RED	6-Drawer Cart, All Red
UTRLU-433369-RED	6-Drawer Cart, All Red
UTRLU-43336P-RED	5-Drawer Cart w/panel, All Red
UXGLA-10PEDS-ACC	10-Drawer Cart, Light Gray Shell/Multi-Colored Drawers

WATERLOO UNI-CART® ER2000™ INSTANT ACCESS EMERGENCY CART

Ergonomic molded top/push handle/rail provides a large unobstructed work area. Removable clear transparent drawer trays includes four 3" and one 6" tray, on full extension roller bearing slides. Additional storage in side door bins. Large yellow instant access breakaway seal mechanism to shear seal and doors open quickly. Shipped assembled and available in red or blue shell. The following accessories come standard and included in price: IV pole, cardiac board and bracket, oxygen tank bracket, plastic numbered seals, ambubag hooks, airway, breathing, circulation labels.

ER2000-B	Instant Access Crash Cart, Blue
ER2000-R	Instant Access Crash Cart, Red

WATERLOO UNI-CART® GENERAL & UTILITY CARTS

Baked red enamel finish. Full access pull-out drawers. Four 4" x 2" plastic casters. Locking tow brakes on two prevent rolling and swivel. Two swivel and two stationary. Pegboard sides allow extra storage. Full-extension drawer.

PC343HC	3-Drawer Cart w/compartments, Red
PC345	5-Drawer Cart, Red

WATERLOO UNI-CART® ISOLATION & INFECTION CONTROL CARTS

Wrap around bumper provides perimeter protection. Ships Assembled. Pre-plated zinc coated steel painted with a superior process of powder coat. Full extension drawers. White ABS plastic top and push handle and key lock standard. Optional push button lock with key override available. Contact Waterloo for more information and part numbers. Twenty brilliant color choices for drawer front colors. Eight color choices for the 9 drawer split drawer cart.

USGKU-499-YEL	3-Drawer Isolation Cart, Grey Shell, Yellow Drawer Front
UTGKU-4699-YEL	4-Drawer Isolation Cart, Grey Shell, Yellow Drawer Front


**LAKESIDE EQUIPMENT STANDS**

Versatile mobile equipment stands come in six different models. Ideal for transporting equipment which needs to be used in many different areas along with the extra supplies in the drawers. Stainless steel construction allows for easy sanitation. One and two drawer models have optional locks available.

356	Stand w/1 Drawer, 3 Shelves, 16 3/4" x 18 3/4" x 30 1/8"H, (4) 3 1/2" Swivel Casters
358	Stand w/2 Drawers, 2 Shelves, 16 3/4" x 18 3/4" x 35"H, (4) 3 1/2" Swivel Casters
359	Stand w/out Drawers, 3 Shelves, 16 3/4" x 18 3/4" x 30 1/8"H, (4) 3 1/2" Swivel Casters
466	Stand w/1 Drawer, 3 Shelves, 16 3/4" x 18 3/4" x 31"H, (4) 4" Swivel Casters
469	Stand w/out Drawers, 3 Shelves, 16 3/4" x 18 3/4" x 31"H, (4) 4" Swivel Casters
474	Stand w/2 Drawers, 2 Shelves, 20 1/8" x 36 3/8" x 35"H, (4) 4" Swivel Casters



339



339S

DETECTO DUAL-READING PHYSICIAN'S SCALE

A scale with pound and kilo readings on both front and back. Capacity is 400 lb x 4 oz (180 kg x 100 grams). Overall height 59". Platform measures 10 1/2" x 14 1/2". Height rod adjusts 30" to 78" (76 to 198 cm). White only on all models. Rapid Ship no longer available.

337	Dual Reading Scale w/o height rod
338	Portable Dual Reading Scale w/height rod & wheels
339	Dual Reading Scale w/height rod
349	Dual Reading Scale w/height rod & handpost
337S	Dual Reading Scale w/out height rod, Stainless Steel
339S	Dual Reading Scale w/height rod, Stainless Steel

Scales

**DETECTO PHYSICIAN'S SCALE**

Heavy duty scale maintains pinpoint accuracy. Capacity is 400 lb x 4 oz or 180 kg x 100 g. Overall height is 59". Platform measures 10 1/2" x 14 1/2". Height rod adjusts 30" to 78". Available in stainless steel and white.

- 2371 Metric Scale w/o height rod
- 2381 Portable Metric Scale w/height rod & wheels
- 2391 Metric Scale w/height rod
- 2491 Metric Scale w/height rod & hand post
- 437 Scale w/out height rod, White
- 438 Portable Scale w/height rod & wheels
- 439 Scale w/height rod
- 448 Scale w/handpost, height rod & wheels
- 449 Scale w/height rod & hand post
- 2371S Metric Scale w/out height rod, Stainless Steel
- 2391S Metric Scale w/height rod, Stainless Steel
- 437S Scale w/out height rod, Stainless Steel
- 439S Scale w/height rod, Stainless Steel



6855




6868

DETECTO PORTABLE HIGH CAPACITY DIGITAL SCALE

Tubular steel, wrap-around handrail. 7" high-contrast LCD display. 6 C size alkaline batteries (included). Capacity: 600 lb x .2lb/270 kg x .1kg. Nonslip platform size: 18"L x 14"W x 1.75"H. White.

- 6855 Portable High Capacity Digital Scale, 600 lb x .2 lb/270 kg x .1 kg
- 6868 Portable High Capacity Digital Scale w/flip-up seat, 800 lb x .2 lb/360 kg x .1kg

**SECA 634 BARIATRIC FLOOR SCALE**

Specially designed for weighing adipose patients. The very high capacity of 300kg and the 60 x 60 centimeter large platform offer a solid base even for larger people. Thanks to the low level platform of only 45 millimeters, the person can effortlessly step onto the nonslip surface. The cable remote display can be placed at eye level or on a table. The graduation of 100g enables recognition of even the smallest weight change, so that the success of diet programs can be precisely determined at any time.

6341321004 Bariatric Floor Scale

**SECA 674 DIGITAL PLATFORM SCALE**

A multi-functional electronic flat scale designed especially for people needing assistance or heavyweights. Thanks to its extremely flat and large platform with built-in ramp, weighing people in wheelchairs or seated during dialysis is no problem at all. Comes with transport castors so it is easy to move. Capacity: 800 lbs/300 kg.

67413210009 Electronic Platform Scale w/high capacity, Can be Transported on Casters

**SECA 674C DIGITAL PLATFORM SCALE**

Represents a step up from the 674 through the incorporation of Smartcare technology. Each unit comes with three holders and cards for use by practitioners in accounting for wheelchair weight. It is no longer necessary to remember the unique weight of each wheelchair, just simply insert the appropriate Smartcard into the scale, weigh the patient and chair, and only the patient's weight appears. Capacity: 800 lbs/300 kg.

67425211139 Electronic Platform Scale w/high capacity, Can be Transported on Casters, Smartcard w/3 holders and cards

Get The Best Seat In The House

Choose from a Full Line of Quality Medical Stools from Clinton

From classic chrome examining stools, contemporary, pneumatic stools with dual nylon casters, extra tall models, to contour seats with comfortable adjustable backrests and hands-free stools with specialty armrests, Clinton has them all. And all are Clinton built for comfort and durability for years of comfortable performance. **Ask about Clinton's full line of quality stools.**



Examine Real Quality™



2145-SP
Hands-Free Stool
with Armrest

- 16" seat
- Height Adjustable
19 $\frac{1}{2}$ " to 24 $\frac{1}{2}$ "

2165-W
Ergonomic Designed
Task Chair

- 18" seat
- Height adjustable
17" to 22 $\frac{1}{2}$ "



2102
Chrome Base Stool

- 14" seat
- Height Adjustable
18 $\frac{1}{2}$ " to 25 $\frac{1}{2}$ "



2135-S
5-Leg Pneumatic Stool

- 16" seat
- Height adjustable
18" to 23"



Detecto Model
439



EYE-LEVEL PHYSICIAN SCALES

Detecto continues to set the pace in the clinical scale industry, backed by over a century of leadership. Now Detecto offers you more weighing ability for your money! The popular weighbeam physician scale line has been increased to a 400 lb / 180 kg capacity. Detecto's model 439 features a height rod in inches and centimeters. The 10.5" x 14.5" platform cover is easily removed for cleaning. The die-cast beam can be easily read from either side.




Detecto Model
6027



WAIST-HIGH DIGITAL SCALES

Detecto's model 6027 waist-high scale features an 11" x 13" platform and is 38" high. Detecto digital medical scales come standard with the enhanced 758C weight display. The 758C display includes a standard serial port for data output to a printer or a PC-based EMR (electronic medical record) software program. With the new time/date feature, the scale can send the weight, time, date, BMI (Body Mass Index), and height directly to a printer or a PC.

Detecto Model
6880



ROLLING CHAIR SCALES

Zero turning radius ability is the hallmark of this lightweight rolling chair scale. Ergonomics dictated this stylish, battery-powered chair scale's contemporary design. This system comes complete with a large LCD display. It has a 440 lb capacity that is sensitive to 0.2 lb. Oversized 3" diameter wheels and sleek styling allow this scale to be extremely portable. The soft, padded armrests fold out of the way to allow patients easy accessibility.

Alko Enterprises

Celebrating 50 Years

of Innovation & Craftsmanship

Hausmann Industries is proud to celebrate our 50th year of manufacturing products for the medical industry. It all started in a small shop in 1955 with the introduction of a basic line of 'H' Brace Tables. Today, with a product offering of more than 400 items, over 90% of the products are manufactured in our high-tech New Jersey based facility featuring the latest in computer operated wood-working machinery. Manufacturing technology has changed over the years, Hausmann's original philosophy and dedication to quality and service has remained steadfast. Hausmann has worked hard to establish a reputation for reliability and performance, and constantly strives for continued improvement before attaching its name to any of its products.



OBESITY IS A DISEASE DECLARES US GOVERNMENT*

The US Government says it will now recognize obesity as a disease. In a statement, Health and Human Services Secretary Tommy Thompson said, "Obesity is a critical public health problem in our country that causes millions of Americans to suffer unnecessary health problems. With this new policy, Medicine will

Are You Equipped?

With the Hausmann's line of Bariatric Equipment designed specifically for obese patients, you will be prepared today for tomorrow's patient.



Model 4795



Model 4074



Model 4086

**POWER
Backrest
Standard!**



Model 4416

*US Department of Health and Human Services Press Release, July 15, 2004.

Preferred by Cardiologists: GREAT PRICE, UNIQUE DESIGN, HAUSMANN QUALITY

POWERMATIC® ECHO-SCAN™ TABLE

features:

- Power height adjustment from 20" wheelchair height to 29" high
- Two direct-drive lift mechanisms
- 10" wide x 12" deep drop-away section on durable piano hinge
- 400 lbs. weight capacity

www.hausmann.com

Model 4790 Drop-Away "Open".
Shown with Option 05 - Headrest.



Alko Enterprises

Nicolet Delivers a Variety of Solutions for all your Doppler Needs!



Nicolet Obstetric Dopplers

- Unmatched accuracy and reliability
- Interchangeable 2 & 3 MHz OB probes
- Crystal-clear sound
- Versatile options: digital display, rechargeable batteries, 2 MHz waterproof probe
- 5-year parts, 1-year labor warranty

Nicolet Vascular Dopplers

- Monitor hard-to-find pulses
- Measure low systolic blood pressures
- Locate vessels for catheterization
- Perform ankle-brachial index (ABI) exam
- Interchangeable 5 MHz probe & 8 MHz pencil style probe
- Reimbursement available with VersaLab bi-directional Doppler systems



VIASYS[®]
HEALTHCARE
Excellence For Life

www.viasyshealthcare.com

Alko Enterprises



SECA 700 MECHANICAL BEAM SCALE

The seca 700 mechanical column scale with its eye-level balance beam sets new standards worldwide - for capacity, graduation, and quality. Each scale has a cast iron base, iron column, and easy-to-use interface. Moreover, the compact packaging helps keep freight costs to a minimum. The included seca 220 measuring rod converts the scale into a complete space-saving measuring and weighing device. Available in kg, lbs or kg/lbs versions.

7001021994	Physician Scale w/height rod & wheels, 160kg Capacity, Kilograms Only
7001121004	Physicians Mechanical Balance Beam Scale w/out height rod, 500lb Capacity, Lbs Display
7001121994	Physician Scale w/height rod & wheels, 500lb Capacity, Lbs Only
7001321994	Physician Scale w/height rod & wheels, 500lb Capacity, Lbs/Kg

SECA 755 MECHANICAL COLUMN SCALE WITH BMI DISPLAY

The seca 755 is one of the first mechanical column scales which, in addition to displaying the weight, also displays the respective nutritional condition of the patient with the BMI function. When the patient mounts the scale, a "BMI disc" of contrasting colour is set in motion below the very precise weight display in the large, clear dial. It is visible in a small window next to which a height scale is printed. Once the weight has been determined, it is easy to read the BMI range which, with contrasting colours, is split up into underweight, normal, and overweight as well as adiposity (abnormally overweight). A flat tread area and practical transport castors are further advantages of the seca 755. It can also be converted into a complete measuring and weighing station by fitting the optional measuring rod seca 220.

7551321994	Mechanical Column Scale with BMI Display
------------	--

SECA 763 DIGITAL WEIGHING AND MEASURING STATION W/AUTOMATIC BMI CALCULATION

With the innovative, digital measuring station seca 763, both the height and weight of a patient can be ascertained in just one step. After measuring and weighing, both values are displayed on two LCD displays integrated into the intuitively operable display. But that is not all; on the basis of these values, the seca 763, at the touch of a button, automatically calculates the BMI of the patient. This value provides information on the general nutritional condition of the patient on the basis of the recognized WHO standard. With its capacity of 550lbs, its low platform and very robust design, the measuring station is also extremely suitable for heavyweight patients.

7631321139	Digital Weighing and Measuring Station
------------	--



SECA 767 DIGITAL COLUMN SCALE

Does not need much energy. With no more than one set of batteries this column scale can weight up to 16,000 times. The big liquid crystal display not only shows the weight but also calculates the patient's Body-Mass-Index. Extend the scale's functionality by fitting it with the optional Seca 200 measuring rod, so the medical practitioner can weigh and measure in just one step. Capacity: 440 lbs/200 kg.

76713210009	Electronic Personal Scale w/large backlit LCD and Low Power Consumption
-------------	---



DETECTO WAIST-HIGH DIAL TYPE SCALES

White. Platform: 11" x 13". Capacity: 300 lb x 1 lb or 150 kg x .5 kg. Dial: single revolution, 7 1/2" reading line diameter, 9" overall diameter.

080	Scales w/out height rod
080KG	Metric Scales w/o height rod

Scales



057



059



6027



6029

DETECTO WAIST-HIGH PERSONAL FITNESS SCALES

White. Platform size: 11" x 13".

- 057 Scale w/o height rod, 340 lb x .5 lb
- 059 Scale w/height rod, 340 lb x .5 lb
- 0591 Metric Scale w/height rod, 150 kg x 200 g
- 6027 Digital Dual Scale w/o height rod, 400 lb x .5 lb/180 kg x .2 kg,
Powered By 6 "C" Cell Batteries or Optional AC Adaptor
- 6027KG Digital Scale w/o height rod, 180 kg x .2 kg, Powered By 6 "C" Cell
Batteries or Optional AC Adaptor
- 6029 Digital Dual Reading Scale w/height rod, 400 lb x .5 lb/180 kg x .2 kg,
Powered By 6 "C" Cell Batteries or Optional AC Adaptor
- 6029KG Digital Scale w/height rod, 180 kg x .2 kg, Powered By 6 "C" Cell
Batteries or Optional AC Adaptor

Accessories

- 4PHTROD-2 Height Rod For Mechanical Waist High Scales)
- 4PHTROD-3 Height Rod For 6029 & 6029KG)

**DETECTO WAIST-HIGH STAND-ON SCALE**

Heavy-duty, yet lightweight. Engineered specifically for weighing obese and unsteady patients and features a 24" square. Low-profile platform with non-skid surface and sturdy handrails with padded grips. Digital indicator with easy-to-use quick weigh keys, plus Body Mass Index calculation. Auto-zero function. Powered by either batteries or AC adapter. Auto-off function.

- 6856 Stand-On Scale, 800 lb x .2 lb/360 kg x .1 kg

**SECA 644 DIGITAL HAND-RAIL SCALE**

Low profile platform makes it easy to use and the curved railing provides a stable anchor for unstable patients. With its locking castors and swiveling control/display panel, it provides the medical staff with a high degree of flexibility. Can be transformed into a complete measuring and weighing station with the addition of Seca 242 or Seca 222 measuring rod options. Capacity: 660 lbs/300 kg.

- 6441321139 Electronic Scale w/standing assistance, High Capacity and Large Platform

**SECA 703 DIGITAL COLUMN SCALE WITH VERY HIGH CAPACITY**

Extremely robust, precise and stable: With a capacity of 550lbs, a large but low platform and the integrated, BMI function, this column scale has the prerequisites for the diagnosis and therapy of heavyweight patients. The durable construction of cast iron can cope with the most extreme loads. The energy-saving battery operation - up to 16,000 weighings with just one set of batteries - means that this scale can be used anywhere as it is not dependent on AC power. The scale can be fitted with the integral measuring rod sec 220, thus facilitating weighing and measuring in one step.

- 7031321009 High Capacity Digital Scale w/height rod and wheels, Model 703

**SECA 780 DIGITAL COLUMN SCALE**

New technology gives longevity without the need for changing batteries or the inconvenience of outlets. Easy to operate, stylish and priced just right. Height rod can be attached. Capacity: 400 lbs/180 kg. Dimensions: 12"W x 35"H x 15 1/4"D.

- 7802321134 Digital Scale w/energy saving battery, 2 Wheels



450



459

DETECTO HEAVY DUTY PEDIATRIC SCALE

Heavy duty lever and beam construction assures long life and durability. Well suited for heavy usage. Overall size is 27"H x 22"W x 21"D. Equipped with sturdy, washable polystyrene tray. White.

450 Pediatric Scale, 130 lb x 1 oz

451 Pediatric Scale, 65 kg x 20 g

459 Dual Reading (Pediatric Scale, 40 lb x 1/2 oz, 17.5 kg x .01 kg)



DETECTO MODEL 8440 INFANT SCALE

Provides accurate, safe weighing of infants and features selectable weight units, lb-oz or kilograms, display hold, and a memory function. Compact, battery powered and portable, perfect for visiting nurses, home or wherever space is limited and performance is paramount. 44 lb x 1/2 oz / 20 kg x 10 g. 4 3/4 digits LCD 1" high. Time and date clock function. Memory for comparison of current and stored weight. 6V DC uses 4 AA batteries or optional AC adapter (not included).

8440 Infant Scale



AP-4KD




PS5-2KD

DETECTO WET DIAPER SCALE/LAP SPONGE SCALE

Fluid output scale has a 4000 gram capacity. Accurate to one gram. For adult incontinence and infant diapers. Red digital 1/2" display. Push buttons simplify accurate readings. Features a roomy, removable stainless steel tray for fast weighing and simple cleanups. Portable, 4.5 lb model. 7"L x 7 1/2"W x 2 3/4"H. AC adapter with 6' power cord.

AP-4KD Diaper Scale, 4000 g x 1 g

PS5-2KD Diaper Scale, 1998 g x 2 g, Enamel Finish



SECA 334 INFANT SCALE

Represents maximum safety and accuracy when weighing babies. Although its weighing tray is designed to be far more generous and comfortable than other models, at only 7 lbs itself, it requires amazingly little space. Has a large, clear LC display and practical handle, which can also be used to hang the scale. Capacity: 44 lbs/20 kg. Dimensions: 25"W x 4"H x 11 3/4"D. Power Supply: battery.

33413210004 Digital Baby Scale w/practical handle

SECA 354 DIGITAL BABY SCALE - 44LB CAPACITY

With their dual function, these baby scales are a small sensation. At the touch of a button, the tray can be detached quickly and easily from the base, leaving you with a robust floor scale for children weighing up to 44lbs. As the two elements are securely connected with one another, it is absolutely safe to lift and transport the scale. The soft, round shape of the tray gives the baby that so very important feeling of safety and security. With the breast-milk-intake function, it is easy to determine the exact amount of milk the baby has taken when being breast-fed.

3541317009 Digital Baby Scale w/detachable tray, Tight graduation

SECA 364 DIGITAL BABY SCALE - 110LB CAPACITY

At the touch of a button, the tray can be detached quickly and easily from the base, leaving you with a robust floor scale for children weighing up to 110lbs. As the two elements are securely connected with one another, it is absolutely safe to lift and transport the scale. The soft, round shape of the tray gives the baby that so very important feeling of safety and security. With the breast-milk-intake function, it is easy to determine the exact amount of milk the baby has taken when being breast-fed.

3641317009 Digital Baby Scale w/detachable tray

SECA 374 DIGITAL BABY SCALE W/EXTRA LARGE WEIGHING TRAY

The shell-shaped design of the seca 374 is consistent with the high demands made on its multiple functions. Suitable for babies and infants weighing up to 40lbs, the large tray provides maximum safety against the baby falling from the scale. The clear controls of the user interface are located on the slightly raised, ergonomically formed side of the scale facing the user. A large LCD display makes it quick and easy to determine the weight. A further advantage of the seca 374 is that it can also be operated using batteries and thus be used anywhere. The scale's extremely robust, hygienic surface of high-quality plastic is easily and quickly cleaned.

3741321009 Shell Shaped Digital Baby Scale

Scales

SECA 725 INFANT SCALE

Classic design and reliable construction make this mechanical baby scale a favorite for over 50 years. Rugged steel construction and cast levers ensure years of durability and accuracy. Removable tray makes transport and cleaning simple and convenient. Capacity: 32 lbs/16 kg. Dimensions: 21 3/4"W x 7 1/2"H x 12"D.

7251121139 Mechanical Baby Scale w/sliding weights

**SECA 727 INFANT SCALE**

A large tray and three-button operation in HOLD, TARE, lbs/kg ensures convenience, precision, comfort and safety. An integrated damping system compensates for movement, and the PC-friendly RS 232 interface eliminates the need for manual input. Capacity: 44 lbs/20kg. Dimensions: 21 1/2"W x 5 1/2"H x 12 1/2"D. Power Supply: 6 AA rechargeable batteries, A/C UL adapter/recharger.

7271321844 Electronic Baby Scale w/digital display and Serial Interface RS232



4250000009



7281321834

SECA 728 INFANT SCALE

Specially designed for situations where the smallest weight fluctuation can be critical. From the integrated digital display to the unique tray design, this baby scale has it all. With the seca 231, the scale becomes a complete weighing and measuring station. Capacity: 44 lbs/20 kg. Dimensions: 21 1/2"W x 6 1/2"H x 12 1/2"D. Power Supply: 6 AA rechargeable batteries, A/C UL adapter/recharger.

7281321834 Electronic Baby Scale w/digital display and Serial Interface RS233

Accessories

4250000009 Carrying Case w/rollers for Seca 725, 727, and 728

**SECA 851 DIAPER SCALE**

Will keep tabs on everything in a baby's diaper down to two grams. Its lightning quick digital readout keeps messy job time to a minimum. And it is as easy to clean as it is to store. Capacity: 4 lbs = 6 oz/ 2 kg. Dimensions: 5 3/4"W x 1 1/2"H x 9 1/2"D. Power Supply: Battery operated (9V included).

8511317009 Electronic Diaper Scale



475



6475K



6875



6876

DETECTO CHAIR SCALE

Ideal for weighing frail, handicapped or elderly patients. Can be used in relatively small spaces on carpet or tile. Lift away arms for easy patient transfer; lift away footrest and oversized wheels for maneuverability. Die-cast beam is easy to read from both sides. 42"H. Black/white. Powered by 6 C cell alkaline batteries or optional AC adaptor.

475 400 lb x 4 oz Chair Scale


4751 Chair Scale Dual Reading Beam, 400 lb x 4 oz (180 kg x 100 g)

6475K Digital Chair Scale, Metric Reading, 180 kg x 1 kg

6875 Euro Chair Scale, Digital Dual Reading, 400 lb x 2 lb/180 kg x .1 kg, Concealed Wheels

6475 Digital Chair Scale, Dual Reading, 400 lb x 2 lb (180 kg x .1 kg)

6876 Flip Seat Euro Chair Scale, 600 lb x 2 lb/270 kg x .1 kg

**DETECTO PORTABLE WHEELCHAIR SCALE**

Fold up scale with four way, dual direction wheels for easier transport. Ambulatory patients hold padded handrails during weighing for added security. Digital display rotates for wide angle viewing. Keypad entry tare, auto power off, auto zero, motion detector, weight lock. 800 lb x 2 lb/360 kg x .1 kg capacity. 28"L x 32"W platform with non-skid mat. Digital indicator with LCD display. Powered by 6 C cell batteries. Optional 115VAC power adapter.


6550 Wheelchair Scale



DETECTO ROLLING CHAIR SCALE

Zero turning radius ability is the hallmark of this lightweight rolling chair scale. Ergonomics dictated this stylishly, battery-powered scale's contemporary design. This system comes complete with a large LCD display and it has a 440 lb capacity that is sensitive to .2 lb. Oversized 3" diameter wheels and sleek styling allow this scale to be extremely portable. The large 0.7" high contrast LCD display is battery powered and is easily readable in the brightest lighting conditions. The soft, padded armrests fold out of the way to allow patients easy accessibility. The 6880's footrests also fold away to promote patient security when getting into or out of the scale. Safety stops on rear wheels for parking in one place without rolling. Contoured, sloping backrest allows the patient to assume a relaxed position.

6880 Rolling Chair Scale



DETECTO SEMI-PORTABLE WHEELCHAIR SCALE

White. Capacity: 1000 lb x 0.5 lb/450 kg x 0.2 kg. Overall dimensions: 32 1/2"L x 43"W x 4 1/2"H. Platform: 32 1/2"L x 30"W x 2 1/2"H. Smooth deck with rubber mat and integral wheel stop. Display: 5-digit, 7 segment, .7" high LCD for wall or table mount. Power source: 6 C batteries or AC power adaptor.

CR-1000D Semi-Portable Wheelchair Scale
CR-500D Wheelchair Scale, Same as CR-1000D Except Capacity is 500 lb x 0.2 lb/225 kg x 0.1 kg



DETECTO STATIONARY GERIATRIC WHEELCHAIR SCALE

Capacity: 1000 lb x .2 lb/450 kg x .1 kg. 7" high-contrast LED display with keypad tare. Powered by AC adaptor. Smooth steel platform covered with black rubberized mat. 24" ramp. 25' cable. Adjustable platform height 2 3/4" to 3 1/4". Blue/black.

FHD-133-II Geri-Chair Scale w/36" x 36" platform
FHD-144-II Geri-Chair Scale w/48" x 48" platform



495



6495

DETECTO STATIONARY WHEELCHAIR SCALES

Platform size: 30"L x 26"W x 2"H with nonskid mat. Padded handrail for ambulatory patients. Platform is designed to accommodate most all sizes of wheelchairs. Ramp provides smooth approach and exit and attached to either side of platform.

495 Mechanical Scale, Capacity 400 lb x 4 oz
4951 Mechanical Scale, Capacity 180 kg x 100 g
6495 Digital Scale, Capacity 400 lb x .2 lb/180 kg x .1 kg, LCD Display w/keypad, Powered By AC Adaptor



DETECTO WHEELCHAIR SCALES

Mechanical eye-level physician scale with 27"W x 30"D ramp, without height rod.

485 Wheelchair Scale, 350 lb x 4 oz Capacity
4851 Wheelchair Scale, 160 kg x 100 g Capacity




SECA 664 DIGITAL WHEELCHAIR SCALE

The high capacity of 800 lbs of this wheelchair scale means that even heavy patients can be weighed with ease. Thanks to the extremely flat platform, the scale is easy for wheelchair-users to mount and the side barriers ensure maximum safety. The integrated pre-TARE function determines the patient's net weight in seconds and a real advantage is its HOLD function, which continues to display weight in large, easy-to-read numbers even after the scale is no longer in use. This scale can be folded up in just a few moves and its transport castors make it easy to move. Capacity: 800 lbs/300 kg. Dimensions: 27"W x 22 1/2"H x 23"D.

6641321004 Electronic Wheelchair Scale, Collapsible, Transport by Wheels

Scales

**SECA 664C DIGITAL WHEELCHAIR SCALE**

All the quality and functionality of the Seca 664 with the added feature of Smartcard technology. Each unit comes with six holders and cards for use by practitioners in accounting for wheelchair weight. It is no longer necessary to remember the unique weight of each wheelchair brand, just simply insert the appropriate Smartcard into the scale, weigh the patient and chair, and only the patient's weight appears.

6642521109 Electronic Wheelchair Scale, Collapsible, Transport by Wheels, Smartcard w/3 holders and cards

**SECA 954 DIGITAL CHAIR SCALE**

Offers traditional design with added mobility, so it can be brought directly to the patient. Compact design and locking wheels make the chair maneuverable. Has swiveling armrests for easy seating and the desirable TARE function for taking out extra weight. Capacity: 440 lbs/200 kg. Dimensions after assembly: 27"W x 22 1/2"H x 23"D.

9541309803 Digital Chair Scale w/4 wheels and swivel armrests, Chair Scale Dismantled, Delivered w/tools



IB400



IB600

DETECTO IN-BED SCALES

Ideal for intensive, coronary, geriatric, surgical and burn care patients. Locking rear casters and base widening stance provide maximum safety. LCD digital indicator powered by 6 "D" cell alkaline batteries has auto sleep mode, auto power shut off, lock/release display, auto/manual zero, motion detector. Thermoplastic Staph Check™ stretcher fabric. Adjustable weighing arm and hydraulic lift mechanism. Optional 115 VAC power adapter.

IB400 In-Bed Scales, 400 lb x .01 lb

IB600 In-Bed Scales, 600 lb x 0.2lb/270 kg x 0.1 kg

Accessories

0046-C007-08 6 ft Adult Stretcher For IB400 Series

0046-C065-08 4 ft Pediatric Stretcher For IB 400 Series

0046-C247-08 6 ft Adult Stretcher For IB600



PL-400



PL-435



PL-475

DETECTO PATIENT LIFT SCALE - DIGITAL

Designed for patient weighing when used with the PL-425 Patient Lift and PL-450 Sling Seat. Keypad controls on/zero, off and lb/kg display select. .75" high-contrast LCD display with annunciator for lb/kg, zero, low battery & display lock. Powered by one 9V battery. Unit measures 4.125"W x 2.125"D x 4"H (center of eye to center of eye; 5" overall).

PL-400	Patient Lift Scale, 400 lb x .2 lb/182 kg x .1 kg
PL-425	Patient Lift Scale w/4.8"H underbed clearance, Hydraulic Lift Mechanism
PL-435	Patient Lift Scale, Sling Seat, Yoke & Suspension
PL-475	Complete 400 Series Patient Lift Scale, Includes PL-400, PL-425 & PL-450
PL-600	Patient Lift Scale, 600 lb x 0.2 lb/270 kg x 0.1 kg


**DETECTO PATIENT LIFT SCALE - MECHANICAL**

To be installed in between the boom and sling of a patient lift. Capacity: 400 lb x 1 lb. Overall dimensions: 24 3/4". Cast iron housing. 7" single dial. Aluminum face with black lettering and single revolution.

11S400H Patient Lift Scale



4340000009



9841321139

SECA 984 BED & DIALYSIS SCALE

Allows bed-ridden patients to be weighed simply, carefully and accurately and is an indispensable constituent of dialysis and intensive care. Integrated lifts allow four measuring platforms to be placed under the castors of the bed. The previously determined weight of the bed is tared using the pre-TARE function and the patient's precise weight accurately determined. The bed scale stores all the previously determined values and retains the information in the event of a power failure using rechargeable batteries. Package contains: Scale w/digital display, four loaded cells w/integrated mechanical lifters and equipment trolley. Capacity: 1100 lbs/500 kg.

9841321139 Electronic Bed and Dialysis Scale w/digital display and 4 Load Cells w/integrated mechanical lifters, Equipment Trolley

Accessories

4340000009 Additional Multifunctional Display for 984 Scale



DETECTO DR400-C PORTABLE DIGITAL SCALE

The perfect scale for those who need a lightweight and portable scale. 400 lb x .5lb/181 kg x .2kg. Platform is 12" square. Powered by 6 AA alkaline batteries or 115V AC adapter is included.

DR400-C Portable Digital Scale



D315



D1130

DETECTO PERSONAL SCALES

A low profile design, heavy duty construction and a non-skid vinyl mat are just a few features of the Detecto D1130. Other features include a white, baked-enamel finish and personal colorful indicators to mark one person's or a whole family's weights. The large diameter dial allows for easy readability, even for tall people. Designed for personal use, fitness centers, and health care facilities.

D315 Personal Scales, 330 lb x 1 lb Capacity, 18 1/2" x 13" x 4 3/4" w/7 3/4" dial, Weigh Tracker™ Indicators Mark Weight Progress, 5 Year Warranty
D1130 Personal Floor Scales, 300 lb x 1 lb Capacity, 11.5"W x 17.25"L x 3"H, w/7 1/4" dial, Weigh Tracker™ Indicators, Five Year Warranty



DETECTO PORTABLE SCALES WITH CARRYING HANDLES

The size of a small briefcase and incredible compact and light, the MD series portable scales are built with sturdy steel construction for durability. 250lb capacity by 1 pound graduation and with a built-in carrying handle, the MD series scales are easy to transport. They feature an easy-to-read dial with shatterproof lens, and the Lock/Hold feature is useful for oversize loads.

MD250PKG Portable Scale w/carrying handles, 13"W x 10.5"D (33cm W x 27cm D), Painted Mild Steel, 2/pkg



SECA 491 BMI CALCULATOR

A simple and quick way to calculate Body Mass Index. Patient data can be entered in pounds, feet and inches or centimeters and kilograms. Also functions as a standard calculator.

491 BMI Calculator



SECA 750 MECHANICAL FLOOR SCALE

Offers both rugged durability and a large, easy-to-read dial. Its design was based on a lower profile and wider platform. Capacity: 320 lbs. Dimensions: 12 1/2"W x 2 1/2"H x 18 1/2"D.

75010190009 Mechanical Personal Scale w/precision weighing mechanism in a classic white housing -KG only

7501119009 Mechanical Personal Scale w/precision weighing mechanism in a classic white housing -LBS only

75013190009 Mechanical Personal Scale w/precision weighing mechanism in a classic white housing -LBS/KG only

Scales

**SECA 754 NOVA SCALE**

Circular dial with nickel plated housing shows the weight in white on black. 150kg capacity.

7541135009 Flat Scale -LBS/KG

**SECA 760 COLORATA SCALE**

Brighten up with a splash of color from a selection of lively shades (white, red, yellow or blue). The housing and the platform with skin-friendly, grooved natural rubber mat are in matching colors. 150 kg capacity.

7601126009 Flat Scale White

**SECA 761 MECHANICAL FLOOR SCALE**

The large easy to read dial displays weight that's accurate to one percent of total body weight. Its rugged all-metal body is enhanced by a heavy-duty lever system. Capacity: 320 lbs. Dimensions: 11 3/4"W x 4 1/4"H x 18"D.

761-11 Mechanical Personal Scale w/large round dial -LBS only
7611119009 Mechanical Personal Scale w/large round dial -LBS only
7611319009 Mechanical Personal Scale w/large round dial -LBS/KG

**SECA 770 HEAVY DUTY ELECTRONIC FLOOR SCALE**

Made for serious traffic. The platform's precision sensor activates the scale and the automatic shut-off adds extra convenience. Battery operation and low power consumption make this the perfect portable scale. Capacity: 440 lbs/ 200 kg. Dimensions: 13"W x 5 1/2"H x 15"D.

7701321004 Electronic Personal Scale w/autohold function, Automatic switch off function, High Capacity

**SECA 841 ELECTRONIC PERSONAL FLOOR SCALE**

8411321009 Electronic Scale w/large digital display

**SECA 843 ELECTRONIC PERSONAL SCALE**

Ergonomically shaped platform and practical handle. Capacity 308 pounds. 14.8"L x 11.8"W x 1.8"H. Battery power supply.

8431321009 Electronic Personal Scale w/ergonomically shaped platform and handle

**SECA 880 ELECTRONIC PERSONAL FLOOR SCALE**

Designed for the private sector. The quick-touch turn on and long-life batteries offer the utmost in operational ease. Carry case #421 available separately. Capacity: 400 lbs/ 180 kg. Dimensions: 12"W x 2 1/2"H x 13"D.

8801321139 Electronic Personal Scale



SECA 881 DIGITAL FLOOR UNISCALE

Ideal if you're looking for weighing tiny tots to grown-ups. Combine the large LC display and activation sensor with the convenient TARE function and you have the seca multipurpose uniscale at your feet. Capacity: 400 lbs/180kg. Dimensions: 12"W x 2 1/2"H x 13"D. Power Supply: 4 AA alkaline batteries.

8811321139 Electronic Personal Scale w/mother-child function



SECA 882 BODY-MASS-INDEX ELECTRONIC FLOOR SCALE

Recognized by the US Department of Health as the standard for normal values of height and weight. The remote head and handle make the BMI easy to operate and move. Carrying case #421 is available separately. Capacity: 400 lbs/180 kg.


8821321134 Electronic Personal Scale w/remote display on cable



SECA 884 PERSONAL ELECTRONIC SCALE

Flat scale equipped with cable remote control, practical carrying handle and battery operation.


8847021099 Certifiable Electronic Personal Scale w/remote display on cable, Additional Functions



DETECTO DIGITAL FLOOR SCALE

Detecto's reliable DR400-758C digital platform scale is lightweight and portable, making it perfect for mobile clinics, home care nurses, or wherever a portable scale with Body Mass Index is needed. The remote digital indicator comes with a wall-mount bracket and easily mounts on a counter top or wall. The indicator features a large .7" / 18mm high-contrast LCD, selectable units, Body Mass Index, pushbutton and keypad tare, and automatic display lock. The 758C quickly calculates Body Mass Index by simply entering the patient's height while weighing. The DR400-758C is battery powered, so you can take the scale anywhere you need it. The unit has a slip-resistant mat to ensure patient safety when getting on and off the scale.

DR400-758C Digital Floor Scale



VET-400



VET-50

DETECTO VETERINARY DIGITAL SCALES

Stainless steel. 7" high-contrast LCD display. 6 C size alkaline batteries (not included).

VET-400 Digital Scales, Capacity: 400 lb x .2 lb/180 kg x .1 kg, 22" x 42" x 1 3/4", 8 ft Coil Cable, Optional Rubber Mat
VET-50 Digital Scales, Capacity: 50 lb x .02 lb/20 kg x .01 kg, 12" x 14" x 3", 10 ft Cable

Accessories

VET-MAT Rubber Mat



SECA 606 ELECTRONIC VETERINARY SCALE

An extremely accurate and dependable scale for weighing pets. The tare function allows a small animal to be weighed in a basket. Simply place a basket on the scale, tare out the basket and place the pet in the basket. Only the pet's weight will show. The dampening function ensures reliable weighing to calm animals, restless animals and very restless animals. The rough surface consists of scratch-resistant plastic and is very easy to clean. The anodized frame construction and the sealed wooden plate at the bottom of the scale make the scale highly moisture resistant when the floor in the vet practice or clinic becomes wet. The intelligent auto-hold function will continue to display the animals weight after it has left the scale. The battery capacity is sufficient for 16,000 weighings. The scale turns off in two minutes when using batteries. 220 lbs capacity. Platform size is 35.4"L x 23.6"W x 3"H.

6061321009 Electronic Veterinary Scale

Scales



285R31



AP-10



AP-4K



AP-6



PS-5A



PS-6A



PS-11



PS-7

DETECTO ELECTRONIC DIGITAL PORTION CONTROL SCALE

Provides accurate, consistent portion control for food preparation and dispensing. Choose a beige, oven-baked enamel base with stainless steel platter or a complete stainless steel unit. Weighs in ounces or grams. Capacity: 70 oz x .1 oz or 1998 grams x 2 grams. Easy-to-read digital display. Operates on four AA batteries. 7"W x 7"D x 2 3/4"H.

AP-10	Portion Control Scale, Capacity: 10 lb - 9.995 lb x .005 lb, Stainless Steel
AP-20	Portion Control Scale, Capacity: 20 lb - 19.99 lb x .01 lb, Stainless Steel
AP-4K	Portion Control Scale, Capacity: 4 kilo - 3999 g x 1 g, Stainless Steel
AP-6	Portion Control Scale, Capacity: 6 lb - 99.95 oz x .05 oz, Stainless Steel
AP-8	Portion Control Scale, Capacity: 8 lb - 7.998 lb x .002 lb, Stainless Steel
PS-5A	Portion Control Dual Scale, Beige, w/AC Adapter
PS-6A	Portion Control Scale, 70 oz x 1/4 oz Capacity Only, Beige, w/AC Adapter
PS-11	Portion Control Scale, Capacity: 11lb x 0.1oz (5000g x 2g), Stainless Steel w/stainless steel commodity tray
PS-7	Portion Control Scale, Capacity: 7lb x 0.1oz (3000g x 1g), Molded Poly Enclosure w/stainless steel commodity tray

Accessories

285R31	Commodity Scoop
--------	-----------------

**DETECTO MECHANICAL 6" DIAL TYPE PORTION SCALES**

6" diameter dial. 5 3/4" square stainless steel platform. All models have temperature compensating springs. 6 1/4"W x 7 1/2"D x 7 1/2"H.

PT-1	Portion Scales, Fixed Dial, Painted Finish, 16 oz x 1/4 oz
PT-1000RK	Portion Scales, Rotating Dial, Painted Finish, 1000 g x 5 g
PT-1000SRK	Portion Scales, Rotating Dial, Stainless Steel, 1000 g x 5 g
PT-2	Portion Scales, Fixed Dial, Painted Finish, 32 oz x 1/4 oz
PT-2-SR	Portion Scales, Rotating Dial, Stainless Steel, 32 oz x 1/4 oz
PT-2R	Portion Scales, Rotating Dial, Painted Finish, 32 oz x 1/4 oz
PT-500RK	Portion Scales, Rotating Dial, Painted Finish, 500 g x 2 g
PT-500SRK	Portion Scales, Rotating Dial, Stainless Steel, 500 g x 2 g

**DETECTO MECHANICAL 8" DIAL TYPE PORTION SCALES**

8" diameter dial. 9" x 9" stainless steel platform. 8 1/2"W x 10 1/4"D x 10 1/2"H. All models have temperature compensating springs. Capacity: 32 oz x 1/8 oz.

T1-KP	Portion Scales, Fixed Dial, Enamel Finish, 1kg x 5g/2.2lb x 1/4 oz
T2-KP	Portion Scales, Fixed Dial, Enamel Finish, 2kg x 10g/4.3lb x 1/2 oz

**DETECTO MODEL 8450 DIGITAL BABY SCALE**

20kg x .01kg/44lb x .02lb/704 oz x .5 oz. Uses a 9V battery. 22"L x 13"W x 2 3/4"D mold tray. Five digit, 1/2" green vacuum fluorescent display.


8450	Digital Baby Scale, Black
------	---------------------------




DETECTO CERTIFIER SCALE - LEGAL FOR TRADE

Special purpose eye-level beam scale with tolerances for certification by weights and measures officials for athletic sporting events.

442 Certifier Scale, Capacity: 350 lb x 4 oz or 140 kg x 100 g, Platform Size: 13 3/4" x 18 1/2", White



GP-400-204



GP-400-758C

DETECTO GENERAL PURPOSE PORTABLE SCALE

Mild steel platform with slip-proof surface, built-in carrying handle.

GP-400-204 18"L x 14"W x 2"H, Digital Scale w/carrying handle, 1" High Contrast LCD, Legal For Trade

GP-400-758C 18" x 14" x 2", Digital Scale, 7" High-Contrast LCD Display, 400 lb Capacity, 6 "C" Cell Alkaline Batteries



3P7030X



728R121



8440-Case



8450-Case



DR-400C-Case

DETECTO SCALE ACCESSORIES

3P7030X	Counterweights, Increases Capacity From 350 lb to 450 lb For Models 439, 438, 437 & 449 Only
3P7039X	Counterweights, Increases Capacity Form 350 lb to 450 lb For Models 337, 338 & 339 Only
3P8030X	Counterweights, Increases Capacity Form 140 kg to 180 kg For Models 2381, 2391 & 2491 Only
3PHNDPST	Handpost (for replacement only)
3PHTROD-1	Height Rod
3PRAMP	Wheelchair Ramp
3PWHL	Wheels For Scale Portability
728R121	115V AC Power Adapter
8440-Case	Carrying Case For 8440
8450-Case	Carrying Case For 8450
DR-400C-Case	Carrying Case For DR-400C



SECA 202 MECHANICAL TELESCOPIC MEASURING ROD

Limited to the essentials, the measuring rod seca 202 is a simpler version of the measuring rod seca 222, but, of course, still suitable for obtaining high precise results.

2021714009 Measuring Rod

Scales

**SECA 207 PORTABLE SLIDER MEASURING DEVICE**

A fast and easy infant measuring rod made of anodized aluminum for long-lasting accuracy and durability. The calipers conveniently fold in for transport or storage. Dimensions: 1 7/8"W x 1 3/8"H x 41"D.

2071814009 Measuring Rod for Infants w/extra wide measuring slide

SECA 210 MEASURE MAT II

The simple and quick way to measure a baby's growth from day one. The mat's soft touch, lightweight rubber is gentle enough for baby, yet portable enough for easy storing and durable enough for wall mounting. It is extremely handy in delivery rooms where newborns are measured upside down. Dimensions: 11 3/4"W x 5 1/2"H x 53 1/2"D.

2101821009 Measuring Mat For Infants & Small Children

**SECA 212 BABY BAND**

Made of durable, non-stretch tape, measure a baby's head circumference in inches or centimeters. Reverse side features a mid-point measure for greater convenience.

2120 Non Stretch Teflon Tape to Measure Head Circumference on Infants and Small Children

SECA 220 TELESCOPIC HEIGHT ROD

The headpiece is secured on both sides for accuracy and folds down for a greater degree of safety. When added to a scale, it becomes a precision weighing and measuring package. Dimensions: 2 1/8"W x 30"H x 1 1/4"D.

2201814004 Telescopic Measuring Rod For Scale Models 766, 767, 780, & 700

**SECA 240 MEASURING ROD (STADIOMETER)**

Mounted onto the wall and the exact measuring result is being displayed in the small viewing window of the wide measuring slide.

2401814009 Wall Mounted Stadiometer

**SECA 404 GRAPH EASE**


A convenient plotting aid for growth and weight percentile charts. Just line up the vertical and horizontal lines on your growth chart and mark the point through a convenient hole at the line intersection. Sheet Size: 8 1/2" x 11" with rounded corners and 3 hole punch. Packaging: 25/pack.

404 Graph Ease, 25/pk


**SECA 416 WEE LENGTHS MECHANICAL INFANTOMETER**

Made of a sturdy, unbreakable hardwood for durability and reliable accuracy. The foot piece is channeled on an aluminum sideboard for smooth wobble-free operation. A handhold enhances transportability. Dimensions: 43 1/2"W x 6 1/2"H x 12"D.

416 Baby Board



405B




405G

SECA GROWTH CHARTS

These charts for boys (blue) and for girls (pink) are the fast, convenient way to track length/height and weight percentiles. Sheet size: 8 1/2" x 11". Packaged 100 sheets per pad.

405B Growth Chart, Boys, Plots Newborns to 36 Months
 405G Growth Chart, Girls, Plots Newborns to 36 Months
 406B Growth Chart, Boys, Plots 2-20 Years
 406G Growth Chart, Girls, Plots 2-20 Years




SECA HEIGHT MEASURING DEVICES

216	Plastic Wall Mounted Height Rod, Lock Knob, 24" - 83" Measuring Range
2221714009	Telescopic Measuring Rod, Wall Mounted w/extra wide slide -cm
2221814009	Telescopic Measuring Rod, Wall Mounted w/extra wide slide -inches
22218149999	Full Length Stadiometer for 644
2251821009	Mobile Measuring Unit on Rollers for Transport
2311717009	Attachable Measuring Rod for Baby Scales #728 and #736
23218170009	Attachable Measuring Rod for Baby Scale #334
2421814009	Electronic Wall Mounted Measuring Rod, Wireless Data Transmission
2421814999	242 Digital Stadiometer for 664

BARJAN SIGNAGE

All of our signage is created from high quality materials to visually present graphics and messages. A.D.A. Indoor/Outdoor. Directional/Way finding. Directory signage created in many sizes, shapes and colors to suit any application or budget. Logos, Subsurface, Vinyl Letters, Braille, Raised Letters, Silk Screening, Engraving, Banners, Flags.


Signage Signage In Several Different Styles



CLINTON EXAM ROOM & OFFICE SIGNS

All animal character exam room signs are available with numbers 1 through 12 and letter signs A through L. Sold individually. Each image is pre-mounted on ridged black plastic that easily mounts to any smooth flat surface with the supplied velcro strips.

EX-21	Exit Sign
EX-23	Doctors Office Sign
EX-24	Sign-In Sign
EX-25	Supplies Sign
EX-26	Waiting Room Sign
EX-27	Lavatory Sign
EX1	Exam Room Number 1 Sign
EX10	Exam Room Number 10 Sign
EX11	Exam Room Number 11 Sign
EX12	Exam Room Number 12 Sign
EX2	Exam Room Number 2 Sign
EX3	Exam Room Number 3 Sign
EX4	Exam Room Number 4 Sign
EX5	Exam Room Number 5 Sign
EX6	Exam Room Number 6 Sign
EX7	Exam Room Number 7 Sign
EX8	Exam Room Number 8 Sign
EX9	Exam Room Number 9 Sign



OMNIMED BEAM® 6 FLAG SIGNAL SYSTEM


Features 6 highly visible flags in 8" or 4" styles. Flags swing 180° to identify room status. Color coded system assures confidentiality and promotes staff efficiency.

291706	8" Flag Signal System w/6 flags
291708	4" Flag Signal System w/6 flags


MJM PORTABLE PRIVACY SCREEN "7000" SERIES

7003	Portable Privacy Screen, 3-Panel 21"W x 51"H
7003-2TW	Portable Privacy Screen, 3-Panel 21"W x 51"H, w/2" Twin Nylon Casters
7003-3TW	Portable Privacy Screen, 3-Panel 21"W x 51"H, w/3" Twin Nylon Casters

Room Accessories



153052



153054

OMNIMED BEAM® BEAMATIC™ SCREEN FRAME

Durable anodized aluminum frame, and unique S shape base provides the screens their fully expandable capability for immediate isolation. Each section measures 68"H x 19"W. Lightweight design makes setup a cinch.

- 153052 Two-Section Screen Frame Only, 2 Panels Needed
- 153053 Three-Section Screen Frame Only, 3 Panels Needed
- 153054 Four-Section Screen Frame Only, 4 Panels Needed
- 153055 Five-Section Screen Frame Only, 5 Panels Needed



153020



153025



153030



153035



153040



153045

OMNIMED BEAM® DESIGNER CLOTH PANELS

Panels are anti-bacterial and meet California flame resistant requirements. Note: King size screens require 2 panels per section.

- 153020 Designer Cloth Screen Panel, Apricot
- 153025 Designer Cloth Screen Panel, Parade
- 153030 Designer Cloth Screen Panel, Freeway
- 153035 Designer Cloth Screen Panel, Norway
- 153040 Designer Cloth Screen Panel, Rouge
- 153045 Designer Cloth Screen Panel, Frost
- 153046 Designer Cloth Screen Panel, Lavender
- 153047 Designer Cloth Screen Panel, Paisley
- 153048 Designer Cloth Screen Panel, Sage

**OMNIMED BEAM® ECONOMY FOLDING SCREEN FRAME**

Available with two, three or four sections. Each section 68"H x 19"W. Anodized aluminum frame. Accordion hinge design allows panels to fold for storage convenience.

- 153092 Two-Section Economy Folding Screen Frame Only
- 153093 Three-Section Economy Folding Screen Frame Only
- 153094 Four-Section Economy Folding Screen Frame Only

OMNIMED BEAM® ECONOMY VINYL PANELS

Can be used with any of our privacy screens. Panels are anti-static, anti-bacterial and meet California flame resistance requirements. Each panel measures 27"W x 52"H and is 8 mil thick. Note: King size screens require 2 panels per section.

- 153010 Economy Vinyl Screen Panel, White
- 153015 Economy Vinyl Screen Panel, Green

**OMNIMED BEAM® KING SIZE SCREEN FRAME**

Features large oversized sections that provide maximum privacy. Strong, yet lightweight anodized aluminum frame. Folds neatly for transporting and storage. Each section measures 73 3/4"H x 32"W. Need 2 panels per section (6 total).

- 153950 Three-Section King-Size Screen Frame Only
- 153961 Three-Section King-Size Screen w/casters Frame Only

OMNIMED BEAM® PRIVACY SCREEN FRAME ACCESSORIES

- 153993 18" Rods For Economy Beamatic Frames
- 153995 31" Rods For King Size Frames

**OMNIMED BEAM® THREE SECTION MOBILE SCREEN FRAME**

Anodized aluminum frame with smooth rolling 2" diameter casters. Excellent for temporary privacy in patient/resident rooms when replacing soiled cubicle curtains. Each section measures 68"H x 19"W.

- 153153 Three Section Folding Screen w/casters Frame Only

BARJAN BLINDS

Select from every name-brand manufacturer: Levolor®, Hunter Douglas®, Graber®, Bali® or Kirsch® plus our own, more economical Barjan brand. Horizontal Blinds: Available as 2" Venetian, 1" mini and 1/2" micro. Mini blinds and Venetians available in aluminum or PVC. 100's of colors available. Vertical Blinds: Huge selection of Barjan Fire Safe patterns and colors. Self-Aligning and lubricating clutch. Removable stems easy to replace. Reversible end controls. Strongest Carrier on the Market.

Blinds Horizontal & Vertical Blinds

BARJAN CUBICLE TRACK SYSTEM

Available in Satin Anodized and White Baked Enamel Finishes. Curves up to 12 inch radius.

CubicleT Cubicle/Privacy Curtain Track System

BARJAN I.V. TRACK SYSTEM

Barjan No. 110 I.V. Track is heavy gauged extruded aluminum. Satin anodized or white baked enamel finish. The track is designed to receive the NC-10 or C-10 carriers. Available with 12" radius curves.

TrackIV I.V. Track System

BARJAN CUBICLE/PRIVACY CURTAINS

Creating a beautiful, dignified environment, promotes the healing process. Pride in Workmanship by Experienced Craftsmen Produce Unsurpassed Quality. Viewed from both patients and visitor's side. Not just a requirement; offers dignity and comfort. Chosen from the largest selection of fire retardant fabrics. Passes NFPA 701. The most visible product in a patient room. Not just decorative; also functional and durable. Available with Safe Guard™ Hanging System. We "FRAME" the mesh on three sides with matching cubicle fabric for strength, durability and beauty. All bottom hems are weighted at corners and joining seams. Barjan Curtains wash better and last longer.

Cubicle Cubicle/Privacy Curtains
Design Designer Series Cubicle/Privacy Curtains

BARJAN DRAPERY

All Barjan inherently fire retardant fabrics meet or exceed Federal, State and Local fire standards. A complete range of drapery applications from window coverings to auditorium/stage drapes and room dividers. Every type of drapery, from basic functional, to the elegant enhancement of a room. Window Drapes by Barjan; Strictly elegant - for the most special occasions. Drapery Room Partitions; One drapery separates two different room functions. Dramatic Drapes; Auditorium Stage Drapes - Box Pleated Drapes, Valances and Side Panels. Installation Available.

Drapes Drapery

BARJAN POPULAR TOP TREATMENT STYLES

Barjan's fabric and style selections will create a finished appearance ranging from warm and casual, to sophisticated and formal. Many styles to choose from.

Top Top Treatment

BARJAN SECURE-GUARD™ DRAPERY SYSTEM


Each uniquely designed for safety and durability. For areas for special use. Perfect for Hotels, Hospitals and Long-term Care Facilities, especially Psychiatric, Developmental Centers and Alzheimer Units. Mountable to ceilings, walls and security screens. Easy to take down and re-install. Inherently fire retardant fabrics. Minimizes down-time, saves time and money. No more drapery pins or broken rods. No more damaged walls or ceilings.

Secure Secure-Guard™ Drapery System

BARJAN PICTURES

Thousands of prints to choose from. Contemporary, Abstracts, Traditional, Landscapes, Still Life, Children's Themes and many more. Framed Pictures: Available under glass or plexiglass with or without matting, Ready to hang, Security locking device available. Laminated Plaques: Beveled Edges, Moisture resistant, Highly durable, Ready to install, Resist scratches.

Pictures Pictures, Framed & Laminated Plaques

**OMNIMED BEAM® ADJUSTABLE KEYBOARD TRAY**

Designed to fit on our 1 3/4" diameter instrument stand columns. Knob provides adjustment for positioning. Stainless steel tray features non-skid deck to hold keyboard or other accessories.

350040 Adjustable Keyboard Tray

**Your representative
has all of the
information to design
a physician office.
Call us for details.**

Specialty Medical Products



With a complete line of high-quality, cost-effective products that perform, Graham makes the clinical setting more comfortable for the patient – and more effective for doctors and staff – with the details that matter to both. Our innovative line includes exam rolls, sheets, towels, gowns, pillows, drapes and more!

STAFF APPAREL



Personal Protection Gown

- Linen-soft, opaque nonwoven for comfort, modesty and warmth
 - Front/back opening and elastic cuff
- | | |
|----------|-----------|
| Item 241 | 50/case |
| Yellow | 30" x 42" |



Protective Gown

- Made of durable spunbond/poly
 - Fluid-repellant to protect wearer
- | | |
|----------|-----------|
| Item 290 | 25/case |
| Blue | 28" x 44" |

PATIENT GOWNS



FabriWear™ Gown

- Made of soft, fabric-like material
 - Comfortable, tear-resistant coverage
- | | |
|----------|-----------|
| Item 234 | 50/case |
| Blue | 30" x 42" |

LAB COATS



LabMates™ Lab Coat

- Snap closure and soft knit cuffs
- Made of fluid-resistant spunbond
- Features raglan sleeves and 3 pockets

Item 254	10/case
White	Sm/Med
Item 255	10/case
White	Lg/X-Lg
Item 46987	10/case
White	XX-Lg

PILLOW CASES



Solacel™ Pillow Case

- Soft, nonwoven material is rugged yet comfortable
 - Innovative material protects pillow from moisture and stains
- | | |
|------------|-----------|
| Item 53157 | 100/case |
| White | 22" x 30" |

TOWELS & WASHCLOTHS



Webtex™ Washcloths & Towels

- Made of strong, soft Webtex
 - Comfortable and absorbent
- | | |
|-----------------|--------------------------------|
| Item 401 - Blue | Washcloth 10" x 13.5" 500/case |
| Item 412 - Blue | Towel 17" x 20" 300/case |



Softsorb® Washcloths

- Made of soft, absorbent nonwoven material
 - Strong enough to withstand rigorous use
- | | |
|------------|-------------|
| Item 49733 | White |
| 500/case | 12" x 13.5" |

**Inspired by innovation,
proven by performance.®**

www.grahammedical.com