

PRECISION DYNAMICS I.V. ARMBOARDS

Sturdy stiffeners and soft foam pads provide needed support and patient comfort during I.V. procedures. Vinyl covers are waterproof, nonallergenic and bacteriostatic. Does not absorb fluids or harbor bacteria.

45001-15-MPN	Adult IV Armboard, Reusable, 3 1/2" x 18", 24/cs
45002-11-MPG	Adult IV Armboard, Disposable, 3" x 17 1/2", 100/cs
45003-11-MPF	Infant IV Armboard, Disposable, 2" x 9", 50/cs
45004-11-MPG	Wrist IV Armboard, Disposable, 3" x 9", 100/cs
45005-15-MPN	Wrist IV Armboard, Reusable, 3 1/2" x 9", 24/cs
45006-15-MPN	Infant IV Armboard, Reusable, 2" x 6", 24/cs
45008-11-MPN	Neonatal IV Armboard, Disposable, 1" x 4", 24/cs
45009-11-MPN	Neonatal IV Armboard, Disposable, 2" x 4", 24/cs
45010-11-MPF	Neonatal IV Armboard, Disposable, 2" x 6", 50/cs
45011-15-MPG	Pediatric IV Armboard, Disposable Foam, 1" x 5", 100/cs
45012-15-MPG	Wrist IV Armboard, Disposable Foam, 3" x 9", 100/cs
45013-15-MPF	Adult IV Armboard, Disposable Foam, 3" x 18", 50/cs
45021-13-MPG	Adult IV Armboard Sleeves, Non-Woven, 4 1/2" x 18", 100/cs

BD I.V. START PAK™ SITE PREP KITS

386119	Kit Includes: Povidone Iodine Prep, Adhesive Bandage, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Impermeable Drape, 25/sp, 8 sp/cs
386120	Kit Includes: Povidone Iodine Prep, Adhesive Bandage, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, 25/sp, 8 sp/cs
386121	Kit Includes: Povidone Iodine Prep, Adhesive Bandage, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Small Exam Gloves, 25/sp, 8 sp/cs
386122	Kit Includes: Povidone Iodine Prep, Adhesive Bandage, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Medium Exam Gloves, 25/sp, 8 sp/cs
386130	I.V. Start Pak w/2" x 3" Ensure-It™ Dressing, 25/bx, 8 bx/cs
386140	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, 25/sp, 8 sp/cs
386141	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Small Exam Gloves, 25/sp, 8 sp/cs
386142	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Medium Exam Gloves, 25/sp, 8 sp/cs
386143	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Large Exam Gloves, 25/sp, 8 sp/cs
386144	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Small Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386145	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Medium Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386149	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Impermeable Drape, 25/sp, 8 sp/cs

BD I.V. START PAK™ SITE PREP KITS (Continued)

386152	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Medium Exam Gloves, 25/sp, 8 sp/cs
386154	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Small Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386155	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Pair Medium Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386159	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, Impermeable Drape, 25/sp, 8 sp/cs
386195	Kit Includes: Persist™ Swabstick, Tegaderm HP Secure Strip Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Pair Medium Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386199	Kit Includes: Persist™ Swabstick, Tegaderm HP Secure Strip Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Impermeable Drape, 25/sp, 8 sp/cs
386702	Kit Includes: Persist™ Swabstick, Tegaderm HP IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, 25/sp, 8 sp/cs
386150	Kit Includes: Povidone Iodine Prep, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, 25/sp, 8 sp/cs
386160	Kit Includes: Povidone Iodine Prep, Bioclusive IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, 25/sp, 8 sp/cs
386162	Kit Includes: Povidone Iodine Prep, Bioclusive IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, (2) Alcohol Wipes, ID Label, (2) Pair Medium Exam Gloves, 25/sp, 8 sp/cs
386170	Kit Includes: Persist™ Swabstick, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, 25/sp, 8 sp/cs
386172	Kit Includes: Persist™ Swabstick, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Pair Medium Exam Gloves, 25/sp, 8 sp/cs
386175	Kit Includes: Persist™ Swabstick, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Pair Medium Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386179	Kit Includes: Persist™ Swabstick, Tegaderm IV Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Impermeable Drape, 25/sp, 8 sp/cs
386180	Kit Includes: Persist™ Swabstick, OpSite IV3000 Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, 25/sp, 8 sp/cs
386182	Kit Includes: Persist™ Swabstick, OpSite IV3000 Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Pair Medium Exam Gloves, 25/sp, 8 sp/cs
386185	Kit Includes: Persist™ Swabstick, OpSite IV3000 Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Pair Medium Exam Gloves, Impermeable Drape, 25/sp, 8 sp/cs
386189	Kit Includes: Persist™ Swabstick, OpSite IV3000 Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Impermeable Drape, 25/sp, 8 sp/cs
386190	Kit Includes: Persist™ Swabstick, Tegaderm HP Secure Strip Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, 25/sp, 8 sp/cs
386192	Kit Includes: Persist™ Swabstick, Tegaderm HP Secure Strip Dressing, (2) 2"x2" Gauze Sponges, 3/4" Roll Transpore® Tape, Tourniquet, Alcohol Wipe, ID Label, Pair Medium Exam Gloves, 25/sp, 8 sp/cs

KIMBERLY-CLARK O.R. ARMBOARD COVER

51300	O.R. Armboard Cover, 27" x 13 1/4", 20/bx, 4 bx/cs
-------	--



KIMBERLY-CLARK IV SUPPORT

29967 IV Support, Infant/Child, 6 1/4" x 2" w/9" straps, 15/bx, 2 bx/cs

BUSSE I.V. START KITS

Latex free. Kits contain: (1) tourniquet, (1) alcohol prep pad, (1) PVP prep pad, (2) 2" x 2" gauze sponges (1) roll tape. Packed in a soft, disposable tray that keeps contents sterile until opened.

818 I.V. Start Kit w/Tegaderm™ Dressing & ChloroPrep® Frepp®, 50/cs
 820 I.V. Start Kit w/Tegaderm™ Dressing, 50/cs

PREMIUM PLASTICS GENT-L-KARE® I.V. START KITS

Kits contain: tourniquet, alcohol prep pad, PVP (ampule or pad), 2" x 2" gauze, dressing (adhesive or transparent), transpore tape and label. Comes in a peel pouch.

2601 IV Start Kit w/adhesive bandage & PVP Prep Pad, 100/cs
 2602 IV Start Kit w/Tegaderm dressing & PVP ointment, 100/cs
 2603 IV Start Kit w/Opsite dressing & PVP Ampule, 100/cs


PREMIUM PLASTICS I.V. START KIT, LATEX FREE

2609 I.V. Start Kit, Latex Free, 100/cs

BD AUTOGUARD PRO™ SHIELDED I.V. CATHETERS

Made of BD Vialon™ biomaterial.

381731 24ga x 0.75" IV Catheter, Notched Needle, 50/sp, 4 sp/cs
 381732 22ga x 1" IV Catheter, Notched Needle, 50/sp, 4 sp/cs
 381733 20ga x 1" IV Catheter, 50/sp, 4 sp/cs
 381734 20ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
 381735 18ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
 381736 16ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
 381737 14ga x 1.16" IV Catheter, 50/sp, 4 sp/cs



BD INSYTE™ AUTOGUARD™ SHIELDED I.V. CATHETERS

IV catheter placement unit with Vialon™ biomaterial. Provides peripheral vascular access in all clinical situations. Clear flash chamber. With the push of a button, AUTOGUARD protects healthcare providers from accidental needlestick injuries associated with I.V. catheter placement.

381412 24ga x .75" IV Catheter, Notched Needle, 50/sp, 4 sp/cs
 381423 22ga x 1" IV Catheter, Notched Needle, 50/sp, 4 sp/cs
 381433 20ga x 1" IV Catheter, 50/sp, 4 sp/cs
 381434 20ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
 381437 20ga x 2" IV Catheter, 50/sp, 4 sp/cs
 381444 18ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
 381447 18ga x 2" IV Catheter, 50/sp, 4 sp/cs
 381454 16ga x 1 1/4" IV Catheter, 50/sp, 4 sp/cs
 381457 16ga x 1 3/4" IV Catheter, 50/sp, 4 sp/cs
 381467 14ga x 1 3/4" IV Catheter, 50/sp, 4 sp/cs
 381512 24ga x .75" IV Winged Catheter, Notched needle, 50/sp, 4 sp/cs
 381523 22ga x 1" IV Winged Catheter, Notched Needle, 50/sp, 4 sp/cs
 381533 20ga x 1" IV Winged Catheter, 50/sp, 4 sp/cs
 381534 20ga x 1.16" IV Winged Catheter, 50/sp, 4 sp/cs
 381537 20ga x 1.88" IV Winged Catheter, 50/sp, 4 sp/cs
 381544 18ga x 1.16" IV Winged Catheter, 50/sp, 4 sp/cs
 381547 18ga x 1.88" IV Winged Catheter, 50/sp, 4 sp/cs
 381554 16ga x 1.16" IV Winged Catheter, 50/sp, 4 sp/cs
 381557 16ga x 1.77" IV Winged Catheter, 50/sp, 4 sp/cs

BD INSYTE-N™ AUTOGUARD™ SHIELDED I.V. CATHETERS

Small vein shielded I.V. catheter made of BD Vialon™ biomaterial.

381411 24ga x .56" IV Catheter, Notched Needle, 50/sp, 4 sp/cs
 381511 24ga x .56" IV Winged Catheter, Notched Needle, 50/sp, 4 sp/cs

**Please ask before
 ordering a product based
 on the pictured item.**



BD INTERLINK® CANNULA & SYRINGE PRODUCTS

303345	Blunt Plastic Cannula For Use w/Baxter Interlink®, Abbott LifeShield® or B.Braun SafeLine® Systems, 100/sp, 10 sp/cs
303346	3mL Syringe w/blunt plastic cannula, For Use w/Interlink® System, 100/spx, 8 sp/cs
303347	5mL Syringe w/blunt plastic cannula, For Use w/Interlink® System, 100/sp, 4 sp/cs
303348	10mL Syringe w/blunt plastic cannula, For Use w/Interlink® System, 100/sp, 4 sp/cs
303365	Y Lock Cannula For Use w/Interlink® System, 100/sp, 10 sp/cs
303367	Vial Access Cannula For Use w/Interlink® System, 100/sp, 10 sp/cs
303369	Threaded Lock Cannula For Use w/Interlink® System, 100/sp, 10 sp/cs
303370	Lever Lock Cannula For Use w/Interlink® System, 100/sp, 10 sp/cs
303390	Twinpak™ Dual Cannula Device For Use w/Baxter Interlink®, Abbott LifeShield® or B.Braun SafeLine® Systems, 100/sp, 10 sp/cs
303391	3mL Syringe w/Twinpak™ Dual Cannula Device, For Use with Interlink®, Abbott LifeShield® or SafeLine® Systems, 100/sp, 8 sp/cs
303392	5mL Syringe w/Twinpak™ Dual Cannula Device, For Use with Interlink®, Abbott LifeShield® or SafeLine® Systems, 100/sp, 4 sp/cs
303393	10mL Syringe w/Twinpak™ Dual Cannula Device, For Use with Interlink®, Abbott LifeShield® or SafeLine® Systems, 100/sp, 4 sp/cs
303401	3mL Syringe w/vial access cannula, For Use w/Interlink® System, 100/sp, 8 sp/cs
303403	5mL Syringe w/vial access cannula, For Use w/Interlink® System, 100/sp, 4 sp/cs
303405	10mL Syringe w/vial access cannula, For Use w/Interlink® System, 100/sp, 4 sp/cs



BD SAF-T-INTIMA™ I.V. CATHETERS

Extension sets with Vialon™ biomaterial. Designed with a built-in safety shield to prevent accidental needlesticks associated with I.V. catheter placement.

383312	24ga x .75" IV Catheter w/wings, PRN & needle shield, 25/sp, 8 sp/cs
383313	24ga x .75" IV Catheter w/wings, Y Adapter & needle shield, 25/sp, 8 sp/cs
383322	22ga x .75" IV Catheter w/wings, PRN & needle shield, 25/sp, 8 sp/cs
383323	22ga x .75" IV Catheter w/wings, Y Adapter & needle shield, 25/sp, 8 sp/cs
383335	20ga x 1" IV Catheter w/wings, PRN & needle shield, 25/sp, 8 sp/cs
383336	20ga x 1" IV Catheter w/wings, Y Adapter & needle shield, 25/sp, 8 sp/cs
383346	18ga x 1" IV Catheter w/wings, Y Adapter & needle shield, 25/sp, 8 sp/cs

BD SAFETY-GARD™ I.V. NEEDLES

305406	18ga Shielded Needle, 100/sp, 10 sp/cs
305407	20ga Shielded Needle, 100/sp, 10 sp/cs

BD ANGIOCATH™ I.V. CATHETERS

For special placement procedures. Made of FEP polymer.

381112	24ga x 3/4" IV Catheter, 50/sp, 4 sp/cs
381137	20ga x 1.88" IV Catheter, 50/sp, 4 sp/cs
381144	18ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
381147	18ga x 1.88" IV Catheter, 50/sp, 4 sp/cs
381154	16ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
381157	16ga x 2" IV Catheter, 50/sp, 4 sp/cs
381164	14ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
381167	14ga x 1.88" IV Catheter, 50/sp, 4 sp/cs
382258	16ga x 3.25" IV Catheter, 10/sp, 5 sp/cs
382259	16ga x 5.25" IV Catheter, 10/sp, 5 sp/cs
382268	14ga x 3.25" IV Catheter, 10/sp, 5 sp/cs
382269	14ga x 5.25" IV Catheter, 10/sp, 5 sp/cs
382277	12ga x 3" IV Catheter, 10/sp, 5 sp/cs
382287	10ga x 3" IV Catheter, 10/sp, 5 sp/cs
381123	22ga x 1.16" IV Catheter, 50/sp, 4 sp/cs
381134	20ga x 1.16" IV Catheter, 50/sp, 4 sp/cs

TERUMO SURFLASH™ I.V. POLYURETHANE CATHETERS

Visual confirmation of catheter placement occurs due to a unique flashback feature. Flashback appears between the specially grooved needle and the catheter, providing medical staff with a higher certainty of proper catheter placement. Made of polyurethane, which is firm at room temperature and softer at body temperatures, patient comfort is enhanced without compromising proper technique.

1SR-FF1451U	14ga x 2" IV Catheter, Orange, 50/bx, 4 bx/cs
1SR-FF1651U	16ga x 2" IV Catheter, Gray, 50/bx, 4 bx/cs
1SR-FF1664U	16ga x 2 1/2" IV Catheter, Gray, 50/bx, 4 bx/cs
1SR-FF1832U	18ga x 1 1/4" IV Catheter w/20ga needle, 50/bx, 4 bx/cs
1SR-FF1851U	18ga x 2" IV Catheter, Green, 50/bx, 4 bx/cs
1SR-FF1864U	18ga x 2 1/2" IV Catheter, Green, 50/bx, 4 bx/cs
1SR-FF2025U	20ga x 1" IV Catheter w/22ga needle, 50/bx, 4 bx/cs
1SR-FF2032U	20ga x 1 1/2" IV Catheter w/22ga needle, 50/bx, 4 bx/cs
1SR-FF2051U	20ga x 2" IV Catheter, Pink, 50/bx, 4 bx/cs
1SR-FF2225U	22ga x 1" IV Catheter w/24ga needle, 50/bx, 4 bx/cs
1SR-FF2419U	24ga x 3/4" IV Catheter w/27ga needle, 50/bx, 4 bx/cs

TERUMO SURFLO® TEFLON I.V. CATHETERS

An ultra-sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible catheter with large inner diameter assures good blood flow. All this plus safe, easy-to-handle, tamper-evident packaging.

1SR-OX1451CA	14ga x 2" IV Catheter w/16ga needle, 50/bx, 4 bx/cs
1SR-OX1632CA	16ga x 1 1/4" IV Catheter w/18ga needle, 50/bx, 4 bx/cs
1SR-OX1651CA	16ga x 2" IV Catheter w/18ga needle, 50/bx, 4 bx/cs
1SR-OX1664CA	16ga x 2 1/2" IV Catheter w/18ga needle, 50/bx, 4 bx/cs
1SR-OX1832CA	18ga x 1 1/4" IV Catheter w/20ga needle, 50/bx, 4 bx/cs
1SR-OX1851CA	18ga x 2" IV Catheter w/20ga needle, 50/bx, 4 bx/cs
1SR-OX1864CA	18ga x 2 1/2" IV Catheter w/20ga needle, 50/bx, 4 bx/cs
1SR-OX2025CA	20ga x 1" IV Catheter w/22ga needle, 50/bx, 4 bx/cs
1SR-OX2032CA	20ga x 1.25" IV Catheter w/22ga needle, 50/bx, 4 bx/cs
1SR-OX2051CA	20ga x 2" IV Catheter w/22ga needle, 50/bx, 4 bx/cs
1SR-OX2225CA	22ga x 1" IV Catheter w/24ga needle, 50/bx, 4 bx/cs
1SR-OX2419CA	24ga x 3/4" IV Catheter w/27ga needle, 50/bx, 4 bx/cs



EXEL SAFELET IV CATHETERS

Latex free. Wingless. Radiopaque teflon.

- 26726 14ga x 1 1/4" IV Catheter, 50/bx, 4 bx/cs
- 26729 14ga x 2" IV Catheter, 50/bx, 4 bx/cs
- 26730 16ga x 1 1/4" IV Catheter, 50/bx, 4 bx/cs
- 26731 16ga x 2" IV Catheter, 50/bx, 4 bx/cs
- 26735 18ga x 1 1/4" IV Catheter, 50/bx, 4 bx/cs
- 26736 18ga x 1 1/2" IV Catheter, 50/bx, 4 bx/cs
- 26738 18ga x 2" IV Catheter, 50/bx, 4 bx/cs
- 26741 20ga x 1" IV Catheter, 50/bx, 4 bx/cs
- 26742 20ga x 1 1/4" IV Catheter, 50/bx, 4 bx/cs
- 26743 20ga x 2" IV Catheter, 50/bx, 4 bx/cs
- 26746 22ga x 1" IV Catheter, 50/bx, 4 bx/cs
- 26751 24ga x 3/4" IV Catheter, 50/bx, 4 bx/cs

MEDEX ACUVANCE® IV SAFETY CATHETER

Simply insert, then thread, remove the introducer - the needle is blunt. Compact size and sleek profile allow for easy one- or two-handed insertions. Smooth, tapered hub for insertion ease and patient comfort. Unique, large-volume chamber allows fast, clear flashback and easy handling. Removable flash plug allows syringe attachment for special procedures.

- 3342 Acuvance® Plus IV Catheter, Straight Hub, 16G x 1 1/4", 50/bx, 4 bx/cs
- 3350 Acuvance® Plus IV Catheter, Straight Hub, 22G x 1", 50/bx, 4 bx/cs
- 3352 Acuvance® Plus IV Catheter, Straight Hub, 16G x 2", 50/bx, 4 bx/cs
- 3353 Acuvance® Plus IV Catheter, Straight Hub, 24G x 5/8", 50/bx, 4 bx/cs
- 3354 Acuvance® Plus IV Catheter, Straight Hub, 18G x 1 3/4", 50/bx, 4 bx/cs
- 3355 Acuvance® Plus IV Catheter, Straight Hub, 18G x 1 1/4", 50/bx, 4 bx/cs
- 3356 Acuvance® Plus IV Catheter, Straight Hub, 20G x 1 1/4", 50/bx, 4 bx/cs
- 3357 Acuvance® Plus IV Catheter, Straight Hub, 20G x 1", 50/bx, 4 bx/cs
- 3358 Acuvance® Plus IV Catheter, Straight Hub, 14G x 2", 50/bx, 4 bx/cs
- 3359 Acuvance® Plus IV Catheter, Straight Hub, 20G x 1 3/4", 50/bx, 4 bx/cs
- 3370 Acuvance® IV Catheter, Winged Hub, 22G x 1", 200/cs
- 3373 Acuvance® IV Catheter, Winged Hub, 24G x 5/8", 200/cs
- 3375 Acuvance® IV Catheter, Winged Hub, 18G x 1 1/4", 200/cs
- 3376 Acuvance® IV Catheter, Winged Hub, 20G x 1 1/4", 200/cs
- 3377 Acuvance® IV Catheter, Winged Hub, 20G x 1", 200/cs

MEDEX CATHLON™ I.V. CATHETERS

Offers all the features and quality of the JELCO™ I.V. Catheter combined with the added feature of a metal hub. FEP Polymer.

- 4412 Non-Radiopaque IV Catheter, 16G x 1 1/4", Grey, 50/bx, 4 bx/cs
- 4418 Non-Radiopaque IV Catheter, 14G x 1 1/4", Orange, 50/bx, 4 bx/cs
- 4424 Non-Radiopaque IV Catheter, 18G x 1 3/4", Green, 50/bx, 4 bx/cs
- 4425 Non-Radiopaque IV Catheter, 18G x 1 1/4", Green, 50/bx, 4 bx/cs
- 4426 Non-Radiopaque IV Catheter, 20G x 1 1/4", Pink, 50/bx, 4 bx/cs
- 4427 Non-Radiopaque IV Catheter, 20G x 1", Pink, 50/bx, 4 bx/cs
- 4442 Radiopaque IV Catheter, 16G x 1 1/4", Grey, 50/bx, 4 bx/cs
- 4448 Radiopaque IV Catheter, 14G x 1 1/4", Orange, 50/bx, 4 bx/cs
- 4454 Radiopaque IV Catheter, 18G x 1 3/4", Green, 50/bx, 4 bx/cs
- 4455 Radiopaque IV Catheter, 18G x 1 1/4", Green, 50/bx, 4 bx/cs
- 4456 Radiopaque IV Catheter, 20G x 1 1/4", Pink, 50/bx, 4 bx/cs
- 4457 Radiopaque IV Catheter, 20G x 1", Pink, 50/bx, 4 bx/cs
- 4461 Striped IV Catheter, 22G x 1", Blue, 50/bx, 4 bx/cs
- 4464 Striped IV Catheter, 18G x 1 3/4", Green, 50/bx, 4 bx/cs
- 4465 Striped IV Catheter, 18G x 1 1/4", Green, 50/bx, 4 bx/cs
- 4466 Striped IV Catheter, 20G x 1 1/4", Pink, 50/bx, 4 bx/cs
- 4467 Striped IV Catheter, 20G x 1", Pink, 50/bx, 4 bx/cs
- 4420 Non-Radiopaque IV Catheter, 22G x 1, 50/bx, 4 bx/cs
- 4472 Radiopaque IV Catheter, 16G x 2", Grey, 50/bx, 4bx/cs
- 4478 Radiopaque IV Catheter, 14G x 2", Orange, 50/bx, 4bx/cs
- 4482 Striped IV Catheter, 16G x 2", Grey, 50/bx, 4bx/cs
- 4488 Striped IV Catheter, 14G x 2", Orange, 50/bx, 4bx/cs
- 4492 Non-Radiopaque IV Catheter, 16G x 2", Grey, 50/bx, 4bx/cs
- 4498 Non-Radiopaque IV Catheter, 14G x 2", Orange, 50/bx, 4bx/cs



4070



4072

MEDEX JELCO™ I.V. CATHETERS

The leading standard short peripheral I.V. catheter made of FEP Polymer. Plastic hub.

- 4042 Radiopaque IV Catheter, 16G x 1 1/4", Grey, 50/bx, 4 bx/cs
- 4048 Radiopaque IV Catheter, 14G x 1 1/4", Orange, 50/bx, 4 bx/cs
- 4050 Radiopaque IV Catheter, 22G x 1", Blue, 50/bx, 4 bx/cs
- 4053 Radiopaque IV Catheter, 24G x 3/4", Yellow, 50/bx, 4 bx/cs
- 4054 Radiopaque IV Catheter, 18G x 1 3/4", Green, 50/bx, 4 bx/cs
- 4055 Radiopaque IV Catheter, 18G x 1 1/4", Green, 50/bx, 4 bx/cs
- 4056 Radiopaque IV Catheter, 20G x 1 1/4", Pink, 50/bx, 4 bx/cs
- 4057 Radiopaque IV Catheter, 20G x 1", Pink, 50/bx, 4 bx/cs
- 4059 Radiopaque IV Catheter, 20G x 1 3/4", Pink, 50/bx, 4 bx/cs
- 4060 Striped Radiopaque IV Catheter, 22G x 1", Blue, 50/bx, 4 bx/cs
- 4063 Striped Radiopaque IV Catheter, 24G x 3/4", Yellow, 50/bx, 4 bx/cs
- 4064 Striped Radiopaque IV Catheter, 18G x 1 3/4", Green, 50/bx, 4 bx/cs
- 4065 Striped Radiopaque IV Catheter, 18G x 1 1/4", Green, 50/bx, 4 bx/cs
- 4066 Striped Radiopaque IV Catheter, 20G x 1 1/4", Pink, 50/bx, 4 bx/cs
- 4067 Striped Radiopaque IV Catheter, 20G x 1", Pink, 50/bx, 4 bx/cs
- 4069 Striped Radiopaque IV Catheter, 20G x 1 3/4", Pink, 50/bx, 4 bx/cs
- 4070 Winged Radiopaque IV Catheter, 22G x 1", Blue, 50/bx, 4 bx/cs
- 4072 Winged Radiopaque IV Catheter, 16G x 1 1/4", Grey, 50/bx, 4 bx/cs
- 4073 Winged Radiopaque IV Catheter, 24G x 5/8", Yellow, 50/bx, 4 bx/cs
- 4074 Winged Radiopaque IV Catheter, 18G x 1 3/4", Green, 50/bx, 4 bx/cs
- 4075 Winged Radiopaque IV Catheter, 18G x 1 1/4", Green, 50/bx, 4 bx/cs
- 4076 Winged Radiopaque IV Catheter, 20G x 1 1/4", Pink, 50/bx, 4 bx/cs
- 4077 Winged Radiopaque IV Catheter, 20G x 1", Pink, 50/bx, 4 bx/cs
- 4078 Winged Radiopaque IV Catheter, 14G x 1 1/4", Orange, 50/bx, 4 bx/cs
- 4082 Radiopaque IV Catheter, 16G x 2", Grey, 50/bx, 4bx/cs
- 4088 Radiopaque IV Catheter, 14G x 2", Orange, 50/bx, 4bx/cs
- 4092 Striped Radiopaque IV Catheter, 16G x 2", Grey, 50/bx, 4bx/cs
- 4098 Striped Radiopaque IV Catheter, 14G x 2", Orange, 50/bx, 4bx/cs

MEDEX JELCO™ INTERMITTENT INJECTION CAP

For use with IV catheters & winged infusion needles.

4600 Intermittent Injection Cap, 50/bx, 4 bx/cs



5042



5050

MEDEX OPTIVA™ I.V. CATHETERS


Ocrilon™ polyurethane is designed to start out firm for a controlled insertion, then soften in the vein to minimize irritation. Plastic hub. Radiopaque.

5042	16G x 1 1/4" IV Catheter, Grey, 50/bx, 4 bx/cs
5048	14G x 1 1/4" IV Catheter, Orange, 50/bx, 4 bx/cs
5050	22G x 1" IV Catheter, Blue, 50/bx, 4 bx/cs
5053	24G x 3/4" IV Catheter, Yellow, 50/bx, 4 bx/cs
5054	18G x 1 3/4" IV Catheter, Green, 50/bx, 4 bx/cs
5055	18G x 1 1/4" IV Catheter, Green, 50/bx, 4 bx/cs
5056	20G x 1 1/4" IV Catheter, Pink, 50/bx, 4 bx/cs
5057	20G x 1" IV Catheter, Pink, 50/bx, 4 bx/cs
5059	20G x 1 3/4" IV Catheter, Pink, 50/bx, 4 bx/cs
5082	16G x 2" IV Catheter, Grey, 50/bx, 4bx/cs
5088	14G x 2" IV Catheter, Orange, 50/bx, 4bx/cs

TERUMO SURFLO® WINGED INFUSION SETS

The soft, flexible interlocking wings readily conform to the body's contour. The needles are bonded to the wings at just under a 3° downward angle, for less risk of vessel damage or occlusion of the lumen. The precisely honed needle assures ease of penetration and continued access with a single insertion.

1SV-18BLK	18ga x 3/4" TW Infusion Set w/12" tubing, 50/bx, 10 bx/cs
1SV-19BI	19ga x 3/4" Intermittent Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-19BLK	19ga x 3/4" TW Infusion Set w/12" tubing, 50/bx, 10 bx/cs
1SV-19BLS	19ga x 3/4" TW Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-21BI	21ga x 3/4" Intermittent Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-21BLK	21ga x 3/4" TW Infusion Set w/12" tubing, 50/bx, 10 bx/cs
1SV-21BLS	21ga x 3/4" TW Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-22BLK	22ga x 3/4" TW Infusion Set w/12" tubing, 50/bx, 10 bx/cs
1SV-23BI	23ga x 3/4" Intermittent Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-23BLK	23ga x 3/4" TW Infusion Set w/12" tubing, 50/bx, 10 bx/cs
1SV-23BLS	23ga x 3/4" TW Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-25BLK	25ga x 3/4" TW Infusion Set w/12" tubing, 50/bx, 10 bx/cs
1SV-25BLS	25ga x 3/4" TW Infusion Set w/3 1/2" tubing, 100/bx, 5 bx/cs
1SV-25EL	25ga x 1/2" TW Infusion Set w/8" tubing, 100/bx, 10 bx/cs
1SV-27EL	27ga x 1/2" TW Infusion Set w/8" tubing, 100/bx, 10 bx/cs

**EXEL ARTERIAL VENOUS FISTULA**

Latex free. 12" tubing.

26781	15ga x 1" Arterial Venous Fistula w/eye, clamp & fixed wing, 50/bx, 10 bx/cs
26782	15ga x 1" Arterial Venous Fistula w/out eye, w/clamp & fixed wing, 50/bx, 10 bx/cs
26786	16ga x 1" Arterial Venous Fistula w/eye, clamp & fixed wing, 50/bx, 10 bx/cs
26791	17ga x 1" Arterial Venous Fistula w/eye, clamp & fixed wing, 50/bx, 10 bx/cs

**EXEL SCALP VEIN (BUTTERFLY) SETS**

Latex free. Manufactured from the highest quality 420 stainless steel. Inner wrap keeps product within the sterile field upon opening. Female luer lock connector.

26702	19ga x 3/4" Scalp Vein Set w/12" tubing, 50/bx
26703	20ga x 3/4" Scalp Vein Set w/12" tubing, 50/bx
26704	21ga x 3/4" Scalp Vein Set w/12" tubing, 50/bx
26706	23ga x 3/4" Scalp Vein Set w/12" tubing, 50/bx
26708	25ga x 3/4" Scalp Vein Set w/12" tubing, 50/bx
26709	27ga x 3/4" Scalp Vein Set w/12" tubing, 50/bx

**EXEL SCALP VEIN SET - SMALL VEIN BLOOD DRAW**


Latex free. With 12" tubing with Multi-Sample Luer Adapter.

26764	21ga x 3/4" Scalp Vein Set, 50/bx, 4 bx/cs
26766	23ga x 3/4" Scalp Vein Set, 50/bx, 4 bx/cs
26768	25ga x 3/4" Scalp Vein Set, 50/bx, 4 bx/cs

EXEL SECURETOUCH SAFETY BUTTERFLY INFUSION SET

Latex free. Safety feature reduces possibility of accidental needle sticks. Sterile. Non-toxic, non-pyrogenic. Siliconized, thin wall needle allows for maximum fluid flow. Flexible, color coded wings allow for easy gauge selection and can be taped or sutured to hold needle in place after venipuncture. Convenient peel pack packaging.

- 27704 Safety Butterfly Infusion Set, 21ga x 3/4", 12" Tube, 50/bx
- 27706 Safety Butterfly Infusion Set, 23ga x 3/4", 12" Tube, 50/bx, 10 bx/cs
- 27708 Safety Butterfly Infusion Set, 25ga x 3/4", 12" Tube, 50/bx



B BRAUN SAFELINE® IV ADMINISTRATION/EXTENSION SETS

- 352434 Filtered Gravity Admin Set, 0.2 Micron Filter, SAFELINE Valve Injection Site 6" Above Distal End, SPIN-LOCK Connector, DEHP-Free, Latex Free, 21.0mL Priming Volume, 91" L, 15 drops/mL, 50/cs
- 375125 Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/SAFELINE Valve Injection Site, Roller Clamp, Slide Clamp, SAFELINE Injection Site 6", SPIN-LOCK Connector, Latex Free, 16.0mL Priming Volume, 88" L, 60 drops/mL, 20/cs
- 375126 Add-On Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/SAFELINE Valve Injection Site SlideClamp, Spike Adapter, Latex Free, 3.0mL Priming Volume, 20" L, 20/cs
- 473993 Filtered Add-On Extension Set, Female Luer Lock, 0.2 Micron Filter, SAFELINE Valve Injection Site 6" Above Distal End, SPIN-LOCK Connector, DEHP-Free, Latex Free, 5.1mL Priming Volume, 16" L, 50/cs
- 473998 Filtered Add-On Extension Set, Female Luer Lock, 1.2 Micron Filter, SAFELINE Valve Injection Site 6" Above Distal End, SPIN-LOCK Connector, DEHP-Free, Latex Free, 5.1mL Priming Volume, 16" L, 50/cs
- NF1150 Admin Set, Control® Clamp, SAFELINE Injection Sites 6" & 40" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 17 mL Priming Volume, 91" Length, 15 Drops/mL, 50/cs
- NF1151 Admin Set, Vented Spike, Control® Clamp, SAFELINE Injection Sites 6" & 38" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 18 mL Priming Volume, 97" Length, 15 Drops/mL, 50/cs
- NF1190 Admin Set, Control® Clamp, SAFELINE Injection Sites 6" & 38" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 18 mL Priming Volume, 96" Length, 60 Drops/mL, 50/cs
- NF1220 Lipid Admin Set, Vented Spike, Non-DEHP Tubing, Control® Clamp, SAFELINE Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 15 mL Priming Volume, 82" Length, 15 Drops/mL, 50/cs
- NF1251 Admin Set, Vented Spike, Check Valve, Control® Clamp, SAFELINE Injection Sites 6", 30" & 88" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 21 mL Priming Volume, 115" Length, 15 Drops/mL, 50/cs
- NF1290 Admin Set, Check Valve, Control® Clamp, SAFELINE Injection Sites 6", 27" & 62" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 17 mL Priming Volume, 88" Length, 60 Drops/mL, 50/cs
- NF1310 Extension Set, SAFELINE Injection Site, Microbore Tubing, Removable Slide Clamp, Male Luer Lock, 0.42 mL Priming Volume, 7" Length, 50/cs
- NF1312 Extension Set, SAFELINE Injection Site, Microbore Tubing, Removable Slide Clamp, Male Luer Slip, 0.4 mL Priming Volume, 7" Length, 50/cs
- NF1318 Y-Type Microbore Extension Set w/attached SAFELINE Injection Site, Female Adapter, 2 Removable Slide Clamps, Male Luer Lock, 0.52 mL Priming Volume, 6" Length, 100/cs
- NF1370 Extension Set, Female Adapter, SAFELINE Injection Sites 6" & 28" Above Distal End, Large Bore Tubing, Slide Clamp, 2-Piece Male Luer Lock, 6 mL Priming Volume, 35" Length, 50/cs
- NF1430 Secondary Set, Control® Clamp, 2-Piece Male Luer Lock, SAFELINE Clip Lock Cannula, 10 mL Priming Volume, 50" Length, 15 Drops/mL, 50/cs

B BRAUN SAFELINE® IV ADMINISTRATION/EXTENSION SETS (Continued)

- NF3447 Admin Set, Check Valve, Control® Clamp, SAFELINE Injection Sites 6", 26" & 88" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 20 mL Priming Volume, 115" Length, 15 Drops/mL, 50/cs
- NF3482 Admin Set, Control® Clamp, SAFELINE Injection Site 26" Above Distal End, 2-Piece Male Luer Lock, 18 mL Priming Volume, 105" Length, 15 Drops/mL, 50/cs
- NF3488 Admin Set, Check Valve, Control® Clamp, SAFELINE Injection Sites 6", 28" & 88" Above Distal End, Slide Clamp, 2-Piece Male Luer Lock, 20 mL Priming Volume, 115" Length, 60 Drops/mL, 50/cs
- NF5140 Admin Set, Y-Type, Drip Chamber w/170µ Blood Filter, 3 Roller Clamps, Slide Clamp, SAFELINE Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 34 mL Priming Volume, 94" Length, 10 Drops/mL, 50/cs
- NF5145 Admin Set, Y-Type, Drip Chamber w/170µ Blood Filter, Pump Chamber, 3 Roller Clamps, Slide Clamp, SAFELINE Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 72 mL Priming Volume, 94" Length, 10 Drops/mL, 50/cs
- NF5300 Extension Set, Female Adapter, Rate Flow® Regulator, SAFELINE Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 3 mL Priming Volume, 18" Length, 50/cs
- NF5322 Admin Set, Universal Adult Spike, 15µ Filter, Roberts Clamp, Rate Flow® Regulator, SAFELINE Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 16 mL Priming Volume, 71" Length, 20 Drops/mL, 50/cs
- NF5932 Admin Set, Universal Adult Spike, 15µ Filter, Check Valve, Rate Flow® Regulator, SAFELINE Injection Sites 6" & 54" Above Distal End, 2-Piece Male Luer Lock, 16 mL Priming Volume, 71" Length, 20 Drops/mL, 50/cs
- NF1210F Filtered Gravity Admin Set, Universal Spike, 0.2 Micron Filter, Backcheck Valve, SAFELINE Valve Injection Sites, SPIN-LOCK Connector, DEHP-Free, Latex Free, 23.0mL Priming Volume, 111" L, 15 drops/mL (352438), 50/cs
- NF1320 Extension Set, SAFELINE Injection Site, Macrobore Kink-Resistant Tubing, Slide Clamp, 2-Piece Male Luer Lock, 0.63 mL Priming Volume, 7" Length, 50/cs
- NF1322 Extension Set, SAFELINE Injection Site, Macrobore Kink-Resistant Tubing, Slide Clamp, 2-Piece Male Luer Slip, 0.58 mL Priming Volume, 7" Length, 50/cs
- NF1323 Extension Set, SAFELINE Injection Site, Microbore Tubing, Adjustable Bridge, Male Luer Slip, 0.32 mL Priming Volume, 4" Length, 50/cs
- NF1330 Extension Set, Female Adapter, SAFELINE Injection Site, Microbore Tubing, Removable Slide Clamp, T-Connector, 2-Piece Male Luer Lock Connection, 0.45 mL Priming Volume, 6" Length, 50/cs
- 352062 Admin Set, Non-Vented Spike, Slide Clamp, Backcheck Valve, Upper Low Profile SAFSITE Valve Injection Site w/roller clamp, Low Profile SAFSITE Valve Injection Site, Distal SPIN-LOCK Connector, Latex Free, DEHP-Free, 106" Length, 50/cs
- 352072 Admin Set, Slide Clamp, Low Profile SAFSITE Needle-Free Y-Site w/backcheck valve, Roller Clamp, Slide Clamp, Low Profile SAFSITE Y-Site, Slide Clamp, Low Profile SAFSITE Low Profile Y-Site, Distal Sliding SPIN-LOCK Connection, 104", 15 Drops/mL, 50/cs
- 352079 Admin Set, Slide Clamp, Low Profile SAFSITE Y-Site w/backcheck valve, Roller Clamp, Slide Clamp, Low Profile SAFSITE Y-Site, Slide Clamp, Low Profile SAFSITE Needle-Free Low Profile Y-Site, Distal Sliding SPIN-LOCK Connection, 104", 60 Drops/mL, 50/cs
- 352207 0.2 Micron SUPOR® Filter Low Volume Extension Set, Air Eliminating In-Line IV Filter w/slide clamp, SAFSITE Latex Free Y-Site, Small Bore Tubing, Distal Luer Lock Connection, DEHP-Free Tubing, 1.4 mL Priming Volume, 9" Length, 50/cs
- 375057 Add-On Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/SAFSITE Valve, Slide Clamp, Spike Adapter, Latex Free, 3.0mL Priming Volume, 120" L, 20/cs
- 375058 Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/SAFSITE Valve, Roller Clamp, Slide Clamp, SAFSITE Valve Injection Site 6" Above Distal End, SPIN-LOCK Connector, Latex Free, 15.3mL Priming Volume, 88" L, 20/cs
- 471905 Small Bore T-Port Extension Set w/low profile SAFSITE Valve, Luer Slip Connection, 0.35 mL Priming Volume, 5" Length, 50/cs
- 471906 Standard Bore T-Port Extension Set w/low profile SAFSITE Valve, Luer Slip Connection, Slide Clamp, 1.0 mL Priming Volume, 5" Length, 100/cs
- 471976 Small Bore Y-Extension Set w/SAFSITE Valve, Luer Slip Connection, 0.5 mL Priming Volume, 6" Length, 100/cs
- 471987 Small Bore T-Port Extension Set w/low profile SAFSITE Valve, SPIN-LOCK® Connection, 0.47 mL Priming Volume, 5" Length, 100/cs
- 471993 Small Bore Y-Extension Set w/SAFSITE Valve, Luer Lock Connection, 0.5 mL Priming Volume, 6" Length, 100/cs
- 471994 Small Bore Y-Extension Set w/SAFSITE Valve, SPIN-LOCK® Connection, 0.5 mL Priming Volume, 6" Length, 100/cs

B BRAUN SAFELINE® IV ADMINISTRATION/EXTENSION SETS

(Continued)

- 472022 Standard Bore Extension Set w/2 low profile SAFSITE Valves at Y-Sites, Distal Luer Slip Connection, 3.4 mL Priming Volume, 23" Length, 100/cs
- 472060 Small Bore Extension Set w/SAFSITE Valve, Luer Lock Connection, 0.42 mL Priming Volume, 6" Length, 100/cs
- 472068 Small Bore T-Port Extension Set w/SAFSITE Valve, SPIN-LOCK® Connection, 0.42 mL Priming Volume, 4" Length, 100/cs
- 472070 Standard Bore Extension Set w/SAFSITE Valve, Luer Slip Connection, 1.27 mL Priming Volume, 6" Length, 100/cs
- 472071 Standard Bore T-Port Extension Set w/SAFSITE Valve, SPIN-LOCK® Connection, 1.2 mL Priming Volume, 7" Length, 100/cs
- 472072 Standard Bore Extension Set w/low profile SAFSITE Valve at Y-Site, In-Line Check Valve, Distal Luer Slip Connection, 1.4 mL Priming Volume, 8" Length, 100/cs
- 473017 Standard Bore Extension Set w/2 low profile SAFSITE Valves at Y-Sites, Distal SPIN-LOCK® Connection, 5.2 mL Priming Volume, 35" Length, 50/cs
- 473027 Small Bore Extension Set w/removable SAFSITE Valve, Luer Lock Connection, 0.6 mL Priming Volume, 13" Length, 100/cs
- 473042 Small Bore Extension Set w/bonded SAFSITE Valve, Luer Lock Connection, 0.45 mL Priming Volume, 13" Length, 100/cs
- 473045 Small Bore Extension Set w/SAFSITE Valve, Luer Slip Connection, 0.42 mL Priming Volume, 6" Length, 100/cs
- 473051 0.2 Micron Charged SUPOR® Filter Extension Set, Endotoxin Retentive, Air Eliminating In-Line IV Filter w/on/off clamp, SAFSITE Latex Free Y-Site, Distal SPIN-LOCK® Connection, DEHP-Free Tubing, 4.2 ml Priming Volume, 14" Length, 50/cs
- 473052 Small Bore T-Port Extension Set w/SAFSITE Valve, Luer Slip Connection, 0.43 mL Priming Volume, 4" Length, 100/cs
- 473149 Standard Bore Y-Extension Set w/low profile SAFSITE Valve, Distal Sliding SPIN-LOCK® Connection, 1.05 mL Priming Volume, 7.5" Length, 50/cs
- 473240 Small Bore Y-Extension Set w/SAFSITE Valve, SPIN-LOCK® Connection, T-Port Connection, 0.6 mL Priming Volume, 6" Length, 100/cs
- 473990 0.2 Micron Filter Extension Set, SAFSITE Injection Site 6" Above Distal End, On/Off Clamp, SPIN-LOCK Connector, DEHP-Free, Latex Free, 4.5mL Priming Volume, 16"L, 50/cs
- 474014 Standard Bore Extension Set w/low profile SAFSITE Valve at Y-Site, Distal Luer Slip Connection, 2.0 mL Priming Volume, 15" Length, 100/cs





EXEL I.V. ADMINISTRATION SETS

Latex free. 78" tubing with Y-injection site. Adelberg III roller clamp plus occlusion clamp (allows user to shut off I.V. without using roller clamp). Sterile. Non toxic. Non pyrogenic fluid path. Priming volume: 23ml/approx. I.D. 3.0 +/- 0.05mm. O.D. 4.1 +0.10/-0.05mm. Vented and non-vented on same set by opening/closing valve.

- 29081 15 Drop/mL IV Administration Set, 50/cs
- 29083 60 Drop/mL IV Administration Set, 50/cs

EXEL MICRON FILTERED IV ADMINISTRATION SETS


- 26085 .2 Micron Filter IV Admin Set w/Y injection site, Vented/Non-Vented, 84" Tube, Two Clamps, 20 Drops/mL, Latex/DEHP Free, 50/cs
- 27085 .2 Micron Filter IV Admin Set w/Y injection site & needle-less, Vented/Non-Vented, 84" Tube, Two Clamps, 20 Drops/mL, Latex/DEHP Free, 50/cs



AMSINO AMSAFE® I.V. ADMINISTRATION SETS

DEHP-free and latex free. Packaged in poly pouch.

- 000606 6" Standardbore Extension Set, One AMSafe® Needleless Injection Site, One Split Septum Injection Site, 50/cs
- 003201 32" Standardbore Extension Set, Two Injection Sites, 50/cs
- 107201 72" IV Admin Set, 10 Drops Per mL, One Injection Site, 50/cs
- 108301 83" IV Admin Set, 10 Drops Per mL, One Injection Site, 50/cs
- 108302 83" IV Admin Set, 10 Drops Per mL, One AMSafe® Needleless & One Split Septum Injection Site, One Back Check Valve, 50/cs
- 108303 83" IV Admin Set, 10 Drops Per mL, One AMSafe® Needleless Injection Site, 50/cs
- 108305 89" IV Admin Set w/pre-attached extension set, 10 Drops Per mL, One AMSafe® Needleless & One Split Septum Injection Site, 50/cs
- 108306 83" IV Admin Set, 10 Drops Per mL, One AMSafe® Needleless & One Split Septum Injection Site, 50/cs
- 109602 96" IV Admin Set, 10 Drops Per mL, Two Injection Sites, Slip Luer w/flashbulb, 50/cs
- 109606Y 96" IV Admin Set, 10 Drops Per mL, 60 Drops Per mL, One AMSafe® Needleless & One Split Septum Injection Sites, 50/cs
- 159602 96" IV Admin Set, 15 Drops Per mL, Two Injection Sites, Slip Luer w/flashbulb, 50/cs
- 607201 72" IV Admin Set, 60 Drops Per mL, One Injection Site, 50/cs
- 608306 83" IV Admin Set, 60 Drops Per mL, One AMSafe® Needleless & One Split Septum Injection Site, 50/cs
- AA3101 78" IV Admin Set, 15 Drops Per mL, One Injection Site, 50/cs
- AC3110 105" IV Admin Set, 10 Drops Per mL, Two AMSafe® Needleless Injection Sites, 50/cs
- AC3112 105" IV Admin Set, 60 Drops Per mL, Two AMSafe® Needleless Injection Sites, 50/cs
- AE3108 8" Standardbore IV Extension Set w/AMSafe® Injection Site, 100/cs
- AS3500 40" IV Secondary Set, 10 Drops Per mL, One Medication Hanger, 50/cs



AMSINO FLOWSAFE™ I.V. ADMINISTRATION SETS

DEHP-free and latex free. Packaged in poly pouch.

- AFS102 100" IV Admin Set, 10 Drops Per mL, One FLOWSafe™ Flow Controller, One AMSafe® Needleless Injection Site, 50/cs
- AFS106 20" Standardbore Extension Set w/FLOWSafe™ Flow Controller, AMSafe® Needleless Injection Site, 50/cs
- AFS107 20" Standardbore Extension Set w/FLOWSafe™ Flow Controller, One Injection Site, 50/cs



B BRAUN ADDITIV® PRIMARY ADMINISTRATION SETS

- V1443 Administration Set, Check Valve, 15 Drops/ml, Injection Sites 6" & 60" Above Distal End, Control® Clamp, 2-Piece Male Luer Lock, 17 mL Priming Volume, 86" Length, 50/cs
- 352401 Filtered Admin Set, 1.2 Micron Filter, 1 Injection Site, Backcheck Valve & SPIN-LOCK Connector, 21mL Priming Volume, 109"CL, 15 drops/mL, 50/cs
- 352424 Filtered Admin Set, 0.2 Micron Filter, Roller Clamp, Slide Clamp, 1 Injection Site & SPIN-LOCK Connector, DEHP-Free, 21mL Priming Volume, 91" L, 15 drops/mL, 50/cs
- 352426 Filtered Admin Set, 0.2 Micron Filter, 2 Injection Sites & SPIN-LOCK Connector, 24mL Priming Volume, 99"CL, 15 drops/mL, 50/cs
- 352427 Filtered Admin Set, 0.2 Micron Filter, 3 Injection Sites, Backcheck Valve & SPIN-LOCK Connector, DEHP-Free, 25mL Priming Volume, 104" L, 15 drops/mL, 50/cs
- 352432 Filtered Admin Set, 0.2 Micron Filter, Non-Vented Spike, 2 Injection Sites, Backcheck Valve & SPIN-LOCK Connector, 15mL Priming Volume, 80" L, 60 drops/mL, 50/cs
- 352436 Filtered Admin Set, 1.2 Micron Filter, 1 Injection Site, SPIN-LOCK Connector, 119mL Priming Volume, 91" L, 15 drops/mL, 50/cs
- 352447 Filtered Admin Set, 0.2 Micron Filter, 3 Injection Sites, Backcheck Valve & SPIN-LOCK Connector, DEHP-Free, 20mL Priming Volume, 83" L, 15 drops/mL, 50/cs
- V1445-44 Administration Set, Check Valve, Injection Sites 6", 12", 48" & 100" Above Distal End, Control® Clamp, 2-Piece Male Luer Lock, 24 mL Priming Volume, 127" Length, 15 Drops/mL, 50/cs
- V1446 Administration Set, Check Valve, Injection Sites 6" & 68" Above Distal End, Control® Clamp, 2-Piece Male Luer Lock, 18 mL Priming Volume, 93" Length, 15 Drops/mL, 50/cs
- V1447 Administration Set, Check Valve, Injection Sites 6", 28" & 90" Above Distal End, Control® Clamp, Slide Clamp, 2-Piece Male Luer Lock, 20 mL Priming Volume, 115" Length, 15 Drops/mL, 50/cs
- V1448 Mini-Drop Admin Set, Check Valve, Injection Sites 6", 28" & 63" Above Distal End, Control® Clamp, 2-Piece Male Luer Lock, 18 mL Priming Volume, 87" Length, 60 Drops/mL, 50/cs
- V1453 Administration Set, Check Valve, Injection Sites 6", 28" & 72" Above Distal End, Control® Clamp, 2-Piece Male Luer Lock, 19 mL Priming Volume, 96" Length, 15 Drops/mL, 50/cs
- V1484 Administration Set, Check Valve, Injection Sites 6" & 58" Above Distal End, Control® Clamp, Slide Clamp, 2-Piece Male Luer Lock, 17 mL Priming Volume, 83" Length, 15 Drops/mL, 50/cs
- V1488 Mini-Drop Admin Set, Check Valve, Injection Sites 6", 28" & 90" Above Distal End, Control® Clamp, Slide Clamp, 2-Piece Male Luer Lock, 20 mL Priming Volume, 115" Length, 60 Drops/mL, 50/cs

B BRAUN ANESTHESIA IV ADMINISTRATION SETS

For use with bags or vented bottles.

- 351510 Admin Set, Universal Spike, Roller Clamp, Y-Site, Luer Lock Connection, Detachable 4-Way Stopcock w/40" extension, Distal SPIN-LOCK Connector, DEHP-Free, Latex Free, 19.0mL Priming Volume, 114" L, 15 drops/mL, 50/cs
- 352088 Admin Set, Non-Vented Spike, Latex Free Upper Y-Site w/Backcheck Valve, Roller Clamp, Removable Extension Set w/one 4-way stopcock, 2 Low Profile Normally Closed Needle-Free Vales, Latex-Free Y-Site & Distal /Luer Slip Connection, 119" Length, 50/cs



B BRAUN BASIC IV ADMINISTRATION SETS

- V1390 Lipid Admin Set, Vented Spike, Non-DEPH Tubing, Control® Clamp, Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 16 mL Priming Volume, 81" Length, 15 Drops/ml, 50/cs
- V1482 Administration Set, Control® Clamp, Slide Clamp, Injection Site 26" Above Distal End, 2- Piece Male Luer Lock, 19 mL Priming Volume, 106" Length, 15 Drops/ml, 50/cs

B BRAUN BASIC Y-TYPE ADMINISTRATION SET

- V1500 Administration Set, Two Roller Clamps, Control® Clamp, Male Luer Slip, 22 mL Priming Volume, 87" Length, 15 Drops/ml, 50/cs

B BRAUN BAXTER COMPATIBLE PUMP ADMINISTRATION SETS

For use with Baxter Flo-Gard® 6200 and 6300 Infusion Pump.


- 375127 Admin Set, Non-Vented Spike, Roller Clamp, Metriset 150mL Burette Chamber, Slide Clamp, Injection Site, SPIN-LOCK Connector, Latex Free, 17mL Priming Volume, 95" L, 60 drops/mL, 20/cs
- NF6413 Admin Set, Check Valve, Control® Clamp, Slide Clamp, SAFELINE® Injection Sites 6", 42" & 88" Above Distal End, 2-Piece Male Luer Lock, 21 mL Priming Volume, 115" Length, 15 Drops/mL, 50/cs
- NF6415 Admin Set, Check Valve, Control® Clamp, SAFELINE® Injection Sites 6", 26" & 72" Above Distal End, 2-Piece Male Luer Lock, 19 mL Priming Volume, 98" Length, 60 Drops/mL, 50/cs
- SS6413 Admin Set, Check Valve, Control® Clamp, Slide Clamp, SAFSITE® Valve 6", 42" & 88" Above Distal End, 2-Piece Male Luer Lock, 21 mL Priming Volume, 115" Length, 15 Drops/mL, 50/cs
- SS6415 Admin Set, Check Valve, Control® Clamp, Slide Clamp, SAFSITE® Valve 6", 26" & 72" Above Distal End, 2-Piece Male Luer Lock, 19 mL Priming Volume, 98" Length, 60 Drops/mL, 50/cs
- US6413 Admin Set, Check Valve, Control® Clamp, Slide Clamp, ULTRASITE® Valve 6" & 88" Above Distal End, 2-Piece Male Luer Lock, 21 mL Priming Volume, 114" Length, 15 Drops/mL, 50/cs
- US6415 Admin Set, Check Valve, Control® Clamp, Slide Clamp, ULTRASITE® Valve 6" & 74" Above Distal End, 2-Piece Male Luer Lock, 19 mL Priming Volume, 100" Length, 60 Drops/mL, 50/cs
- V6400 Admin Set, Control® Clamp, Injection Site 26" Above Distal End, 2-Piece Male Luer Lock, 18 mL Priming Volume, 101" Length, 15 Drops/mL, 50/cs
- V6401 Admin Set, Vented Spike, Control® Clamp, Injection Site 26" Above Distal End, 2-Piece Male Luer Lock, 19 mL Priming Volume, 101" Length, 15 Drops/mL, 50/cs
- V6405 Admin Set, Control® Clamp, Injection Site 26" Above Distal End, 2-Piece Male Luer Lock, 18 mL Priming Volume, 101" Length, 60 Drops/mL, 50/cs
- V6410 Admin Set, Check Valve, Control® Clamp, Injection Sites 6" & 76" Above Distal End, 2-Piece Male Luer Lock, 19 mL Priming Volume, 103" Length, 15 Drops/mL, 50/cs
- V6413 Admin Set, Check Valve, Slide Clamp, Control® Clamp, Injection Sites 6", 42" & 88" Above Distal End, 2-Piece Male Luer Lock, 21 mL Priming Volume, 113" Length, 15 Drops/mL, 50/cs
- V6415 Admin Set, Check Valve, Control® Clamp, Injection Sites 6", 26" & 72" Above Distal End, 2-Piece Male Luer Lock, 19 mL Priming Volume, 98" Length, 60 Drops/mL, 50/cs

B BRAUN FILTERED EXTENSION SETS

- 352205 0.2 Micron SUPOR® Filter Low Volume Extension Set, Air Eliminating IV Filter w/slide clamp, Latex Free Injection Site, Small Bore Tubing & Distal Luer Lock Connection, DEHP-Free Tubing, 1.4 mL Priming Volume, 9" Length, 50/cs
- 352212 1.2 Micron SUPOR® Filter Low Volume Extension Set, Air Eliminating IV Filter w/slide clamp, Latex Free Injection Site, Small Bore Tubing & Distal Luer Lock Connection, DEHP-Free Tubing, 1.4 mL Priming Volume, 9" Length, 50/cs
- 473036 0.2 Micron SUPOR® Low Volume Air Eliminating IV Filter w/luer lock connections, 0.8 mL Priming Volume, 50/cs
- 473040 5.0 Micron Filter Extension Set, In-Line IV Filter w/slide clamp, Latex Free Injection Site & Distal Luer Slip Connection, 2.7 mL Priming Volume, 12" Length, 100/cs
- 473043 0.2 Micron Charged SUPOR® Filter Extension Set, Endotoxin Retentive, Air Eliminating In-Line IV Filter w/on/off clamp, Latex Free Injection Site & Distal SPIN-LOCK® Connection, DEHP-Free Tubing, 4.0 mL Priming Volume, 14" Length, 50/cs
- 473061 0.2 Micron SUPOR® Micro Volume Filter Extension Set, Air Eliminating IV Filter w/slide clamp, Latex Free Injection Site, Small Bore Tubing & Distal SPIN-LOCK® Connection, DEHP-Free Tubing, 0.6 mL Priming Volume, 7" Length, 50/cs
- 473989 Extension Set, Female Luer Connector, 0.2 Micron Filter, 1 Injection Site, On/Off Clamp & SPIN-LOCK Connector, DEHP-Free, Latex Free, 5.1mL Priming Volume, 16"L, 50/cs
- 473994 Extension Set, 1.2 Micron Filter, Luer Lock Connector, DEHP-Free, 4.3mL Priming Volume, 10"L, 50/cs
- 473995 Extension Set, 1.2 Micron Filter, 1 Injection Site & SPIN-LOCK Connector, DEHP-Free, 5.1mL Priming Volume, 16"L, 50/cs
- N2023 5.0 PMF™ Particulate Matter Filter, 200/cs

B BRAUN HYPODERMOCLYSIS ADMINISTRATION SET FOR SUBCUTANEOUS INFUSION

- V1527 Admin Set, 2 Roller Clamps, 2 Flashtubes, 2 Male Luer Slips, 2 Unattached 22G Needles, 24 mL Priming Volume, 68" Length, 15 Drops/mL, 50/cs



B BRAUN IV ACCESSORIES


- 654499 Luer Slip Adaptor, 100/cs
- N2030-100 Extended IV Catheter Plug, Male Luer Lock, 100/cs
- N2044-100 IV Catheter Plug, Male Luer Lock, 100/cs
- N2046-100 IV Catheter Plug, 100/cs
- V5600 Dual Injection Site, 50/cs

B BRAUN METRASET® MEASURED VOLUME SOLUTION ADMINISTRATION SETS

- 375059 Add-On Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/ULTRASITE Valve, Slide Clamp, Spike Adapter, DEHP-Free, Latex-Free, 3.0mL Priming Volume, 20"L, 20/cs
- 375112 Add-On Admin Set w/o automatic shutoff, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Slide Clamp, Spike Adapter, Latex Free, 3mL Priming Volume, 20"L, 20/cs
- 375113 Mini-Drop Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Roller Clamp, Injection Site 6" Above Distal End, SPIN-LOCK Connector, Latex Free, 15mL Priming Volume, 87"L, 60 drops/mL, 20/cs
- 375115 Add-On Admin Set w/automatic shutoff, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Slide Clamp, Spike Adapter, Latex Free, 3mL Priming Volume, 20"L, 20/cs
- 375123 Blood Admin Set, Drip Chamber w/170 micron blood filter, Roller Clamp, 150mL Chamber, Roller Clamp, Injection Site 6" Above Distal End, SPIN-LOCK Connector, 33mL Priming Volume, 89"⊘L, 20 drops/mL, 50/cs
- 375124 Mini-Drop Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Roller Clamp, Filterflow® 0.2 Micron Air-Eliminating Filter, 2 Injection Sites, SPIN-LOCK Connector, Latex Free, 20mL Priming Volume, 112"L, 60 drops/mL, 20/cs
- 375128 Mini-Drop Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Roller Clamp, Injection Sites 26" & 6" Above Distal End, SPIN-LOCK Connector, 18mL Priming Volume, 109"⊘L, 60 drops/mL, 20/cs
- 375134 Mini-Drop Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Roller Clamp, Injection Site 6" Above Distal End, SPIN-LOCK Connector, Latex Free, 17mL Priming Volume, 90"L, 60 drops/mL, 20/cs
- 375137 Admin Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/ULTRASITE Valve, Roller Clamp, Slide Clamp, ULTRASITE Injection site 6", SPIN-LOCK Connector, DEHP-Free, Latex Free, 16.0mL Priming Volume, 88"L, 60 drops/mL, 20/cs

B BRAUN PRIMARY ADDITIV® MINI-DROP ADMINISTRATION SETS WITH CHECK VALVE

- V1449-01 Administration Set, Check Valve, Injection Sites 6", 28" & 63" Above Distal End, Controll® Clamp, 2-Piece Male Luer Lock, 16 mL Priming Volume, 87" Length, 60 Drops/mL, 50/cs



B BRAUN RATE FLOW™ REGULATOR ADMINISTRATION SETS

- US5300 Admin Set, Female Adapter, Rate Flow Regulator, ULTRASITE Y-Site 6" Above Distal End, SPIN-LOCK Connection, Latex Free, 2.7 mL Priming Volume, 18" Length, 50/cs
- V5200 Admin/Extension Set, Female Adaptor, Rate Flow™ Regular, Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 2.5 mL Priming Volume, 19"⊘ Length, 50/cs
- V5922 Admin/Adult IV Set, Universal Adult Spike, 15µ Filter, Roberts Clamp, Rate Flow™ Regulator, Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 14 mL Priming Volume, 84" Length, 20 Drops/mL, 50/cs
- V5926 Admin/Pediatric IV Set, Universal Mini-Drop Spike, 15µ Filter, Roberts Clamp, Rate Flow™ Regulator, Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 14 mL Priming Volume, 84" Length, 60 Drops/mL, 50/cs
- V5932 Admin Set, Universal Adult Spike, Check Valve, 15µ Filter, Roberts Clamp, Rate Flow™ Regulator, Injection Sites 6" & 54" Above Distal End, 2-Piece Male Luer Lock, 16 mL Priming Volume, 90" Length, 20 Drops/mL, 50/cs

B BRAUN SAFELINE® ACCESSORIES

N2319-100	Protected Needle, 19G Thin Wall, 100/cs
NF9100	SAFELINE Valve Injection Site: Split Septum Access w/SPIN-LOCK Connector, 0.25mL Priming Volume, 400/cs
NF9200	Clip Lock™ Cannula, Locking Cannula w/female luer, 400/cs
NF9210	Blunt Cannula w/female luer, 400/cs
NF9285	Medic Plastic Anti-Stick Needle/Connector, 500/cs

B BRAUN STOPCOCK EXTENSION SETS

456027	Extension Set, 4-Way, 2 Female Luer Lock Ports, 20" Extension Tubing w/Luer Slip Adapter, Port Covers, 1.8 mL Priming Volume, 22" Length, 50/cs
456030	Extension Set, 4-Way, 2 Female Luer Lock Ports, 33" Extension Tubing w/SPIN-LOCK® Adapter, Port Covers, 2.6 mL Priming Volume, 35" Length, 50/cs
456032	Extension Set, 4-Way, 2 Female Luer Lock Ports, 33" Extension Tubing w/Luer Slip Adapter, Port Covers, 2.6 mL Priming Volume, 35" Length, 50/cs
471961	Extension Set, 4-Way, 2 Female Luer Lock Ports, 41" Extension Tubing, One Latex Free Y-Site w/Luer Slip Male Adapter, Port Covers, 6.3 mL Priming Volume, 42" Length, 100/cs
473054	Extension Set, Three 4-Way, 4 Female Luer Lock Ports, 40" Extension Tubing w/SPIN-LOCK® Adapter, Port Covers, 6.8 mL Priming Volume, 43" Length, 100/cs
473543	Extension Set, High-Flow, 4-Way, 2 Female Luer Lock Ports, 33" Extension Tubing w/SPIN-LOCK® Adapter, Port Covers, 5.5 mL Priming Volume, 35" Length, 50/cs
456025	Extension Set, 4-Way, 2 Female Luer Lock Ports, 20" Extension Tubing w/SPIN-LOCK® Adapter, Port Covers, 1.8 mL Priming Volume, 22" Length, 50/cs

B BRAUN ANESTHESIA STOPCOCK EXTENSION SETS

473547	Extension Set, 3-Way Anesthesia Stopcock w/injection site, 2 Full-Threaded Female Ports, 20" Extension Tubing w/SPIN-LOCK® Adapter, Port Covers, 3.3 mL Priming Volume, 22" Length, Latex Free, 50/cs
473549	Extension Set, 3-Way Anesthesia Stopcock w/injection site, 2 Full-Threaded Female Ports, 33" Extension Tubing w/SPIN-LOCK® Adapter, Port Covers, 5.2 mL Priming Volume, 35" Length, Latex Free, 50/cs
473551	Extension Set, 3-Way Anesthesia Stopcock, 3 Full-Threaded Female Ports, 20" Extension Tubing w/Luer Slip Adapter, Port Covers, 3.3 mL Priming Volume, 22" Length, Latex Free, 50/cs
473552	Extension Set, 3-Way Anesthesia Stopcock, 3 Full-Threaded Female Ports, 33" Extension Tubing w/Luer Slip Adapter, Port Covers, 5.0 mL Priming Volume, 35" Length, Latex Free, 50/cs

BD VALVES & ACCESSORIES

513094	470V Two Way Automatic Valve-Female Luer-Lok to Male Luer-Lok w/side hose end for 1/8" I.D. plastic tubing (Non-Return Type)
513095	470V-1 Modified Two Way Automatic Valve-AT Shaped w/Female Luer-Lok & Female Luer-Lok Side Connection (Non-Return Type)
513397	470W-Black Buna-N Rubber Washer For 470V & 3410V, 12/cs
513576	470RVS Set-Four Buna-N Rubber Valve Tacks & 2 Stainless Steel Springs

MEDEX ADAPTERS & CONNECTORS

MX453L	T-Connector w/male luer lock, 50/cs
MX492	Injection Adapter w/male luer lock, 7/8", Sterile, Latex Free, 100/cs

MEDEX STOPCOCKS

Sterile. Single. Latex free.

MX5311L	3-Way Stopcock w/male luer lock, Gamma/ETO, PVC Free, Phthalate Free, 50/cs
---------	---

B BRAUN VISTA® PUMP & GRAVITY ADMINISTRATION SETS

352229	Basic Set, Universal Spike, Free-Flow Protector Clip, Injection Site 72" Above Distal Connector, Roller Clamp, SPIN-LOCK® Connector, Latex Free, DEHP-Free, 113" Length, 15 Drops/mL, 50/cs
352230	Macro Drop Set, Universal Spike, Slide Clamp, Backcheck Valve, Free-Flow Protector Clip, Injection Sites 83", 58" & 6" Above Distal Connector, Roller Clamp, SPIN-LOCK® Connector, Latex Free, DEHP-Free, 111" Length, 15 Drops/mL, 50/cs
352232	Nitro Set, Universal Spike, Free-Flow Protector Clip, Roller Clamp, SPIN-LOCK Connector, Latex Free, 125" Length, 60 Drops/mL, 50/cs
352233	Half Set, Female Luer Lock Connector, Free-Flow Protector Slip, Roller Clamp, Injection Site 6" Above Distal Connector, SPIN-LOCK Connector, Latex Free, DEHP-Free, 87" Length, 50/cs
352234	Micro Drop Set, Non-Vented Spike, Slide Clamp, Backcheck Valve, Free-Flow Protector Clip, Injection Sites 83", 58" & 6" Above Distal Connector, Roller Clamp, SPIN-LOCK® Connector, Latex Free, 111" Length, 60 Drops/mL, 50/cs
352337	Pump Set, Universal Spike, Free Flow Protector Clip, ULTRASITE Injection Site 72" Above Distal End, Roller Clamp, SPIN-LOCK Connection, DEHP-Free, Latex Free 113" Length, 15 Drops/mL, 50/cs
352338	Pump Set, Universal Spike, Backcheck Valve, Free Flow Protector Clip, ULTRASITE Y-Sites 6" & 83" Above Distal End, Roller Clamp, SPIN-LOCK Connection, DEHP-Free, 109" Length, 15 Drops/mL, 50/cs
352339	Pump Set, Non-Vented Spike, Backcheck Valve, Free Flow Protector Clip, ULTRASITE Y-Sites 6" & 83" Above Distal End, Roller Clamp, SPIN-LOCK Connection, DEHP-Free, 109" Length, 60 Drops/mL, 50/cs
352355	Pump Set, Universal Spike, Backcheck Valve, Free Flow Protector Clip, ULTRASITE Y-Sites 6", 58" & 83" Above Distal End, Roller Clamp, SPIN-LOCK Connection, DEHP-Free, 111" Length, 15 Drops/mL, 50/cs
352403	Pump Set, Universal Spike, Backcheck Valve, Injection Sites 87" & 6" Above Distal End, Free-Flow Protector Clip, Roller Clamp, 0.2 Micron Filter, Slide Clamp, SPIN-LOCK Connector, DEHP-Free, 23.7mL Priming Volume, 111" L, 15 drops/mL, 50/cs
352437	Pump Set, Universal Spike, Free-Flow Protector Clip, Roller Clamp, 0.2 Micron Filter, ULTRASITE Injection Sites 6" & 87" Above Distal End, SPIN-LOCK Connector, Latex Free, DEHP-Free, 41" Length, 15 drops/mL, 50/cs
375114	Pump Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Free-Flow Protector Clip, Roller Clamp, Injection Sites 28" & 6" Above Distal End, Slide Clamp, SPIN-LOCK connector, Latex Free, 22mL Priming Volume, 112" L, 20 drops/mL, 20/cs
375129	Pump Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber, Free-Flow Protector Clip, Roller Clamp, Injection Sites 28" & 6" Above Distal End, Slide Clamp, SPIN-LOCK connector, Latex Free, 22mL Priming Volume, 112" L, 60 drops/mL, 20/cs
375130	Pump Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/ULTRASITE Valve, Free-Flow Protector Clip, Roller Clamp, 2 Injection Sites, Slide Clamp, SPIN-LOCK Connector, Latex Free, 22mL Priming Volume, 112" L, 20 drops/mL, 20/cs
375131	Pump Set, Non-Vented Spike, Roller Clamp, 150mL Burette Chamber w/ULTRASITE Valve, Free-Flow Protector Clip, Roller Clamp, 2 Injection Sites, Slide Clamp, SPIN-LOCK Connector, Latex Free, 22mL Priming Volume, 112" L, 60 drops/mL, 20/cs
US8812F	Pump Set, Universal Spike, Free-Flow Protector Clip, 1.2 Micron Filter, Roller Clamp, Slide Clamp, ULTRASITE Injection Site 6" Above Distal End, SPIN-LOCK Connector, DEHP-Free, 21.0mL Priming Volume, 109" Length, 15 drops/mL (352400), 50/cs
US8850	Blood Set, Y-Type, 2 Upper Roller Clamps, Drip Chamber w/170µ Blood Filter, Free Flow Protector Clip, ULTRASITE Y-Sites 6" Above Distal End, Roller Clamp, SPIN-LOCK Connection, 112" Length, 10 Drops/mL (352342), 50/cs
US8870	Pump Set, Female Adapter, Free Flow Protector Clip, Roller Clamp, ULTRASITE Y-Site 6" Above Distal End, SPIN-LOCK Connection, 112" Length (352343), 50/cs
V8850	Y-Type Blood Set, Non-Vented Spikes, 2 Upper Roller Clamps, Drip Chamber w/170 Micron Blood Filter, Free-Flow Protector Clip, Roller Clamp, Injection Site 6" Above Distal Connector, SPIN-LOCK Connector, Latex Free, 112" Length, 10 Drops/mL (352231), 50/cs

B BRAUN ADDITIVE CAPS


315000	Big Tab® Additive Cap, X-Large, Ribbed Finger Grip w/tamper-proof feature, Fits 28mm Bottle Tops, Sterile, 100/cs
--------	---

B BRAUN ADMIXTURE ACCESSORIES

- 415060 Backcheck Valve w/17G x 1" Needle For The Injection of Multiple Additives to an IV Bag Injection Port, Latex Free, DEHP-Free, 0.26 mL Priming Volume, 100/cs
- 415125 Transofix® Transfer Spike, Vented Double Spike For Rapid Transfer of Medication or Diluent Between Rubber Stoppered Containers, 100/cs
- 418105 Vented Spike Adapter To Permit Connection of a Non-Vented IV Spike to an IV Bottle Requiring a Vent, Latex Free, 20/cs
- N2050 Non-Vented Spike Adapter, Permits Connection of a Vented IV Spike to a Semi-Rigid Plastic IV Container, Latex Free, 24/cs
- N7890 Add-A-Vial® 20 mm Binary Connector, 200/cs
- N7893 Add-A-Via® 13 mm Binary Connector, 200/cs
- N7990 Add Base™ 20 mm Binary Connector, 200/cs
- N7993 Add Base™ 13 mm Binary Connector, 200/cs

B BRAUN ANCILLARY DISPOSABLE PRODUCTS

- 418020 Injection Cap, Distal Luer Locking Connection For Attachment To Female Luer Lock or Slip Connection, 0.2 mL Priming Volume, 1 1/8" Length, 100/cs
- 418025 Injection Cap, Distal Luer Connection For Attachment To Female Luer Lock or Slip Connection, 0.2 mL Priming Volume, 1 1/8" Length, 100/cs
- 418030 Injection Cap, Distal Luer Locking Connection For Attachment To Female Luer Lock or Slip Connection, 0.1 mL Priming Volume, 3/4" Length, 100/cs
- 456080 Fistula Needle Recirculation Connector, 100/cs

**B BRAUN BLOOD ADMINISTRATION SETS**

- V2203 Plasma Transfer Admin Set, Roller Clamp, Male Luer Slip, 27" Length 50/cs
- V2400 Admin Set, Drip Chamber w/170µ Blood Filter, Roller Clamp, Flashtube, Male Luer Slip, 31 mL Priming Volume, 76" Length, 10 Drops/mL, 50/cs
- V2410 Admin Set, Drip Chamber w/170µ Blood Filter, Roller Clamp, Injection Site 6" Above Distal End, Flashtube, Male Luer Slip, 30 mL Priming Volume, 75" Length, 10 Drops/mL, 50/cs
- V2460 Admin Set, Roller Clamp, Drip Chamber w/170µ Blood Filter, Pump Chamber, Injection Site 6" Above Distal End, Male Luer Slip, 70 mL Priming Volume, 85" Length, 10 Drops/mL, 50/cs
- V2500 Admin Set, Y-Type, 3 Roller Clamps, Drip Chamber w/170µ Blood Filter, Flashtube, Male Luer Slip, 34 mL Priming Volume, 87" Length, 10 Drops/mL, 50/cs
- V2510 Admin Set, Y-Type, 3 Roller Clamps, Drip Chamber w/170µ Blood Filter, Injection Site 6" Above Distal End, Flashtube, Male Luer Slip, 33 mL Priming Volume, 78" Length, 10 Drops/mL, 50/cs
- V2511 Admin Set, Y-Type, 2 Roller Clamps, Drip Chamber w/170µ Blood Filter, Controll® Clamp, Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 34mL Priming Volume, 98" Length, 10 Drops/mL, 50/cs
- V2520 Admin Set, Y-Type, 3 Roller Clamps, Drip Chamber w/170µ Blood Filter, Injection Sites 6" & 28" Above Distal End, 2-Piece Male Luer Lock, 36 mL Priming Volume, 96" Length, 10 Drops/mL, 50/cs
- V2560 Admin Set, Y-Type, 3 Roller Clamps, Drip Chamber w/170µ Blood Filter, Pump Chamber, Injection Site 6" Above Distal End, Male Luer Slip, 69 mL Priming Volume, 85" Length, 10 Drops/mL, 50/cs
- V2560-37 Admin Set, Y-Type, 3 Roller Clamps, Drip Chamber w/170µ Blood Filter, Pump Chamber, Injection Sites 6" & 41" Above Distal End, Male Luer Slip, 71 mL Priming Volume, 98" Length, 10 Drops/mL, 50/cs
- V2560-39 Admin Set, Y-Type, 3 Roller Clamps, Drip Chamber w/170µ Blood Filter, Pump Chamber, Injection Site 6" Above Distal End, 2-Piece Male Luer Lock, 80 mL Priming Volume, 148" Length, 10 Drops/mL, 50/cs
- V2950 Add-On Blood Filter Set, Drip Chamber w/170µ Blood Filter, Slide Clamp, Spike Adapter, 19 mL Priming Volume, 9" Length, 50/cs

B BRAUN CHEMO DISPENSING PIN

- 412014 Chemo Dispensing Pin w/Distal MINI-SPIKE®, Proximal Female Luer Lock, Built-In 0.2µ Hydrophobic Air-Venting Filter, 50/cs

B BRAUN FILTERED MEDICATION TRANSFER DEVICES

- 414995 1.2µ SUPOR® Aspiration/Injection Disc Filter, Proximal & Distal Luer Lock Connection Designed For Fungal & Particulate Filtration of Medication, Used For Aspiration or Injection, 0.3 mL Fluid Retention, Clear, 50/cs
- 415002 0.2µ SUPOR® Aspiration/Injection Disc Filter, Proximal & Distal Luer Lock Connections Designed For Bacteria Retentive Filtration of Medication, Used For Aspiration or Injection, 0.3 mL Fluid Retention, Green, 50/cs
- 415003 0.2µ PTFE Aspiration/Injection Disc Filter, Hydrophobic Membrane, Proximal Luer Lock & Distal Luer Slip Connections Designed Chemically Resistant Filter Designed For Sterile Filtration of Aggressive Solutions, For Aspiration or Injection, White, 50/cs
- 415005 Large Volume 5µ Conical Final Filter w/distal luer locking connection, Proximal Female Luer Lock Connects to Any Standard IV Device in Pharmacy, 50/cs
- 415020 FILTER STRAW®, 5µ Filter, 4" Flexible Straw For Fluid Aspiration From Glass Ampules, 100/cs
- 415021 FILTER STRAW®, 5µ Filter, 1 3/4" Flexible Straw For Fluid Aspiration From Glass Ampules, 100/cs
- 415025 Filter Needle, 5µ Filter in Female Luer Lock Connection, 20G x 1 1/2" Thinwall Needle For Withdrawal or Injection of Medication From Rubber-Stoppered Vial, 100/cs
- 415030 Filter Needle, 5µ Filter in Female Luer Lock Connection, 19G x 1 1/2" Thinwall Needle For Withdrawal or Injection of Medication From Rubber-Stoppered Vial, 100/cs
- 415035 Filter Needle, 5µ Filter in Female Luer Lock Connection, 20G x 1" Thinwall Needle For Withdrawal or Injection of Medication From Rubber-Stoppered Vial, 100/cs
- 415040 Filter Needle, 5µ Filter in Female Luer Lock Connection, 19G x 1" Thinwall Needle For Withdrawal or Injection of Medication From Rubber-Stoppered Vial, 100/cs
- 415041 High Flow Filter Needle, 19G x 7/8" Needle w/5µ Filter, Recommended For Use w/viscous medications, 100/cs
- 415042 Filter Needle II, Removable 5µ Filter, 19G x 1" Thinwall Needle For Withdrawal or Injection of Medication From Rubber-Stoppered Vial, 100/cs
- 418021 Filter Hub, 5µ Filter in Male/Female Luer Lock Hub For Use w/standard luer connections & SAFSITE System, 100/cs

B BRAUN FLUID DISPENSING CONNECTORS

- 415080 Fluid Dispensing Connector For Transfer of Blood to Vacutainer, 100/cs
- 415081 Fluid Dispensing Connector, Proximal Female Luer Lock & Distal Female Luer Slip Connections, Used For Transferring Medication From Syringe to Syringe or Between a Syringe & a Master Container w/male luer fitting, Blue, 100/cs

B BRAUN FLUID TRANSFER SETS

356043	Fluid Transfer Set, Vented Proximal Spike, Wide Bore Tubing, Distal Non-Vented Spike w/on/off clamp, 24" Length, 50/cs
356048	Fluid Transfer Set, Vented Proximal Spike, Wide Bore Tubing, On/Off Clamp, Distal Luer Slip, 30" Length, 50/cs
356049	Fluid Transfer Set, Proximal & Distal Non-Vented Spikes, Wide Bore Tubing w/on/off clamp, 24" Length, 50/cs
356050	Fluid transfer Set, Vented Spike w/gripping flange, Large Bore Tubing 23", Plastic On/Off Clamp, 17G Needle w/gripping flange, Used For Direct Transfer of Fluids Between Large Volume Parenterals, 23" Length, 50/cs
356052	Fluid Transfer Set, Same Features as 356050 But w/Built-In Latex Free Y-Site Medicinal Entry Point For Multiple Injections, 23" Length, 50/cs
356053	Fluid Transfer Set, Same Features as 356052 But w/30" Length & Distal Luer Lock Connection w/o needle, 30" Length, 50/cs
356054	Fluid Transfer Set, Same Features as 356048 But 22" Length w/Built-In Latex Free Y-Site Medicinal Entry Point & Distal Luer Slip w/o needle, 22" Length, 50/cs
356055	Fluid Transfer Set w/Built-In Large Volume 5µ Conical Filter, Vented Spike w/gripping flange, Large Bore Tubing 25", On/Off Clamp, Distal 17G Needle w/gripping flange, 25" Length, 50/cs
356057	Fluid Transfer Set, Same Features as 356055 But w/Latex Free Y-Site Medicinal Entry Above The 5µ Conical Filter, 25" Length, 50/cs
358000	Fluid Transfer Set, Vented Spike w/gripping flange, 0.2µ SUPOR® Filter, On/Off Clamp, Distal Luer Slip w/o needle, 35" Length, 50/cs
V1904	Fluid Transfer Set, Roller Clamp, 17G Unattached Needle, 25" Length, 50/cs

B BRAUN MICRO-COMP® ADMIXTURE COMPOUNDING SYSTEMS

473060	Standard Bore Connector Set, Male-To-Male Luer Lock Extension Set For Connecting 2 MICRO-COMP Units, 6" Length, 50/cs
--------	---

B BRAUN MULTI-AD® FLUID DISPENSING SYSTEMS

513500	Fluid Dispensing System, 10cc Luer Lock Syringe w/same features as 513000, 10/cs
513501	Fluid Dispensing System, 10cc Luer Lock Syringe w/same features as 513500 But w/30" length, Wide Diameter Transfer Set w/Weighted Stainless Sinkers & Tubing Clip, 10/cs
513502	Fluid Dispensing System, Hand-Held Style, 10cc Luer Lock Syringe w/automatic spring return, 2-Way Valve, Syringe Volume Adjusting Device & 41" Vented Transfer Set, For Reconstitution & Dispensing of IV Additives & Unit Dose Fluids, 10/cs
513506	Fluid Dispensing System, Hand-Held Style, 10cc Luer Lock Syringe w/automatic spring return, 2-Way Valve, Syringe Volume Adjusting Device Located Below Handle & 41" Vented Transfer Set, For Reconstitution & Dispensing IV Additives & Unit Dose Fluids, 10/cs
513508	Fluid Dispensing System, Hand-Held Style w/ergonomically designed hand grip, Same Features as 513502 Except Volume Adjustment Device Located Below Handle, 10/cs
513540	Fluid Dispensing System, 10cc Luer Lock Syringe w/same features as 513500, But w/a 40" Extension Set Affixed to the 2-Way Valve, 10/cs
513553	Fluid Dispensing System, 50cc Luer Lock Syringe w/2-way valve, 41" Vented Transfer Set, Distal Aspiration Valve, 3-Way Stopcock, Re-Circulates Air Into Reservoir Container From Cassette During Filling Procedure, 10/cs
513554	Fluid Dispensing System, 50cc Luer Lock Syringe w/2-Way Valve, 41" Vented Transfer Set, Permanently Connected 40" Luer Lock Extension Set, 10/cs

B BRAUN NEEDLE-FREE DISPENSING PINS

412000	MINI-SPIKE® Dispensing Pin w/SAFSITE® Valve, Automatic 2-Way Valve For Aspiration or Injection of Medication From Inverted Multi-Dose Vials, Utilizes a Bacterial Retentive Air-Venting Filter, 50/cs
412005	Standard Spike Dispensing Pin w/SAFSITE® Valve, Automatic 2-Way Valve For Aspiration or Injection of Medication From Inverted Multi-Dose Vials, Utilizes a Bacterial Retentive Air-Venting Filter, 50/cs
412006	Standard Spike Dispensing Pin w/SAFSITE® Valve, Same Features as 412005 Except w/luer lock clear shroud, 50/cs
412007	Standard Spike Dispensing Pin w/24" Extension Line & SAFSITE® Valve For Aspiration of Medication From Inverted Multi-Dose Vials, Utilizes a Bacterial Retentive Air-Venting Filter, 50/cs
412012	MINI-SPIKE® IV Additive Dispensing Pin For Preparing & Dispensing Diluent or Additive From Multi-Dose Rubber-Stoppered Vials, Utilizes a Bacterial Retentive Air-Venting Filter, 50/cs
412013	Standard Spike Dispensing Pin For Preparing & Dispensing Diluent or Additive From Multi-Dose Rubber-Stoppered Vials, Utilizes a Bacterial Retentive Air-Venting Filter, 50/cs
412022	Non-Vented Dispensing Pin w/1-Way Valve For Aspiration of Medication From Inverted Bags or Semi-Rigid Plastic Containers, Luer Lock Connection, 50/cs
412023	Non-Venting Dispensing Pin w/one-way valve & ULTRASITE® Valve, For Aspiration of Medication From Inverted Bags or Semi-Rigid Plastic Containers, DEHP-Free, Latex Free, 50/cs
412026	MINI-SPIKE Dispensing Pin w/ULTRASITE Valve For Rubber-Stoppered Vials, 2-Way Valve For Aspiration or Injection, Utilizes a Bacterial Retentive Air-Venting Filter, DEHP-Free, Latex Free, 50/cs
413500	Non-Vented Dispensing Pin w/SAFSITE Valve, Automatic 2-Way Valve For Aspiration or Injection of Medication From Inverted Bags or Semi-Rigid Plastic Containers, 50/cs
413501	Non-Vented Dispensing Pin w/SAFSITE Valve, Luer Lock Clear Shroud, Automatic 2-Way Valve For Aspiration or Injection of Medication From Inverted Bags or Semi-Rigid Plastic Containers, 50/cs
413503	Non-Vented MINI-SPIKE® Dispensing Pin For Withdrawal or Injection of Medication From Multi-Dose Rubber-Stoppered Vials, 100/cs
413504	Mini-Spike Dispensing Pin w/ULTRASITE Valve & Security Clip For Attachment to Standard 20mm Rubber-Stoppered Vials, Two-Way Valve For Aspiration & Injection, DEHP-Free, Latex Free, 50/cs
415019	MICRO-PIN™ Non-Vented, Single Use Dispensing Pin For Withdrawal or Injection of Medication From Rubber-Stoppered Vials, 100/cs

B BRAUN PLASMA VOLUME IV MANAGEMENT

L6511	500 mL HESpan®, 6% Hetastarch in 0.9% Sodium Chloride Injection, EXCEL® Container, 12/cs
-------	--

B BRAUN PVC-FREE FLUID PATH MICROBORE EXTENSION SETS

V6203	Extension Set w/PVC-Free Fluid Path, 0.03" ID, Male & Female Luer Locks, Use w/higher flow rate & viscous fluids, 0.6 mL Priming Volume, 36" Length, 50/cs
-------	--

B BRAUN SAFSITE® ACCESS PINS

415075	SAFSITE Access Pin, Shielded Plastic Pin For Accessing Y-Injection Site, 100/cs
415077	SAFSITE Access Pin w/valve, Shielded Plastic Pin For Accessing Y-Injection Site w/pre-attached SAFSITE valve, 100/cs
415079	SAFSITE Access Pin w/valve, Shielded Plastic Pin For Accessing Y-Injection Site w/pre-attached SAFSITE valve, 100/cs
415082	SAFSITE Access Pin, Shielded Plastic Pin For Accessing Y-Injection Site, 100/cs
415084	SAFSITE Access Pin, Shielded Plastic Pin For Accessing Y-Injection Site, 100/cs
415085	SAFSITE Access Pin w/valve, Shielded Plastic Pin For Accessing Y-Injection Site w/pre-attached SAFSITE valve, 100/cs

B BRAUN SAFSITE® REPLACEMENT CAPS

418017	Blue Catheter Cap, Non-Injectable, Male Luer Lock, 100/cs
418027	Blue Replacement Cap, 100/bx, 10 bx/cs
474900	White Replacement Cap, Female Luer Lock, 100/cs
654491	Blue Cap, Dual Function Blue Cap w/male & female end, 100/bx, 10 bx/cs
654496	RED CAP™, Dual Function Red Cap w/male & female end, 100/bx, 10 bx/cs
BW1000	Replacement Caps, Blue Male Luer Lock Cap & White Female Luer Lock Cap (418016), 100/cs

B BRAUN SAFSITE® VALVES

415062	Normally Closed Check Valve Requires Syringe Pressure to Open, Luer Lock, 0.12 mL Priming Volume, 100/cs
415067	SAFSITE Valve Allows Aspiration, Injection or Gravity Flow of Fluids, Luer Taper Operated, Normally Closed, Includes Extra B1000 Replacement Cap, 0.12 mL Priming Volume, 100/cs
415068	SAFSITE Valve Allows Aspiration, Injection or Gravity Flow of Fluids, Luer Taper Operated, Normally Closed, 0.12 mL Priming Volume, 100/cs
415069	SAFSITE Valve w/19G x 1" Needle Pre-Attached, Not Bonded, 0.26 mL Priming Volume, 100/cs
415078	Luer Adapter For Pre-Filled Rubber-Stoppered Cartridges, 100/cs

B BRAUN STIMUPLEX® NERVE STIMULATOR PRODUCTS

333542	STIMUPLEX Kit Containing STIMUPLEX-DIG Nerve Stimulator, STIMUPLEX Electrode Connecting Cable, 9V Battery & Operating Manual
333545	CONTIPLEX Connecting Cable
333680	20G x 6" Insulated Needle w/extension set, For Use w/STIMUPLEX Nerve Stimulator, 25/cs
333686	21G x 4" Insulated Needle w/extension set, For Use w/STIMUPLEX Nerve Stimulator, 25/cs
333690	22G x 2" Insulated Needle w/extension set, For Use w/STIMUPLEX Nerve Stimulator, 25/cs
333694	24G x 1" Insulated Needle w/extension set, For Use w/STIMUPLEX Nerve Stimulator, 25/cs

B BRAUN SYRINGE PUMP EXTENSION SETS

V5450	Extension Set, Male Luer Locks at Each End & 1 Slide Clamp, 0.7 mL Priming Volume, 31" Length, 50/cs
V5456	Extension Set, Male Luer Locks at Each End & 1 Slide Clamp, 1.3 mL Priming Volume, 63" Length, 50/cs
V6200	Extension Set w/PVC-Free Fluid Path, 0.02" ID, Male & Female Luer Locks, 0.3 mL Priming Volume, 36" Length, 50/cs
V6212	Extension Set w/PVC-Free Fluid Path, 0.02" ID, Male & Female Luer Locks, 0.4 mL Priming Volume, 60" Length, 50/cs
V6212-22	Extension Set w/PVC-Free Fluid Path & 22G Protected Needle, 0.02" ID, Male & Female Luer Locks, 0.5 mL Priming Volume, 60" Length, 50/cs
V6213	Extension Set w/PVC-Free Fluid Path, 0.03" ID, Male & Female Luer Locks, Use w/higher flow rates & viscous fluids, 0.8 mL Priming Volume, 60" Length, 50/cs
V6215	Extension Set w/PVC-Free Fluid Path & 0.2 Micron Filter, 0.03" ID, Male & Female Luer Locks, 1.3 mL Priming Volume, 60" Length, 50/cs
V6510	Extension Set, Microbore Tubing, 0.023" ID, Male & Female Luer Locks, 0.3 mL Priming Volume, 36" Length, 100/cs
V6512	Extension Set, Microbore Tubing, 0.023" ID, Male & Female Luer Locks, 0.4 mL Priming Volume, 60" Length, 100/cs
V6516	Extension Set, Microbore Tubing, 0.050" ID, Male & Female Luer Locks, Use w/higher flow rates & viscous fluids, 2.5 mL Priming Volume, 60" Length, 100/cs

B BRAUN TRANSFER DEVICES VIA DIRECT VACUUM

Used for transfer of medication from rubber-stoppered vial to vacuumized IV container.

415017	Filtered Double-Ended Transfer Needle, Proximal, 17G x 3/8" Needle Opposite 17G x 1" Needle, 100/cs
--------	---

B BRAUN ULTRASITE® NEEDLE-FREE ACCESS PIN

For use on in-line Y-injection site.

415093	ULTRASITE On 19G Protected Needle For Use On In-Line Y-Site, 100/cs
415095	Access Pin w/ULTRASITE Valve, 100/cs

B BRAUN ULTRASITE® NEEDLE-FREE VALVE


For aspiration, injection or gravity flow of fluid upon insertion of a male luer fitting

415110	ULTRASITE Capless Positive Pressure Valve, 100/cs
--------	---

**EXEL INJECTION PLUG**

Latex free. Used with any I.V. system for rapid injection. Sterile, non-toxic and non-pyrogenic. Ready for use.

26539	Injection Plug w/cap, 50/bx, 4 bx/cs
-------	--------------------------------------

**HYOGUARD COMFORT INFUSION SETS**

Sleek design, easy-to-hide under clothing or a swimsuit. Allergy tested dressing, eliminates most allergy problems for the user. Flexibility, 20 to 45 degree insertion angle according to individual people's needs.

999000	Comfort 43" Infusion Set (standard), 5 sets/bx, 20 bx/cs
999001	Comfort 23" Infusion Set (standard), 5 sets/bx, 20 bx/cs
999008	Comfort 43" Infusion Set (hospital), 10 sets/bx, 20 bx/cs
999009	Comfort 23" Infusion Set (hospital), 10 sets/bx, 20 bx/cs

KENDALL MONOJECT SYRINGE TIP CAPS

Flexible. Integral luer plug fits luer tip and luer lock syringes. Caps positioned upright in tray wells. Autoclavable. Sterile.

8881682085	Syringe Tip Cap, 100/bg, 10 bg/ctn
8881682101	Syringe Tip Cap, 10/tray, 100 tray/ctn
8881682119	Syringe Tip Cap, 25/tray, 40 tray/ctn

TERUMO SURFLO® INJECTION PLUGS

Large target area allows for easy needle penetration. Tamper resistant seal and easy to open packaging.

3SR-IP1	Luer Slip Injection Plug, 100/bx, 10 bx/cs
3SR-IP2	Luer Lock Injection Plug, 100/bx, 10 bx/cs



309570

309585

BD 3 ML SYRINGES & NEEDLES

- 301073 3mL Syringe Only w/Luer-Lok™ tip, Non-Sterile, Bulk, 1600/cs
- 309570 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 25ga x 5/8", 100/sp, 8 sp/cs
- 309571 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 23ga x 1", 100/sp, 8 sp/cs
- 309572 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 1", 100/sp, 8 sp/cs
- 309573 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 1 1/4", 100/sp, 8 sp/cs
- 309574 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 1 1/2", 100/sp, 8 sp/cs
- 309575 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga x 1", 100/sp, 8 sp/cs
- 309577 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga x 1 1/2", 100/sp, 8 sp/cs
- 309578 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 20ga x 1", 100/sp, 8 sp/cs
- 309579 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 20ga x 1 1/2", 100/sp, 8 sp/cs
- 309581 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 25ga x 1", 100/sp, 8 sp/cs
- 309582 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 25ga x 1 1/2", 100/sp, 8 sp/cs
- 309585 3mL Syringe Only w/Luer-Lok™ tip, 100/bx, 8 bx/cs
- 309586 3mL Syringe Only w/slip tip, 100/bx, 8 bx/cs
- 301077 3mL Syringe Only w/slip tip, Non-Sterile, Bulk, 1600/cs
- 305060 3mL Syringe w/blunt fill needle & Luer-Lok™ tip, 18ga x 1 1/2", 100/sp, 8 sp/cs
- 309569 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 3/4", 100/sp, 8 sp/cs
- 309580 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 18ga x 1 1/2", 100/sp, 8 sp/cs
- 309587 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 26ga Sub-Q x 5/8", 100/sp, 8 sp/cs
- 309588 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 23ga IM x 1", 100/sp, 8 sp/cs
- 309589 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 23ga IM x 1 1/2", 100/sp, 8 sp/cs
- 309598 3 mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga IV x 1", 100/sp, 8 sp/cs
- 309599 3mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga IV x 1 1/2", 100/sp, 8 sp/cs
- 309702 3mL Syringe w/luer slip tip, Sterile Convenience Tray Pak, 25/tray, 8 trays/cs



301603

309630

BD 5 ML SYRINGES & NEEDLES

- 301603 5mL Syringe Only w/slip tip, 100/sp, 4 sp/cs
- 309603 5mL Syringe Only w/Luer-Lok™ tip, 100/sp, 4 sp/cs
- 309630 5mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 1", 100/sp, 4 sp/cs
- 309631 5mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 1 1/2", 100/sp, 4 sp/cs
- 309632 5mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga x 1", 100/sp, 4 sp/cs
- 309633 5mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga x 1 1/2", 100/sp, 4 sp/cs
- 309634 5mL Syringe/Needle Combination w/Luer-Lok™ tip, 20ga x 1", 100/sp, 4 sp/cs
- 309635 5mL Syringe/Needle Combination w/Luer-Lok™ tip, 20ga x 1 1/2", 100/sp, 4 sp/cs
- 301027 5mL Syringe Only w/Luer-Lok™ tip, Non-Sterile, Bulk, 1400/cs
- 301028 5mL Syringe Only w/slip tip, Non-Sterile, Bulk, 1400/cs
- 305062 5mL Syringe w/blunt fill needle & Luer-Lok™ tip, 18ga x 1 1/2", 100/sp, 8 sp/cs
- 309703 5mL Syringe w/luer slip tip, Sterile Convenience Tray Pak, 25/tray, 8 trays/cs




301604

309640

BD 10 ML SYRINGES & NEEDLES

- 301604 10mL Syringe Only w/slip tip, 100/sp, 4 sp/cs
- 309604 10mL Syringe Only w/Luer-Lok™ tip, 100/sp, 4 sp/cs
- 309640 10mL Syringe/Needle Combination w/Luer-Lok™ tip, 22ga x 1", 100/sp, 4 sp/cs
- 309642 10mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga x 1", 100/sp, 4 sp/cs
- 309643 10mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga x 1 1/2", 100/sp, 4 sp/cs
- 309644 10mL Syringe/Needle Combination w/Luer-Lok™ tip, 20ga x 1", 100/sp, 4 sp/cs
- 309645 10mL Syringe/Needle Combination w/Luer-Lok™ tip, 20ga x 1 1/2", 100/sp, 4 sp/cs
- 301029 10mL Syringe Only w/Luer-Lok™ tip, Non-Sterile, Bulk, 850/cs
- 301030 10mL Syringe Only w/slip tip, Non-Sterile, Bulk, 850/cs
- 305064 10mL Syringe w/blunt fill needle & Luer-Lok™ tip, 18ga x 1 1/2", 100/sp, 4 sp/cs
- 305462 10mL Lab Syringe w/eccentric tip, Non-Sterile, 100/sp, 4 sp/cs
- 305482 10mL Syringe Only w/eccentric tip, 100/sp, 4 sp/cs
- 309605 10mL Syringe w/Luer-Lok™ tip, Sterile Convenience Pak, 20/tray, 12 trays/cs
- 309639 10mL Syringe/Needle Combination w/Luer-Lok™ tip, 21ga IV x 1 1/2" TW, 100/sp, 4 sp/cs
- 309695 10mL Luer-Lok Tip Control Syringe, 25/sp, 4 sp/cs




BD 20 ML SYRINGES

- 301625 20mL Syringe Only w/slip tip, 40/sp, 4 sp/cs
- 309661 20mL Syringe Only w/Luer-Lok™ tip, 40/sp, 4 sp/cs
- 300613 20mL Syringe Only w/eccentric tip, 120/sp, 4 sp/cs
- 301031 20mL Syringe Only w/Luer-Lok™ tip, Non-Sterile, Bulk, 325/cs
- 301032 20mL Syringe Only w/slip tip, Non-Sterile, Bulk, 325/cs
- 305617 20mL Syringe w/Luer-Lok™ tip, Sterile Convenience Pak Tray, 10/tray, 12 trays/cs

BD 30 ML SYRINGES

- 309651 30mL Syringe Only w/slip tip, 40/sp, 4 sp/cs
- 301033 30mL Syringe Only w/Luer-Lok™ tip, Non-Sterile, Bulk, 225/cs
- 301034 30mL Syringe Only w/slip tip, Non-Sterile, Bulk, 225/cs
- 305618 30mL Syringe w/Luer-Lok™ tip, Sterile Convenience Pak Tray, 10/tray, 12 trays/cs
- 309650 30mL Syringe Only w/Luer-Lok™ tip, 40/sp, 4 sp/cs



BD 60 ML SYRINGES

- 309654 60mL Syringe Only w/slip tip, 40/sp, 4 sp/cs
- 300866 60mL Syringe Only w/eccentric tip, 60/sp, 4 sp/cs
- 301035 60mL Syringe Only w/Luer-Lok™ tip, Non-Sterile, Bulk, 125/cs
- 301036 60mL Syringe Only w/slip tip, Non-Sterile, Bulk, 125/cs
- 309653 60mL Syringe Only w/Luer-Lok™ tip, 40/sp, 4 sp/cs
- 309680 60mL Syringe w/Luer-Lok™ tip, Sterile Convenience Pak Tray, 20/tray, 6 trays/cs




BD ECLIPSE™ NEEDLES & SYRINGES

- 305761 25ga x 1" Needle, 100/sp, 12 sp/cs
- 305762 23ga x 1" Needle, 100/sp, 12 sp/cs
- 305763 22ga x 1 1/2" Needle, 100/sp, 12 sp/cs
- 305764 21ga x 1" Needle, 100/sp, 12 sp/cs
- 305765 21ga x 1 1/2" Needle, 100/sp, 12 sp/cs
- 305766 18ga x 1 1/2" Needle, 100/sp, 12 sp/cs
- 305780 1mL 25ga x 5/8" Syringe w/needle, 50/sp, 6 sp/cs
- 305781 3mL 25ga x 5/8" Syringe w/needle, 50/sp, 6 sp/cs
- 305782 3mL 23ga x 1" Syringe w/needle, 50/sp, 6 sp/cs
- 305783 3mL 22ga x 1 1/2" Syringe w/needle, 50/sp, 6 sp/cs
- 305785 5mL 22ga x 1 1/2" Syringe w/needle, 50/sp, 6 sp/cs
- 305786 10mL 22ga x 1 1/2" Syringe w/needle, 50/sp, 6 sp/cs
- 305789 1mL 27ga x 1/2" Syringe w/needle, 50/sp, 6 sp/cs
- 301799 Bifurcated Needle (Non-CE), 50/sp, 6 sp/cs
- 305767 25ga x 1 1/2" Needle, 100/sp, 12 sp/cs
- 305768 22ga x 1" Needle, 100/sp, 12 sp/cs
- 305778 1mL 30ga x 1/2" Syringe w/needle, 50/sp, 6 sp/cs
- 305779 3mL 21ga x 1" Syringe w/needle, 50/sp, 6 sp/cs
- 305784 3mL 21ga x 1 1/2" Syringe w/needle, 50/sp, 6 sp/cs
- 305787 3mL 25ga x 1" Syringe w/needle, 50/sp, 6 sp/cs
- 305788 3mL 22ga x 1" Syringe w/needle, 50/sp, 6 sp/cs



BD INTEGRA™ SYRINGES

- 305269 3mL Syringe w/25ga x 5/8" Detachable Needle, 100/sp, 4 sp/cs
- 305270 3mL Syringe w/25ga x 1" Detachable Needle, 100/sp, 4 sp/cs
- 305271 3mL Syringe w/23ga x 1" Detachable Needle, 100/sp, 4 sp/cs
- 305272 3mL Syringe w/22ga x 1 1/2" Detachable Needle, 100/sp, 4 sp/cs
- 305273 3mL Syringe w/21ga x 1" Detachable Needle, 100/sp, 4 sp/cs
- 305274 3mL Syringe w/21ga x 1 1/2" Detachable Needle, 100/sp, 4 sp/cs
- 305278 1mL Syringe w/25ga x 1" Permanently Attached Needle, 100/sp, 8 sp/cs
- 305279 1mL TB Syringe w/26ga x 3/8" Permanently Attached Needle, 100/sp, 8 sp/cs
- 305282 1mL Insulin Syringe w/29ga x 1/2" Permanently Attached Needle, 100/spk, 8 sp/cs
- 305283 3mL Syringe Only, 100/sp, 4 sp/cs
- 305277 1mL Syringe w/25ga x 5/8" Permanently Attached Needle, 100/sp, 8 sp/cs
- 305298 1mL TB Syringe w/27ga x 1/2" Permanently Attached Needle, 100/sp, 8 sp/cs



BD SAFETY-LOK™ SAFETY SYRINGES & NEEDLES

- 329464 1mL Insulin Syringe, 29ga x 1/2" U-100 Ultra Fine, 100/sp, 5 sp/cs
- 305553 1mL TB Syringe, 27ga x 1/2" w/permanently attached needle, 100/sp, 5 sp/cs
- 305554 1mL TB Syringe, 25ga x 5/8" w/permanently attached needle, 100/sp, 5 sp/cs
- 305558 5mL Syringe Only w/Luer-Lok™ tip, 50/sp, 8 sp/cs
- 305559 10mL Syringe Only w/Luer-Lok™ tip, 50/sp, 8 sp/cs
- 305561 5mL Syringe, 21ga x 1 1/2" w/detachable needle, 50/sp, 8 sp/cs
- 305564 10mL Syringe, 21ga x 1 1/2" w/detachable needle, 50/sp, 8 sp/cs
- 309592 3mL Syringe, 25ga x 5/8" w/detachable needle, 100/sp, 8 sp/cs
- 309593 3mL Syringe, 22ga x 1 1/2" w/detachable needle, 100/sp, 8 sp/cs
- 309594 3mL Syringe, 23ga x 1" w/detachable needle, 100/sp, 8 sp/cs
- 309595 3mL Syringe, 21ga x 1 1/2" w/detachable needle, 100/sp, 8 sp/cs
- 309596 3mL Syringe, 22ga x 1" w/detachable needle, 100/sp, 8 sp/cs
- 309606 3mL Syringe Only w/Luer-Lok™ Tip, 100/sp, 8 sp/cs

TERUMO SURSAVER™ SYRINGES

Designed with a permanently attached needle which eliminates the needle hub. This design greatly reduces dead space as compared to standard syringes. The needle hub on standard syringes traps drugs that are not expelled after injection. Over the course of one year, these accumulated "wasted" drugs can represent a significant cost to your facility.

- 3SS01D2313 1cc Syringe w/23ga x 1/2" Fixed Needle, 100/bx, 5 bx/cs
- 3SS01D2516 1cc Syringe w/25ga x 5/8" Fixed Needle, 100/bx, 5 bx/cs
- 3SS01D2713 1cc Syringe w/27ga x 1/2" Fixed Needle, 100/bx, 5 bx/cs

TERUMO HYPODERMIC NEEDLE & SYRINGE

Clear syringe barrel. Bold, precise scale markings. Translucent needle hub for quick observation of flashback. Plastic blister package is tamper evident, for assurance of sterility. Peel-packs provide waste disposal cost savings. Syringes with needles offer maximum patient comfort delivered through ultimate sharpness and smoothness.

3SS-10LT	Pharmacy Tray, 10cc Luer Lock Tip w/out needle, 0.2cc Graduations, 25/tray, 8 tray/cs
3SS-20LT	Pharmacy Tray, 20cc Luer Lock Tip w/out needle, 1.0cc Graduations, 25/tray, 8 tray/cs
3SS-30LT	Pharmacy Tray, 30cc Luer Lock Tip w/out needle, 1.0cc Graduations, 10/tray, 12 tray/cs
3SS-60LT	Pharmacy Tray, 60cc Luer Lock Tip w/out needle, 1.0cc Graduations, 10/tray, 12 tray/cs
3SS-03L2025	3cc 20ga x 1" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2038	3cc 20ga x 1 1/2" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2125	3cc 21ga x 1" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2138	3cc 21ga x 1 1/2" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2225	3cc 22ga x 1" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2238	3cc 22ga x 1 1/2" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2325	3cc 23ga x 1" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2338	3cc 23ga x 1 1/2" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2516	3cc 25ga x 5/8" Syringe w/needle, 100/bx, 10 bx/cs
3SS-03L2525	3cc 25ga x 1" Syringe w/needle, 100/bx, 10 bx/cs
3SS-05L2025	5cc 20ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
3SS-05L2038	5cc 20ga x 1 1/2" Syringe w/needle, 100/bx, 6 bx/cs
3SS-05L2125	5cc 21ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
3SS-05L2138	5cc 21ga x 1 1/2" Syringe w/needle, 100/bx, 6 bx/cs
3SS-05L2225	5cc 22ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
3SS-05L2238	5cc 22ga x 1 1/2" Syringe w/needle, 100/bx, 6 bx/cs
3SS-10L2025	10cc 20ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
3SS-10L2038	10cc 20ga x 1 1/2" Syringe w/needle, 100/bx, 6 bx/cs
3SS-10L2125	10cc 21ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
3SS-10L2138	10cc 21ga x 1 1/2" Syringe w/needle, 100/bx, 6 bx/cs
3SS-10L2238	10cc 22ga x 1 1/2" Syringe w/needle, 100/bx, 6 bx/cs

TERUMO HYPODERMIC SYRINGE WITHOUT NEEDLE

Clear syringe barrel. Bold, precise scale markings. Latex free.

3SS-03L	3cc Syringe w/out needle, Luer Lock, 100/bx, 10 bx/cs
3SS-03S	3cc Syringe w/out needle, Luer Taper, 100/bx, 10 bx/cs
3SS-05L	5cc Syringe w/out needle, Luer Lock, 100/bx, 6 bx/cs
3SS-05S	5cc Syringe w/out needle, Luer Taper, 100/bx, 6 bx/cs
3SS-10ES	10cc Syringe w/out needle, Eccentric Luer Tip, 100/bx, 6 bx/cs
3SS-10L	10cc Syringe w/out needle, Luer Lock, 100/bx, 6 bx/cs
3SS-10S	10cc Syringe w/out needle, Luer Taper, 100/bx, 6 bx/cs
3SS-20ES	20cc Syringe w/out needle, Eccentric Luer Tip, 25/bx, 8 bx/cs
3SS-20L	20cc Syringe w/out needle, Luer Lock, 25/bx, 8 bx/cs
3SS-20S	20cc Syringe w/out needle, Luer Taper, 25/bx, 8 bx/cs
3SS-30ES	30cc Syringe w/out needle, Eccentric Luer Tip, 25/bx, 8 bx/cs
3SS-30L	30cc Syringe w/out needle, Luer Lock, 25/bx, 8 bx/cs
3SS-30S	30cc Syringe w/out needle, Luer Taper, 25/bx, 8 bx/cs
3SS-60C	60cc Syringe w/out needle, Catheter Tip, 25/bx, 4 bx/cs
3SS-60ES	60cc Syringe w/out needle, Eccentric Luer Tip, 25/bx, 4 bx/cs
3SS-60L	60cc Syringe w/out needle, Luer Lock, 25/bx, 4 bx/cs

TERUMO SURGUARD™ SAFETY NEEDLE PROTECTION SYSTEM

Single-handed needle protection is just a click away. Unique locking mechanism helps prevent needlestick injuries. One-handed activation meets key criteria for safety and ease of use. Compatible with all types of syringes.

3SG05M2813	.5cc 28ga x 1/2" Insulin Syringe w/permanently attached needle, 100/bx, 6 bx/cs
3SG05M2913	.5cc 29ga x 1/2" Insulin Syringe w/permanently attached needle, 100/bx, 6 bx/cs
3SG10M2813	1cc 28ga x 1/2" Insulin Syringe w/permanently attached needle, 100/bx, 6 bx/cs
3SG10M2913	1cc 29ga x 1/2" Insulin Syringe w/permanently attached needle, 100/bx, 6 bx/cs
3SG1825	18ga x 1" Safety Needle, 100/bx, 8 bx/cs
3SG1925	19ga x 1" Safety Needle, 100/bx, 8 bx/cs
3SG30M2913	.3cc 29ga x 1/2" Insulin Syringe w/permanently attached needle, 100/bx, 6 bx/cs



TERUMO SURGUARD2™ 3CC SYRINGES W/SAFETY NEEDLE

3SG2-03L2025	3cc Syringe w/20ga x 1" Safety Needle, 100/bx, 4 bx/cs
3SG2-03L2125	3cc Syringe w/21ga x 1" Safety Needle, 100/bx, 4 bx/cs
3SG2-03L2225	3cc Syringe w/22ga x 1" Safety Needle, 100/bx, 4 bx/cs
3SG2-03L2238	3cc Syringe w/22ga x 1.5" Safety Needle, 100/bx, 4 bx/cs
3SG2-03L2325	3cc Syringe w/23ga x 1" Safety Needle, 100/bx, 4 bx/cs
SG2-03L2138	21ga x 1.5" Safety Needle, 100/bx 4 bx/cs
SG2-03L2516	25ga x 5/8" Safety Needle, 100/bx, 4 bx/cs
SG2-03L2525	25ga x 1", 100/bx, 4 bx/cs



TERUMO SURGUARD2 5CC SYRINGE W/SAFETY NEEDLE

SG2-05L2025	20ga x 1" Safety Needle, 100/bx, 4 bx/cs
SG2-05L2038	20ga x 1.5" Safety Needle, 100/bx, 4 bx/cs
SG2-05L2138	21ga x 1.5" Safety Needle, 100/bx, 4 bx/cs



TERUMO SURGUARD2 10CC SYRINGE W/SAFETY NEEDLE

SG2-10L2025	20ga x 1" Safety Needle, 100/bx, 4 bx/cs
SG2-10L2038	20ga x 1.5" Safety Needle, 100/bx, 4 bx/cs
SG2-10L2125	21ga x 1.5" Safety Needle, 100/bx, 4 bx/cs



EXEL LUER LOCK SYRINGE WITH NEEDLE

Latex free.

26100	3cc Luer Lock Syringe w/needle, 25ga x 5/8", 100/bx, 10 bx/cs
26101	3cc Luer Lock Syringe w/needle, 23ga x 1", 100/bx, 10 bx/cs
26102	3cc Luer Lock Syringe w/needle, 22ga x 1", 100/bx, 10 bx/cs
26103	3cc Luer Lock Syringe w/needle, 22ga x 1 1/4", 100/bx, 10 bx/cs
26104	3cc Luer Lock Syringe w/needle, 22ga x 1 1/2", 100/bx, 10 bx/cs
26105	3cc Luer Lock Syringe w/needle, 21ga x 1", 100/bx, 10 bx/cs
26107	3cc Luer Lock Syringe w/needle, 21ga x 1 1/2", 100/bx, 10 bx/cs
26108	3cc Luer Lock Syringe w/needle, 20ga x 1", 100/bx, 10 bx/cs
26109	3cc Luer Lock Syringe w/needle, 20ga x 1 1/2", 100/bx, 10 bx/cs
26110	3cc Luer Lock Syringe w/needle, 18ga x 1 1/2", 100/bx, 10 bx/cs
26111	3cc Luer Lock Syringe w/needle, 25ga x 1", 100/bx, 10 bx/cs
26112	3cc Luer Lock Syringe w/needle, 25ga x 1 1/2", 100/bx, 10 bx/cs
26115	3cc Luer Lock Syringe w/needle, 22ga x 3/4", 100/bx, 10 bx/cs
26127	3cc Luer Lock Syringe w/needle, 21ga x 1 1/4", 100/bx, 10 bx/cs
26129	3cc Luer Lock Syringe w/needle, 27ga x 1 1/4", 100/bx, 10 bx/cs
26210	5cc Luer Lock Syringe w/needle, 22ga x 1", 100/bx, 8 bx/cs
26211	5cc Luer Lock Syringe w/needle, 22ga x 1 1/2", 100/bx, 8 bx/cs
26212	5cc Luer Lock Syringe w/needle, 21ga x 1", 100/bx, 8 bx/cs
26213	5cc Luer Lock Syringe w/needle, 21ga x 1 1/2", 100/bx, 8 bx/cs
26214	5cc Luer Lock Syringe w/needle, 20ga x 1", 100/bx, 8 bx/cs
26218	5cc Luer Lock Syringe w/needle, 21ga x 1 1/4", 100/bx, 8 bx/cs
26219	5cc Luer Lock Syringe w/needle, 22ga x 1 1/4", 100/bx, 8 bx/cs
26250	10cc Luer Lock Syringe w/needle, 22ga x 1", 100/bx, 8 bx/cs
26251	10cc Luer Lock Syringe w/needle, 22ga x 1 1/2", 100/bx, 8 bx/cs
26252	10cc Luer Lock Syringe w/needle, 21ga x 1", 100/bx, 8 bx/cs
26253	10cc Luer Lock Syringe w/needle, 21ga x 1 1/2", 100/bx, 8 bx/cs
26254	10cc Luer Lock Syringe w/needle, 20ga x 1", 100/bx, 8 bx/cs
26255	10cc Luer Lock Syringe w/needle, 20ga x 1 1/2", 100/bx, 8 bx/cs
26257	10cc Luer Lock Syringe w/needle, 18ga x 1", 100/bx, 8bx/cs



EXEL LUER SLIP SYRINGE WITH NEEDLE

Latex free.

26616	3cc Luer Slip Syringe w/needle, 26ga x 3/8", 100/bx, 10 bx/cs
26601	3cc Luer Slip Syringe w/needle, 22ga x 3/4", 100/bx, 10 bx/cs
26602	3cc Luer Slip Syringe w/needle, 25ga x 5/8", 100/bx, 10 bx/cs
26603	3cc Luer Slip Syringe w/needle, 23ga x 1", 100/bx, 10 bx/cs
26605	3cc Luer Slip Syringe w/needle, 22ga x 1", 100/bx, 10 bx/cs
26608	3cc Luer Slip Syringe w/needle, 21ga x 1", 100/bx, 10 bx/cs
26611	3cc Luer Slip Syringe w/needle, 20ga x 1", 100/bx, 10 bx/cs



EXEL LUER LOCK SYRINGES

Latex free.

26200	3cc Luer Lock Syringe, 100/bx, 10 bx/cs
26230	5cc Luer Lock Syringe, 100/bx, 8 bx/cs
26265	10cc Luer Lock Syringe, 100/bx, 8 bx/cs
26280	20cc Luer Lock Syringe, 50/bx, 10 bx/cs
26290	30cc Luer Lock Syringe, 50/bx, 5 bx/cs
26300	50-60cc Luer Lock Syringe, 25/bx, 6 bx/cs



EXEL LUER SLIP SYRINGES

Latex free.

26301	50-60cc Luer Slip Syringe, Eccentric, 25/bx, 6 bx/cs
26201	3cc Luer Slip Syringe, 100/bx, 10 bx/cs
26231	5cc Luer Slip Syringe, 100/bx, 8 bx/cs
26266	10cc Luer Slip Syringe, 100/bx, 8 bx/cs
26281	20cc Luer Slip Syringe, Eccentric, 50/bx, 10 bx/cs
26291	30cc Luer Slip Syringe, Eccentric, 50/bx, 5 bx/cs



HYO GUARD FUTURA™ SAFETY SYRINGE

The patented, state-of-the-art needle retraction design requires little change to injection technique, and virtually eliminates exposure time. One-handed operation keeps hands behind needle during procedures. No recapping via accessory shields or sheaths required. Barrel harbors the needle post-retraction. Safety mechanism disables device. Cannot be reused. Accurate dose administration. Sterile, non-toxic, and non-pyrogenic. Latex free.

103224	3cc 22ga x 1.5" Safety Syringe, 50/bx, 10 bx/cs
103233	3cc 23ga x 1" Safety Syringe, 50/bx, 10 bx/cs
103252	3cc 25ga x 5/8" Safety Syringe, 50/bx, 10 bx/cs



KENDALL MONOJECT 3CC SAFETY SYRINGES

Latex free. Safety shield covers entire needle length. True transport position. Restrictive end collar. Bold graduations: 0.1cc. Ultra-sharp, tri-beveled, anti-coring needle. Luer lock syringe tip. Rigid pack and Soft Pack packaging. Sterile.

8881033056	3cc Safety Syringe, 20 x 1 1/2 SP, 80/bx, 6 bx/cs
8881033130	3cc Safety Syringe, 21 x 1 SP, 80/bx, 6 bx/cs
8881033239	3cc Safety Syringe, 22 x 1 SP, 80/bx, 6 bx/cs
8881033254	3cc Safety Syringe, 22 x 1 1/2 SP, 80/bx, 6 bx/cs
8881033338	3cc Safety Syringe, 23 x 1 SP, 80/bx, 6 bx/cs
8881033908	3cc Safety Syringe, Luer Lock Tip, SP, 80/bx, 6 bx/cs
8881533056	3cc Safety Syringe, 20 x 1 1/2 RP, 100/bx, 5 bx/cs
8881533130	3cc, Safety Syringe, 21 x 1 RP, 100/bx, 5 bx/cs
8881533155	3cc Safety Syringe, 21 x 1 1/2 RP, 100/bx, 5 bx/cs
8881533239	3cc, Safety Syringe, 22 x 1 RP, 100/bx, 5 bx/cs
8881533254	3cc, Safety Syringe, 22 x 1 1/2 RP, 100/bx, 5 bx/cs
8881533510	3cc, Safety Syringe, 25 x 5/8 RP, 100/bx, 5 bx/cs
8881533908	3cc, Safety Syringe Luer Lock Tip RP, 100/bx, 5 bx/cs



KENDALL MONOJECT SOFTPACK 6CC SYRINGES

Latex free. Bold graduations: 0.2cc. Ultra-sharp, tri-beveled, anti-coring needle. Luer lock syringe tip on combination units. Color-coded translucent needle hub. Tip cap included on syringe only units. Color-coded gauge stripe and bar code on individual package. Shippable cartons. Sterile.

1180600555	6cc Syringe, Regular Tip, 100/bx, 4 bx/cs
1180600777	6cc, Syringe Luer Lock Tip, 100/bx, 4bx/cs
1181620112	6cc Syringe, 20 x 1 1/2, 100/bx, 4 bx/cs
1181621100	6cc, Syringe, 21 x 1, 100/bx, 4 bx/cs
1181621112	6cc Syringe, 21 x 1 1/2, 100/bx, 4 bx/cs
1181622112	6cc, Syringe, 22 x 1 1/2, 100/bx, 4 bx/cs



KENDALL MONOJECT SOFTPACK 12CC SYRINGES

0.2cc bold graduations. Ultra-sharp, tri-beveled, anti-coring needle. Luer lock syringe tip on combination units. Color-coded translucent needle hub. Tip cap included on syringe only units. Color-coded gauge stripe and bar code on individual package. Sterile.

1181200555	12cc Syringe, Regular Luer Tip, 80/bx, 6 bx/cs
1181200777	12cc Syringe Luer Lock Tip, 80/bx, 6 bx/cs
1181218100	12cc, Syringe, 18 x 1, 80/bx, 6 bx/cs
1181220112	12cc Syringe, 20 x 1 1/2, 80/bx, 6 bx/cs
1181221100	12cc Syringe, 21 x 1, 80/bx, 6 bx/cs
1181221112	12cc Syringe, 21 x 1 1/2, 80/bx, 6 bx/cs
1181222112	12cc Syringe, 22 x 1 1/2, 80/bx, 6 bx/cs



KENDALL MONOJECT 12CC SAFETY SYRINGES

Latex free. Safety shield covers entire needle length. True transport position. Restrictive end collar. Bold graduations: 0.2cc. Ultra-sharp, tri-beveled, anti-coring needle. Luer lock syringe tip. Rigid pack packaging. Sterile.

8881522018	12cc, Safety Syringe, 20 x 1 1/2, 50/bx, 6 bx/cs
------------	--

SAFEPRO PLUS SAFETY SYRINGES & NEEDLES

7032025-6	3cc 20ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
7032238-6	3cc 22ga x 1.5" Syringe w/needle, 100/bx, 6 bx/cs
7032325-6	3cc 23ga x 1" Syringe w/needle, 100/bx, 6 bx/cs
7032516-6	3cc 25ga x 5/8" Syringe w/needle, 100/bx, 6 bx/cs

SAFEPRO SAFETY SYRINGES & NEEDLES

6012516TB	1cc 25ga x 5/8" TB Syringe w/fixed needle, 100/bx, 10 bx/cs
6012610A	1cc 26ga x 3/8" Allergy Syringe w/fixed needle, 25/bx, 40 bx/cs
6012610TB	1cc 26ga x 3/8" TB Syringe w/fixed needle, 100/bx, 10 bx/cs
6012713A	1cc 27ga x 1/2" Allergy Syringe w/fixed needle, 25/bx, 40 bx/cs
6012713TB	1cc 27ga x 1/2" TB Syringe w/fixed needle, 100/bx, 10 bx/cs
6012813INS	1cc 28ga x 1/2" Insulin Syringe w/fixed needle, 100/bx, 10 bx/cs
6012913INS	1cc 29ga x 1/2" Insulin Syringe w/fixed needle, 100/bx, 10 bx/cs
6032025-6	3cc 20ga x 1" Syringe w/fixed needle, 100/bx, 6 bx/cs
6032125-6	3cc 21ga x 1" Syringe w/fixed needle, 100/bx, 6 bx/cs
6032225-6	3cc 22ga x 1" Syringe w/fixed needle, 100/bx, 6 bx/cs
6032238-6	3cc 22ga x 1.5" Syringe w/fixed needle, 100/bx, 6 bx/cs
6032325-6	3cc 23ga x 1" Syringe w/fixed needle, 100/bx, 6 bx/cs
6032516-6	3cc 25ga x 5/8" Syringe w/fixed needle, 100/bx, 6 bx/cs
6032525-6	3cc 25ga x 1" Syringe w/fixed needle, 100/bx, 6 bx/cs



BD INSULIN SYRINGES & NEEDLES

329622	1mL Insulin Syringe w/detachable needle, 25ga x 1", U-100, Orange/Blue, 100/sp, 8 sp/cs
329650	1mL Insulin Syringe Only, U-100 Slip Tip, 100/sp, 8 sp/cs
329651	1mL Insulin Syringe w/detachable needle, 25ga x 5/8", U-100, Orange/Blue, 100/sp, 8 sp/cs
329652	1mL Insulin Syringe w/detachable needle, 26ga x 1/2", U-100, Orange/Tan, 100/sp, 8 sp/cs
309300	3/10mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 28ga x 1/2", Self Contained, U-100 Micro-Fine™ IV, Orange, 200/sp, 8 sp/cs
309301	3/10mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 29ga x 1/2", Self Contained, U-100 Ultra-Fine™, 200/sp, 8 sp/cs
309305	1/2mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 28ga x 1/2", Self Contained, U-100 Micro-Fine™ IV, Orange, 200/sp, 8 sp/cs
309306	1/2mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 29ga x 1/2", Self Contained, U-100 Ultra-Fine™, 200/sp, 8 sp/cs
309309	1mL Insulin Syringe w/permanently attached needle, 28ga x 1/2", Self Contained, U-100 Micro-Fine™ IV, Orange, 200/sp, 8 sp/cs
309310	1mL Insulin Syringe w/permanently attached needle, 27ga x 5/8", Self Contained, U-100 Micro-Fine™ IV, Orange, 200/sp, 8 sp/cs
309311	1mL Insulin Syringe w/permanently attached needle, 29ga x 1/2", Self Contained, U-100 Ultra-Fine™, 200/sp, 8 sp/cs
309330	1/2mL 27ga x 1/2" Syringe, 10/bg, 20 bg/bx, 8 bx/cs
309331	1mL 27ga x 1/2" Syringe, 10/bg, 20 bg/bx, 8 bx/cs
328418	1mL Insulin Syringe w/permanently attached needle, 31ga x 5/16", Self-Contained, U-100 Ultra-Fine™ Short, 100/sp, 5 sp/cs
328438	3/10mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 31ga x 5/16", Self Contained, U-100 Ultra-Fine™ Short, 100/sp, 5 sp/cs
328440	3/10mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 31ga x 5/16", Self Contained, U-100 Ultra-Fine™ Short 1/2 Unit, 100/sp, 5 sp/cs
328468	1/2mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 31ga x 5/16", Self Contained, U-100 Ultra-Fine™ Short, 100/sp, 5 sp/cs
329410	1mL Insulin Syringe w/permanently attached needle, 28ga x 1/2", Self-Contained, U-100 Micro-Fine™ IV, Orange, 100/sp, 5 sp/cs
329411	1mL Insulin Syringe w/permanently attached needle, 29ga x 1/2", Self-Contained, U-100 Ultra-Fine™, 100/sp, 5 sp/cs
329412	1mL Insulin Syringe w/permanently attached needle, 27ga x 5/8", Self-Contained, U-100 Micro-Fine™ IV, Orange, 100/sp, 5 sp/cs
329420	1mL Insulin Syringe w/permanently attached needle, 28ga x 1/2", Blister Package, U-100 Micro-Fine™ IV, Single Unit Graduation, 100/sp, 5 sp/cs
329424	1mL Insulin Syringe w/permanently attached needle, 28ga x 1/2", Blister Package, U-100 Micro-Fine™ IV, Orange, 100/sp, 5 sp/cs
329461	1/2mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 28ga x 1/2", Blister Pkg, U-100 Micro-Fine™ IV, Orange, 100/sp, 5 sp/cs
329465	1/2mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 28ga x 1/2", Self-Contained, U-100 Ultra-Fine™ IV, Orange, 100/sp, 5 sp/cs
329466	1/2mL Lo-Dose™ Insulin Syringe w/permanently attached needle, 29ga x 1/2", Self-Contained, U-100 Ultra-Fine™ 100/sp, 5 sp/cs

BD SAFETYGLIDE™ INSULIN SYRINGES

Trays with permanently attached SafetyGlide™ needle.

305930	1mL Insulin Syringe w/28ga x 1/2" permanently attached needle, 100/sp, 4 sp/cs
305932	1/2mL Insulin Syringe w/28ga x 1/2" permanently attached needle, 100/sp, 4 sp/cs
305934	1/2mL Insulin Syringe w/30ga TW x 5/16" permanently attached needle, 100/sp, 4 sp/cs
305935	3/10mL Insulin Syringe w/29ga x 1/2" permanently attached needle, 100/sp, 4 sp/cs
305937	3/10mL Insulin Syringe w/31ga TW x 5/16" permanently attached needle, 100/sp, 4 sp/cs

TERUMO INSULIN SYRINGES

Low dead space design, for use with U-100 insulin only, with attached needle. Based on extreme needle sharpness and smooth plunger movement, they are designed for maximum comfort. Clear, bold scale markings assure dosage control and accuracy.

3SS05M2713	.5cc 27ga x 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs
3SS05M2813	.5cc 28ga x 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs
3SS05M2913	.5cc 29ga x 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs
3SS10M2713	1cc 27ga x 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs
3SS10M2813	1cc 28ga x 1 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs
3SS10M2913	1cc 29ga x 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs
3SS30M2913	.3cc 29ga x 1/2" Insulin Syringe w/attached needle, 100/bx, 5 bx/cs

EXEL COMFORT POINT INSULIN PEN NEEDLES

Optimum treatment with maximum simplicity. Sterile, non-toxic, non-pyrogenic. Double point system. Patient point with triple beveling for easy and atraumatic penetration. Cartridge point sharpened for a perfect punch with no fragmentation. Compatible with all available pen models. Easy flow of medication without any impurities.

26003	29ga x 1/2" Insulin Needle, 100/bx, 12 bx/cs
26004	31ga x 1/4" Insulin Needle, 100/bx, 12 bx/cs
26005	31ga x 1/3" Insulin Needle, 100/bx, 12 bx/cs



EXEL INSULIN SYRINGE WITH NEEDLE

Latex free.

26014	30ga x 5/16", 3/10cc Insulin Syringe w/needle, 10/bg, 10 bg/bx, 5 bx/cs
26015	30ga x 5/16", 1/2cc LO-DOSE Insulin Syringe w/needle, 10/bg, 10 bg/bx, 5 bx/cs
26016	30ga x 5/16", 1cc Insulin Syringe w/needle, 10/bg, 10 bg/bx, 5 bx/cs
26018	29ga x 1/2", 3/10cc Insulin Syringe w/needle, 10/bg, 10 bg/bx, 5 bx/cs
26026	28ga x 1/2", 1/2cc Insulin Syringe w/needle, Retail Pack, 10/bg, 10bg/bx, 5bx/cs
26027	28ga x 1/2", 1cc Insulin Syringe w/needle, Retail Pack, 10/bg, 10bg/bx, 5bx/cs
26028	29ga x 1/2", 1/2cc Insulin Syringe w/needle, Retail Pack, 10/bg, 10bg/bx, 5bx/cs
26029	29ga x 1/2", 1cc Insulin Syringe w/needle, Retail Pack, 10/bg, 10bg/bx, 5bx/cs
26032	28ga x 1/2", 1/2cc Insulin Syringe w/needle, Nursing Home Pack, 10/bg, 50 bg/cs
26034	28ga x 1/2", 1cc Insulin Syringe w/needle, Nursing Home Pack, 10/bg, 50 bg/cs



KENDALL MONOJECT INSULIN SAFETY SYRINGES

Latex free. Safety shield covers entire needle length. True transport position. Ultra-comfort, anti-coring, permanently attached needle. Accu-tip flat plunger tip. Low dead space. For U-100 insulin. Sterile.

- 888151110 1cc Insulin Safety Syringe, 29ga x 1/2", 100/bx, 5 bx/cs
- 888151136 1/2cc Insulin Safety Syringe, 29ga x 1/2", 100/bx, 5 bx/cs



1188127012



8881601358

KENDALL MONOJECT SOFTPACK INSULIN SYRINGES

Latex free. Ultra-comfort, anti-coring needle. Accu-tip flat plunger tip. Accu-scale graduations. Clear barrel construction and bolder scale. Low dead space permanently attached needle. For U-100 insulin. Sterile.

- 1188127012 1cc, Insulin Syringe, 27 x 1/2 Det, 100/bx, 5 bx/cs
- 1188128012 1cc, Insulin Syringe, 28 x 1/2, 100/bx, 5 bx/cs
- 8881600350 1/2cc Insulin Safety Syringe 29ga x 1/2", 100/bx, 3 bx/cs
- 8881601358 1cc Insulin Safety Syringe, 29ga x 1/2", 100/bx, 3 bx/cs

BD 1 ML SYRINGES & NEEDLES

- 309628 1mL Luer-Lok™ Syringe Only, 1/100 mL, 100/bx, 8 sp/cs
- 301025 1mL Syringe Only w/slip tip, Non-Sterile, Bulk, 3000/cs
- 309597 1mL Syringe/Needle Combination, 26ga SUB-Q x 5/8", Slip Tip, 1/100mL, 100/sp, 8 sp/cs
- 309637 1mL Luer-Lok™ Syringe/Needle Combination, 20ga x 1", 100/sp, 8 sp/cs
- 309701 1mL Syringe, Luer Slip Tip, Sterile Convenience Tray Pak, 25/tray, 8 trays/cs

BD SAFETYGLIDE™ TUBERCULIN SYRINGES

- 305945 1mL TB Syringe w/27ga x 1/2" Regular Bevel Permanently Attached Needle, 100/sp, 4 sp/cs
- 305946 1mL TB Syringe w/26ga x 3/8" Intradermal Bevel Permanently Attached Needle, 100/sp, 4 sp/cs

BD TUBERCULIN SYRINGES & NEEDLES

- 309602 1mL Tuberculin Syringe Only, Slip Tip, 100/sp, 8 sp/cs
- 309623 1mL Tuberculin Syringe w/detachable needle, Slip Tip, 27ga x 1/2", 100/sp, 8 sp/cs
- 309624 1mL Tuberculin Syringe w/detachable needle, Slip Tip, 21ga x 1", 100/sp, 8 sp/cs
- 309625 1mL Tuberculin Syringe w/detachable needle, Slip Tip, 26ga x 3/8", Intradermal Bevel, 100/sp, 8 sp/cs
- 309626 1mL Tuberculin Syringe w/detachable needle, Slip Tip, 25ga x 5/8", 100/sp, 8 sp/cs
- 305620 1/2mL Tuberculin Syringe w/permanently attached needle, 27ga x 1/2", Regular Bevel, 100/sp, 5 sp/cs



EXEL TB TUBERCULIN SYRINGES

Latex free. Zero Dead Space, P.A. needle.

- 26046 Tuberculin Syringe w/needle, 25ga x 5/8", Zero Dead Space, 10/bg, 10 bg/bx, 5 bx/cs



EXEL TB TUBERCULIN SYRINGES WITH LUER SLIP

Latex free. 1 cc.

- 26040 Tuberculin Syringe w/needle, 27ga x 1/2", 100/bx, 10 bx/cs
- 26043 Tuberculin Syringe w/needle, 26ga x 1/2", 100/bx, 10 bx/cs
- 26044 Tuberculin Syringe w/needle, 25ga x 5/8", 100/bx, 10 bx/cs
- 26048 Syringe Only, Slip Tip, 100/bx, 10 bx/cs



KENDALL MONOJECT SOFTPACK TUBERCULIN SYRINGES

Latex free. ULTRA-COMFORM, anti-coring needle. ACCU-TIP flat plunger tip. ACCU-SCALE graduations. Low dead space permanently attached needle. Shippable cartons. Sterile.

- 1180100555 1cc, TB Syringe Regular Tip, 100/bx, 5 bx/cs
- 1180125058 1cc, TB Syringe w/25 x 5/8 Det Needle, 100/bx, 5 bx/cs
- 1180125158 1cc, TB Syringe w/25 x 5/8 Needle, 100/bx, 5 bx/cs
- 1180126038 1cc, TB Syringe w/26 x 3/8 Det Needle, 100/bx, 5 bx/cs
- 1180127012 1cc, TB Syringe w/27 x 1/2 Det Needle, 100/bx, 5 bx/cs
- 1180128012 1cc, TB Syringe w/28 x 1/2 Needle, 100/bx, 5 bx/cs
- 1180528012 1/2cc, TB Syringe with 25 x 1/2 Needle, 100/bx, 5 bx/cs



KENDALL MONOJECT TUBERCULIN SAFETY SYRINGES

Latex free. Safety shield covers entire needle length. True transport position. ULTRA-COMFORT, anti-coring, permanently attached needle. ACCU-TIP flat plunger tip. ACCU-SCALE graduations: 0.01cc. Low dead space. Sterile.

- 8881511201 1cc, TB Safety Syringe, 28 x 1/2, 100/bx, 5 bx/cs
- 8881511235 1cc, TB Safety Syringe, 25 x 5/8, 100/bx 5 bx/cs



KENDALL MONOJECT TUBERCULIN SYRINGES

Latex free. Ultra-comfort, anti-coring needle. Accu-scale graduations: 0.01cc. Accu-tip flat plunger tip. Sterile.

- 8881500105 1/2cc, TB Syringe, 28 x 1/2, 100/bx, 5 bx/cs
- 8881501160 1cc TB Syringe, 25ga x 5/8", 100/bx, 5 bx/cs
- 8881501178 1cc, TB Syringe, 26 x 3/8, 100/bx, 5 bx/cs
- 8881501368 1cc, TB Syringe, 27 x 1/2, 100/bx, 5 bx/cs
- 8881501400 1cc TB Syringe Regular Luer Tip, 100/bx, 5 bx/cs
- 8881501608 1cc, TB Syringe, 28 x 1/2", Brown, 100/bx, 5 bx/cs



TERUMO SURGUARD2™ 1CC TB SYRINGE W/SAFETY NEEDLE

- 3SG2-01T2609 1cc Syringe w/26ga x 3/8" Safety Needle, 100/bx, 4 bx/cs
- 3SG2-01T2516 1cc Syringe w/25ga x 5/8" Safety Needle, 100/bx, 4 bx/cs
- 3SG2-01T2713 1cc Syringe w/27ga x 1/2" Safety Needle, 100/bx, 4 bx/cs

TERUMO TB SYRINGES

Provides clear, bold scale markings to assure dosage control and accuracy. The package is tamper evident, for greater protection of the sterile contents.

- 3SS-01T 1cc TB Syringe Only, 100/bx, 10 bx/cs
- 3SS-01T2516 1cc 25ga x 5/8" TB Syringe w/removable needle, 100/bx, 10 bx/cs
- 3SS-01T2609 1cc 26ga x 3/8" TB Syringe w/removable needle, 100/bx, 10 bx/cs
- 3SS-01T2713 1cc 27ga x 1/2" TB Syringe w/removable needle, 100/bx, 10 bx/cs

BD ALLERGIST TRAYS

Trays with permanently attached SafetyGlide™ needle.

- 305535 1/2mL Allergist Tray w/permanently attached needle, 27ga x 1/2", Regular Bevel, 25/tray, 40 trays/cs
- 305536 1/2mL Allergist Tray w/permanently attached needle, 27ga x 3/8", Intradermal Bevel, 25/tray, 40 trays/cs
- 305537 1mL Allergist Tray w/permanently attached needle, 26ga x 1/2", Regular Bevel, 25/tray, 40 trays/cs
- 305538 1mL Allergist Tray w/permanently attached needle, 26ga x 1/2", Intradermal Bevel, 25/tray, 40 trays/cs
- 305539 1mL Allergist Tray w/permanently attached needle, 26ga x 3/8", Intradermal Bevel, 25/tray, 40 trays/cs
- 305540 1mL Allergist Tray w/permanently attached needle, 27ga x 1/2", Regular Bevel, 25/tray, 40 trays/cs
- 305541 1mL Allergist Tray w/permanently attached needle, 27ga x 3/8", Intradermal Bevel, 25/tray, 40 trays/cs
- 305542 1mL Allergist Tray w/permanently attached needle, 27ga x 3/8", Regular Bevel, 25/tray, 40 trays/cs

BD SAFETYGLIDE™ ALLERGIST TRAYS

- 305950 Allergist Tray, 1mL w/27ga x 1/2" Regular Bevel, 25/tray, 40 trays/cs
- 305951 Allergist Tray, 1mL w/26ga x 3/8" Intradermal Bevel, 25/tray, 40 trays/cs

BD ALLERGY SYRINGE

- 305500 1mL Allergy Syringe w/permanently attached needle, 28ga x 1/2", Regular Bevel, 100/sp, 5 sp/cs

TERUMO ALLERGY SYRINGES & TRAYS

Low dead space for maximum cost savings. The barrel is clear and easy to read. Sharp needles minimize patient discomfort.

- 3SS-01A2713 27ga x 1/2" Allergy Syringe w/regular bevel needle, 100/bx, 5 bx/cs
- 3SS01A2609T 26ga x 3/8" Allergy Syringe w/regular bevel needle, 25/tray, 40 tray/cs
- 3SS01A2713T 27ga x 1/2" Allergy Syringe w/regular bevel needle, 25/tray, 40 tray/cs
- 3SS10A2709IDT 27ga x 3/8" Allergy Syringe w/intradermal bevel needle, 25/tray, 40 tray/cs
- 3SS10A2713IDT 27ga x 1/2" Allergy Syringe w/intradermal bevel needle, 25/tray, 40 tray/cs

TERUMO SURGUARD SAFETY ALLERGY TRAYS

- SG01A2613T 1cc 26ga x 1/2" Allergy Syringe Tray w/regular bevel needle, 25/tray, 40 tray/cs
- SG01A2713T 1cc 27ga x 1/2" Allergy Syringe Tray w/regular bevel needle, 25/tray, 40 tray/cs
- SG10A2609IDT 1cc 26ga x 3/8" Allergy Syringe Tray w/intradermal bevel needle, 25/tray, 40 tray/cs
- SG10A2709IDT 1cc 27ga x 3/8" Allergy Syringe Tray w/intradermal bevel needle, 25/tray, 40 tray/cs

BD EPILOR™ LOSS OF RESISTANCE SYRINGE

- 405291 7cc Luer-Lok™ Plastic Loss Of Resistance Syringe, 10/sp, 5 sp/cs
- 405292 7cc Luer-Slip Plastic Loss Of Resistance Syringe, 10/sp, 5 sp/cs

BD ORAL SYRINGE WITH TIP CAP

- 305207 1mL Amber Syringe, 100/sp, 5 sp/cs
- 305208 5mL Amber Syringe, 100/sp, 5 sp/cs
- 305209 10mL Amber Syringe, 100/sp, 5 sp/cs
- 305210 3mL Amber Syringe, 100/sp, 5 sp/cs
- 305217 1mL Clear Syringe, 100/sp, 5 sp/cs
- 305218 5mL Clear Syringe, 100/sp, 5 sp/cs
- 305219 10mL Clear Syringe, 100/sp, 5 sp/cs
- 305220 3mL Clear Syringe, 100/sp, 5 sp/cs

BD SYRINGE TIP CAP

- 305819 Syringe Tip Cap, Sterile, 200/sp, 10 sp/cs
- 308341 Luer Tip Cap, 10/tray, 200 trays/cs

BD MULTIFIT™ REUSABLE SYRINGES

- 512130 2ML Syringe Zone 1
- 512131 5ML Syringe Zone 1
- 512132 10ML Syringe Zone 1
- 512133 20ML Syringe Zone 1
- 512134 30ML Syringe Zone 1
- 512135 50ML Syringe Zone 1
- 512136 3ML Syringe Zone 1
- 512140 2ML Syringe Zone 2
- 512141 5ML Syringe Zone 2
- 512142 10ML Syringe Zone 2
- 512143 20ML Syringe Zone 2
- 512144 30ML Syringe Zone 2
- 512145 50ML Syringe Zone 2
- 512146 3ML Syringe Zone 2
- 512150 2ML Syringe Zone 3
- 512151 5ML Syringe Zone 3
- 512152 10ML Syringe Zone 3
- 512153 20ML Syringe Zone 3
- 512154 30ML Syringe Zone 3
- 512155 50ML Syringe Zone 3
- 512156 3ML Syringe Zone 3
- 512430 B2ML Syringe Zone 1, 36/cs
- 512432 B10ML Syringe Zone 1, 36/cs
- 512433 B20ML Syringe Zone 1, 12/cs
- 512434 B30ML Syringe Zone 1, 6/cs
- 512435 B50ML Syringe Zone 1, 6/cs
- 512440 B2ML Syringe Zone 2, 36/cs
- 512441 B5ML Syringe Zone 2, 36/cs
- 512442 B10ML Syringe Zone 2, 36/cs
- 512443 B20ML Syringe Zone 2, 12/cs
- 512444 B30ML Syringe Zone 2, 6/cs
- 512445 B50ML Syringe Zone 2, 6/cs
- 512450 B2ML Syringe Zone 3, 36/cs
- 512452 B10ML Syringe Zone 3, 36/cs
- 512453 B20ML Syringe Zone 3, 12/cs
- 512454 B30ML Syringe Zone 3, 6/cs
- 512455 B50ML Syringe Zone 3, 6/cs
- 512463 B20MM Syringe Zone 1, 12/cs
- 512464 B30MM Syringe Zone 1, 6/cs
- 512465 B50MM Syringe Zone 1, 6/cs
- 512470 B2MM Syringe Zone 2, 36/cs
- 512471 B5MM Syringe Zone 2, 36/cs
- 512472 B10MM Syringe Zone 2, 36/cs
- 512473 B20MM Syringe Zone 2, 12/cs
- 512474 B30MM Syringe Zone 2, 6/cs
- 512475 B50MM Syringe Zone 2, 6/cs
- 512480 B2MM Syringe Zone 3, 36/cs
- 512481 B5MM Syringe Zone 3, 36/cs
- 512482 B10MM Syringe Zone 3, 36/cs
- 512484 B30MM Syringe Zone 3, 6/cs
- 512485 B50MM Syringe Zone 3, 6/cs
- 512900 3MLC 3mL Control Syringe w/Luer-Lok™ tip, Zone 1
- 512901 3MLC 3mL Control Syringe w/Luer-Lok™ tip, Zone 2
- 512902 3MLC 3mL Control Syringe w/Luer-Lok™ tip, Zone 3
- 512903 5MLC 5mL Control Syringe w/Luer-Lok™ tip, Zone 1
- 512904 5MLC 5mL Control Syringe w/Luer-Lok™ tip, Zone 2
- 512905 5MLC 5mL Control Syringe w/Luer-Lok™ tip, Zone 3
- 512906 10MLC 10mL Control Syringe w/Luer-Lok™ tip, Zone 1
- 512907 10MLC 10mL Control Syringe w/Luer-Lok™ tip, Zone 2
- 512908 10MLC 10mL Control Syringe w/Luer-Lok™ tip, Zone 3

BD YALE™ SYRINGES

- 512300 2Y 2mL Syringe
- 512301 5Y 5mL Syringe
- 512302 10Y 10mL Syringe
- 512303 20Y 20mL Syringe
- 512304 30Y 30mL Syringe
- 512305 50Y 50mL Syringe
- 512306 3Y 3mL Syringe
- 512310 2YL 2mL Syringe
- 512311 5YL 5mL Syringe
- 512312 10YL 10mL Syringe
- 512313 20YL 20mL Syringe
- 512314 30YL 30mL Syringe
- 512315 50YL 50mL Syringe
- 512316 3YL 3mL Syringe
- 512317 100YL 100mL Syringe
- 512001 1/4YT 1/4mL in 1/100mL & 1/2 min. TB Syringe, Precision w/glass luer tip
- 512002 1/2YT 1/2mL in 1/100mL & 1/2 min. TB Syringe, Precision w/glass luer tip
- 512004 1YT 1mL in 1/100 mL & 1/2 min. TB Syringe, Precision w/glass luer tip, 6/cs
- 512019 2YB 2mL in 1/20 mL TB Syringe, Precision w/glass luer tip
- 512027 1YTL 1mL Luer-Lok Tip TB Syringe, Precision w/glass luer tip

BUSSE KRATON™ IRRIGATION SYRINGES

Reusable, red ear/ulcer syringe. Sturdy, flexible, and latex free allows increased suction, better control and a better grip when wet. Packed sterile in a peel-back blister pouch.

- 140 2 oz. Syringe, Clean, 200/cs
- 147 3 oz. Syringe, Clean, 100/cs
- 148 2 oz. Syringe, Sterile, 24/cs
- 149 3 oz. Syringe, Sterile, 12/cs



EXEL CATHETER TIP SYRINGES

Latex free.

- 26292 30cc Catheter Tip Syringe, Eccentric, 50/bx, 5 bx/cs
- 26302 50-60cc Catheter Tip Syringe, Eccentric, 25/bx, 6 bx/cs
- 26302BN 60cc Catheter Tip Syringe, Non-Sterile, Bulk, 400/cs



KENDALL MONOJECT CURVED TIP SYRINGE

Ungraduated, 12ml capacity curved tip syringe. Reusable plastic finger grip extender. Ideal for pressure irrigation.

- 8881412012 Curved Tip Syringe, 50/bx, 10 bx/cs



KENDALL MONOJECT ORAL MEDICATION SYRINGES

Permanent, heat-etched dual graduations. Tip will not accept hypodermic needle. Polypropylene. Ribbed tip caps packaged separately.

8881901006	Amber 1mL Syringe, 100/bx, 5 bx/cs
8881901014	Clear 1mL Syringe, 100/bx, 5 bx/cs
8881903002	Clear 3mL Syringe, 100/bx, 5 bx/cs
8881903010	Amber 3mL Syringe, 100/bx, 5 bx/cs
8881906005	Amber 6mL Syringe, 100/bx, 5 bx/cs
8881906104	Clear 6mL Syringe, 100/bx, 5 bx/cs
8881907003	Amber 10mL Syringe, 100/bx, 5 bx/cs
8881907102	Clear 10mL Syringe, 100/bx, 5 bx/cs

B BRAUN ORAL SYRINGES

115000	1 cc Clear Polypropylene Syringe w/sealing tip cap, 100/cs
115015	1 cc Amber Light-Resistant Polypropylene Syringe w/sealing tip cap, 500/cs
133000	3 cc Clear Polypropylene Syringe w/sealing tip cap, 100/cs
133015	3 cc Amber Light-Resistant Polypropylene Syringe w/sealing tip cap, 500/cs
170010	5 cc Clear Polypropylene Syringe w/sealing tip cap, 100/cs
170025	5 cc Amber Light-Resistant Polypropylene Syringe w/sealing tip cap, 500/cs
196015	10 cc Amber Light-Resistant Polypropylene Syringe w/sealing tip cap, 500/cs

WELCON FLAT TOP PISTON SYRINGE


Bulk packaged.

W018	Piston Syringe, 150/cs
W034	Piston Syringe, 30/cs
W035	Piston Syringe, 300/cs
W045	Piston Syringe, 30/bg, 10 bg/cs

WELCON THUMB CONTROL RING PISTON SYRINGE

Bulk packaged.

W016	Piston Syringe, 150/cs
W031	Piston Syringe, 30/cs
W032	Piston Syringe, 300/cs
W042	Piston Syringe, 30/bg, 10 bg/cs



AS00501



AS00501S

AMSINO AMSURE® EAR/ULCER SYRINGE

Latex free.

AS00501	1 oz Vinyl Ear/Ulcer Syringe, Non-Sterile, 300/cs
AS00501S	1 oz Ear/Ulcer Syringe, Form Fill Seal Package, Sterile, 50/cs
AS00502	2 oz Vinyl Ear/Ulcer Syringe, Non-Sterile, 200/cs
AS00502S	2 oz Ear/Ulcer Syringe, Form Fill Seal Package, Sterile, 50/cs
AS00503	3 oz Vinyl Ear/Ulcer Syringe, Non-Sterile, 150/cs
AS00503S	1 oz Ear/Ulcer Syringe, Form Fill Seal Package, Sterile, 50/cs

B BRAUN PERIFIX® PLASTIC LOSS-OF-RESISTANCE SYRINGES

332150	10 cc Plastic Luer Slip Loss-of-Resistance Syringe, 50/cs
332152	8 cc Plastic Luer Lock Loss-of-Resistance Syringe, 50/cs

B BRAUN SYRINGE CAPS

418004	Tamper-Evident Syringe Cap, Used For Capping Luer Lock Tip IV Syringes, Red, 10/pk, 50 pk/cs
418010	MULTI-AD® Syringe Cap, Used For Capping Luer Tip IV Syringes, Green, 10/pk, 50 pk/cs
418013	MULTI-AD® Luer Lock Syringe Cap, Used For Capping Luer Slip or Luer Lock IV Syringes, Blue, 10/pk, 50 pk/cs
418200	Flexible Syringe Cap, Used For Capping Luer Tip IV Syringes, White, 30/pk, 20 pk/cs
418202	Flexible Syringe Cap, Used For Capping Luer Tip IV Syringes, Blue, 30/pk, 20 pk/cs
418203	Flexible Syringe Cap, Used For Capping Luer Tip IV Syringes, Yellow, 30/pk, 20 pk/cs
418204	Flexible Syringe Cap, Used For Capping Luer Tip IV Syringes, White, 10/pk, 50 pk/cs
418205	Flexible Syringe Cap, Used For Capping Luer Tip IV Syringes, Blue, 10/pk, 50 pk/cs
N7820	ADD-PRO® Syringe Caps, Used For Capping Luer Tip IV Syringes, Beige, 20/pk, 25 pk/cs
418012	MULTI-AD® Luer Lock Syringe Cap, Used For Capping Luer Slip or Luer Lock IV Syringes, Red, 10/pk, 50 pk/cs

BD CORNWALL™ DISPOSABLE SYRINGE SYSTEM

305224	10 mL Fluid Dispensing Syringe, 10/cs
305225	Syringe Tip Connector, 50/cs

BD CORNWALL™ SYRINGES

512190	1240S 1 mL Control Syringes w/Luer-Lok™ tip
512191	1250S 2 mL Control Syringes w/Luer-Lok™ tip
512192	1260S 5 mL Control Syringes w/Luer-Lok™ tip
512193	1270S 10 mL Control Syringes w/Luer-Lok™ tip

BD POSIFLUSH™ HEPARIN LOCK FLUSH SYRINGES

306509	Heparin Lock Flush Syringes, 10 Units/mL, 6mL Fill in 10mL Syringe, 30/sp, 4 sp/cs
306510	Heparin Lock Flush Syringes, 10 Units/mL, 5mL Fill in 10mL Syringe, 30/sp, 4 sp/cs
306511	Heparin Lock Flush Syringes, 10 Units/mL, 5mL Fill in 5mL Syringe, 30/sp, 4 sp/cs
306512	Heparin Lock Flush Syringes, 10 Units/mL, 3mL Fill in 3mL Syringe, 30/sp, 4 sp/cs
306513	Heparin Lock Flush Syringes, 100 Units/mL, 5mL Fill in 10mL Syringe, 30/sp, 4 sp/cs
306514	Heparin Lock Flush Syringes, 100 Units/mL, 3mL Fill in 10mL Syringe, 30/sp, 4 sp/cs
306515	Heparin Lock Flush Syringes, 100 Units/mL, 5mL Fill in 5mL Syringe, 30/sp, 4 sp/cs
306516	Heparin Lock Flush Syringes, 100 Units/mL, 3mL Fill in 5mL Syringe, 30/sp, 4 sp/cs
306517	Heparin Lock Flush Syringes, 100 Units/mL, 3mL Fill in 3mL Syringe, 30/sp, 4 sp/cs
306521	Heparin Lock Flush Syringes, 10 Units/mL, 3mL Fill in 10mL Syringe, 30/sp, 4 sp/cs
306525	Heparin Lock Flush Syringes, 10 Units/mL, 5mL Fill in 10mL Syringe w/blunt plastic cannula, 30/sp, 4 sp/cs
306527	Heparin Lock Flush Syringes, 10 Units/mL, 2mL Fill in 3mL Syringe, 30/sp, 4 sp/cs
306528	Heparin Lock Flush Syringes, 100 Units/mL, 2mL Fill in 3mL Syringe, 30/sp, 4 sp/cs
306531	Heparin Lock Flush Syringes, 100 Units/mL, 5mL Fill in 10mL Syringe w/blunt plastic cannula, 30/sp, 4 sp/cs

BD SYRINGE PARTS & ACCESSORIES

513560	1220FO Filling Outfit, Includes 470V Valve, 20" 470T Tubing, 330S Metal Sinker
511789	1250NR 14ga x 4 Laboratory Cannula
513092	470C Adapter Catheter Tip
513328	470F Fistula Tip
513329	470P Pomeroy Ear Tip



305109, 305115

BD REGULAR BEVEL NEEDLES

303005	18ga x 1 1/2" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
303006	19ga x 1 1/2" Needle, Non-Sterile, Regular Bevel, TW, 2500/bg, 2 bg/cs
303007	20ga x 1 1/2" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
303008	21ga x 1 1/2" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
303009	22ga x 1 1/2" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
303010	25ga x 5/8" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
303011	23ga x 1" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
303012	25ga x 1" Needle, Non-Sterile, Regular Bevel, 2500/bg, 2 bg/cs
305109	27ga x 1/2" Needle, Regular Bevel, 100/sp, 10 sp/cs
305110	26ga x 3/8" Needle, Intradermal Bevel, 100/sp, 10 sp/cs
305111	26ga x 1/2" Needle, Regular Bevel, 100/sp, 10 sp/cs
305115	26ga Sub-Q x 5/8" Needle, Regular Bevel, 100/sp, 10 sp/cs
305122	25ga x 5/8" Needle, Regular Bevel, 100/sp, 10 sp/cs
305125	25ga x 1" Needle, Regular Bevel, 100/sp, 10 sp/cs
305127	25ga x 1 1/2" Needle, Regular Bevel, 00/sp, 10 sp/cs
305143	23ga x 3/4" Needle, Regular Bevel, 100/sp, 10 sp/cs
305145	23ga x 1" Needle, Regular Bevel, 100/sp, 10 sp/cs
305155	22ga x 1" Needle, Regular Bevel, 100/sp, 10 sp/cs
305156	22ga x 1 1/2" Needle, Regular Bevel, 100/sp, 10 sp/cs
305165	21ga x 1" Needle, Regular Bevel, 100/sp, 10 sp/cs
305167	21ga x 1 1/2" Needle, Regular Bevel, 100/sp, 10 sp/cs
305175	20ga x 1" Needle, Regular Bevel, 100/sp, 10 sp/cs
305176	20ga x 1 1/2" Needle, Regular Bevel, 100/sp, 10 sp/cs
305177	21ga IV x 1" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305185	18ga FILL x 1 1/2" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305186	19ga x 1" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305187	19ga x 1 1/2" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305190	21ga IV x 1 1/2" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305193	23ga IM x 1" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305194	23ga IM x 1 1/2" TW Needle, Regular Bevel, 100/sp, 10 sp/cs
305195	18ga x 1" Needle, Regular Bevel, 100/sp, 10 sp/cs
305196	18ga x 1 1/2" Needle, Regular Bevel, 100/sp, 10 sp/cs
511000	LNR 27ga x 1/4" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511002	LNR 27ga x 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511005	LNR 27ga x 1" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511012	LNR 26ga x 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511017	LNR 26ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511019	LNRI 26ga x 3/8" Needle, Intradermal Bevel 12/cs
511023	LNR 25ga x 5/8" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs

511025	LNR 25ga x 1" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511028	LNR 25ga x 2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511034	LNR 24ga x 3/4" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511042	LNR 23ga x 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511047	LNR 23ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511048	LNR 23ga x 2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511055	LNR 22ga x 1" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511057	LNR 22ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511059	LNR 22ga x 3" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511067	LNR 21ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511077	LNR 20ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511078	LNR 20ga x 2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511089	LNR 19ga x 3" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511097	LNR 18ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511098	LNR 18ga x 2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511101	LNR 17ga x 3 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511108	LNR 15ga x 3 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511110	LNR 13ga x 3 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511252	LNR 30ga x 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs
511258	LNR 30ga x 1 1/2" Needle, Regular Bevel Metal Hub, Nickel Chrome Plated, 12/cs

BD SHORT BEVEL NEEDLES

303013	18ga x 1 1/2" Needle, Non-Sterile, Short Bevel, 2500/bg, 2 bg/cs
303015	22ga x 1 1/2" Needle, Non-Sterile, Short Bevel, 2500/bg, 2 bg/cs
305159	22ga x 1 1/2" Needle, Short Bevel, 100/sp, 10 sp/cs
305178	20ga x 1" Needle, Short Bevel, 100/sp, 10 sp/cs
305179	20ga x 1 1/2" Needle, Short Bevel, 100/sp, 10 sp/cs
305188	19ga x 1" TW Needle, Short Bevel, 100/sp, 10 sp/cs
305189	19ga x 1 1/2" TW Needle, Short Bevel, 100/sp, 10 sp/cs
305199	18ga x 1 1/2" Needle, Short Bevel, 100/sp, 10 sp/cs
511104	LNRS 16ga x 2" Needle, Short Bevel 12/cs
511107	LNRS 15ga x 2" Needle, Short Bevel 12/cs
511109	LNRS 13ga x 2" Needle, Short Bevel, 12/cs



BD SAFETYGLIDE™ NEEDLES & SYRINGES

305902	23ga x 1" Needle, 50/sp, 10 sp/cs
305903	1mL Syringe w/needle 25ga x 5/8", 50/sp, 8 sp/cs
305904	3mL Syringe w/needle 25ga x 5/8", 50/sp, 8 sp/cs
305909	3mL Syringe w/needle 21ga x 1 1/2", 50/sp, 8 sp/cs
305916	25ga x 1" Needle, 50/sp, 10 sp/cs
305917	21ga x 1 1/2" Needle, 50/sp, 10 sp/cs
305918	18ga x 1 1/2" Needle, 50/sp, 10 sp/cs
305900	22ga x 1 1/2" Needle, 50/sp, 10 sp/cs
305901	25ga x 5/8" Needle, 50/sp, 10 sp/cs
305905	3mL Syringe w/needle 23ga x 1", 50/sp, 8 sp/cs
305906	3mL Syringe w/needle 22ga x 1 1/2", 50/sp, 8 sp/cs
305907	5mL Syringe w/needle 22ga x 1 1/2", 50/sp, 8 sp/cs
305908	10mL Syringe w/needle 22ga x 1 1/2", 50/sp, 8 sp/cs
305915	21ga x 1" Needle, 50/sp, 10 sp/cs

TERUMO HYPODERMIC NEEDLES

Translucent hub and regular bevel.

3NN-1825R	18ga x 1" TW Needle, 100/bx, 10 bx/cs
3NN-1838R	18ga x 1 1/2" TW Needle, 100/bx, 10 bx/cs
3NN-1925R	19ga x 1" TW Needle, 100/bx, 10 bx/cs
3NN-1938R	19ga x 1 1/2" TW Needle, 100/bx, 10 bx/cs
3NN-2025R	20ga x 1" UTW Needle, 100/bx, 10 bx/cs
3NN-2038R	20ga x 1 1/2" UTW Needle, 100/bx, 10 bx/cs
3NN-2125R	21ga x 1" UTW Needle, 100/bx, 10 bx/cs
3NN-2138R	21ga x 1 1/2" UTW Needle, 100/bx, 10 bx/cs
3NN-2225R	22ga x 1" UTW Needle, 100/bx, 10 bx/cs
3NN-2238R	22ga x 1 1/2" UTW Needle, 100/bx, 10 bx/cs
3NN-2325R	23ga x 1" UTW Needle, 100/bx, 10 bx/cs
3NN-2338R	23ga x 1 1/2" UTW Needle, 100/bx, 10 bx/cs
3NN-2516R	25ga x 5/8" TW Needle, 100/bx, 10 bx/cs
3NN-2525R	25ga x 1" RW Needle, 100/bx, 10 bx/cs
3NN-2538R	25ga x 1 1/2" TW Needle, 100/bx, 10 bx/cs
3NN-2613R	26ga x 1/2" RW Needle, 100/bx, 10 bx/cs
3NN-2713R	27ga x 1/2" RW Needle, 100/bx, 10 bx/cs
3NN-3013R	30ga x 1/2" RW Needle, 100/bx, 10 bx/cs

EXEL ALUMINUM HUB NEEDLES

Designed for veterinary use. Sterile, Non-toxic, Pyrogen free and autoclavable. Needles and boxes are color-coded for easy identification. Each needle is individually packaged and sealed with a tamper evident heat stake.

26450	14ga x 1", 100/bx, 10 bx/cs
26451	14ga x .75", 100/bx, 10 bx/cs
26452	14ga x 1.5", 100/bx, 10 bx/cs
26454	14ga x 2", 100/bx, 10 bx/cs
26455	16ga x 5/8", 100/bx, 10 bx/cs
26456	16ga x .75", 100/bx, 10 bx/cs
26457	16ga x 1", 100/bx, 10 bx/cs
26458	16ga x 1.5", 100/bx, 10 bx/cs
26461	18ga x .75", 100/bx, 10 bx/cs
26462	18ga x 1", 100/bx, 10 bx/cs
26464	18ga x 1.5", 100/bx, 10 bx/cs
26486	Radiopaque 16ga x .75", 100/bx, 10 bx/cs
26487	Radiopaque 16ga x 1", 100/bx, 10 bx/cs
26488	Radiopaque 18ga x 1", 100/bx, 10 bx/cs
26489	Radiopaque 18ga x 1.5", 100/bx, 10 bx/cs



EXEL HYPODERMIC NEEDLES

Latex free.

26400	27ga x 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26402	26ga x 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26403	25ga x 5/8" Hypodermic Needle, 100/bx, 20 bx/cs
26405	25ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26406	25ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26407	23ga x 3/4" Hypodermic Needle, 100/bx, 20 bx/cs
26408	23ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26410	22ga x 3/4" Hypodermic Needle, 100/bx, 20 bx/cs
26411	22ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26412	22ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26414	21ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26416	21ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26417	20ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26418	20ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26419	18ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26420	18ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26424	19ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs
26426	27ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26427	27ga x 1 1/4" Hypodermic Needle, 100/bx, 20 bx/cs
26437	30ga x 1/2" Specialty Use Needle, 100/bx, 20 bx/cs
26438	19ga x 1 1/2" Hypodermic Needle, 100/bx, 20 bx/cs
26439	30ga x 1" Hypodermic Needle, 100/bx, 20 bx/cs

B BRAUN HUBER POINT NEEDLES

Straight, metal hub.

322015	20G x 1 1/2" Huber Needle, 12/cs
322215	22G x 1 1/2" Huber Needle, 12/cs


B BRAUN SURECAN® SAFETY HUBER INFUSION NEEDLES

Female luer connector, on/off clamp, 90° Surecan Safety Huber needle. DEHP-free, latex free. Priming volume: 0.2 mL. 6" length.

471741	19G x 3/4" Safety Huber Infusion Needle, 20/cs
471742	19G x 1" Safety Huber Infusion Needle, 20/cs
471743	19G x 1 1/2" Safety Huber Infusion Needle, 20/cs
471744	20G x 3/4" Safety Huber Infusion Needle, 20/cs
471745	20G x 1" Safety Huber Infusion Needle, 20/cs
471746	22G x 1/2" Safety Huber Infusion Needle, 20/cs
471747	22G x 3/4" Safety Huber Infusion Needle, 20/cs
471748	22G x 1" Safety Huber Infusion Needle, 20/cs
471749	19G x 3/4" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471750	19G x 1" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471751	20G x 1/2" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471752	20G x 3/4" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471753	20G x 1" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471754	20G x 1 1/2" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471755	22G x 3/4" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs
471756	22G x 1" Safety Huber Infusion Needle w/Ultrasite® Valve, 20/cs

B BRAUN WHIN® WINGED HUBER INFUSION NEEDLES


471682	19G x 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471685	22G x 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471686	22G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471687	22G x 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471688	20G x 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471691	20G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471692	20G x 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471693	20G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471694	20G x 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471695	19G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471696	19G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471697	19G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471698	19G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471707	19G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471709	22G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471714	22G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471716	20G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471720	20G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471722	19G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471725	20G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471726	22G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471727	22G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, Latex Free Y-Injection Site, On/Off Clamp, 20/cs
471729	19G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471730	19G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471731	20G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471732	20G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471733	22G x 3/4" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471734	22G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, On/Off Clamp, 20/cs
471735	19G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471736	20G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471737	20G x 1 1/2" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
471738	22G x 1" Huber Needle 90° Bend, 6" DEHP-Free Tubing, SAFSITE® Y-Site, On/Off Clamp, 20/cs
472104	22G x 3/4" Huber Needle, 90° Bend, DEHP-Free Tubing, ULTRASITE Y-Site, On/Off Clamp, 6" Length, 20/cs
472105	22G x 1" Huber Needle, 90° Bend, DEHP-Free Tubing, ULTRASITE Y-Site, On/Off Clamp, 6" Length, 20/cs



EXEL HUBER INFUSION SETS

Latex free. Non-coring, 90° 316L stainless steel needle. Y injection site available. Polypropylene hubs, color-coded according to gauge for easy identification. A rotating wing ensures reduced stress on both the tissue and the port. Also available DEHP Free. DEHP is a bonding agent that is used to make tubing soft and pliable. DEHP Free are required when administering some medications. Used on patients who are receiving drug therapy through an implanted port. A Huber needle is used to prevent damage to the port. Packaging: 50/box.

26946R	20ga x 1 1/2" Huber Infusion Set w/Y" Site, Yellow, 50/box
27940R	19ga x 1" Huber Infusion Set, Non DEHP w/Y" Site, Brown, 50/box
27941R	19ga x 3/4" Huber Infusion Set, Non DEHP w/Y" Site, Brown, 50/box
27944R	20G x 1" Huber Infusion Set, Non DEHP w/Y" Site, Yellow, 50/box
27945R	20ga x 3/4" Huber Infusion Set, Non DEHP w/Y" Site, Yellow, 50/box
27946R	20ga x 1 1/2" Huber Infusion Set, Non DEHP w/Y" Site, Yellow, 50/box
27948R	22ga x 1" Huber Infusion Set, Non DEHP w/Y" Site, Black, 50/box
27949R	22ga x 3/4" Huber Infusion Set, Non DEHP w/Y" Site, Black, 50/box
27950R	22ga x 1 1/2" Huber Infusion Set, Non DEHP w/Y" Site, Black, 50/box
27954R	20ga x 1" Huber Infusion Set, Non DEHP w/o "Y" Site, Yellow, 50/box
27955R	20ga x 3/4" Huber Infusion Set, Non DEHP w/o "Y" Site, Yellow, 50/box
27959R	22ga x 3/4" Huber Infusion Set, Non DEHP w/o "Y" Site, Black, 50/box
27958R	22ga x 1" Huber Infusion Set, Non DEHP w/o "Y" Site, Black, 50/box



EXEL HUBER NEEDLES

Latex free. Non-coring huber needles are used on patients who have an implanted port and are receiving drug therapy. Huber needles are used to prevent damage to the port.

26901	20ga x 3/4" Huber Needle, Right Angle, 50/box
26902	22ga x 3/4" Huber Needle, Right Angle, 50/box
26903	22ga x 1 1/4" Huber Needle, Right Angle, 50/box
26904	22ga x 1" Huber Needle, Right Angle, 50/box
26905	19ga x 3/4" Huber Needle, Right Angle, 50/box
26906	22ga x 1 1/2" Huber Needle, Right Angle, 50/box
26907	19ga x 1" Huber Needle, Right Angle, 50/box
26908	20ga x 1" Huber Needle, Right Angle, 50/box
26909	20ga x 1 1/2" Huber Needle, Right Angle, 50/box
26911	19ga x 1 1/2" Huber Needle, Right Angle, 50/box
26921	22ga x 1" Huber Needle, Straight, 50/box
26922	22ga x 1 1/2" Huber Needle, Straight, 50/box
26923	20ga x 1 1/2" Huber Needle, Straight, 50/box
26924	20ga x 1" Huber Needle, Straight, 50/box
26925	20ga x 1 1/4" Huber Needle, Straight, 50/box

EXEL SECURETOUCH SAFETY HUBER INFUSION SETS

Non-coring infusion sets with 8" DEHP free tubing.

27980S	19ga x 1" Huber Infusion Set w/Y site, 20/bx
27981S	19ga x 3/4" Huber Infusion Set w/Y site, 20/bx
27984S	20ga x 1" Huber Infusion Set w/Y site, 20/bx
27985S	20ga x 3/4" Huber Infusion Set w/Y site, 20/bx
27986S	20ga x 1 1/2" Huber Infusion Set w/Y site, 20/bx
27988S	22ga x 1" Huber Infusion Set w/Y site, 20/bx
27989S	22ga x 3/4" Huber Infusion Set w/Y site, 20/bx
27990S	22ga x 1 1/2" Huber Infusion Set w/Y site, 20/bx
27994S	20ga x 1" Huber Infusion Set, 20/bx
27995S	20ga x 3/4" Huber Infusion Set, 20/bx
27998S	22ga x 1" Huber Infusion Set, 20/bx
27999S	22ga x 3/4" Huber Infusion Set, 20/bx

EXEL DENTAL NEEDLES

Latex free. Plastic hub with a red dot to indicate bevel position.

26555	Dental Needle, 30ga Short, 100/bx, 10 bx/cs
26557	Dental Needle, 27ga Long, 100/bx, 10 bx/cs
26558	Dental Needle, 27ga Short, 100/bx, 10 bx/cs
26562	Dental Needle, 30ga Long, 100/bx, 10 bx/cs
26556	Dental Needle 30ga X-Short, 100/bx, 10bx/cs

KENDALL 400 PLASTIC HUB DENTAL NEEDLE

Rigid autoclavable polypropylene hard pack. Tamper-evident heat stake. Self-threading, self aligning universal fit hub design. Ultra sharp, tri-beveled, anti-coring lancet point. Anti-coring butt end point design.

8881400074	400 Plastic Hub Dental Needle, 30ga Short, 3/4" (21mm) L, Blue, 100/bx, 10 bx/cs
------------	--

TERUMO DENTAL NEEDLES

For conduction anesthesia and infiltration anesthesia. The ultra-sharp lancet point minimizes pain and discomfort.

1DN-2530K	25ga x 1 3/16" Dental Needle, Long, 100/bx, 20 bx/cs
1DN-2721K	27ga x 13/16" Dental Needle, Short, 100/bx, 20 bx/cs
1DN-2730K	27ga x 1 3/16" Dental Needle, Long, 100/bx, 20 bx/cs
1DN-3021K	30ga x 13/16" Dental Needle, Short, 100/bx, 20 bx/cs

**We have the right
products to meet
your needs.**

**Please call for help
locating certain items.**

**BD SPINAL NEEDLES**

Sterile, disposable needles have tinted translucent Luer-Lok hub for better visualization of cerebrospinal fluid. Quick point spreads, rather than cuts, longitudinal dural fibers.

405073	25ga x 1" Neonatal Spinal Needle, Blue, 10/sp, 5 sp/cs
405078	25ga x 2" Neonatal Spinal Needle, Blue, 25/sp, 4 sp/cs
405081	27ga x 3 1/2" Spinal Needle, Gray, 25/sp, 4 sp/cs
405148	22ga x 5" Spinal Needle, Black, 10/sp, 5 sp/cs
405149	22ga x 7" Spinal Needle, Black, 10/sp, 5 sp/cs
405160	20ga x 1 1/4" Spinal Introducer, Yellow, 25/sp, 4 sp/cs
405161	22ga x 1 1/2" Spinal Needle, Black, 25/sp, 4 sp/cs
405162	20ga x 1 1/2" Spinal Needle, Yellow, 25/sp, 4 sp/cs
405164	26ga x 3 1/2" Spinal Needle, Tan, 25/sp, 4 sp/cs
405171	22ga x 3" Spinal Needle, Black, 25/sp, 4 sp/cs
405180	25ga x 3 1/2" Spinal Needle, Blue, 25/sp, 4 sp/cs
405181	22ga x 3 1/2" Spinal Needle, Black, 25/sp, 4 sp/cs
405182	20ga x 3 1/2" Spinal Needle, Yellow, 25/sp, 4 sp/cs
405211	20ga x 6" Spinal Needle, Yellow, 10/sp, 5 sp/cs
405234	25ga x 4 11/16" Spinal Needle, Blue, 10/sp, 5 sp/cs
400106	23ga x 3 1/2" Spinal Needle, Turquoise, 25/sp, 4 sp/cs
405071	20ga x 2 1/2" Spinal Needle, Yellow, 25/sp, 4 sp/cs
405074	22ga x 2 1/2" Spinal Needle, Black, 25/sp, 4 sp/cs
405170	25ga x 3" Spinal Needle, Blue, 25/sp, 4 sp/cs
405172	20ga x 3" Spinal Needle, Yellow, 25/sp, 4 sp/cs
405174	18ga x 3" Spinal Needle, Pink, 25/sp, 4 sp/cs
405184	18ga x 3 1/2" Spinal Needle, Pink, 25/sp, 4 sp/cs
408360	18ga x 6" Spinal Needle, Pink, 10/sp, 5 sp/cs
511293	18ga x 6" TW Spinal Needle w/fitted stylet, T462LNR
511310	14ga x 2" TW Spinal Needle w/fitted stylet, T462LNR

BD WHITACRE PENCIL POINT SPINAL NEEDLES

Sterile.

405010	22ga x 3 1/2" Spinal Needle, Black, 10/sp, 5 sp/cs
405133	24ga x 3 1/2" Spinal Needle, High Flow, Purple, 10/sp, 5 sp/cs
405138	25ga x 3 1/2" Spinal Needle, High Flow, Blue, 10/sp, 5 sp/cs
405140	25ga x 5" Spinal Needle, High Flow, Blue, 10/sp, 5 sp/cs
405144	27ga x 5" Spinal Needle, High Flow, Grey, 10/sp, 5 sp/cs
405146	25ga x 1 1/2" Spinal Needle, Blue, 10/sp, 5 sp/cs
409442	25ga x 4 11/16" Spinal Needle, Blue, 10/sp, 5 sp/cs
409443	27ga x 4 11/16" Spinal Needle, Gray, 10/sp, 5 sp/cs
511336	22ga x 3 1/2" Spinal Needle, W462LNR
405079	27ga x 3 1/2" Spinal Needle, High Flow, Grey, 10/sp, 5 sp/cs

BD YALE™ REUSABLE SPINAL NEEDLES

Quincke.

511114	18ga x 3 1/2" Spinal Needle, 462LNR
511120	22ga x 2" Spinal Needle, 462LNR
511136	22ga x 3 1/2" Spinal Needle, 462LNR
511139	20ga x 3 1/2" Spinal Needle, 462LNR
511140	20ga x 4" Spinal Needle, 462LNR
511142	19ga x 3 1/2" Spinal Needle, 462LNR

B BRAUN PENCAN® PENCIL POINT SPINAL NEEDLES

333853	25G x 3 1/2" PENCAN Pencil Point Spinal Needle, 25/cs
333866	22G x 3 1/2" PENCAN Pencil Point Spinal Needle, 25/cs
333873	27G x 3 1/2" PENCAN Pencil Point Spinal Needle, 25/cs
333875	25G x 5" PENCAN Pencil Point Spinal Needle, 25/cs
333876	27G x 5" PENCAN Pencil Point Spinal Needle, 25/cs
333877	25G x 6" PENCAN Pencil Point Spinal Needle, 25/cs

B BRAUN SPINOCAN® SPECIAL PROCEDURAL SPINAL NEEDLES

333308	25G x 4 3/4" Spinal Needle w/clear hub, 50/cs
333315	22G x 4 3/4" Spinal Needle w/clear hub (333318), 50/cs

B BRAUN SPINOCAN® SPINAL NEEDLE WITH QUINCKE BEVEL

333290	27G x 3 1/4" Thinwall Spinal Needle w/clear hub, 50/cs
333300	26G x 3 1/2" Spinal Needle w/clear hub, 50/cs
333320	22G x 3 1/2" Spinal Needle w/clear hub, 50/cs
333325	22G x 1 1/2" Spinal Needle w/clear hub, 50/cs
333335	20G x 3 1/2" Spinal Needle w/clear hub, 50/cs
333350	18G x 3 1/2" Spinal Needle w/clear hub, 50/cs
333313	25G x 3 1/2" Spinal Needle w/clear hub (333310), 50/cs

B BRAUN SPROTTE® PENCIL POINT SPINAL NEEDLES

333721	22G x 3 1/2" SPROTTE Pencil Point Spinal Needle w/introducer, 25/cs
333741	24G x 3 1/2" SPROTTE Pencil Point Spinal Needle w/introducer, 25/cs
333751	25G x 3 1/2" SPROTTE Pencil Point Spinal Needle, 25/cs
333769	27G x 4 7/8" SPROTTE Pencil Point Spinal Needle, 25/cs
333731	25G x 4 7/8" SPROTTE Pencil Point Spinal Needle (333749), 25/cs

**EXEL SPINAL NEEDLES**

Latex free. Plastic luer-lock hub. Sterile. Disposable. Color coded to size.

26960	18ga x 3 1/2" Spinal Needle, 20/bx, 5 bx/cs
26964	20ga x 3 1/2" Spinal Needle, 50/bx, 2 bx/cs
26967	22ga x 3 1/2" Spinal Needle, 50/bx, 2 bx/cs
26970	25ga x 3 1/2" Spinal Needle, 50/bx, 2 bx/cs

TERUMO HYPODERMIC SPINAL NEEDLES

The tip is beveled for smoother penetration. However, it has just the right amount of penetration resistance to give you a fingertip feeling of exactly where the needle point is. The blister pack's peel-apart feature allows you to dispense the needle into a sterile field, and the protector keeps it sterile until the very moment of use. Color-coding is used on the stylet-hub, blister pack, and shelf package for rapid gauge identification at every stage.

1SN-1890	18ga x 3 1/2" Spinal Needle, 20/bx, 10 bx/cs
1SN-2070	20ga x 2 3/4" Spinal Needle, 20/bx, 10 bx/cs
1SN-2090	20ga x 3 1/2" Spinal Needle, 20/bx, 10 bx/cs
1SN-2270	22ga x 2 3/4" Spinal Needle, 20/bx, 10 bx/cs
1SN-2290	22ga x 3 1/2" Spinal Needle, 20/bx, 10 bx/cs
1SN-2590	25ga x 3 1/2" Spinal Needle, 20/bx, 10 bx/cs

B BRAUN PERIFIX® EPIDURAL NEEDLES

332160	18G x 6" Tuohy Epidural Needle w/clear hub, 25/cs
332164	17G x 3 1/2" Tuohy Needle w/Weiss wings, 25/cs
332165	18G x 3 1/2" Hustead Needle w/clear hub, 25/cs
332166	18G x 3 1/2" Tuohy Needle w/clear hub, 25/cs
332167	20G x 3 1/2" Non-Winged Tuohy Needle, 25/cs
332168	20G x 3 1/2" Winged Tuohy Epidural Needle, 25/cs
332171	22G x 3" Winged Tuohy Epidural Needle, 25/cs
332180	17G x 3 1/2" Tuohy Needle w/clear hub (332162), 25/cs
332181	17G x 4 1/2" Winged Tuohy Epidural Needle (332169), 25/cs

RUSCH SPROTTE® SPECIAL EPIDURAL NEEDLES

For continuous epidural technique or combined technique. Sterile. Minimized the risk of PLPH/PDPH after an unintentional puncture of the dura. Lateral opening is located far up in the conical point. Catheter channel has been rounded to avoid the risk of a catheter being severed.

00111514	18G x 3 1/2" (90mm) Epidural Needle, 25/bx
00211514	17G x 3 1/2" (90mm) Epidural Needle, 25/bx

RUSCH TUOHY EPIDURAL NEEDLES

1400	Extra Length Epidural Needle, 4 1/2", 17Ga, Sterile, 10/bx
1401	Extra Length Epidural Needle, 6", 17Ga, Sterile, 10/bx

BD CHIBA FINE NEEDLE ASPIRATION BIOPSY

408268	Chiba Fine Needle Aspiration Biopsy/Transhepatic Cholangiography Needle, 22ga x 6" TW, 10/sp, 5 sp/cs
408269	Chiba Fine Needle Aspiration Biopsy/Transhepatic Cholangiography Needle, 22ga x 8" TW, 10/sp, 5 sp/cs

BD LIVER, KIDNEY & SOFT TISSUE NEEDLES

511282	469LNR 14ga x 3 3/8" Silverman
511283	468LRT 13ga x 3 1/4" Aspiring Trocar
511403	01-0109 16ga x 4" Klatskin TW
511547	01-0008 14ga x 4 1/2" Silverman-B-D, Modified Franklin, TW
511548	01-0008 14ga x 3 3/8" Silverman-B-D, Modified Franklin, TW

BD WESTCOTT FINE NEEDLE ASPIRATION BIOPSY

408262	Westcott Fine Needle Aspiration Biopsy, 20ga x 3 1/2", 10/spx, 5 sp/cs
408263	Westcott Fine Needle Aspiration Biopsy, 20ga x 5 1/2", 10/spx, 5 sp/cs
408264	Westcott Fine Needle Aspiration Biopsy, 20ga x 7", 10/spx, 5 sp/cs
408265	Westcott Fine Needle Aspiration Biopsy, 22ga x 3 1/2", 10/spx, 5 sp/cs
408266	Westcott Fine Needle Aspiration Biopsy, 22ga x 5 1/2", 10/spx, 5 sp/cs
408267	Westcott Fine Needle Aspiration Biopsy, 22ga x 7", 10/spx, 5 sp/cs

BUSSE ENHANCED J-STYLE MARROW BIOPSY/ASPIRATION NEEDLE

The easily gripped handle of the 11g x 4" J-Style needle provides the physician with greater torque than similar instruments. Combined with the ultra-sharp trocar tip on the stylet, Busse's exclusive design offers much easier penetration of bone tissue. Once the stylet is withdrawn, the cannula tip presents four sharp, beveled planes with pointed serrations at their junctures, and produces tightly sheared marrow samples as the needle is advanced. A blunt sample extractor is provided for easy removal of sample. The handle allows for both Luer lock and Luer slip connections to accommodate either style of syringes.

630 11ga x 4" J-Style Bone Marrow/Aspiration Needle, 10/cs

BUSSE I-STYLE STERNAL-ILIAC BONE MARROW BIOPSY/ASPIRATION NEEDLE

Busse's 15g I-style needle allows variable depth adjustment from 3/8" to 1-7/8". Luer lock allows for easy connection to syringe. The handle allows for both Luer lock and Luer slip connections to accommodate both style of syringes.

631 15ga x 38" to 1-7/8" Adjustable, I-Style Sternal-Iliac Aspiration Needle, 10/cs

KENDALL MONOJECT BONE MARROW ASPIRATION NEEDLES

Plastic hub. Traditional bone marrow aspiration needle design. Precision ground and fitted needle point and stylet tip. Sterile.

8881240009 Bone Marrow Needle, 16 x 2", Rosenthal Type, Sterile, 20/cs

8881240017 Bone Marrow Needle, 16 x 3/4", Rosenthal Type, Sterile, 20/cs

BD ANGIOGRAPHY NEEDLES

- 408270 Potts-Courmand Needle, TW w/obturator & cannula extender, 18ga x 2 3/4" For .032"-.038" Guide Wire, 10/sp, 5 sp/cs
- 408273 Arterial Needle Outer Blunt Cannula, 18ga x 2 7/8" For .032"-.038" Guide Wire, 10/sp, 5 sp/cs
- 408275 Amplatz Needle, FEP Polymer Catheter Over Inner Needle w/fitted stylet, Sterile Single Use, 18ga x 2 1/2" For .032"-.038" Guide Wire, 10/sp, 5 sp/cs
- 408295 Single Wall Puncture (M/P) Procedure Needle, 18ga x 2 3/4" One Part Needle w/wings For .032"-.038" Guide Wire, 10/sp, 5 sp/cs
- 408300 Single Wall Puncture (M/P) Procedure Needle, 18ga x 2 3/4" One Part Needle For .032"-.038" Guide Wire, 10/sp, 5 sp/cs

BD BONE MARROW, ASPIRATION & TRANSFUSION NEEDLE

- 511341 460LNR 18ga x 1" Rosenthal
- 511346 460LNR 16ga x 1 5/16" Rosenthal
- 511348 458LNR 18ga x 1" Osgood
- 511349 458LNR 16ga x 1 5/16" Osgood
- 511420 01-0081 10ga x 3" Large Silverman
- 511573 01-0004 10ga x 4 1/4" Westerman-Jensen

BD BONE MARROW BIOPSY NEEDLE

408627 11ga x 4" Bone Marrow Biopsy Needle, 10/cs

BD BREAST LESION LOCALIZATION NEEDLES

- 409407 PrecisionGuide™ Breast Lesion Localization Needles, 20ga x 3 1/2" w/hook wire, 10/sp, 5 sp/cs
- 409408 PrecisionGuide™ Breast Lesion Localization Needles, 20ga x 2 1/2" w/hook wire, 10/sp, 5 sp/cs
- 409409 PrecisionGuide™ Breast Lesion Localization Needles, 20ga x 1 1/2" w/hook wire, 10/sp, 5 sp/cs

BD FILTER NEEDLES

- 305200 19ga x 1 1/2" TW Needle, 5 Micron, 100/sp, 10 sp/cs
- 305201 18ga x 1 1/2" TW Needle w/Nokor™ point, 5 Micron, 10/sp, 10 sp/cs



BD INTEGRA™ NEEDLES

- 305275 18ga x 1 1/2" Blunt Fill Needle, 100/sp, 10 sp/cs
- 305310 25ga x 5/8" Retracting Needle, 100/sp, 10 sp/cs
- 305311 25ga x 1" Retracting Needle, 100/sp, 10 sp/cs
- 305312 23ga x 1" Retracting Needle, 100/sp, 10 sp/cs
- 305313 22ga x 1 1/2" Retracting Needle, 100/sp, 10 sp/cs
- 305833 18ga x 1 1/2" Blunt Filter Needle, 50/sp, 10 sp/cs

BD MENGHINI TECHNIQUE ASPIRATING NEEDLE

- 511485 Menghini Technique Aspirating, 16ga x 4 3/4" TW
- 511491 Menghini Technique Aspirating, 16ga x 2 3/4" TW

BD MENGHINI BIOPSY RETAINING PIN

511591 Menghini Biopsy Retaining Pin, 16ga, 12/cs

MILTEX MENGHINI BIOPSY NEEDLES

Improved model with built-in stop plug that never requires replacement.

- 13-130 1.2mm Biopsy Needle, 18G For Liver Biopsy, Adult Size, 70mm Long
- 13-132 1.4mm Biopsy Needle, 17G For Liver Biopsy, Adult Size, 70mm Long
- 13-156 1.8mm Biopsy Needle, 16G For Kidney Biopsy, 120mm Long

BD NOKOR™ ADMIX NEEDLES

- 305215 18ga x 1" TW Non-Coring Needle, 10/sp, 10 sp/cs
- 305216 16ga x 1" TW Non-Coring Needle, 10/sp, 10 sp/cs

BD NOKOR™ VENTED NEEDLES

- 305213 16ga x 1" TW Non-Coring Needle, 10/sp, 10 sp/cs
- 305214 18ga x 1" TW Non-Coring Needle, 10/sp, 10 sp/cs



BD SPECIALTY USE NEEDLES

305106	30ga x 1/2" Specialty Use Needle, 100/sp, 10 sp/cs
305120	23ga x 1 1/4" Specialty Use Needle, 100/sp, 10 sp/cs
305124	25ga x 7/8" Specialty Use Needle, 100/sp, 10 sp/cs
305128	30ga x 1" Specialty Use Needle, 100/sp, 10 sp/cs
305129	21ga x 2" Specialty Use Needle, 100/sp, 10 sp/cs
305136	27ga x 1 1/4" Specialty Use Needle, 100/sp, 10 sp/cs
305197	16ga x 1" Specialty Use Needle, 100/sp, 10 sp/cs
305198	16ga x 1 1/2" Specialty Use Needle, 100/sp, 10 sp/cs

B BRAUN CONTIPLEX® CONTINUOUS NERVE BLOCK SET

Contents: catheter connector, threading assist guide, 20G polyamide radiopaque catheter with open tip, 18G x 1 3/4" catheter over a 20G short bevel introducer needle.

331692	Continuous Nerve Block Set, 12/cs
--------	-----------------------------------

B BRAUN CONTIPLEX® CONTINUOUS NERVE BLOCK TRAY

Contents: prep solution well, (3) prep sponges, fenestrated drape, (4) 4" x 4" gauze sponges, paper towel, 4" filter straw, 3 mL plastic measuring syringe with 18G x 1 1/2" needle attached, 50 mL medicine well, 5 mL plastic syringe, 20 mL plastic syringe, 25G x 5/8" skin wheal needle, 20G radiopaque polyamide catheter, threading assist guide, catheter, 18G x 1 3/4" catheter over a 20G short bevel introducer needle.

331700	Continuous Nerve Block Tray, 10/cs
--------	------------------------------------

B BRAUN CONTIPLEX® NEEDLE CONTINUOUS NERVE BLOCK SETS

Contents: 18G insulated tuohy needle, hemostasis valve with sideport, 20G polyamide nylon catheter with closed tip, catheter connector.

331693	Continuous Nerve Block Set w/18G x 3 1/2" Insulated Tuohy Needle, 12/cs
331694	Continuous Nerve Block Set w/18G x 6" Insulated Tuohy Needle, 12/cs
331695	Continuous Nerve Block Set w/18G x 2" Insulated Tuohy Needle, 12/cs

B BRAUN VENTED NEEDLES

415070	Vented Needle, Vented Piercing Pin For Drug Reconstitution, 100/cs
415072	Lateral Flow Vented Needle, Vented Piercing Pin, Specifically Designed to Minimize Drug Foaming During the Reconstitution Process, 100/cs

BD BLUNT FILL NEEDLE

305180	18ga x 1 1/2" Blunt Fill Needle, 100/sp, 10 sp/cs
--------	---

BD BLUNT FILTER NEEDLES

305211	18ga x 1 1/2" TW Needle w/blunt fill tip, 5 Micron, 100/sp, 10 sp/cs
--------	--

BD DISPENSING PINS

305263	Sure-Med™ Dispensing Pin, .45 Micron Filter, 50/sp, 4 sp/cs
305264	Safety-Med™ Chemotherapy Dispensing Pin, .20 Micron Filter, 50/sp, 5 sp/cs

BD ECHOGENIC PROCEDURE NEEDLES

For use with ultrasound.

400727	EchoGenic Procedure Needles, Ultra-vue 22ga x 3 1/2" Spinal Type Point, 10/sp, 5 sp/cs
400730	EchoGenic Procedure Needles, Ultra-vue 20ga x 3 1/2" Spinal Type Point, 10/sp, 5 sp/cs

BD EPIDURAL ANESTHESIA NEEDLES

Sterile.

405018	18ga x 3 1/2" Tuohy Needle, No Wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405028	20ga x 3 1/2" Tuohy Needle, No Wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405092	17ga X 3 1/2" Tuohy Needle, No Wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405185	20ga x 2" Weiss Needle w/fixed wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405186	20ga x 3 1/2" Weiss Needle w/fixed wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405188	18ga x 3 1/2" Weiss Needle w/fixed wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405190	18ga x 5" Weiss Needle w/fixed wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405191	17ga x 3 1/2" Weiss Needle w/fixed wings, Modified Tuohy Point, 10/sp, 5 sp/cs
405193	17ga x 5" Weiss Needle w/fixed wings, Modified Tuohy Point, 10/sp, 5 sp/cs
511219	17ga x 3" Tuohy TW Needle, T466LNRH
511285	17ga x 7" Tuohy TW Needle, T466LNRH
405528	17ga x 3 1/2" Weiss Needle, No Wings, Tactile Point, 10/sp, 5 sp/cs
405529	18ga x 3 1/2" Weiss Needle, No Wings, Tactile Point, 10/sp, 5 sp/cs
405530	17ga x 3 1/2" Tuohy Needle, No Wings, Tactile Point, 10/sp, 5 sp/cs
405531	18ga x 3 1/2" Tuohy Needle, No Wings, Tactile Point, 10/sp, 5 sp/cs

BD LOCAL ANESTHESIA SECURITY BEAD NEEDLES

511359	22ga x 4" Security Bead Needle, 45LNR, 12/sp
511362	22ga x 2" Security Bead Needle, 45LNR, 12/sp
511364	20ga x 6" Security Bead Needle, 45LNR
511369	20ga x 4" Security Bead Needle, 45LNR, 12/sp

BD PERIPHERAL BLOCK NEEDLES

Sterile. Single use.

404204	22ga x 2 1/8", Stimex™ Insulated Needle, FEP Polymer-Coated, Block Bevel, 10/sp, 5 sp/cs
404247	25ga x 2 3/16" Insulated Needle Sets, 15° Bevel (electrical leads & tubing, polyester coating), 10/sp, 5 sp/cs
404249	22ga x 2" Insulated Needle Sets, 15° Bevel (electrical leads & tubing, polyester coating), 10/sp, 5 sp/cs
404250	22ga x 3 3/16" Insulated Needle Sets, 15° Bevel (electrical leads & tubing, polyester coating), 10/sp, 5 sp/cs
404251	22ga x 4 3/4" Insulated Needle Sets, 15° Bevel (electrical leads & tubing, polyester coating), 10/sp, 5 sp/cs
408348	22ga x 1 1/2" Regional Block Needle, Security Bead, Block Bevel, 10/sp, 5 sp/cs

BD PUDENDAL & PARASYMPATHETIC BLOCK NEEDLES

511317	20ga x 5" Sympathetic Nerve Needle, 462LNRS
511318	19ga x 5" Sympathetic Nerve Needle, 462LNRS
511319	20ga x 6" Parasympathetic Nerve Needle, 462LNRS
511320	18ga x 6" Sympathetic Nerve Needle, 462LNRS

BD RETROBULBAR ANESTHESIA NEEDLE

511583 23ga x 1 3/8" Atkinson Anesthesia Needle, 01-0001

BD TONSIL NEEDLES

511118 LLX Extension
 511168 23ga x 1/2" Tonsil Needle (Curved), 465LNR
 511169 23ga x 3/4" Tonsil Needle (Straight), 465LNR

BD VISITEC™ IRRIGATION & ASPIRATION CANNULAS

581624 1.2mm TW Outer Coaxial I&A Cannula, .80mm Inner Cannula w/.30mm Side Aspiration Port (Pearce), 5/sp
 585047 .60mm/.60mm TW Coaxial I&A Cannula, .30mm Aspiration Port (Simcoe), 5/sp
 585048 .60mm/.60mm TW Coaxial I&A Cannula, .30mm Aspiration Port (Reverse) (Simcoe), 5/sp
 585086 .60mm/.60mm TW Coaxial I&A Cannula, .30mm Aspiration Port, Metal Hub (Simcoe), 3/sp
 585104 .60mm/.60mm TW Coaxial I&A Cannula, .40mm Aspiration Port, Metal Hub, Micro Etched Tips (Dahan/Simcoe), 3/sp
 585110 .60mm/.60mm TW Coaxial I&A Cannula, .40mm Aspiration Port, Metal Hub (Simcoe), 3/sp
 585222 .60 x 17mm Bimanual I&A Cannula, Curved 45°, .40mm Top Aspiration Port, 10/sp
 585223 .80 x 17mm Bimanual I&A Cannula, TW, Curved 45°, .51mm End Port, 10/sp

BD VISITEC™ LACRIMAL CANNULAS

581276 .45 x 28mm Cannula, 10/sp
 581615 .80 x 35mm Cannula, Curved, 10/sp
 585016 1.1 x 20mm Shaft Attached to .60 x 12mm Malleable Cannula, 5/sp

BD VISITEC™ LASIK CANNULAS

585216 1.6 x 12mm Cannula, Two .50 x 13mm Tip Extensions At 45° Angle (Vidaum), 5/sp
 585234 .60 x 22mm Cannula, Angled 45°, 4mm From Spatulated Tip (Hersh), 20/sp
 585250 .50 x 22mm Cannula, Three .25mm Ports On 7mm Extension, Angled 45° (Lindstrom), 10/sp
 585252 .60 x 22mm Cannula, Angled 45°, 9mm From Spatulated Tip (Hersh), 20/sp
 585259 .80 x 22mm Cannula, Angled 45°, 12mm From End, Flattend .30mm, 5mm From Tip (Daya), 10/sp
 585266 .50 x 22mm Cannula, Angled 45°, 4mm From Spatulated Tip (Hersh), 20/sp
 585267 .40 x 22mm Cannula, Angled 45°, 4mm From Spatulated Tip (Hersh), 20/sp
 585268 .50 x 22mm Cannula, Angled 45°, 9mm From Spatulated Tip (Hersh), 20/sp
 585269 .40 x 22mm Cannula, Angled 45°, 9mm From Spatulated Tip (Hersh), 20/sp
 585271 .60 x 20mm Cannula, Angled 45°, 10mm From End, Flattend .30mm, 1mm From Tip (Daya), 10/sp
 585278 .40 x 22mm Cannula, Angled 25°, Formed 7mm From End, 10/sp
 585279 .60 x 22mm Cannula, Flattened Tip On Horizontal Plane (Donnenfeld), 10/sp
 585280 .50 x 22mm Cannula, TW, Angled Tip w/three .25mm Ports On 8mm Extension (Solomon), 10/sp
 585283 23ga LASIK Irrigating Cannula (Lindstrom), 10/sp

BD VISITEC™ VITREORETINAL CANNULAS

585002 High Viscosity Injector, 1.6 x 9mm, 20° Beveled Up (Peyman), 5/sp
 585173 High Viscosity Injector, .90 x 4mm, Ultra TW Cannula, 45° Beveled Tip, 10cm Tubing (Hammer), 5/sp
 585184 High Viscosity Injector, .90 x 6mm, Ultra TW Cannula, 45° Beveled Tip, 10cm Tubing (Hammer), 5/sp
 585205 .90 x 28mm Cannula, 10/sp
 585208 High Viscosity Injector, .90 x 4mm, Ultra TW Cannula, 45° Beveled Tip, 16cm Tubing (Hammer), 5/sp
 585228 High Viscosity Injector, .90 x 6mm, Ultra TW Cannula, 45° Beveled Tip, 16cm Tubing (Hammer), 5/sp

MADAJET XL DENTAL NEEDLE-FREE INJECTOR

A proven means of reducing fear and trauma associated with traditional needle and syringe injections. Injects local anesthetics into palatal and soft tissues with greater ease and virtually no pain. Provides instant local anesthesia or regional blocking for many dental and maxillofacial procedures. Provides adequate anesthesia for deep needle insertion. Has interchangeable Extenda Tips® for easy sterilization. Consistent depth of penetration of 2-2.5mm below epithelium and makes a wheal at the base of injection of 5-6mm in diameter. Consistent volume of 0.1cc per injection intradermally. Permits approximately 38 injections with single loading (to 4.0cc). May be sterilized by autoclaving or cold sterilization. Can be used for all age groups. Full one year warranty.

301 MadaJet Instrument w/Angle Head, (2) #300 Extenda Tips 1/2", Madacide-FD (32 oz), ME30 Wrench, 307 Holder, (3) Pyrex Fill Chambers, (3) Stylets & Custom Case
 302 MadaJet Instrument Same as 301 w/1 Extra Head Assembly (Angle Head)
 303 Head Assembly w/Extenda Tip & Sheath
 305-C Conversion MadaJet Dental of Old Style to New Style (XL) w/2 #300 Extenda Tips w/Sheath 1/2", MadaCide-FD (32 oz), ME30 Wrench, (3) Stylets, 307 Holder, (3) Pyrex Fill Chambers & Custom Case (Note: must be returned to MADA for conversion)
 305-C2 Same as 305-C w/Extra Head Assembly

Accessories

300 Extenda Tip, 1/2" w/sheath
 307 Holder
 308 Sheaths For Extenda Tip, 1/2", 30/bx
 400-1 1/2D Extenda Tip, 1 1/2" w/sheath
 602 Pyrex Fill Chambers, 5/pkg
 604 Stylets, 5/pkg
 606 Instructional Video Tape
 ME30 Wrench

We have the right product to meet your need.

MADAJET XL MEDICAL NEEDLE-FREE INJECTOR

Designed for patient comfort to reduce fear and trauma associated with traditional needle injections. Injects local anesthetics and other medicaments with greater ease and virtually no pain. Not only functions as an injection device to inject skin tissue, but is the only jet injector that allows physicians to inject through all body orifices for internal tissue infiltration (i.e., oral cavity & laryngeal area, vagina & cervix, rectum, ear & nose, open wounds & surgical openings). Provides instant local anesthesia for many procedures as well as adequate anesthesia for deep needle insertion. Due to the small volume of local anesthetics used by jet injection, tissue distention (ballooning) experienced with needle injections are virtually eliminated. May be used with local anesthetics, steroids and other medicaments. Has interchangeable Extenda Tips and sheaths for easy sterilization. Consistent depth of penetration of 4-4.5mm below the epithelium and makes a wheal at the base of the injection of 5-6mm in diameter. Consistent volume of 0.1cc per injection intradermally. Permits approximately 9 injections per 1cc, loading capacity of fill chamber of 4cc. Can be autoclaved or gas sterilized or disinfected with MadaCide solution. Full one year warranty.

- 401 MadaJet Instrument w/Angle Head, (2) #400 Extenda Tips 3/4", MadaCide-FD (32 oz), ME30 Wrench, 407 Holder, (3) Pyrex Fill Chambers, (3) Stylets & Custom Case
- 401-6T MadaJet Instrument w/Angle Head, (6) Extenda Tips (3/4", 1 1/2", 2 1/2", 3 1/2", 4 1/2" 6 1/2") MadaCide-FD (32 oz), ME30 Wrench, 407 Holder, (3) Pyrex Fill Chambers, (3) Stylets & Custom Case
- 402 MadaJet Instrument Same as #401 w/Extra Head Assembly (Angle Head)
- 402-6T MadaJet Instrument Same as #401-6T w/Extra Head Assembly (Angle Head)
- 405-C Conversion MadaJet Medical of Old Style to New Style (XL) w/2 #400 Extenda Tips w/sheath 3/4", MadaCide-FD (32 oz), ME30 Wrench, (3) Stylets, 407 Holder, (3) Pyrex Fill Chambers & Custom Case (Note: must be returned to MADA for conversion)
- 405-C2 Same as 405-C w/extra head assembly

MADAJET XL MEDICAL NEEDLE-FREE INJECTOR ACCESSORIES

- 400-1 1/2 Extenda Tip, 1 1/2" w/sheath
- 400-3 1/2 Extenda Tip, 3 1/2" w/sheath
- 400-4 1/2 Extenda Tip, 4 1/2" w/sheath
- 400-6 1/2 Extenda Tip, 6 1/2" w/sheath
- 403 Head Assembly w/Extenda Tip & Sheath (Angle)
- 407 Holder
- 408 Sheath 3/4" For Extenda Tip, 30/bx
- 409 Sheath 1 1/2" For Extenda Tip, 30/bx
- 410 Sheath 2 1/2" For Extenda Tip, 30/bx
- 412 Sheath 4 1/2" For Extenda Tip, 30/bx
- 414 Instructional Video Tape

MADAJET XL PODIATRY NEEDLE-FREE INJECTOR

Easy to use and virtually painless compared to needles. Provides instant local anesthesia or regional blocking for most podiatric procedures. Minimizes tissue trauma and facilitates suturing as there is no tissue distention around the injection site. Provides adequate anesthesia for deep needle insertion. May be used with anesthetics, steroids and other medicaments. Has interchangeable Extenda Tips for easy sterilization between patients. Can be used on all age groups. Consistent depth of penetration of 4-5.5mm below the epithelium and makes a wheal at the base of the injection of 5-6mm in diameter. Consistent volume of 0.1cc per injection intradermally. Permits approximately 38 injections with single loading (to 4.0cc). May be sterilized by autoclaving or your usual sterilization process - do not use dry heat. Full one year warranty.

- 501 MadaJet Instrument w/503 Head Assembly (Straight), (2) #400 Extenda Tips w/sheath 3/4", MadaCide-FD (32 oz), ME30 Wrench, 507 Holder, (3) Pyrex Fill Chambers, (3) Stylets & Custom Case
- 502 MadaJet Instrument Same as #501 w/Extra Head Assembly (Straight) w/#400 Extenda Tip on Head 3/4"
- 503 Head Assembly, Straight w/3/4" Extenda Tip
- 505-C Conversion MadaJet Podiatry of Old Style to New Style (XL) w/2 #400 Extenda Tips w/sheath 3/4", MadaCide-FD (32 oz), ME30 Wrench, (3) Stylets, 507 Holder, (3) Pyrex Fill Chambers & Custom Case (Note: must be returned to MADA for conversion)
- 505-C2 Same as 505-C w/Extra Head Assembly
- 507 Holder
- 509 Instructional Video Tape

OWEN MUMFORD AUTOJECT® 2 INJECTION AID DEVICE

Specially designed for those who have difficulty injecting themselves manually with a syringe. At the touch of a button, it automatically inserts the needle and contents of the syringe into the skin without the user having to see the needle. Available in two models, for either fixed or non-fixed needle syringes, features adjustable penetration depth, a safety lock to prevent accidental firing and a dose completion indicator. The one-handed operation makes it possible to use a wider range of injection sites than manually injecting with only a syringe.

- AJ1300 Autoject® 2, Supplied w/wallet, Depth Adjusters & Instructions, Fixed Needle
- AJ1311 Autoject® 2, Supplied w/wallet, Depth Adjusters & Instructions, Non-Fixed Needle

**Your sales representative
can assist you with
product information
and current pricing.**